

DISCHARGE MONITORING REPORT - QUALITY ASSURANCE (DMR-QA) STUDY 18

To Participating Chemistry Laboratories:

In March of 1998, the U.S. Environmental Protection Agency (USEPA) Office of Enforcement and Compliance Assurance sent announcement letters to designated National Pollutant Discharge Elimination System (NPDES) permittees, informing them that the USEPA and related state agencies were continuing the NPDES Discharge Monitoring Report - Quality Assurance (DMR-QA) studies. These announcements contained the chemistry sample ordering form that you subsequently submitted. The purpose of these studies is to evaluate the analytical and reporting ability of laboratories routinely performing chemistry and whole-effluent toxicity self-monitoring analyses required in the designated NPDES permits. For further information, please refer to the Chemistry Laboratory Requirements on page 1 inside.

Participation of NPDES permittees in this program, including proper analyses, reporting and record retention, is mandatory based on the authority of Section 308(a) of the Clean Water Act. The Agency's legal opinions, dated August 11, 1977 and January 9, 1989, reaffirm this authority.

If you have questions regarding the DMR-QA, please contact your state DMR-QA coordinator, whose name address and telephone number are listed in the attached Chemistry Laboratory Instructions. These offices will play an important role in providing any subsequent study follow-up action or guidance. Please include the complete name, address, USEPA Lab Code and telephone number of your laboratory in all correspondence.

Thank you for your cooperation in this national program to improve the quality of NPDES self-monitoring data.

Enclosures (3):

- 1) Chemistry Laboratory Instructions
- 2) NPDES Chemistry Laboratory Data Reporting Form
- 3) Chemistry Samples (per your request)

DMR-QA CHEMISTRY LABORATORY INSTRUCTIONS

TABLE OF CONTENTS

	Starts on Page:
Chemistry Laboratory Requirements.....	1
Chemistry Sample/Information Replacement Order Form.....	2
Ampul Disposal Information.....	5
Reporting Chemistry Data:	
Completing the NPDES Chemistry Laboratory Data Report Form.....	6
What to Do With the Form After Completion.....	7
Sample Preparation Instructions.....	9
Method Code Reference Guide.....	25
Material Safety Data Sheets (MSDS).....	39
List of DMR-QA State Coordinators.....	Back Page

CHEMISTRY LABORATORY REQUIREMENTS

1. Some tests you normally perform for your permittee may not be included in the DMR-QA Study. You are only required to perform analyses which have been requested by a permittee, provided they are included in the study. Participating laboratories are required to perform only one of each type of analysis, regardless of the number of permittees that request it. **Review the analytes listed below to determine which analyses must be performed to satisfy the DMR-QA reporting requirements for each NPDES Permittee you support.** To the extent feasible, analyze each DMR-QA sample the same way you would analyze a routine sample, i.e., usual analyst, usual analytical system, usual procedures, etc.
 - **Chemistry Analytes:** aluminum, arsenic, cadmium, (total) chromium, cobalt, copper, iron, lead, manganese, mercury, nickel, selenium, vanadium, zinc, pH, non-filterable residue (total suspended solids), oil and grease, ammonia nitrogen, nitrate nitrogen, orthophosphate, total Kjeldahl nitrogen, total phosphorus, 5-day biochemical oxygen demand, chemical oxygen demand, total organic carbon, 5-day carbonaceous biochemical oxygen demand, total cyanide, total phenolics (by the 4AAP method), and total residual chlorine.
2. **Ensure that you have received the appropriate study samples, reporting forms and instructions:**
 1. **Replacement samples, reporting forms or instructions** may be acquired by returning the completed Chemistry Sample / Information Replacement Form on the next page. If you have analytical problems, the available concentrate volume in each ampul is usually sufficient to prepare a second sample¹.
 2. **Unused samples** may be handled in one of three ways:
 - Properly dispose of unused samples ensuring compliance with all federal, state and local regulations governing waste management. For more information regarding waste disposal, see the Ampul Disposal Information on page 3.
 - Return unopened ampuls to the USEPA Contractor, ManTech Environmental, in accordance with 49 CFR 173.4, and using the return address on the sample shipment (4907 S. Alston Ave., Durham, NC 27713).
 - Retain samples for internal quality control. If stored at room temperature and out of direct light, unopened samples are stable for at least one year. A complete list of "true"/calculated values may be obtained after the results of the study have been distributed to permittees, which is expected to occur in January, 1999. Your permittees may request a copy by contacting their Regional/State Coordinator. Users of a Personal Computer w/modem may use the USEPA Office of Research and Development Electronic Bulletin Board, modem number (513) 569-7700 or (513) 569-7610. If you need assistance with use of the Bulletin Board, call (513) 569-7272 (do not call this number for general information regarding PE Studies).
3. **Perform the analyses requested by your permittees;** refer to the enclosed NPDES Chemistry Laboratory Data Report Form for a convenient listing of the analyses available in the study.
4. **Report the results of your analyses and additional required information** as specified on page 6, Completing the NPDES Chemistry Laboratory Data Report Form.
5. **Submit your reporting form** to the requesting permittees in accordance with the instructions on page 7, Chemistry Laboratory Reporting Procedure.
6. **After official evaluation of the study data by the USEPA, permittees will receive Performance Evaluation Reports and laboratories will be asked by the supported NPDES permittee(s) to explain reported values that were not within the study acceptance limits.** Permittees will then provide this explanation to the appropriate USEPA or State Agency.

¹Not applicable for the pH, demand and residue samples.

CHEMISTRY SAMPLE / INFORMATION REPLACEMENT ORDER FORM

If replacement samples, reporting form or instructions are needed, copy this form, complete the information and return it as indicated below.

1. **USEPA Chemistry Laboratory Code:** If you are unsure of your USEPA Lab Code, it can be found on the label used to deliver your original shipment or the label on page 1 of your NPDES Chemistry Laboratory Data Report Form.

--	--	--	--	--	--	--

(State Code + 5 digits, for example NC12345)

2. **Shipping Address:**

Laboratory Name _____ Phone (____) _____

Contact Name _____ Title _____

Laboratory Shipping Address (No PO Boxes) _____

City _____ State _____ Zip Code _____

3. **Samples: Check the appropriate block(s).** You may order only one of each ampul (*The same analytical result for each analyte may be reported to all your permittees to satisfy their reporting requirements*).

☐ Demand

☐ Total Cyanide

☐ Residue

☐ Total Phenolics

☐ Nutrients #1

☐ pH

☐ Nutrients #2

☐ Total Residual Chlorine

☐ Oil & Grease

☐ NPDES Chemistry Laboratory Data Report Form

☐ Trace Metals #1

☐ DMR-QA Chemistry Laboratory Instructions

4. **Reason for Request** [check the appropriate block(s)]:

☐ Broken in shipment

☐ Missing in shipment

☐ Lab accident

☐ Other _____

5. **Fax completed form to the USEPA Contractor: ManTech Environmental (919) 406-2246 attn: Terry Bundy.**

Reminder: Maintain a copy of completed order form for your records.

AMPUL DISPOSAL INFORMATION

U.S. EPA WATER POLLUTION PERFORMANCE EVALUATION SAMPLES


WP Inorganic		Other Contents	Chemical Formula	Max. Conc. mg/L Ampul
Ampul	Solvent			
Trace Metals conc-1	5%HNO ₃ / H ₂ O	aluminum	Al	410
		arsenic trioxide	As ₂ O ₃	150
		cadmium	Cd	110
		chromium trioxide	CrO ₃	210
		cobalt	Co	110
		copper	Cu	110
		iron	Fe	410
		lead nitrate	Pb(NO ₃) ₂	500
		manganese	Mn	410
		mercuric oxide	HgO	3.2
		methyl mercury chloride	CH ₃ HgCl	3.8
		nickel	Ni	310
		selenium dioxide	SeO ₂	300
		vanadium pentoxide	V ₂ O ₅	2000
		zinc oxide	ZnO	260
Trace Metals conc-2	5%HNO ₃ / H ₂ O	beryllium	Be	110
		molybdenum trioxide	MoO ₃	150
		potassium antimony tartrate	KSbC ₄ H ₄ O ₇	270
		silver nitrate	AgNO ₃	160
		strontium nitrate	Sr(NO ₃) ₂	97
		thallium nitrate	TlNO ₃	130
		titanium	Ti	40
Minerals conc-1A	H ₂ O	calcium chloride	CaCl ₂ +2H ₂ O	41000
		hydrochloric acid	HCl	150
		magnesium chloride	MgCl ₂ +6H ₂ O	34000
Minerals conc-1B	H ₂ O	potassium chloride	KCl	7800
		potassium sulfate	K ₂ SO ₄	9100
		sodium bicarbonate	NaHCO ₃	20000
		sodium chloride	NaCl	11000
		sodium fluoride	NaF	910
		sodium sulfate	Na ₂ SO ₄	19000
pH-Minerals conc-2	H ₂ O	hydrochloric acid	HCl	1700
		potassium dihydrogen phosphate	KH ₂ PO ₄	6800
		potassium hydrogen phthalate	KHC ₈ H ₄ O ₄	10000
		sodium hydroxide	NaOH	1900
		sodium borate....."borax"	B ₄ Na ₂ O ₇ +10H ₂ O	4600
		tris(hydroxymethyl) aminomethane....."THAM"	H ₂ NC(CH ₂ OH) ₃	6100
Nutrient-1	H ₂ O	ammonium chloride	NH ₄ Cl	7600
		potassium dihydrogen phosphate	KH ₂ PO ₄	2400
		potassium nitrate	KNO ₃	30000
Nutrient-2	H ₂ O	glycine	NH ₂ CH ₂ COOH	19000
		sodium B-glycerolphosphate	(HOCH ₂) ₂ CHOPO(ONa) ₂ +5H ₂ O	11000
Phenolics	H ₂ O	2-chlorophenol	C ₆ H ₄ OH	210
		2,4-dichlorophenol	Cl ₂ H ₃ OH	270
		2,4-dinitrophenol	(NO ₂) ₂ C ₆ H ₃ OH	300
		phenol	C ₆ H ₅ OH	310
Demands	H ₂ O	glucose	C ₆ H ₁₁ O ₅ CHO	6800
		glutamic acid	HO ₂ C ₂ H ₄ NH ₂ (CH ₂) ₂ CO ₂ H	6800
Cyanide	H ₂ O	potassium ferricyanide	K ₃ Fe(CN) ₆	230
		potassium hydroxide	KOH	29
		sodium hydroxide	NaOH	29
Residual Chlorine	H ₂ O	sodium hypochlorite....."Clorox"	NaClO	440
Residue (TSS-nonfilt)	earth+salt solid mix solid mix	infusorial earth (calcined powder)	(a natural clay)	Maximum percent of mix 50
		sodium chloride	NaCl	99

COMPLETING THE NPDES CHEMISTRY LABORATORY DATA REPORT FORM

1. **USEPA Lab Code:** Enter your USEPA Lab Code in the spaces provided at the top of page 1. The USEPA Laboratory Code begins with your state abbreviation and is followed by five numbers. Do not submit a permit number, state laboratory code or certification ID in lieu of the USEPA Laboratory Code. If you are unsure of your USEPA Lab Code, it can be found on the label in the middle of page 1.
2. **Enter the appropriate data into each field.** The headings listed below relate to information on pages 3 & 4 of the Report Form.
 - a. **Analyte Name and No.:** Identifies each analyte available and their corresponding analyte numbers (do not alter).
 - b. **Sample No.:** Identifies the sample number which was analyzed to produce each result (do not alter).
 - c. **USEPA Laboratory Code:** Enter your USEPA Code next to each analytical result produced by your laboratory and on any additional documentation that you may submit.
 - d. **Voluntary Analyte:** Chemistry laboratories should not enter anything in this column.
 - e. **Method Codes:** Using the enclosed list of USEPA approved methods, identify which method was used to produce each result being reported. A method code of "99" will require additional documentation describing the method used.
 - f. **Usage of "<" or ">":** All analytes are present at quantifiable concentrations, so use of this field should not be necessary. However, if you are not capable of quantifying the sample concentration, you must enter the "<" or ">" portion of your analytical result in this space.
 - g. **Quantity:**
 - Report only one value for each analyte, the result of a single analysis; do not report the mean of multiple analytical results.
 - Each box utilized must contain only a number or a decimal (no dashes, letters or symbols for any reason). **A decimal point requires its own box.**
 - Report your data to three significant digits, i.e., 0.XXX, X.XX, XX.X or XXX, if method allows.
 - "Right justify" each result in the spaces provided.
 - Report all data in the units specified on the data report form.
3. **Complete the Checklist and Certification Statement** on page 1 of the Chemistry Laboratory Data Report Form. The Certification Statement must be signed by the Director or an authorized representative of the laboratory. Most importantly, the laboratory is responsible for certifying the truth, accuracy, and completeness of the data and that **all work has been performed by the personnel that routinely conduct the self-monitoring required for each permittee they report to.** Reference to this section may be found in 40 CFR Part 122.22.

CHEMISTRY LABORATORY REPORTING PROCEDURE

- 1) **Make photo copies** of the completed NPDES Chemistry Laboratory Data Report Form (pink colored form).
- 2) **Retain the original form for your files.**
- 3) **Contract labs should send 1 copy to each NPDES Permittee (your clients) who requested your participation, by September 4, 1998. *This date should allow each permittee enough time to meet their reporting deadline of September 30, 1998.***
- 4) **In-house labs may disregard the NPDES Chemistry Laboratory Data Report Form and report directly on the NPDES Permittee Data Report Form (orange colored form). This form will be sent in a separate mailing to each NPDES Permittee during the summer of 1998).**
- 5) **Do not submit NPDES Chemistry Laboratory Report Forms (pink colored forms) to the USEPA Contractor, USEPA or to a State Agency; NPDES Permittees must report using NPDES Permittee Data Report Forms (orange colored forms).**



DMR-QA STUDY 18

GENERAL INSTRUCTIONS

- Immediately upon receipt of samples, check for breakage and correctness of sets shipped; see list of requested sets on the bottom of the report form label. To request replacements for broken and/or missing ampuls, FAX your request to **Terry Bundy at ManTech (919) 406-2246 (your EPA assigned lab code must be included on all sample requests)** using the form attached to the cover letter for this mailing; to obtain answers to questions, call your State or regional coordinator (see list included in these instructions). Sample replacement takes at least 2 to 3 weeks.
- Read instructions carefully before opening ampuls.
- Treat each PE study sample the same way you would treat a routine sample, i.e., usual analyst, usual analytical system, usual procedures, etc.
- For each analyte you report, provide results of a single analysis, as would be done with a routine sample.
- For best results, ampul solutions and reagent water should be at room temperature (near 20° Centigrade) when used.
- Amber-colored total residual chlorine ampuls as well as the clear cyanide ampuls should be stored away from light until the analyst is ready to use them.
- Open ampuls cautiously to avoid cuts or flying glass. A suggested procedure is to wrap the ampul in cloth or a thick paper towel and snap off the ampul top at the narrow part of the neck.
- Unless otherwise stated, samples need not be filtered.
- Use only Class A volumetrics.
- "Reagent water" is equivalent to Type II Reagent Water as specified in ANSI/ASTM Standard D 1193-91.
- Please Note: There is only one concentration per analyte.

TRACE METAL SAMPLE

An ampul containing 5% nitric acid (v/v) is provided for preparation of the study sample.

SAMPLE PREPARATION

To prepare sample 1, fill a 1-L volumetric flask approximately 3/4 full with reagent water, pipette 10.0 mL (which contains 0.5 mL of concentrated nitric acid) of ampul 1 and 1.0 mL of trace-metal-grade concentrated nitric acid into the flask, dilute to volume with reagent water, and mix well. The resulting samples should be analyzed for the following analytes at levels within the stated concentration ranges.

Analyte	Concentration, micrograms per liter (µg/L) after dilution
Aluminum	<5000
Arsenic	<2000
Cadmium	<2000
Chromium	<2000
Cobalt	<2000
Copper	<2000
Iron	<4000
Lead	<4000
Manganese	<4000
Mercury	<60
Nickel	<4000
Selenium	<4000
Vanadium	<11000
Zinc	<4000

Notes:

- Samples, standards and reagent water blanks should all have the same acid matrix. The solution resulting from the above dilution will contain 0.15% v/v nitric acid; additional or different acids may be used as appropriate for analytical procedures or to match standards.
- **Mercury is present as a mixture of organic and inorganic forms. If a cold vapor method is used, an aliquot of the sample must be digested before analysis for mercury. However, digestion of this sample for any other metal analysis is unnecessary.**
- Chromium should be determined only as total chromium. When chromium is analyzed by direct air-acetylene AA, use of matrix modifiers or the less sensitive oxidizing flame may be necessary to control interfering ions. Preferably, use a fuel-rich nitrous oxide-acetylene flame for best results.
- Some analytes may occur at concentrations that are too low for determination with simple direct AA and must be determined with ICP, furnace AA, or AA after chelation extraction (where approved) to produce reportable results.

DMR-QA STUDY 18

- With some methods (such as colorimetric chromium method 3500D from Reference 74), preparatory extraction may be required to remove relatively high concentrations of interferences before analysis.
- Because many analysts will use atomic absorption or inductively coupled plasma spectrometry, multiple aspirations or measurements may be obtained to determine an analyte concentration, as with routine samples.
- Report results as micrograms per liter with three figures.
- If AA is used, specify which method (i.e., air-acetylene, nitrous oxide-acetylene, platform or non-platform) was used.

**Ampuls contain 5% nitric acid.
CORROSIVE/OXIDIZER - TOXIC - SEVERE IRRITANT
Material Safety Data Sheets are attached
immediately following the sample preparation instructions.**

END OF INSTRUCTIONS FOR TRACE METALS

DMR-QA STUDY 18

pH SAMPLE

SAMPLE PREPARATION

To prepare the pH sample, pour the contents of the pH ampul directly into a ≤ 50 -mL beaker; do not dilute before determining pH.

NOTES:

- We acknowledge that, ideally: 1) the strong memory effects of the calibrating buffers are removed by allowing the electrode(s) to soak in some of the sample prior to analysis, and 2) laboratories should be able to use their entire routine analytical procedures in performance evaluation studies. However, there are practical limits to the amount of pH sample we can provide to you in our studies.

Therefore, we propose that you soak your electrode(s) in distilled water to remove the memory effect of the calibration buffers before you analyze our pH sample as directed in the study instructions. We believe this adjustment to your normal procedure is reasonable, effective and should not adversely affect your chances of successful performance.

- If you cannot produce a valid pH measurement from your analytical system using a sample of 20 mL (or 40 mL using the contents of an additional ampul which is available upon request), your analytical system is beyond the scope of this study. For example, systems for continuous pH measurement may require a minimum sample size which is incompatible with the current study design. If your pH system is incompatible with the study design, you should discuss this with your State DMR-QA Coordinator (see list at back of this set of instructions).

END OF INSTRUCTIONS FOR pH

NUTRIENT SAMPLES

Two ampuls containing reagent-water solutions are provided to prepare two separate samples. Sample 1 will contain inorganic nitrogen and phosphorus and may be analyzed for ammonia nitrogen, nitrate nitrogen, and orthophosphate phosphorus. Sample 2 will contain organically-bound forms and may be analyzed for total Kjeldahl nitrogen and total phosphorus. When prepared according to instructions, samples will contain <50 mg/L of each nutrient.

SAMPLE PREPARATION

To prepare sample 1, fill a 1-L volumetric flask approximately 3/4 full with ammonia-free, reagent water, pipette 10.0 mL of ampul 1 solution into the flask, dilute to volume with reagent water, and mix well. To prepare sample 2, repeat the procedure with ampul 2.

NOTES:

- Solutions in ampuls were preserved by autoclaving the sealed ampuls. Because this preservative treatment is not effective after ampuls are opened, open ampuls just before analysis.
- If a solution is not to be analyzed within two days after preparation, preserve a sample aliquot with 0.2 mL of concentrated sulfuric acid per 100 mL of sample (pH <2) and store at 4°C until just before analysis. Analyze the preserved sample as soon as possible within the maximum 28-day holding time.
- Report results as milligrams per liter with three figures.

END OF INSTRUCTIONS FOR NUTRIENTS

DEMAND SAMPLE

One ampul containing reagent-water solution is provided to prepare one sample. The sample may be analyzed to determine 5-day biochemical oxygen demand (BOD), 5-day carbonaceous BOD (CBOD), chemical oxygen demand (COD), and total organic carbon (TOC).

SAMPLE PREPARATION

To prepare the sample, fill a 1-L volumetric flask approximately 3/4 full with reagent water, pipette 20.0 mL from the ampul into the flask, dilute to volume with reagent water, and mix well. The sample is ready for analysis to determine 5-day BOD, 5-day CBOD, COD and TOC. Each measured value should be <300 mg/L.

NOTES:

- The ampul contains about 21 mL of solution; if clean dry pipettes are used, this is sufficient to prepare the sample. More ampuls may be requested if both 5-day BOD and 5-day CBOD are to be determined. If an additional sample is needed, complete the "Request for Sample Replacement" sheet attached to the cover memorandum and submit as directed.
- The solution in this ampul was preserved by autoclaving the sealed ampul. Because this preservative treatment is not effective after the ampul is opened, open the ampul just before analysis.
- Begin the BOD and/or CBOD as soon as possible after sample preparation. If cooled to 4°C, the maximum holding time for these analyses is 48 hours. If COD and/or TOC cannot be analyzed within your usual holding time after sample preparation, an aliquot should be preserved (pH <2) using 0.2 mL of concentrated sulfuric acid per 100 mL (or an equivalent amount of another acid which is compatible with your chosen methods). This acid preserved portion should be stored at 4°C and analyzed as soon as possible within the 28-day maximum holding time.
- 5-day BOD and 5-day CBOD: For proper oxygen depletion, the sample should be tested in a dilution series such as 4%, 8%, and 16%. Use a natural surface water, diluted domestic sewage, or commercially available seeds for seeded dilution water and avoid entrainment of air. For a 5-day CBOD, the use of a nitrifying inhibitor is recommended.
- Be sure to determine the "blank" BOD of your seeded dilution water so that the proper seed correction can be made (See STANDARD METHODS, 17th edition, Page 5-6).
- Report results as milligrams per liter with three figures.

END OF INSTRUCTIONS FOR DEMAND

CYANIDE SAMPLE

One ampul containing a reagent-water solution is provided to prepare one sample.

SAMPLE PREPARATION

To prepare the sample, fill a 1-L volumetric flask approximately 3/4 full with reagent water and add 1 mL of 1N NaOH. Pipette 10.0 mL of the solution from the ampul into the flask, dilute to volume with reagent water, and mix well.

NOTES:

- **CAUTION:** The sample should be analyzed as soon as possible after opening and diluting.
- The cyanide concentration will be <2 mg/L.
- To prevent a concentration change, the ampul and, even more importantly, the prepared sample should not be exposed to sunlight or other ultraviolet light any longer than necessary. Therefore, perform the analysis immediately after sample preparation.
- The specified addition of 1 mL of 1 normal sodium hydroxide is only a recommended minimum to prevent the loss of cyanide. If another amount, or even another base (such as KOH) is part of the analyst's "normal routine", it is his/her option to change the recommendation.
- Since all the cyanide is present in soluble forms, filtration of the sample is **not** necessary, however, the sample will require distillation.
- The presence of a precipitate in the ampul is **not** a problem, unless it is blue in color.
- Report results as milligrams per liter with three figures.

END OF INSTRUCTIONS FOR CYANIDE

NON-FILTERABLE RESIDUE (TOTAL SUSPENDED SOLIDS) SAMPLE

One vial containing homogeneously mixed solids is provided to prepare one sample.

SAMPLE PREPARATION

To prepare the sample, fill a 1-L volumetric flask about 3/4 full with reagent water. Soon after opening and without drying, weigh 500 mg of solids from the vial and quantitatively transfer to the volumetric flask using a glass funnel and wash bottle filled with reagent water. After rinsing the glass funnel thoroughly and removing it, dilute flask contents to volume with additional reagent water (see notes) and mix by inverting the flask 10 times. Let the sample flask stand for 10 min and invert again 10 times. Immediately after the second set of inversions, rapidly measure a 500-mL aliquot into a graduated cylinder by repetitively mixing and pouring so as not to allow solids to settle in the flask before the transfer is complete. Filter this 500-mL aliquot and determine the non-filterable residue (total suspended solids) according to your chosen routine method.

NOTES:

- Do not dry the vial contents prior to weighing the 500 mg of mixed solids. Weigh the solids as received without any preparatory action.
- If a significant portion of the weighed (500 mg) solids fails to disperse evenly and stubbornly floats on the surface, add a few drops of diluted detergent to the partially filled sample preparation flask. For example, dissolve approximately 10 mL of detergent in 100 mL of reagent water and mix. Then add 2 drops of this diluted detergent to the sample preparation flask, mix, observe, and repeat (if necessary) the 2 drop additions until the suspended sample appears to be dispersed evenly. Choose a detergent which will rinse clean from the Gooch crucible and filter pad without leaving a measurable residue requiring a blank subtraction. Possible conveniently accessible choices might be Triton-X-100, or any house-hold liquid dish-washing detergent void of insoluble ingredients. If your detergent fails to disperse the solids evenly and it was not Triton-X-100, try Triton-X-100. If all else fails, request a new ampul.
- CAUTION: The analysis should be performed soon after the 500 mL aliquot is obtained.
- Report results as milligrams per liter with three figures.

**Ampuls contain NaCl.
IRRITANT**

**Material Safety Data Sheets are attached
immediately following the sample preparation instructions.**

END OF INSTRUCTIONS FOR NON-FILTERABLE RESIDUE

OIL AND GREASE SAMPLE

One ampul containing a solution of oil and grease in n-propanol and glycerol is provided to prepare one sample.

SAMPLE PREPARATION

To prepare the sample, add 994 mL of reagent water and 5 mL of 1:1 hydrochloric acid into a 2-L separatory funnel and mix. Pipette 1.0 mL of the ampul solution into the separatory funnel and mix well. Analyze the sample using whatever analytical method you most often use on routine samples. Note that data produced using hexane for extraction must be reported in a different place on the report form than data produced using Freon for extraction.

NOTES:

- To avoid a concentration change, the sample must be analyzed soon after the ampul has been opened.
- The sample should be analyzed with an approved gravimetric method, because it is not suitable for an infrared method.
- For calculations, assume that initial sample volume was 1000 mL.
- Report results as milligrams per liter with three figures.

**Ampuls contain n-propanol (flammable liquid)
and glycerol (irritant)
Material Safety Data Sheets are attached
immediately following the sample preparation instructions.**

END OF INSTRUCTIONS FOR OIL AND GREASE

TOTAL PHENOLICS SAMPLE

One ampul containing a reagent-water solution is provided to prepare one sample.

SAMPLE PREPARATION

To prepare the sample, fill a 1-L volumetric flask approximately 3/4 full with reagent water. Pipette 10.0 mL from the ampul into the flask, dilute to volume with reagent water, and mix well. Analyze the sample using a 4AAP method ONLY.

Notes:

- In each sample, the resulting total phenolic concentration will be <6 mg/L.
- The solution in the ampul was preserved by autoclaving the sealed ampul. Because this preservative treatment is not effective after the ampul is opened, open the ampul just before analysis.
- Do not report phenolic results determined with GC/MS, because GC/MS data are not equivalent to data obtained with the 4AAP methods.
- Report results as milligrams per liter with three figures.

END OF INSTRUCTIONS FOR TOTAL PHENOLICS

DMR-QA STUDY 18**TOTAL RESIDUAL CHLORINE SAMPLE**

One ampul containing a reagent-water solution is provided to prepare one sample.

SAMPLE PREPARATION

To prepare the sample, fill a 1-L volumetric flask approximately 3/4 full with chlorine-free, reagent water. Pipette 10.0 mL from the ampul into the volumetric flask, dilute to volume with reagent water, and mix well.

NOTES:

- In this sample, the resulting concentration of total residual chlorine will be <6 mg/L.
- To prevent a concentration change, the ampul and prepared sample should not be exposed to sunlight or other strong light or agitation any longer than necessary. Therefore, perform the analysis immediately after sample preparation.
- Report results as milligrams per liter with three figures.

END OF INSTRUCTIONS FOR TOTAL RESIDUAL CHLORINE

ABD00060322
DMR-QA STUDY 18

Method Code Reference Guide

Analyte	Method Code	Ref.-Method	Description
Trace Metals:			
Aluminum (ICP)	13	R19 - 200.7	Inductively Coupled Plasma
	14	R24 - 202.1	Atomic Absorption (AA): Direct aspiration unless otherwise noted
	15	R24 - 202.2	AA; Furnace
	26	R89 - 3111D	AA
	27	R89 - 3113B	AA; Furnace
	28	R89 - 3500D	Eriochrome Cyanine R
	29	R89 - 3120B	ICP
	31	R25 - D4190-82(88)	Direct Current Plasma
	43	R79 - I-3051-85	AA
	56	R80 - AES0029	DCP
	99	Other	
Antimony	13	R24 - 204.1	AA
	14	R24 - 204.2	AA; Furnace
	15	R19 - 200.7	ICP
	23	R89 - 3111B	AA
	24	R89 - 3113B	AA; Furnace
	25	R89 - 3120B	ICP
	99	Other	
Arsenic	14	R19 - 200.7	ICP
	15	R24 - 206.2	AA; Furnace
	16	R24 - 206.3	AA; Gaseous hydride
	17	R24 - 206.4	Silver Diethyl- dithiocarbamate
	20	R89 - 3114B4d	AA; Gaseous hydride
	21	R89 - 3113B	AA; Furnace
	22	R89 - 3120B	ICP
	29	R89 - 3500C	Silver Diethyl- dithiocarbamate
	31	R25 - D2972-88A	Silver Diethyl- dithiocarbamate
	32	R25 - D2972-88B	AA; Gaseous hydride
	33	R25 - D2972-88C	AA; Furnace
	45	R79 - I-3062-85	AA; Gaseous hydride
	46	R79 - I-3060-85	Silver Diethyl- dithiocarbamate
	99	Other	

ABD00060323
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Cadmium	13	R19 - 200.7	ICP
	14	R24 - 213.1	AA
	15	R24 - 213.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	23	R89 - 3120B	ICP
	24	R89 - 3500D	Dithizone
	35	R25 - D3557-90A	AA
	36	R25 - D3557-90B	AA; Chelation-extraction
	37	R25 - D4190-82(88)	DCP
	38	R25 - D3557-90C	Voltametry
	39	R25 - D3557-90D	AA; Furnace
	44	R79 - I-3135-85	AA
	45	R79 - I-3136-85	AA
	46	R79 - I-1472-85	ICP
	52	R8 - p. 37	AA
	55	R81 - 974.27	AA
	56	R80 - AES0029	DCP
	99	Other	
Chromium	13	R19 - 200.7	ICP
	14	R24 - 218.1	AA
	15	R24 - 218.2	AA; Furnace
	16	R24 - 218.3	AA; Chelation-extraction
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	23	R89 - 3120B	ICP
	24	R89 - 3500D	Diphenylcarbazide
	35	R25 - D1687-92B	AA
	36	R25 - D4190-82(88)	DCP
	38	R25 - D1687-92C	AA; Furnace
	44	R79 - I-3236-85	AA
	55	R81 - 974.27	AA
	56	R80 - AES0029	DCP
	99	Other	

ABD00060324
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Cobalt	13	R19 - 200.7	ICP
	14	R24 - 219.1	AA
	15	R24 - 219.2	AA; Furnace
	25	R89 - 3111B	AA
	26	R89 - 3111C	AA; Chelation-extraction
	27	R89 - 3113B	AA; Furnace
	28	R89 - 3120B	ICP
	34	R25 - D3558-90A	AA
	35	R25 - D3558-90B	AA; Chelation-extraction
	36	R25 - D4190-82(88)	DCP
	37	R25 - D3558-90C	AA; Furnace
	44	R79 - I-3239-85	AA
	51	R8 - p. 37	AA
	56	R80 - AES0029	DCP
	99	Other	
Copper	13	R19 - 200.7	ICP
	14	R24 - 220.1	AA
	15	R24 - 220.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	27	R89 - 3120B	ICP
	28	R89 - 3500D	Neocuproine
	29	R89 - 3500E	Bicinchoninate
	30	R25 - D1688-90C	AA; Furnace
	35	R25 - D1688-84A(88)	Neocuproine
	36	R25 - D1688-90A	AA
	37	R25 - D1688-90B	AA; Chelation-extraction
	39	R25 - D4190-82(88)	DCP
	44	R79 - I-3271-85	AA
	45	R79 - I-3270-85	AA
	52	R8 - p. 37	AA
	54	R41 - 8506	Bicinchoninate
	56	R81 - 974.27	AA
	57	R80 - AES0029	DCP
	99	Other	

ABD00060325
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Iron	13	R19 - 200.7	ICP
	14	R24 - 236.1	AA
	15	R24 - 236.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	27	R89 - 3120B	ICP
	28	R89 - 3500D	Phenanthroline
	30	R25 - D1068-90C	AA; Furnace
	36	R25 - D1068-90A	AA
	37	R25 - D1068-90B	AA; Chelation-extraction
	38	R25 - D1068-90D	Phenanthroline
	39	R25 - D4190-82(88)	DCP
	43	R79 - I-3381-85	AA
	53	R44 - 8008	Phenanthroline
	56	R81 - 974.27	AA
	57	R80 - AES0029	DCP
	99	Other	
Lead	13	R19 - 200.7	ICP
	14	R24 - 239.1	AA
	15	R24 - 239.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	23	R89 - 3120B	ICP
	24	R89 - 3500D	Dithizone
	35	R25 - D3559-90A	AA
	36	R25 - D3559-90B	AA; Chelation-extraction
	37	R25 - D4190-82(88)	DCP
	38	R25 - D3559-90C	Voltametry
	39	R25 - D3559-90D	AA; Furnace
	43	R79 - I-3399-85	AA
	56	R81 - 974.27	AA
	57	R80 - AES0029	DCP
	99	Other	

ABD00060326
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Manganese	13	R19 - 200.7	ICP
	14	R24 - 243.1	AA
	15	R24 - 243.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	28	R89 - 3120B	ICP
	29	R89 - 3500D	Persulfate
	35	R25 - D858-90A	AA
	36	R25 - D858-90B	AA; Chelation-extraction
	37	R25 - D4190-82(88)	DCP
	38	R25 - D858-90C	AA; Furnace
	43	R79 - I-3454-85	AA
	54	R41 - 8034	Periodate (pp. 2-113 and 2-117)
	56	R81 - 974.27	AA
	57	R81 - 920.203	Persulfate
	58	R80 - AES0029	DCP
	99	Other	
Mercury	13	R24 - 245.1	Manual Cold Vapor
	14	R24 - 245.2	Automated Cold Vapor
	23	R89 - 3112B	Manual Cold Vapor
	33	R25 - D3223-91	Manual Cold Vapor
	43	R79 - I-3462-85	Manual Cold Vapor
	51	R81 - 977.22	Manual Cold Vapor
	99	Other	
Molybdenum	13	R19 - 200.7	ICP
	14	R24 - 246.1	AA
	15	R24 - 246.2	AA; Furnace
	23	R89 - 3111D	AA
	24	R89 - 3113B	AA; Furnace
	25	R89 - 3120B	ICP
	42	R79 - I-3490-85	AA
	52	R80 - AES0029	DCP
	99	Other	
Nickel	13	R19 - 200.7	ICP
	14	R24 - 249.1	AA
	15	R24 - 249.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	27	R89 - 3120B	ICP
	28	R89 - 3500D	Heptoxime
	34	R25 - D1886-90A	AA
	35	R25 - D1886-90B	AA; Chelation-extraction
	36	R25 - D4190-82(88)	DCP
	37	R25 - D1886-90C	AA; Furnace
	43	R79 - I-3499-85	AA
	52	R80 - AES0029	DCP
	99	Other	

ABD00060327
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Selenium	12	R19 - 200.7	ICP
	13	R24 - 270.2	AA; Furnace
	25	R89 - 3113B	AA; Furnace
	26	R89 - 3114B	AA; Gaseous hydride
	27	R89 - 3120B	ICP
	32	R25 - D3859-88A	AA; Gaseous hydride
	43	R79 - I-3667-85	AA; Gaseous hydride
	99	Other	
Silver	13	R19 - 200.7	ICP
	14	R24 - 272.1	AA
	15	R24 - 272.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3113B	AA; Furnace
	23	R89 - 3120B	ICP
	43	R79 - I-3720-85	AA
	52	R8 - p. 37	AA
	56	R81 - 974.27	AA
	57	R80 - AES0029	DCP
	99	Other	
Strontium	99	Other	
Thallium	13	R19 - 200.7	ICP
	14	R24 - 279.1	AA
	15	R24 - 279.2	AA; Furnace
	23	R89 - 3111B	AA
	24	R89 - 3120B	ICP
	99	Other	
Titanium	13	R24 - 283.1	AA
	14	R24 - 283.2	AA; Furnace
	23	R89 - 3111D	AA
	57	R80 - AES0029	DCP
	99	Other	
Vanadium	13	R19 - 200.7	ICP
	14	R24 - 286.1	AA
	15	R24 - 286.2	AA; Furnace
	26	R89 - 3111D	AA
	27	R89 - 3120B	ICP
	28	R89 - 3500D	Gallic Acid
	33	R25 - D4190-82(88)	DCP
	57	R80 - AES0029	DCP
	99	Other	

ABD00060328
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Zinc	13	R19 - 200.7	ICP
	14	R24 - 289.1	AA
	15	R24 - 289.2	AA; Furnace
	20	R89 - 3111B	AA
	21	R89 - 3111C	AA; Chelation-extraction
	22	R89 - 3500F	Zincon
	27	R89 - 3120B	ICP
	28	R89 - 3500E	Dithizone
	34	R25 - D1691-90B	AA; Chelation-extraction
	35	R25 - D1691-90A	AA
	36	R25 - D4190-82(88)	DCP
	43	R79 - I-3900-85	AA
	52	R8 - p. 37	AA
	54	R41 - 8009	Zincon (pp. 2-231 and 2-333)
	55	R81 - 974.27	AA
	57	R80 - AES0029	DCP
	99	Other	

Nutrients:

Ammonia, as Nitrogen	14	R24 - 350.2	Titration
	15	R24 - 350.3	Electrode
	16	R24 - 350.1	Automated Phenate
	17	R24 - 350.2	Nesslerization
	20	R89 - 4500C	Nesslerization
	21	R89 - 4500E	Titration
	22	R89 - 4500H	Automated Phenate
	23	R89 - 4500G	Electrode
	24	R89 - 4500F	Electrode
	32	R25 - D1426-89A	Nesslerization
	33	R25 - D1426-89B	Electrode
	44	R79 - I-3520-85	Nesslerization
	45	R79 - I-4523-85	Automated Phenate
	53	R40 - 379-75WE	Automated Electrode
	54	R81 - 973.49	Nesslerization
	99	Other	

ABD00060329
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Nitrate, as Nitrogen	15	R24 - 352.1	Brucine Sulfate
	16	R24 - 353.1*	Automated Hydrazine Reduction
	17	R24 - 353.2*	Automated Cadmium Reduction
	18	R24 - 353.3*	Cadmium Reduction
	20	R89 - 4500E*	Cadmium Reduction
	21	R89 - 4500F*	Automated Cadmium Reduction
	22	R89 - 4500H*	Automated Hydrazine Reduction
	27	R3 - 419D	Brucine Sulfate
	35	R25 - D3867-90A*	Automated Cadmium Reduction
	36	R25 - D3867-90B*	Cadmium Reduction
	43	R79 - I-4545-85*	Automated Camium Reduction
	52	R8 - p. 28	Brucine Sulfate
	54	R81 - 973.50	Brucine Sulfate
	99	Other	
Orthophosphate, as P	13	R24 - 365.1	Automated Ascorbic Acid Reduction
	14	R24 - 365.2	Manual Ascorbic Acid Reduction
	15	R24 - 365.3	Manual Two Reagent
	25	R89 - 4500F	Automated Ascorbic Acid Reduction
	26	R89 - 4500E	Manual Ascorbic Acid Reduction
	33	R25 - D515-88A	Manual Ascorbic Acid Reduction
	43	R79 - I-4601-85	Automated Ascorbic Acid Reduction
	54	R81 - 973.56	Automated Ascorbic Acid Reduction
	55	R81 - 973.55	Manual Ascorbic Acid Reduction
	99	Other	

* Normally, use of these methods for nitrate would require a prior analysis of nitrite and subtraction of the nitrite results to get an estimate for nitrate, however, WP study samples do not contain any nitrite.

ABD00060330
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Kjeldahl Nitrogen	14	R24 - 351.3	Nesslerization
	15	R24 - 351.4	Potentiometric
	16	R24 - 351.2	Semi-automated Block Digestion
	17	R24 - 351.3	Titration
	18	R24 - 351.3	Electrode
	19	R24 - 351.1	Automated Phenate
	34	R25 - D3590-89A	Titration
	35	R25 - D3590-89A	Nesslerization
	36	R25 - D3590-89A	Potentiometric
	37	R25 - D3590-89B	Semi-automated Block Digestion
	42	R26 - I-4551-78	Automated Phenate
	53	R81 - 973.48	Titration
	70	R89 - 4500H	Automated Phenate
	71	R89 - 4500B+E	Macro Digestion + Titration
	72	R89 - 4500B+C	Macro Digestion + Nesslerization
	73	R89 - 4500B+F	Macro Digestion + Electrode
	74	R89 - 4500B+G	Macro Digestion + Electrode
	75	R89 - 4500C+E	Semi-micro Digestion + Titration
	76	R89 - 4500C	Semi-micro Digestion + Nesslerization
	77	R89 - 4500C+F	Semi-micro Digestion + Electrode
	78	R89 - 4500C+G	Semi-micro Digestion + Electrode
	99	Other	
Total Phosphorus, as P	13	R24 - 365.2	Manual Ascorbic Acid Reduction
	14	R24 - 365.3	Man. Ascorbic Acid Reduction (2 reagents)
	15	R24 - 365.1	Automated Ascorbic Acid Reduction
	16	R24 - 365.4	Semi-automated Block Digestion
	26	R89 - 4500E	Manual Ascorbic Acid Reduction
	27	R89 - 4500F	Automated Ascorbic Acid Reduction
	33	R25 - D515-88A	Manual Ascorbic Acid Reduction
	34	R25 - D515-88B	Semi-automated Block Digestion
	43	R79 - I-4600-85	Automated Ascorbic Acid Reduction
	54	R81 - 973.56	Automated Ascorbic Acid Reduction
	55	R81 - 973.55	Manual Ascorbic Acid Reduction
	99	Other	

ABD00060331
DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Demands:			
Chemical	12	R24 - 410.1	Dichromate Reflux
Oxygen	13	R24 - 410.2	Dichromate Reflux
Demand (COD)	14	R24 - 410.3	Dichromate Reflux
	15	R24 - 410.4	Spectrophotometric
	23	R89 - 5220C	Dichromate Reflux
	24	R89 - 5220D	Spectrophotometric
	33	R25 - D1252-88A	Dichromate Reflux
	34	R25 - D1252-88B	Spectrophotometric
	45	R79 - I-3560-85	Dichromate Reflux
	46	R79 - I-3562-85	Dichromate Reflux
	47	R79 - I-3561-85	Spectrophotometric
	52	R8 - p. 17	Dichromate Reflux
	54	R41 - 8000	Spectrophotometric
	55	R42 - p. 1	Spectrophotometric
	56	R81 - 973.46	Dichromate Reflux
	99	Other	
Total	12	R24 - 415.1	Combustion or oxidation
Organic	23	R89 - 5310B	Combustion or oxidation
Carbon (TOC)	24	R89 - 5310C	Combustion or oxidation
	25	R89 - 5310D	Combustion or oxidation
	32	R25 - D2579-85A	Combustion or oxidation
	33	R25 - D2579-85B	Combustion or oxidation
	42	R82 - p. 14	Combustion/Infrared
	52	R81 - 973.47	Combustion or oxidation
	99	Other	
Biochemical	11	R24 - 405.1	Winkler (Azide Modification)
Oxygen Demand	12	R24 - 405.1	Electrode
(5 day BOD)	24	R89 - 5210B	Winkler (Azide Modification)
	25	R89 - 5210B	Electrode
	42	R26 - I-1578-78	Winkler (Azide Modification)
	51	R8 - p. 17	Winkler (Azide Modification)
	53	R81 - 973.44	Winkler (Azide Modification)
	99	Other	
Carbonaceous BOD	22	R89 - 5210B	Carbonaceous BOD
	99	Other	

DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
Total Cyanide	13	R24 - 335.2	Manual Spectrophotometric
	14	R24 - 335.3	Automated Spectro- photometric
	25	R89 - 4500D	Titrimetric
	26	R89 - 4500E	Manual Spectrophotometric
	34	R25 - D2036-91A	Manual Spectrophotometric
	43	R79 - I-3300-85	Manual Spectrophotometric
	51	R8 - p. 22	Titration or Colorimetric
	99	Other	
Non-filterable Residue (Total Suspended Solids)	12	R24 - 160.2	Glass fiber 103° to 105°C
	24	R89 - 2540D	Glass fiber 103° to 105°C
	42	R79 - I-3765-85	Glass fiber 103° to 105°C
	99	Other	
Oil and Grease	12	R24 - 413.1	Freon Extraction/Gravimetric
	13	R102 - 1664	Hexane Extraction/ Gravimetric
	23	R89 - 5520B	Freon Extraction/Gravimetric
	99	Other	

DMR-QA STUDY 18

Analyte	Method Code	Ref.-Method	Description
pH	12	R24 - 150.1	Electrode
	23	R89 - 4500B	Electrode
	32	R25 - D1293-84A(90)	Electrode
	33	R25 - D1293-84B(90)	Electrode
	43	R79 - I-1586-85	Electrode
	52	R43 - p. 1	Automated Electrode
	54	R81 - 973.41	Electrode
	99	Other	
Total Phenolics	11	R24 - 420.1	Manual 4-AAP with distillation
	12	R24 - 420.2	Automated 4-AAP with distillation (with the manual distillation in 420.1)
	21	R3 - 510A + 510B	Manual 4-AAP with distillation and chloroform extraction
	22	R3 - 510A + 510C	Manual 4-AAP with distillation
	99	Other	
Total Residual Chlorine	11	R24 - 330.1	Titrimetric - amperometric
	12	R24 - 330.2	Back-starch end point
	13	R24 - 330.3	Iodometric
	14	R24 - 330.4	DPD-FAS
	15	R24 - 330.5	Spectrophotometric, DPD
	16	R24 - 330.2	Back-amperometric end point
	20	R89 - 4500B	Iodometric
	26	R89 - 4500C	Back-starch end point
	27	R89 - 4500D	Titrimetric - amperometric
	28	R89 - 4500F	DPD - FAS
	29	R89 - 4500G	Spectrophotometric, DPD
	31	R25 - D1253-86(92)	Titrimetric - amperometric
	51	R49 - p. 1	Electrode
	99	Other	

DMR-QA STUDY 18

METHOD REFERENCES

- R3. "Standard Methods for the Examination of Water and Waste Water," 14th ed., 1976. Available from the American Public Health Association, 1015 18th Street, N.W., Washington, DC 20036.
- R8. American National Standards on Photographic Processing Effluents, April 2, 1975. Available from ANSI, 1430 Broadway, New York, NY 10018.
- R19. "Inductively Coupled Plasma-Atomic Emission Spectrometric Method for Trace Element Analysis of Water and Wastes," Appendix C, 40 CFR, Part 136.
- R24. "Methods of Chemical Analysis of Water and Wastes." U.S. Environmental Protection Agency, NERL, Cincinnati, OH 45268, March, 1983 (EPA-600/4-79-020). Available from ORD Publications, CERL, EPA, Cincinnati, OH 45268.
- R25. Annual Book of ASTM Standards, Volumes 11.01 and 11.02 for Water, American Society for Testing and Materials, 1916 Race Street, Philadelphia, PA 19103.
- R26. Techniques of Water-Resources Investigation of the U.S. Geological Survey, Chapter A-1, "Methods for Determination of Inorganic Substances in Water and Fluvial Sediments," Book 5, 1979, Stock #024-001-03177-9. Available from Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402.
- R40. Ammonia, Automated Electrode Method, Industrial Method Number 379-75WE, dated February 19, 1976, Bran & Luebbe (Technicon) AutoAnalyzer II, Bran & Luebbe Analyzing Technologies, Inc., Elmsford, NY 10523.
- R41. Hach Handbook of Water Analysis, 1979, Hach Chemical Company, P.O. Box 389, Loveland, CO 80537.
- R42. COD Method, Oceanography International Corporation, 512 West Loop, P.O. Box 2980, College Station, TX 77840
- R43. Hydrogen Ion (pH) Automated Electrode Method, Industrial Method Number 378-75WA, October, 1976, Bran & Luebbe (Technicon) AutoAnalyzer II, Bran & Luebbe Analyzing Technologies, Inc., Elmsford, NY 10523.
- R44. Iron, 1,10-Phenanthroline Method, Method 8008, 1980, Hach Chemical Company, P.O. Box 389, Loveland, CO 80537.

DMR-QA STUDY 18

- R49. ORION Research Instruction Manual, Residual Chlorine Electrode Model 97-70, 1977, ORION Research, Inc., 840 Memorial Drive, Cambridge, MA 02138.
- R79. Fishman, M. J., et al. "Methods for Analysis of Inorganic Substances in Water and Fluvial Sediments," U.S. Department of the Interior, Techniques of Water Resources Investigations of the U.S. Geological Survey, Denver, CO, revised 1989.
- R80. "Direct Current Plasma (DCP) Optical Emission Spectrometric Method for Trace Elemental Analysis of Water and Wastes, Method AES0029," 1986 - Revised 1991, Applied Research Laboratories, Inc., 24911 Stanford Avenue, Valencia, CA 91355.
- R81. "Official Methods of Analysis of the Association of Official Analytical Chemists," Methods manual, 15th ed. (1990). Available from the AOAC, 2200 Wilson Blvd., Suite 400, Arlington, VA 22201-3301. (703) 522-3032.
- R82. Wershaw, R.L., et al, "Methods for Analysis of Organic Substances in Water and Fluvial Sediments," Techniques of Water-Resources Investigation of the U.S. Geological Survey, Book 5, Chapter A3, (1972 Revised 1987).
- R89. Standard Methods for Examination of Water and Wastewater, 18th ed., 1992. Available from the American Public Health Association, 1015 18th Street, N.W., Washington, DC 20036.
- R102. "Method 1664; N-Hexane Extractable Material (HEM) and Silica GelTreated N-Hexane Extractable Mataerial (SGT-HEM) by Extraction and Gravimetry (Oil and Grease and Total Petroleum Hydrocarbons)," EPA-821-B-94-004b, April, 1995. Available from the Water Resources Center, Mail Code RC-4100, 401 M Street SW, Washington, DC 20460.

J. T. BAKER CHEMICAL CO. 222 RED SCHOOL LANE, PHILLIPSBURG, NJ 08865
H A T E R I A L S A F E T Y D A T A S H E E T
24-HOUR EMERGENCY TELEPHONE -- (201) 859-2151
CHEMTREC # (800) 424-9300 -- NATIONAL RESPONSE CENTER # (800) 424-8802

N3660 -01 NITRIC ACID PAGE: 1
EFFECTIVE: 10/01/85 ISSUED: 01/23/86

SECTION I - PRODUCT IDENTIFICATION

PRODUCT NAME: NITRIC ACID
FORMULA: HNO3
FORMULA WT: 63.01
CAS NO.: 07697-37-2
NOSH/RTCS NO.: Q05775000
COMMON SYNONYMS: HYDROGEN NITRATE
PRODUCT CODES: 4801, 9602, 9598, 9606, 5371, 9597, 9600, 5113, 9616, 9605, 9601
PRECAUTIONARY LABELLING
BAKER SAF-T-DATA(TH) SYSTEM

HEALTH - 3 (POISON)
FLAMMABILITY - 0
REACTIVITY - 3 (OXIDIZER)
CONTACT - 4 (CORROSIVE)

LABORATORY PROTECTIVE EQUIPMENT

GOGGLES & SHIELD: LAB COAT & APRON: VENT HOOD: PROPER GLOVES

PRECAUTIONARY LABEL STATEMENTS

POISON DANGER
STRONG OXIDIZER - CONTACT WITH OTHER MATERIAL MAY CAUSE FIRE
LIQUID AND VAPOR CAUSE SEVERE BURNS - MAY BE FATAL IF SWALLOWED
HARMFUL IF INHALED AND MAY CAUSE DELAYED LUNG INJURY
SPILLAGE MAY CAUSE FIRE OR LIBERATE DANGEROUS GAS
KEEP FROM CONTACT WITH CLOTHING AND OTHER COMBUSTIBLE MATERIALS. DO NOT
STORE NEAR COMBUSTIBLE MATERIALS. DO NOT GET IN EYES, ON SKIN, OR CLOTHING.
DO NOT BREATHE VAPOR. KEEP IN TIGHTLY CLOSED CONTAINER. USE WITH ADEQUATE
VENTILATION. WASH THOROUGHLY AFTER HANDLING. IN CASE OF FIRE, FLOOD WITH
WATER. FLUSH SPILL AREA WITH WATER SPRAY.

SECTION II - HAZARDOUS COMPONENTS
COMPONENT % CAS NO.

NITRIC ACID 65-75 7697-37-2
SECTION III - PHYSICAL DATA
BOILING POINT: 120 C (248 F) VAPOR PRESSURE(MM HG): 2.9
RELING POINT: -42 C (-44 F) VAPOR DENSITY(AIR=1): 2.5
SPECIFIC GRAVITY: 1.50 EVAPORATION RATE: N/A
(H2O=1) (BUTYL ACETATE=1)

CONTINUED ON PAGE: 2

J. T. BAKER CHEMICAL CO. 222 RED SCHOOL LANE, PHILLIPSBURG, NJ 08865
H A T E R I A L S A F E T Y D A T A S H E E T
24-HOUR EMERGENCY TELEPHONE -- (201) 859-2151
CHEMTREC # (800) 424-9300 -- NATIONAL RESPONSE CENTER # (800) 424-8802

N3660 -01 NITRIC ACID PAGE: 2
EFFECTIVE: 10/01/85 ISSUED: 01/23/86

SECTION III - PHYSICAL DATA (CONTINUED)

SOLUBILITY(H2O): COMPLETE (IN ALL PROPORTIONS) & VOLATILES BY VOLUME: 100
APPEARANCE & ODOR: COLORLESS LIQUID, WITH CHOKING ODOR.
SECTION IV - FIRE AND EXPLOSION HAZARD DATA
FLASH POINT: N/A NEPA 704M RATING: 3-0-0 OXY
FIRE EXTINGUISHING MEDIA
USE WATER SPRAY.

SPECIAL FIRE-FIGHTING PROCEDURES

FIREFIGHTERS SHOULD WEAR PROPER PROTECTIVE EQUIPMENT AND SELF-CONTAINED
BREATHING APPARATUS WITH FULL FACEPIECE OPERATED IN POSITIVE PRESSURE MODE
MOVE CONTAINERS FROM FIRE AREA IF IT CAN BE DONE WITHOUT RISK. USE WATER
TO KEEP FIRE-EXPOSED CONTAINERS COOL.

UNUSUAL FIRE & EXPLOSION HAZARDS

STRONG OXIDIZER. CONTACT WITH OTHER MATERIAL MAY CAUSE FIRE.

TOXIC GASES PRODUCED
NITROGEN OXIDES

SECTION V - HEALTH HAZARD DATA

THRESHOLD LIMIT VALUE (TLV/TWA): 5 MG/M3 (2 PPM)

SHORT-TERM EXPOSURE LIMIT (STEL): 10 MG/M3 (4 PPM)

EFFECTS OF OVEREXPOSURE

LIQUID MAY CAUSE SEVERE BURNS TO SKIN AND EYES.
INHALATION OF VAPORS MAY CAUSE SEVERE IRRITATION OF THE RESPIRATORY SYSTEM
OR UNCONSCIOUSNESS.
INGESTION MAY BE FATAL.

EMERGENCY AND FIRST AID PROCEDURES

CALL A PHYSICIAN.
IF SWALLOWED, DO NOT INDUCE VOMITING; IF CONSCIOUS, GIVE WATER, MILK, OR
MILK OF MAGNESIA.
IF INHALED, REMOVE TO FRESH AIR. IF NOT BREATHING, GIVE ARTIFICIAL
RESPIRATION. IF BREATHING IS DIFFICULT, GIVE OXYGEN.
IN CASE OF CONTACT, IMMEDIATELY FLUSH EYES OR SKIN WITH PLENTY OF WATER FO
AT LEAST 15 MINUTES WHILE REMOVING CONTAMINATED CLOTHING AND SHOES.
WASH CLOTHING BEFORE RE-USE.

CONTINUED ON PAGE: 3

J. T. BAKER CHEMICAL CO. 222 RED SCHOOL LANE, PHILLIPSBURG, NJ 08865
M A T E R I A L S A F E T Y D A T A S H E E T
24-HOUR EMERGENCY TELEPHONE -- (201) 859-2151
CHEMTREC # (800) 424-9300 -- NATIONAL RESPONSE CENTER # (800) 424-8802

N3660 -01
EFFECTIVE: 10/01/85
NITRIC ACID
PAGE: 3
ISSUED: 01/23/86

SECTION VI - REACTIVITY DATA

STABILITY: STABLE
HAZARDOUS POLYMERIZATION: WILL NOT OCCUR

CONDITIONS TO AVOID: HEAT, LIGHT

INCOMPATIBILITIES: STRONG BASES, COMBUSTIBLE MATERIALS,
STRONG REDUCING AGENTS

DECOMPOSITION PRODUCTS: OXIDES OF NITROGEN

SECTION VII - SPILL AND DISPOSAL PROCEDURES

STEPS TO BE TAKEN IN THE EVENT OF A SPILL OR DISCHARGE
WEAR SELF-CONTAINED BREATHING APPARATUS AND FULL PROTECTIVE CLOTHING. STO
LEAK IF YOU CAN DO SO WITHOUT RISK. VENTILATE AREA. NEUTRALIZE SPILL WITH
SODA ASH OR LIME. WITH CLEAN SHOVEL, CAREFULLY PLACE MATERIAL INTO CLEAN,
DRY CONTAINER AND COVER; REMOVE FROM AREA. FLUSH SPILL AREA WITH WATER.
KEEP COMBUSTIBLES (WOOD, PAPER, OIL, ETC.) AWAY FROM SPILLED MATERIAL.

J. T. BAKER NEUTRASORB(R) OR NEUTRASOL(R) "LOW NA+" ACID NEUTRALIZERS
ARE RECOMMENDED FOR SPILLS OF THIS PRODUCT.

DISPOSAL PROCEDURE
DISPOSE IN ACCORDANCE WITH ALL APPLICABLE FEDERAL, STATE, AND LOCAL
ENVIRONMENTAL REGULATIONS.

EPA HAZARDOUS WASTE NUMBER: 0002, 0003 (CORROSIVE, REACTIVE WASTE)

SECTION VIII - PROTECTIVE EQUIPMENT

USE GENERAL OR LOCAL EXHAUST VENTILATION TO MEET
TLV REQUIREMENTS.

RESPIRATORY PROTECTION: RESPIRATORY PROTECTION REQUIRED IF AIRBORNE
CONCENTRATION EXCEEDS TLV. AT CONCENTRATIONS UP
TO 100 PPM, A CHEMICAL CARTRIDGE RESPIRATOR WITH
ACID CARTRIDGE IS RECOMMENDED. ABOVE THIS LEVEL,
A SELF-CONTAINED BREATHING APPARATUS IS ADVISED.

EYE/SKIN PROTECTION: SAFETY GOGGLES AND FACE SHIELD, UNIFORM,
PROTECTIVE SUIT, ACID-RESISTANT GLOVES ARE
RECOMMENDED.

SECTION IX - STORAGE AND HANDLING PRECAUTIONS

"AF"-T-DATA(TM) STORAGE COLOR CODE: YELLOW

SPECIAL PRECAUTIONS
KEEP CONTAINER TIGHTLY CLOSED. STORE SEPARATELY AND AWAY FROM FLAMMABLE

CONTINUED ON PAGE: 4

J. T. BAKER CHEMICAL CO. 222 RED SCHOOL LANE, PHILLIPSBURG, NJ 08865
M A T E R I A L S A F E T Y D A T A S H E E T
24-HOUR EMERGENCY TELEPHONE -- (201) 859-2151
CHEMTREC # (800) 424-9300 -- NATIONAL RESPONSE CENTER # (800) 424-8802

N3660 -01
EFFECTIVE: 10/01/85
NITRIC ACID
PAGE: 4
ISSUED: 01/23/86

SECTION IX - STORAGE AND HANDLING PRECAUTIONS (CONTINUED)

AND COMBUSTIBLE MATERIALS.

SECTION X - TRANSPORTATION DATA AND ADDITIONAL INFORMATION

DOMESTIC (O.O.T.)

PROPER SHIPPING NAME
HAZARD CLASS
UN/NA
LABELS
REPORTABLE QUANTITY
NITRIC ACID (OVER 40%)
OXIDIZER
UN2031
OXIDIZER, CORROSIVE
1000 LBS.

INTERNATIONAL (I.M.O.)

PROPER SHIPPING NAME
HAZARD CLASS
UN/NA
LABELS
NITRIC ACID
8
UN2031
CORROSIVE

(TH) AND (R) DESIGNATE TRADEMARKS.
N/A = NOT APPLICABLE OR NOT AVAILABLE

THE INFORMATION PUBLISHED IN THIS MATERIAL SAFETY DATA SHEET HAS BEEN COMPILED
FROM OUR EXPERIENCE AND DATA PRESENTED IN VARIOUS TECHNICAL PUBLICATIONS. IT IS
THE USER'S RESPONSIBILITY TO DETERMINE THE SUITABILITY OF THIS INFORMATION FOR
THE ADOPTION OF NECESSARY SAFETY PRECAUTIONS. WE RESERVE THE RIGHT TO REVISE
MATERIAL SAFETY DATA SHEETS PERIODICALLY AS NEW INFORMATION BECOMES AVAILABLE.

-- LAST PAGE --



SODIUM CHLORIDE

Material Safety Data Sheet

Emergency Telephone Number
314-982-5000

Mallinckrodt Inc.
Science Products Division
P.O. Box M
Paris, Kentucky 40361

Effective Date: 08-07-85

PRODUCT IDENTIFICATION:

Synonyms: Salt

Formula CAS No.: 7647-14-5

Molecular Weight: 58.44

Chemical Formula: NaCl

Hazardous Ingredients:
None.PRECAUTIONARY MEASURES

As part of good industrial and personal hygiene and safety procedure, avoid all unnecessary exposure to the chemical substance and insure prompt removal from skin, eyes and clothing.

SEE SECTION 5.

DOT Hazard Class: Not Regulated

Physical Data

Appearance: White crystalline.

Odor: Odorless.

Solubility: 36g/100cc water @ 20°C (68°F)

Boiling Point: 1413°C (2575°F)

Vapor Density (Air=1): No Information
found.

Melting Point: 801°C (1474°F)

Vapor Pressure (mm Hg): 1.0 @ 845°C
(1553°F)

Specific Gravity: 2.16

Evaporation Rate: No Information found

-2-

Fire and Explosion InformationSECTION 2

Fire: Not considered to be a fire hazard.

Explosion: Not considered to be an explosion hazard.

Fire Extinguishing Media:

Use any means suitable for extinguishing surrounding fire.

Special Information:

In the event of a fire, wear full protective clothing and NIOSH-approved self-contained breathing apparatus with full facepiece operated in the pressure demand or other positive pressure mode.

Reactivity DataSECTION 3

Stability:

Stable under ordinary conditions of use and storage.

Hazardous Decomposition Products:

When heated to above 801°C (1474°F) it emits toxic fumes of chloride and sodium oxide

Hazardous Polymerization:

Will not occur.

Incompatibilities:

Lithium, bromide, trifluoride.

Leak/Spill Disposal InformationSECTION 4

Spills: Sweep up and containerize for reclamation or disposal. Vacuuming or wet sweeping may be used to avoid dust dispersal.

Disposal: Whatever cannot be saved for reclamation may be delivered to an approved waste disposal facility, or if local ordinances allow, can be dissolved in sufficient amounts of water to meet water quality standards, and flushed down a sewer drain.

Ensure compliance with local, state and federal regulations.

-3-

Health Hazard Information

A. Exposure/Health Effects

SECTION 3

Inhalation:

Inhalation of dust may cause mild irritation to mucous membranes, nose and throat. Symptoms may include coughing, dryness, and sore throat.

Ingestion:

Very large doses can cause vomiting, diarrhea, and prostration. Dehydration and congestion occur in most internal organs. Hypertonic salt solutions can produce violent inflammatory reactions in the gastrointestinal tract.

Skin Contact:

Not expected to be a health hazard.

Eye Contact:

May cause irritation.

Chronic Exposure:

No information found.

Aggravation of Pre-existing Conditions:

No information found.

B. FIRST AID

Inhalation:

Remove to fresh air. Get medical attention for any breathing difficulty.

Ingestion:

If large amounts were swallowed, get medical advice.

Skin Exposure:

Wash exposed area with soap and water. Get medical advice if irritation develops.

Eye Exposure:

Wash thoroughly with running water. Get medical advice if irritation develops.

C. TOXICITY DATA (NTPCS, 1982)

Oral rat LD50: 3000 mg/kg.
Reproductive effects cited.

-4-

Occupational Control Measures

SECTION 4

Airborne Exposure Limits:

None established

Ventilation Systems:

In general, dilution ventilation is a satisfactory health hazard control for this substance. However, if conditions of use create discomfort to the worker, a local exhaust system should be considered.

Personal Respirators: (NIOSH Approved)

For conditions of use where exposure to the dust is apparent, a dust/mist respirator may be worn. For emergencies, a self-contained breathing apparatus may be necessary.

Skin Protection:

Wear protective gloves and clean body-covering clothing.

Eye Protection:

Use chemical safety goggles. Contact lenses should not be worn when working with this material.

Maintain eye wash fountain and quick-drench facilities in work area.

Storage and Special Information

SECTION 7

Keep in a tightly closed container, stored in a cool, dry, ventilated area. Protect against physical damage.

The information contained herein is provided in good faith and is believed to be correct as of the date hereof. However, Mallinckrodt, Inc. makes no representation as to the comprehensiveness or accuracy of the information. It is expected that individuals receiving the information will exercise their independent judgment in determining its appropriateness for a particular purpose. Accordingly, Mallinckrodt, Inc. will not be responsible for damages of any kind resulting from the use of or reliance upon such information. NO REPRESENTATIONS, OR WARRANTIES, EITHER EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR OF ANY OTHER NATURE ARE MADE HEREUNDER WITH RESPECT TO THE INFORMATION SET FORTH HEREIN OR TO THE PRODUCT TO WHICH THE INFORMATION REFERS.



MATERIAL SAFETY DATA SHEET

EM SCIENCE
A Division of EM Industries
P.O. Box 70
490 Democrat Rd.
Gibbstown, N.J. 08027

BIOMETRICS CORP.
16 TRANCE DRIVE
LINCINNATI, OH
43146

Preparation Date: OCT 27, '88
Data Sent to Customer: NOV 5, '88
Information Phone Number:
(609) 354-5200
Emergency Number:
(800) 424-9300

NFPA HAZARD RATINGS

Health: 0
Reactivity: 0
Special Hazards: N/A
Flammability: 3
Catalog Number(s): SECTION I - GENERAL INFORMATION
PX1815
PX1824

Chemical Name: Propyl Alcohol
Trade Name: 71-259
Chemical Family: Alcohol
Formula: CH₃CH₂CH₂OH
Molecular Weight: 60.10
Boiling Point: 97.2°C
DOT Shipping Name: N/A

SECTION II - HAZARDOUS INGREDIENTS

N/A

SECTION III - PHYSICAL DATA

Boiling Point (C, 760 mm Hg): 97.2°C
Melting Point (C): -126°C
Vapor Pressure (mm Hg): 150
Vapor Density (Air = 1): 2.0
Relative Volatility (Air = 1): 1.1
Vapor Density (Air = 1): 2.0
Vapor Density (Air = 1): 2.0
Reactivity: 0
Stability: Stable
Solubility: Soluble
Specific Gravity: 0.81
Appearance: Clear, colorless liquid
Odor: Faint, alcoholic

SECTION IV - FIRE & EXPLOSION HAZARD DATA

Flash Point (F): 74F (CC)
Flammable Limits (LFL, UFL): 3.1, 13.1
Explosion Limits (LFL, UFL): 3.1, 13.1
Explosion Index (MEL): 13.1
Fire Fighting Media: "Alcohol" foam, Water spray
Fire Fighting Precautions: Wear self-contained breathing apparatus, protective clothing
Fire & Explosion Hazards: Vapor can travel distances to ignition source and flash back

SECTION V - HEALTH HAZARD DATA (ACUTE AND CHRONIC)

ACUTE TOXICITY (LD50): 5700 mg/kg
Toxicity Data: orl-ret LD50: 1870 mg/kg
Symptoms of Exposure: Irritation to eyes, nose, and throat
Vapor irritating to eyes, nose, and throat
Effects from Inhalation or Ingestion: Headache, nausea, dizziness, irritation from ingestion of high concentrations
Inhalation: Irritation of respiratory tract
Ingestion: Irritation of gastrointestinal tract
Routes of Entry: Inhalation, ingestion, skin contact
Cardiovascular: No data available
Carcinogenicity: No data available
Reproductive: No data available
Emergency First Aid: Wash thoroughly with soap/water
Skin: Wash thoroughly with soap/water
Eyes: Wash thoroughly with water for at least 15 minutes
Inhalation: Remove to fresh air; give artificial respiration if breathing has stopped
Ingestion: If conscious, induce vomiting
GET MEDICAL ASSISTANCE FOR ALL CASES OF OVEREXPOSURE

SECTION VI - REACTIVITY DATA

MSDS-PX1815 Page 8 : 01
The statements contained herein are offered for informational purposes only and are based upon technical data that EM SCIENCE believes to be accurate. It is intended for use only by persons having the necessary technical skill and a thorough understanding of the chemical and physical properties of the material. EM SCIENCE makes NO WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS OR OTHERWISE.



Preparation Date: OCT 27, '88
Data Sent to Customer: NOV 5, '88
Information Phone Number:
(609) 354-5200
Emergency Number:
(800) 424-9300

NFPA HAZARD RATINGS

Health: 0
Reactivity: 0
Special Hazards: N/A
Flammability: 3
Catalog Number(s): SECTION I - GENERAL INFORMATION
PX1815
PX1824

Chemical Name: Propyl Alcohol
Trade Name: 71-259
Chemical Family: Alcohol
Formula: CH₃CH₂CH₂OH
Molecular Weight: 60.10
Boiling Point: 97.2°C
DOT Shipping Name: N/A

SECTION II - HAZARDOUS INGREDIENTS

N/A

SECTION III - PHYSICAL DATA

Boiling Point (C, 760 mm Hg): 97.2°C
Melting Point (C): -126°C
Vapor Pressure (mm Hg): 150
Vapor Density (Air = 1): 2.0
Relative Volatility (Air = 1): 1.1
Vapor Density (Air = 1): 2.0
Vapor Density (Air = 1): 2.0
Reactivity: 0
Stability: Stable
Solubility: Soluble
Specific Gravity: 0.81
Appearance: Clear, colorless liquid
Odor: Faint, alcoholic

SECTION IV - FIRE & EXPLOSION HAZARD DATA

Flash Point (F): 74F (CC)
Flammable Limits (LFL, UFL): 3.1, 13.1
Explosion Limits (LFL, UFL): 3.1, 13.1
Explosion Index (MEL): 13.1
Fire Fighting Media: "Alcohol" foam, Water spray
Fire Fighting Precautions: Wear self-contained breathing apparatus, protective clothing
Fire & Explosion Hazards: Vapor can travel distances to ignition source and flash back

SECTION V - HEALTH HAZARD DATA (ACUTE AND CHRONIC)

ACUTE TOXICITY (LD50): 5700 mg/kg
Toxicity Data: orl-ret LD50: 1870 mg/kg
Symptoms of Exposure: Irritation to eyes, nose, and throat
Vapor irritating to eyes, nose, and throat
Effects from Inhalation or Ingestion: Headache, nausea, dizziness, irritation from ingestion of high concentrations
Inhalation: Irritation of respiratory tract
Ingestion: Irritation of gastrointestinal tract
Routes of Entry: Inhalation, ingestion, skin contact
Cardiovascular: No data available
Carcinogenicity: No data available
Reproductive: No data available
Emergency First Aid: Wash thoroughly with soap/water
Skin: Wash thoroughly with water for at least 15 minutes
Eyes: Wash thoroughly with water for at least 15 minutes
Inhalation: Remove to fresh air; give artificial respiration if breathing has stopped
Ingestion: If conscious, induce vomiting
GET MEDICAL ASSISTANCE FOR ALL CASES OF OVEREXPOSURE

SECTION VI - REACTIVITY DATA

MSDS-PX1815 Page 8 : 01
The statements contained herein are offered for informational purposes only and are based upon technical data that EM SCIENCE believes to be accurate. It is intended for use only by persons having the necessary technical skill and a thorough understanding of the chemical and physical properties of the material. EM SCIENCE makes NO WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS OR OTHERWISE.

Stability: Stable
Reactivity: 0
Special Hazards: N/A
Flammability: 3
Catalog Number(s): SECTION I - GENERAL INFORMATION
PX1815
PX1824

Stability: Stable
Reactivity: 0
Special Hazards: N/A
Flammability: 3
Catalog Number(s): SECTION I - GENERAL INFORMATION
PX1815
PX1824

SECTION VII - ENVIRONMENTAL PROTECTION PROCEDURES

Spill Response:
- Eliminate ignition sources; collect on absorbent
- Flush area with water
- Waste disposal must be performed in compliance with all current local, state and federal regulations.

SECTION VIII - SPECIAL PROTECTION INFORMATION

Ventilation, respiratory protection, protective clothing, eye protection:
- Material should be handled or transferred in an approved fume hood
or with adequate ventilation
- Protective gloves (nitrile or equivalent) should be worn
- Safety glasses with side shields must be worn at all times

SECTION IX - SPECIAL PROTECTION INFORMATION

Handling & Storage:
- Keep container closed
- Store in well-ventilated area away from ignition sources
- Do not breathe vapor
- Do not get in eyes; avoid skin contact
- Electrically grounded all equipment when handling this product
- Residual residue may make empty containers hazardous; use caution
- Wash thoroughly after handling. Do not eat, drink, or smoke
- Internally, eye wash and safety equipment should be readily available.

SECTION X - OTHER INFORMATION

Comments: Tests on laboratory animals indicate material may produce adverse mutagenic effects and cause tumors
- Tests on laboratory animals indicate material may produce adverse mutagenic effects and cause tumors

Revision History: OCT 27, '87
N/A = Not available

Rev. 6/86

MSDS-PX1815 Page 8 : 02
The statements contained herein are offered for informational purposes only and are based upon technical data that EM SCIENCE believes to be accurate. It is intended for use only by persons having the necessary technical skill and a thorough understanding of the chemical and physical properties of the material. EM SCIENCE makes NO WARRANTY, EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS OR OTHERWISE.



J. T. Baker Inc.

222 Reg School Lane
Phillipsburg, N.J. 08865
24 Hour Emergency Telephone: (201) 859-2151
Chemtec # (800) 424-9300
National Response Center # (800) 424-8802

MATERIAL SAFETY DATA SHEET

Page: 1
Effective: 04/15/87
Issued: 11/09/89

SECTION I - PRODUCT IDENTIFICATION

Product Name: Glycerol, Anhydrous
Formula: $\text{HOCH}_2\text{CHOHCH}_2\text{OH}$
Formula Wt: 92.10
CAS No.: 56-81-5
HIOER/RECS No.: H48050000
Common Synonyms: 1,2,3-Propanetriol; Glycerin; Glycyl Alcohol
Product Codes: 2141, 2135, M778, 2136

PRECAUTIONARY LABELLING

BAKER SAF-T-DATA System

HEALTH	FLAMMABILITY	REACTIVITY	CONTACT
1	1	0	1
SLIGHT	SLIGHT	NONE	SLIGHT

Laboratory Protective Equipment



Precautionary Label Statements

CAUTION!
MAY CAUSE IRRITATION
MAY BE HARMFUL IF SWALLOWED
During use avoid contact with eyes, skin, clothing. Wash thoroughly after handling. When not in use keep in tightly closed container.

SAF-T-DATA Storage Color Code: Orange (general storage)

SECTION II - HAZARDOUS COMPONENTS

Component	CAS No.
Glycerol, Anhydrous	90-100 56-81-5

SECTION III - PHYSICAL DATA

Boiling Point: 290°C (554°F) Vapor Pressure (mmHg): <0.1

Continued on Page: 2

0007



J. T. Baker Inc.

222 Reg School Lane
Phillipsburg, N.J. 08865
24 Hour Emergency Telephone: (201) 859-2151
Chemtec # (800) 424-9300
National Response Center # (800) 424-8802

MATERIAL SAFETY DATA SHEET

Page: 2
Effective: 04/15/87
Issued: 11/09/89

SECTION III - PHYSICAL DATA (Continued)

Melting Point: 18°C (64°F) Vapor Density (air=1): 3.1
Specific Gravity: 1.26 Evaporation Rate: N/A
($\text{H}_2\text{O}=1$) (Butyl Acetate=1)
Solubility (H_2O): Complete (in all proportions) & Volatiles by Volume: 0

Appearance & Odor: Clear, colorless, odorless syrupy liquid.

SECTION IV - FIRE AND EXPLOSION HAZARD DATA

Flash Point (Closed Cup): 199°C (390°F) NFPA 704M Rating: 1-1-0

Flammable Limits: Upper - N/A & Lower - 0.9 %

Fire Extinguishing Media

Use extinguishing media appropriate for surrounding fire.

Special Fire-Fighting Procedures

Firefighters should wear proper protective equipment and self-contained breathing apparatus with full facepiece operated in positive pressure mode

Unusual Fire & Explosion Hazards

Contact with strong oxidizers may cause fire or explosion.

Toxic Gases Produced: carbon monoxide, carbon dioxide, acrolein

SECTION V - HEALTH HAZARD DATA

Threshold Limit Value (TLV/TWA): 10 mg/m^3 (ppm)

Toxicity:	LD ₅₀ (oral-rat) (g/kg)	- 12.6
	LD ₅₀ (ipr-mouse) (mg/kg)	- 63
	LD ₅₀ (acu-rat) (mg/kg)	- 100
	LD ₅₀ (iu-rat) (mg/kg)	- 5566

Carcinogenicity: NTP: No IARC: No 2 List: No OSHA reg: No

Effects of Overexposure

Inhalation of vapors may cause coughing and difficult breathing.
Contact with skin or eyes may cause irritation.
Prolonged skin contact may cause dermatitis.

Continued on Page: 3

0007

ABD00060341



MATERIAL SAFETY DATA SHEET

J. T. Baker Inc.
222 Red School Lane
Phillipsburg, N.J. 08865
24 Hour Emergency Telephone - (201) 859-2151
Chemtrec # (800) 424-9300
National Response Center # (800) 424-8802

Glycerol, Anhydrous

C4836 -05 Page: 3
Effective: 04/15/87 Issued: 11/09/88

SECTION U - HEALTH HAZARD DATA (Continued)

Ingestion may cause headache, nausea, vomiting, gastrointestinal irritation, convulsions, and unconsciousness.

Target Organs

eyes, skin, mucous membranes

Medical Conditions Generally Aggravated By Exposure
None Identified

Routes Of Entry

eye contact, skin contact, inhalation, ingestion

Emergency and First Aid Procedures

INGESTION: If swallowed and the person is conscious, immediately give large amounts of water. Get medical attention.

INHALATION: If a person breathes in large amounts, move the exposed person to fresh air. Get medical attention.

EYE CONTACT: Immediately flush with plenty of water for at least 15 minutes. Get medical attention.

SKIN CONTACT: Immediately wash with plenty of soap and water for at least 15 minutes.

SECTION UI - REACTIVITY DATA

Stability: Stable Hazardous Polymerization: Will not occur

Incompatibles: strong oxidizing agents, strong acids, strong bases

Decomposition Products: carbon monoxide, carbon dioxide, acrolein

SECTION UII - SPILL AND DISPOSAL PROCEDURES

Steps to be taken in the event of a spill or discharge

Wear suitable protective clothing. Take up with sand or other noncombustible absorbent material and place into container for later disposal. Flush spill area with water.

Disposal Procedures

Dispose in accordance with all applicable federal, state, and local environmental regulations.

SECTION UIII - INDUSTRIAL PROTECTIVE EQUIPMENT

Ventilation: Use general or local exhaust ventilation to meet TLU requirements.

Respiratory Protection: Respiratory protection required if airborne concentration exceeds TLU. At concentrations above 10 mg/m³, a dust/mist respirator is

Continued on Page: 4

0007



J. T. Baker Inc.
222 Red School Lane
Phillipsburg, N.J. 08865
24 Hour Emergency Telephone - (201) 859-2151
Chemtrec # (800) 424-9300
National Response Center # (800) 424-8802

Glycerol, Anhydrous

C4836 -05 Page: 4
Effective: 04/15/87 Issued: 11/09/88

SECTION UIII - INDUSTRIAL PROTECTIVE EQUIPMENT (Continued)

recommended.

Eye/Skin Protection: Safety goggles, uniform, apron, proper gloves are recommended.

SECTION IX - STORAGE AND HANDLING PRECAUTIONS

SAF-T-DATA* Storage Color Code: Orange (general storage)

Special Precautions

Keep container tightly closed. Suitable for any general chemical storage area.

Isolate from incompatible materials.

SECTION X - TRANSPORTATION DATA AND ADDITIONAL INFORMATION

DOMESTIC (D.O.T.)

Proper Shipping Name Chemicals, n.o.s. (Non-regulated)

INTERNATIONAL (I.M.O.)

Proper Shipping Name Chemicals, n.o.s. (Non-regulated)

N/A - Not Applicable or Not Available

The information published in this Material Safety Data Sheet has been compiled from our experience and data presented in various technical publications. It is the user's responsibility to determine the suitability of this information for the adoption of necessary safety precautions. We reserve the right to revise Material Safety Data Sheets periodically as new information becomes available. J. T. Baker Inc. makes no warranty or representation about the accuracy or completeness nor fitness for purpose of the information contained herein.

COPYRIGHT 1988 J.T. BAKER INC.
* TRADEMARKS OF J.T. BAKER INC.

-- LAST PAGE --

0007

DMR-QA STATE COORDINATORS

*** Alabama**

E. John Williford
AL Dept. of Env. Mgt.
1890 A. Cong. W.L. Dickinson Dr.
Montgomery, AL 36109
(334) 260-2700

Alaska

Karen Hoover
US EPA, Region X
OW-133
1200 Sixth Ave.
Seattle, WA 98101
(206) 553-1213

American Samoa

Sheila Wiegman
American Samoa EPA
Environmental Quality Commission
Pago Pago, AS 96799
(684) 633-2304

Arizona

Gary Brussels
AZ Dept. of Health Services
Environmental Lab Licensure
3443 North Central Ave., Suite 810
Phoenix, AZ 85012
(602) 255-3454

*** Arkansas**

Dick Cassat
AR Dept of Poll. Ctrl. & Ecology
8001 National Drive
Little Rock, AR 72209
(501) 682-0744

California

Bill Ray
CA State Water Res. Control Board
PO Box 944213
Sacramento, CA 95814-0100
(916) 657-1123

Colorado

Carla Lenkey
CO Dept. of Public Division
Field Support Section
4300 S. Cherry Creek Dr.
Denver, CO 80222-1530
(303) 692-3603

Connecticut

Donald Gonyea
CT Dept. of Env. Protection
Water Management Bureau
PO Box 5066
Hartford, CT 06106-5066
(860) 424-3827

Delaware

Joe Mulrooney
DNR & EC
Water Pollution Branch
PO Box 1401
Dover, DE 19903
(302) 739-5731

District of Columbia

William Ruby
Environmental Regulation Admin.
Water Resources Management Div.
2100 M.L. King Jr. Ave, SE
Washington, DC 20020
(202) 645-6601 ext. 3032

*** Florida**

Ed Sims
US EPA, Region IV
Water Management Div. (CWAES)
100 Alabama St., SW - Atlanta Federal
Center
Atlanta, GA 30303-3104
(404) 562-9768

Carlos Boueres
FL Dept. of Env. Protection
Quality Assurance Section
2600 Blear Stone Road, MS 6505
Tallahassee, FL 32399
(904) 488-2796

Georgia

Jeff Larson
GA Dept of Natural Resources
Env. Protection Div.
4244 Int. Pkwy, Suite 110 - Atlanta
Tradeport
Atlanta, GA 30354
(404) 362-2680

Guam

Jesus Salas
Guam EPA
PO Box 22439-GMF
Tiyen Barrigada, GU 96921
(671) 472-8863

*** Hawaii**

Rendy Chow
HI Dept. of Health
State Laboratories Div.
2725 Waimano Home Rd.
Pearl City, HI 96782
(808) 453-6684

Idaho

Karen Hoover
US EPA, Region X
OW-133
1200 Sixth Ave.
Seattle, WA 98101
(206) 553-1213

*** Illinois**

Erin Rednour
Illinois EPA
Bureau of Water
PO Box 19276
Springfield, IL 62702
(217) 782-9720

Indiana

Steve Kim
IN Dept. of Env. Mgt., OWM
Oper. Assist. & Training Section
PO Box 6015
Indianapolis, IN 46204
(317) 232-8793

Iowa

Charles Furrey
IA Dept. of Natural Resources
Henry A. Wallace Bldg.
900 E. Grand
Des Moines, IA 50319
(515) 281-4067

Kansas

Jack McKenzie
Kansas Dept. of Health & Env.
Laboratory Services & Research
Forbes Field Bldg. #740
Topeka, KS 66620-0001
(913) 296-1639

Kentucky

Donna Drury
KY Dept. for Env. Protection
Division of Water
14 Reilly Rd. Ft. Boone Plaza
Frankfort, KY 40601
(502) 564-3410 ext. 461

Louisiana

Elaine Sorbet
LDEQ Water Lab
8618 GSRI
Baton Rouge, LA 70808
(504) 765-2406

Maine

David Dodge
ME Dept. of Env. Protection
Div. of Water Resource Regulation
State House, Station 17
Augusta, ME 04333
(207) 287-7659

*** Maryland**

Marlene Patillo
MD Dept. of the Environment
Div. of Municipal Compliance
2500 Broening Highway
Baltimore, MD 21224
(410) 631-3646

Melvin Knott
MD Dept. of the Environment
Industrial Wastewater Program
2500 Broening Highway
Baltimore, MD 21224
(410) 631-3906

Massachusetts

Ping Yee
MA Dept. of Env. Protection
Div. of Water Pollution Control
Training Center - Route 20
Milbury, MA 01527
(508) 756-7281

*** Michigan**

Clyde Marion
US EPA, Region V
(WC-15J)
77 West Jackson Blvd.
Chicago, IL 60604
(312) 353-5966

Minnesota

Kim Sandrock
MN Pollution Control Agency
520 Lafayette Road
St. Paul, MN 55155
(612) 296-7387

Mississippi

Phillip Bass
MS Dept of Env. Quality
Office of Pollution Control
PO Box 10385
Jackson, MS 39204
(601) 961-5143

Missouri

Jack Pate
MO Dept. of Natural Resources
Water Pollution Control Program
205 Jefferson Street
Jefferson City, MO 65101
(573) 751-1399

Montana

Mike Pasichnyk
MT Dept of Environmental Quality
Water Protection Bureau
PO Box 200901
Helena, MT 59620-0901
(406) 444-5326

Nebraska

Chris Helms
Dept of Env. Quality
Water Quality Division
PO Box 98922
Lincoln, NE 68509-8922
(402) 471-2186

Nevada

Wendall McCurry
NV Div. of Environmental Protection
Bureau of Water Quality Planning
333 W. Nye Lane, Capitol Complex
Carson City, NV 89701
(702) 687-4670 ext. 3098

New Hampshire

Stephanie Larson
NH Dept. of Env. Services
Water Supply & Poll. Ctrl. Div.
P.O. Box 95
Concord, NH 03301
(603) 271-1493

*** New Jersey**

Linda Mauel
US EPA, Region II
2890 Woodbridge Ave. (MS-220)
Edison, NJ 08837-3679
(732) 321-6766

New Mexico

Patrick Hanson
NM Environmental Dept.
Rennels Bldg., Rm. N 2050
PO Box 26110
Santa Fe, NM 87501
(505) 827-2799

ABD00060344
DMR-QA STATE COORDINATORS

*** New York**

Linda Maue
US EPA, Region II
2890 Woodbridge Ave. (MS-220)
Edison, NJ 08837-3679
(732) 321-6766

North Carolina

Jim Meyer
NC DEM/Laboratory Section
NC Dept. of Env. Health & Nat. Res.
4405 Reedy Creek Rd.
Raleigh, NC 27607-4405
(919) 733-3908

North Dakota

Jean Pfeifer
ND Dept of Health
Div. of Water Quality, Missouri Office Bldg.
1200 Missouri Avenue
Bismarck, ND 58505-5520
(701) 328-5228

Northern Islands

John I. Castro Jr.
Div. of Env. Quality
Mariana Island
PO Box 1304
Saipan, CM 96950
(670) 234-6950

*** Ohio**

Susan Plank
Ohio EPA
1571 Perry Street
Columbus, OH 43201
(614) 644-4240

Oklahoma

Aaron Milligan
OK Dept of Environmental Qlty.
1000 NE Tenth, 10th Floor
Oklahoma City, OK 73117-1212
(405) 271-5205

*** Oregon**

Judy Johndahl
OR Dept. of Environmental Quality
Executive Bldg.
811 SW Sixth Ave.
Portland, OR 97204
(503) 229-6896

Pennsylvania

R. Laurie Wyrick
PA DEP - Bureau of Water Quality Mgt.
Div. of Permits & Compliance
PO Box 8465
Harrisburg, PA 17105-8465
(717) 783-2940

Cary Pesek

PA DEP
Northwest Region, Water Mgt. Program (DMR-QA)
101 South Mercer Street
New Castle, PA 16101
(412) 656-3267

David Long

PA DEP
North-central Region, Water Mgt. Program (DMR-QA)
208 W. Third St., Suite 101
Williamsport, PA 17701-6448
(717) 327-3781

Randy King

PA DEP
South-central Region, Water Mgt. Program (DMR-QA)
One Ararat Boulevard
Harrisburg, PA 17110
(717) 657-4671

Thomas Sherk

PA DEP
Northeast Region, Water Mgt. Program (DMR-QA)
2 Public Square
Wilkes-Barre, PA 18711-0790
(717) 826-2533

David Burke

PA DEP
Southeast Region, Water Mgt. Program (DMR-QA)
555 N. Lane, Lee Park, Suite 6010
Conshohocken, PA 19428
(610) 832-6106

Charles Brethauer

PA DEP
Southwest Region, Water Mgt. Program (DMR-QA)
400 Water Front Drive
Pittsburgh, PA 15222-4745
(412) 442-4328

Puerto Rico

Linda Maue
US EPA, Region II
2890 Woodbridge Ave. (MS-220)
Edison, NJ 08837-3679
(732) 321-6766

Rhode Island

Benjamin Lovesky
RI Dept. of Env. Mgt.
Div. of Water Resources
235 Promenade Street
Providence, RI 02908-5767
(401) 277-3961 ext. 7268

South Carolina

Wayne Davis
SC Dept of Health & Env. Control
Laboratory Certification
PO Box 72
State Park, SC 29147
(803) 935-6856

South Dakota

Brian Zinda
SD Dept. of Env. & Natural Res.
Point Source Control Program
523 E. Capital, Joe Foss Bldg.
Pierre, SD 57501-3181
(605) 773-3351

Tennessee

Pamela Townsend
TN Dept. of Env. & Conservation
Div. of Water Pollution Control
401 Church Street
Nashville, TN 37243-1234
(615) 532-0677

Texas

Mary Stordal
TNRCC
5144 E. Sam Houston Pkwy. N.
Houston, TX 77015
(281) 457-5229

Trust Territory

Patricia Mack
US EPA, Region IX
Laboratory
1337 South 46th St., Bldg. 201
Richmond, CA 94804-4698
(510) 412-2333

Carolyn Tambwekar

US EPA, Region IX
1337 South 46th Street, Bldg 201
Richmond, CA 94804-4698
(510) 412-2383

*** Utah**

Mike Herkimer
UT Dept. of Env. Quality
Div. of Water Quality
288 North 1460 West
Salt Lake City, UT 84114-4870
(801) 538-6146

Vermont

Andrew Fish
VT Dept. of Env. Conservation
103 S. Main Street, Sewing Bldg.
Waterbury, VT 05676
(802) 241-3739

Virgin Islands

Linda Maue
US EPA, Region II
2890 Woodbridge Ave. (MS-220)
Edison, NJ 08837-3679
(732) 321-6766

Virginia

Elizabeth Ziomek
VA Dept. of Env. Qlty.
Water Division
PO Box 10009
Richmond, VA 23240
(804) 698-4181

Washington

Stewart Lombard
WA State Dept. of Ecology
Quality Assurance Section
PO Box 488
Manchester, WA 98353
(360) 895-4649

West Virginia

Don E. Caldwell
State of West Virginia
Dept. of Natural Resources
1201 Greenbrier Street
Charleston, WV 25305
(304) 558-2108

Wisconsin

Mike Kvitrud
WI Dept of Natural Resources
PO Box 7921
Madison, WI 53715
(608) 621-8459

Tom Basher

WI Dept. of Natural Resources
North Central District
PO Box 818
Rhineland, WI 54501
(715) 369-8964

Don Domincick

WI Dept. of Natural Resources
Southeast District
PO Box 12436
Milwaukee, WI 53212
(414) 263-8717

Colleen Higgins

WI Dept. of Natural Resources
Western District
PO Box 4001
Eau Claire, WI 54702
(715) 839-1603

Janet LaRose

WI Dept. of Natural Resources
Northwest District
810 W. Maple St.
Spooner, WI 54801-1255
(715) 635-4067

Roy Lemke

WI Dept. of Natural Resources
Southern District
3911 Fish Hatchery Road
Fitchburg, WI 53711
(608) 275-3283

Linda Vogen

WI Dept. of Natural Resources
Lake Michigan District
PO Box 10448
Green Bay, WI 54307-0448
(414) 492-5876

Wyoming

Edward Mock
WY Dept. of Env. Quality
Water Quality Division
122 W. 25th St.
Cheyenne, WY 82002
(307) 777-7317