

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPTS. 2,3,4,5,7,
8,20, 44,
45,46,48,
56,60
METHOD 602.009 (WC)
Rev. 7/19/71
Page 1 of 3

MATERIAL

Product SM, AN, MeOH, Bz, Toluene, EB and all other clear and colorless materials requiring an APHA color measurement, - DOP, DIDP, -711- Other Alcohols, Caustic.

ANALYSIS REQUIRED

Color

SAMPLE

One quart bottle

APPARATUS

Nessler (short) tubes 100 ml marked at 150 mm from inside bottom.
Lumetron Model 450 with color filter #B-420.

REAGENTS

Distilled water XXX
APHA 500 Standard

PROCEDURE

1. With switch off, set instrument to 0% transmission by means of the lever in the back.
2. Fill a tube to 150 mm mark with clear distilled water. Insert tube in machine and turn switch on.
3. Adjust the instrument to 100% transmission by means of the adjustment level at the top. Turn switch off and remove tube.
4. Fill a tube with sample and place in instrument. Turn switch on and read % transmission from the scale.

7/20/71

RSV 0014954

PREPARATION OF REAGENTS

None

STANDARDIZATION AND/OR CALIBRATION

Prepare standards by diluting APHA 500 stock at rate of 0, 1, 2, 3, 4, & 5 ml to 100 ml with distilled water. The standards represent 0, 5, 10, 15, 20 & 25 APHA colors. Read the % transmission of these vs. the distilled water used. Plot the %-transmissions thus obtained vs. the corresponding APHA colors. Draw in the best representative line. Use this curve for routine work or prepare a chart from it showing %T vs. APHA color.

SAFETY

Use normal precautions in handling chemicals, glassware and electrical equipment.

ANALYTICAL TIME

5 minutes

PRECISION OF RESULTS AND SIGNIFICANT FIGURES TO REPORT

95% confidence limits are ± 2 at a level of 10 APHA. Report results to the nearest 1 APHA.

SCOPE

The color of the sample is compared to the APHA scale.


REFERENCES

Method 31-16-6

/jw
7/20/71

RSV 0014955

LIGHT ADJUSTMENT ON LUMETRON



-----STEPS IN ADJUSTING LIGHT SOURCE-----

1. When bulb* burns out; remove light housing and, replace bulb.
2. Adjust iris diaphragm to 50% transmittance.
3. Adjust Bar A for maximum deflection to the right.
4. Loosen thumb screw B and adjust housing sideways for maximum deflection to the right and retighten.
5. Loosen thumb screw C and adjust backward or forward for maximum deflection to the right.
6. Adjust dark current to zero.
7. The instrument is ready to use.
8. Follow usual procedures for adjusting light source on a water blank for 100% transmittance.

* The bulb used has a double contact on the base and is coded GE 1130.

/jw
7/20/71

RSV 0014956

MONSANTO CHEMICAL COMPANY
TEXAS CITY, TEXAS

DEPTS. 10, 60
METHOD 602.016 (GC)
Rev. 1 7/20/71
Page 1 of 4

MATERIAL

10-D7-OH, Gulf C₂H₄, Humble C₂H₄, Dept. 2 & and Dept. 3
Meter Runs, Dept. 2 & 3 C₂H₄ Composite, Dept. 23 C₂H₄
Composite, Dept. 41 C₂H₄

ANALYSIS REQUIRED

CH₄, C₂H₂, C₂H₄, C₂H₆

SAMPLE

1.7 liter oxygen bomb with positive sample pressure.

APPARATUS

Perkin-Elmer Vapor Fractometer, Model 154B/or equivalent
Gas charging valve with 1 cc sample loop
1/4" O.D. aluminum or stainless steel tubing

REAGENTS

30-60 mesh silica gel
Methane, Phillips Research Grade
Ethane, Phillips Research Grade
Ethylene, Phillips Research Grade
Acetylene, Commercial Grade - water scrubbed and dried.

PROCEDURE

A. Column Preparation

1. Pack a 1.5 meters column of 1/4" tubing with 30-60 mesh silica gel.
2. Install the column in the Fractometer. Adjust the temperature to 50°C., determine optimum detector current and set at this value, and adjust the pressure such that the ethylene peak appears in 3-4 minutes after charging the sample. Allow the instrument to reach equilibrium at these conditions.

RSV 0014957

B. Analysis

1. Turn recorder on and set recorder range to 1.
Recorder trace should show no drift for five minutes
before Fractometer is ready for use.
2. Purge 1 cc gas sample loop with sample to be analyzed,
shut off valve from bomb, and enter 1 cc sample
into Fractometer.
3. Record peaks for air, methane, ethane, ethylene, and
acetylene or allow enough time for acetylene to
appear as determined by retention times from a
previous chart run under the same conditions.

CALCULATIONS

1. Draw base lines for all peaks except air.
2. Use rule with 50 divisions per inch and measure peak
height and width at one-half peak height for all peaks.
Measurements should be made to the nearest one-tenth of
a division.
3. Calculate corrected area for all peaks as follows:

$$\text{Corrected area} = h \times w \times r \times F$$

Where: h = peak height
w = width at one-half peak height
r = recorder range
F = calibration factor for that component

4. Calculate mol percentage of each component in the
sample as follows:

$$\text{Mol \%} = \frac{A_{\text{corr.}} \times 100}{\sum (A_{\text{corr.}})}$$

Where: $A_{\text{corr.}}$ = corrected area from Step 3

$\sum (A_{\text{corr.}})$ = total of all corrected areas
calculated in Step 3.

PREPARATION OF REAGENT SOLUTIONS

None necessary.

CALIBRATION

1. Prepare at least two synthetics containing methane, ethane, and acetylene in approximately 95% ethylene. Allow to stand overnight to insure thorough mixing.
2. Analyze as in the procedure under Analysis.
3. Measure peaks as in Steps 1 and 2 under Calculations.
4. Calculate area of each peak as follows:

$$A = h \times w \times r$$

Where: A = area

h = peak height

w = width at one-half peak height

r = recorder range

5. Calculate factors for each component according to the following equation:

$$F = \frac{\%}{A}$$

Where: F = factor

% = mol % of component in synthetic

A = area of component peak calculated in Step 4 above

6. All factors should be normalized such that the factor for ethylene equals 1.000.

SAFETY PRECAUTIONS

No special precautions need be observed other than those normally observed when handling gaseous hydrocarbon samples.

ANALYTICAL TIME

The average time of analysis using this method should be 0.4 hours.

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES FOR REPORTING RESULTS

95% confidence limits using this method should be $\pm$ or - 0.5%. Report results to nearest 0.01%.

SCOPE

This method is designed for samples containing only methane and the C₂'s. Samples containing propane, but not propylene, may also be analyzed using this method as propane has a retention time slightly less than acetylene. If "acetylene" appears in a sample normally known not to contain this component, the sample should be analyzed by another method; eg., infrared, to determine if it is acetylene or propane.

/jw
7/20/71

RSV 0014960

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPTS. 10, 60
METHOD 602.040 (GC)
Rev. 1 7/20/71
Page 1 of 4

MATERIAL

High Purity Ethylene

ANALYSIS REQUIRED

Carbon Dioxide

SAMPLE

Minimum of liter stainless steel bomb with positive sample pressure.

APPARATUS

Perkin-Elmer Vapor Fractometer, Model 154B/or equivalent.

Gas charging valve with 5 cc sample loop.

1/4" O.D. aluminum or stainless steel tubing.

REAGENTS

Activated Charcoal (Absorbent-A, used as purchased from Burrell Corp.)

Di-2-ethyl hexyl sebacate (Octoil "S" or Narcoil 20)

High purity ethylene - CO₂ free

CO₂

PROCEDURE

A. Column Preparation

1. Pack 6 feet of 1/4" aluminum or stainless steel tubing with activated charcoal modified with 4 wt.% di-2-ethyl hexyl sebacate.
2. Install the column in the chromatograph. Adjust the temperature to 70°C, determine optimum detector current and set at this value, and adjust the pressure such that the CO₂ peak appears in approximately 5 minutes after charging the sample. Allow the instrument to reach equilibrium at these conditions.

PROCEDURE (cont'd)

B. Analysis

1. Turn recorder on and set recorder range to 1. Recorder trace should show no drift for five minutes before chromatograph is ready for use.
2. Purge 5 cc gas sample loop with sample to be analyzed, shut off valve from bomb, and enter 5 cc sample into chromatograph.
3. Record all peaks up to and including CO₂ peak or allow enough time for CO₂ to appear as determined by retention time from a previous chart run under the same conditions.

CALCULATIONS

1. Draw base line for CO₂ peak only.
2. Use rule with 50 divisions per inch and measure peak height and width at one-half peak height for the CO₂ peak. Measurements should be made to the nearest one-tenth of a division.
3. Calculate ppm of CO₂ as follows:
$$\text{ppm CO}_2 = h \times w \times r \times F$$

Where h = peak height
w = width at one-half peak height
r = recorder range
F = calibration factor for CO₂, expressed as ppm/unit area

PREPARATION OF REAGENT SOLUTIONS

None necessary.

CALIBRATION

1. Prepare at least two synthetics containing CO₂ in high purity ethylene. These synthetics should be in the range of 10 - 50 ppm CO₂. Allow to stand overnight to insure thorough mixing.
2. Analyze as in the procedure under Analysis.
3. Measure the CO₂ peak as in Steps 1 and 2 under Calculations.
4. Calculate area of CO₂ peak as follows:

$$A = h \times w \times r$$

Where A = Area

h = peak height

w = width at one-half peak height

r = recorder range

5. Calculate CO₂ factor according to the following equation:

$$F = \frac{C}{A}$$

Where F = Factor

C = Concentration of CO₂ in synthetic expressed in ppm

A = Area of CO₂ peak calculated in Step 4 above.

6. Calibrations should be checked daily.

SAFETY PRECAUTIONS

No special precautions need be observed other than those normally observed when handling gaseous hydrocarbon samples.

ANALYTICAL TIME

The average time of analysis using this method should be 0.2 hour.

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES FOR REPORTING RESULTS

95% confidence limits using this method should be $\pm$ or ± 2 ppm.
Report results to nearest ppm.

SCOPE

This method is designed for trace CO₂ analysis of high purity ethylene, although it may be used for the determination of CO₂ in any inert or hydrocarbon gas mixture.

If used for analysis of streams containing larger concentrations of CO₂, it may be desirable to use a smaller sample size. If the sample size is changed, the factor should be adjusted for the new sample size. Samples containing large amounts of butanes and heavier will eventually "over-modify" the charcoal and cause a decrease in the retention time of CO₂. If this occurs, it may be necessary to replace column.

/jw
7/20/71

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPTS. 10, 60
METHOD 602.041 (GC)
REV. 7/20/71
Page 1 of 4

MATERIAL

All gas samples, hydrocarbon and inert samples.

ANALYSIS REQUIRED

Oxygen, Nitrogen and Carbon Monoxide

SAMPLE

Minimum of one liter stainless steel bomb with positive sample pressure.

APPARATUS

Perkin-Elmer Vapor Fractometer, Model 154B/ or equivalent.

Gas charging valve with 25 cc sample loop.

1/4" O.D. aluminum or stainless steel tubing.

REAGENTS

5 Angstrom Molecular Sieve/or 13 Angstrom Molecular Sieve

Air - scrubbed with caustic and dried.

CO

Helium

PROCEDURE

A. Column Preparation

1. Pack a 2 meter length of 1/4" aluminum or stainless steel tubing with 30 - 60 mesh molecular sieve which has been dried in an oven overnight at 110°C.
2. Install the column in the chromatograph. If a 5 Angstrom sieve is used, adjust temperature to 70°C; if a 13 Angstrom sieve is used, adjust temperature to 50°C. Determine optimum detector current and set at this value, and adjust the pressure such that the oxygen peak appears in one minute after charging the sample. Allow the instrument to reach equilibrium at these conditions.

PROCEDURE (Cont'd)

B. Analysis

1. Turn recorder on and set recorder range to 1. Recorder trace should show no drift for five minutes before chromatograph is ready for use.
2. Purge 25 cc gas sample loop with sample to be analyzed, shut off valve from bomb, and enter 25 cc sample into chromatograph.
3. Record peaks for oxygen, nitrogen, methane, and carbon monoxide or allow enough time for carbon monoxide to appear as determined by retention time from a previous chart run under the same conditions.

CALCULATIONS

1. Draw base lines for all peaks except methane.
2. Use rule with 50 divisions per inch and measure peak height and width at one-half peak height for all peaks. Measurements should be made to the nearest one-tenth of a division.
3. Calculate concentration of all components, except methane, in ppm as follows:

$$\text{ppm} = h \times w \times r \times F$$

Where h = peak height

w = width at one-half peak height

r = recorder range

F = calibration factor for that component
expressed as ppm/unit area.

PREPARATION OF REAGENT SOLUTIONS

None necessary.

CALIBRATION

1. Prepare at least two synthetics containing air, which has been caustic-scrubbed and dried, and carbon monoxide in helium. 2.46 mm of air per 100 psia synthetic will be 100 ppm oxygen. Synthetics should be in the range of 25-100 ppm oxygen and carbon monoxide. Allow to stand overnight to insure thorough mixing.

CALIBRATION (Cont'd)

2. Analyze as in the procedure under Analysis.
3. Measure peaks as in Steps 1 and 2 under Calculations.
4. Calculate area of each peak as follows:

$$A = h \times w \times r$$

Where A = Area
h = peak height
w = width at one-half peak height
r = recorder range

5. Calculate factors for each component according to the following equation:

$$F = \frac{C}{A}$$

Where F = Factor
C = concentration of component in synthetic
expressed in ppm
A = area of component peak calculated in
Step 4 above.

6. Calibrations should be checked daily.

SAFETY PRECAUTIONS

No special precautions need be observed other than those normally observed when handling gaseous hydrocarbon samples.

ANALYTICAL TIME

The average time of analysis using this method should be 0.3 hour for a 5 Angstrom column and 0.2 hour for a 13 Angstrom column.

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES FOR REPORTING RESULTS

95% confidence limits using this method should be ± 4 ppm.
Report results to nearest ppm.

SCOPE

This method is designed for high purity ethylene samples; however, it may be used on any inert or gaseous hydrocarbon mixture. If argon is present in the sample, it will be eluted with oxygen. Periodically, depending on usage, the column may need to be regenerated. This is done by purging with nitrogen or helium for 2 hours while the column is maintained at 500° to 550°F. If sample volume is changed, all factors should be adjusted to new sample volume.

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 10, 60
METHOD 602.042 (GC)
Rev. 7/20/71
Page 1 of 5

MATERIAL

High Purity Ethylene

ANALYSIS REQUIRED

Hydrogen

SAMPLE

Minimum of one liter stainless steel sample bomb with positive-sample pressure.

APPARATUS

Loenco, Model 15-B - or equivalent.
Gas charging valve with a 25 cc sample loop.
1/4" O.D. aluminum or stainless steel tubing.

REAGENTS

13 Angstrom molecular sieve, 20-60 mesh.
Nitrogen (to be used as carrier gas)
Hydrogen, electrolytic
Ethylene, Phillips Research Grade

PROCEDURE

A. Column Preparation

1. Activate the molecular sieve by drying in an oven at 400 °F for 3 hours, cool in a desiccator and store in screw-cap bottles. Pack 12 feet of 1/4" aluminum or stainless steel tubing with the activated sieve.
2. Install the column in the chromatograph. With the instrument at room temperature, determine the optimum detector current and set at this value. Adjust the pressure such that the hydrogen peak appears in approximately 1.5 - 2.0 minutes after charging the sample. Allow the instrument to reach equilibrium at these conditions.

RSV 0014969

B. Analysis

1. Turn on recorder and line out equipment.
2. Purge 25 cc sample loop with sample for 2 to 3 minutes or until all air has been removed.
3. Stop purge and after pressure in sample loop has stabilized, enter ethylene sample into instrument.
4. The recorder pen will deflect off scale when sample is entered but should return to base-line point after one minute.
5. Record hydrogen peak.

CALCULATIONS

1. Read peak height and calculate ppm hydrogen in sample.

$$\text{ppm H}_2 = \text{Peak height} \times \text{attenuation} \times \text{factor}$$

PREPARATION OF REAGENT SOLUTIONS

None necessary

CALIBRATION

1. Prepare standards samples of hydrogen in ethylene in the 1 to 50 ppm range.
2. Analyze these standard samples following the steps under procedure.
3. Read peak heights and calculate calibration factor.

$$F = \frac{\text{Concentration of H}_2}{\text{Peak height} \times \text{attenuation}}$$

SAFETY PRECAUTIONS

No special precautions need be observed other than those normally observed when handling gaseous hydrocarbon samples.

ANALYTICAL TIME

The average time of analysis using this method should be
0.2 hour.


PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES FOR REPORTING RESULTS

95% confidence limits using this method should be + or - .5 /ppm
in the 0-50 ppm range. Report results to nearest 1 ppm.

SCOPE

This method is applicable to any gaseous mixture. Some changes
may be desirable as to sample size depending on the concentration
of hydrogen in the sample.

LYS:ff
10/19/64







FIGURE 1
CONCENTRATION & SAMPLING EQUIPMENT

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.043 (G.C.)
Rev. 7/20/71
Page 1 of 4

MATERIAL

High Purity Ethylene

ANALYSIS REQUIRED

Propylene and other olefins.

SAMPLE

Minimum of one liter stainless steel bomb with positive-sample pressure.

APPARATUS

Perkin-Elmer Vapor Fractometer, Model 154B/or equivalent.

Gas charging valve with 1 cc. sample loop (modified for evacuating the sample loop)

$\frac{1}{8}$ " O.D. aluminum or stainless steel tubing.

Vacuum pump capable of producing 50 microns vacuum or better, with a free air capacity of at least 21 liters per minute.

Mercury manometer.

REAGENTS

Celite, 30-80 mesh; a satisfactory source is the Celite Division, Johns Manville Company, New York.

Diethyl Ether, for column preparation.

Helium.

Hexamethylphosphoramide, available from Monsanto Chemical Company, Organic Chemicals Division, 800 N. 12th Street, St. Louis, Missouri.

Propylene, Phillips Research Grade.

RSV 0014974

PROCEDURE

A. Column Preparation

1. Pack 25 feet of $\frac{1}{4}$ " aluminum or stainless steel tubing with 30-80 mesh Celite aggregate containing 0.2 g of hexamethylphosphoramide per gram of aggregate and shape column to fit into chromatograph.
2. Install the column in the chromatograph and purge with helium (approx. 20 ml/min.) for at least 24 hours before using.
3. With the instrument at room temperature, determine the optimum detector current and set at this value, and adjust the pressure such that the propylene peak appears in approximately 7 to 8 minutes after charging the sample. Allow the instrument to reach equilibrium at these conditions.

B. Analysis

1. Turn recorder on and set recorder range to 1. Recorder trace should show no drift for five minutes before chromatograph is ready for use.
2. Purge 1 cc gas sample loop with sample to be analyzed, shut off valve from bomb, and enter 1 cc sample at one atmosphere pressure into chromatograph.
3. Record all peaks up to and including propylene on recorder attenuation of 1, or allow enough time for propylene peak to appear as determined by retention time from a previous chart run under the same conditions.

CALCULATIONS

1. Draw base line for propylene peak only.
2. Use rule with 50 divisions per inch and measure the peak height of the propylene peak. Measurements should be made to nearest one-tenth of a division.

CALCULATIONS (cont'd)

3. Calculate concentration of propylene as follows:

$$C = \frac{D}{S}$$

Where: C = concentration of propylene in sample
expressed in mol per cent.
D = deflection of propylene peak in sample
chromatogram measured in step 2 above.
S = sensitivity of propylene, expressed as
deflection per mole per cent as deter-
mined from the calibration curve.

PREPARATION OF REAGENT SOLUTIONS

None necessary.

CALIBRATION

1. Prepare a sensitivity curve for propylene by adding pure propylene at 150, 100, 50, and 25 mm of pressure in the 1 cc sample loop.
2. Plot the peak height (scale of 50 divisions per inch) versus pressure and extrapolate the curve to zero pressure.
3. Calculate the sensitivity in scale divisions (50 divisions per inch) of deflection per mol per cent for a 1 cc sample from the slope of the curve below 20 mm pressure.
4. Before analyzing each set of samples, check the sensitivity by measuring a secondary standard and make suitable corrections to the calibrations.

SAFETY PRECAUTIONS

No special precautions need be observed other than those normally observed when handling gaseous hydrocarbon samples.

ANALYTICAL TIME

The average time of analysis using this method should be 0.3 hour.

PRECISION OF ANALYSIS AND SIGNIFICANT
FIGURES FOR REPORTING RESULTS

95% confidence limits using this method should be
/ or - 0.0003 mol per cent in the 0-0.025 mol per cent
range. Report results to nearest .001 mol per cent.

SCOPE

This method is for the determination of propylene and
other olefins in high purity ethylene. It may be used
on other gaseous hydrocarbon mixtures with possible
expansion of calibration curves.

/jw
7/21/71

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPTS. 2,3,20, 33,60
METHOD 602.063 (WC)
July 21, 1971
Page 1 of 2

MATERIAL

NaOH receipts, Circulating Caustic, Waste Caustic

ANALYSIS REQUIRED

NaOH content

SAMPLE SIZE

Pint bottle

APPARATUS

50 ml Burette
Oleum bulbs
300 ml Erlenmeyer flask

REAGENTS

Methyl red indicator
0.5 N HCl

PROCEDURE

Place about 1 gram of sample into a tared oleum bulb. Seal and reweigh bulb. Slide the oleum bulb into a 300 ml Erlenmeyer flask containing about 25 ml of distilled water. Cautiously, break the oleum bulb with a stirring rod which has been flattened on one end. Be sure that oleum bulb stem is also crushed with stirring rod. Add about 25 ml more of distilled water. Add a few drops of Methyl red indicator, swirl and titrate to a red end point with 0.5 N HCl. Analyze receipts (only) in duplicate. Single analysis for other caustic samples.

CALCULATIONS

$$\frac{(\text{Vol}_{\text{HCl}}) (\text{N}_{\text{HCl}}) (4.00)}{\text{Sample Weight}} = \% \text{ NaOH}$$

On caustic receipts also report % Na₂O
(% NaOH)
(Average of) (1.0195 = % Na₂O (@ 76%
(Duplicates)

RSV 0014978

PREPARATION OF REAGENTS

Refer to solutions manual

STANDARDIZATION AND CALIBRATION

Refer to solutions manual for standardization of 0.5 N HCl.

SAFETY PRECAUTIONS

Caustic receipts are approximately 50% NaOH. Caution should be used to prevent contact with skin, eyes or clothing. If 50% NaOH does come in contact with the body, wash affected part with copious quantities of water.

ANALYTICAL TIME

0.3 Hours

PRECISION AND SIGNIFICANT FIGURES

Report to the nearest .01%. Duplicate results should agree within 0.5%. Precision is $\leq$ 0.4% at 50% caustic level.

SCOPE

Method is applicable to any range of caustic by varying sample size and normality of HCl.

/jw
7/21/71

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPTS. 2,3,10,20,
33,60
METHOD 602.212 (WC)
7/21/71
Page 1 of 2

MATERIAL

Specification Acrylonitrile and Product Methanol
Muriatic Acid - Caustic Receipts

ANALYSIS REQUIRED

Appearance -

SAMPLE

Quart Bottle

APPARATUS

APHA Color Tube

REAGENTS

None

PROCEDURE

Observe the sample in a clean dry tube of size large enough to hold about 100 ml of sample. Sample should be clear and free of suspended matter. Report results accordingly.

CALCULATIONS

None

STANDARDIZATION AND/OR CALIBRATION

None

SAFETY PRECAUTIONS

Avoid skin contact with AN, Caustic, HCl, MeOH.
Avoid breathing AN and MeOH, HCl, Muriatic Acid and caustic soda vapors.
Keep out of eyes.

RSV 0014980

ANALYTICAL TIME

0.1 hour

PRECISION OF ANALYSIS AND SIGNIFICANT
FIGURES FOR REPORTING RESULTS

Report appearance as clear and free of suspended matter,
cloudy or suspended matter present.

SCOPE

This method is capable of detecting any foreign matter of
large enough particle size to be observed visually.

REFERENCES

Method 20-29-3, MLO, 7-21-51, 10-22-51

/jw
7/21/71

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.244
Rev. 5/31/71
Page 1 of 6

MATERIAL

Acrylamide

ANALYSIS REQUIRED

Melting range

APPARATUS

- (a) Bath - JFQ Standard Apparatus Drawing 5
(b) Bath Liquid - Use appropriate liquid as follows:

<u>Bath Temperature</u>	<u>Liquid</u>
40-85°C	Water
80-200°C	Dow Corning No. 550 Silicone Fluid
190-270°C	Dow Corning No. 550 Silicone Fluid

- (c) Thermometer - Precision Thermometer and Instrument Company, Philadelphia, Pa., round front precision type thermometer, or equivalent, calibrated against a standard or platinum resistance thermometer, selected as follows:

<u>Bath Temperature</u>	
40-85°C	40-85°C range, 3" immersion 0.2°C division, 12" length
80-200°C	80-200°C range, 3" immersion 0.2°C division, 18" length
190-270°C	190-270°C range, 4" immersion 0.2°C division, 20" length

RSV 0014982

APPARATUS (cont'd)

- (d) Magnifier - Surth Number UG-49-D without base, Safety Incorporated, St. Louis, Missouri, type magnifier containing a 2" x 4" lens with a 3" focal length.
- (e) Melting Point Tubing - Friedrich and Dimmack, Inc., P.O. Box 230, Millville, New Jersey.

7" length, 1.4-1.6 mm OD purchased with both ends sealed. Score with a chip of carborundum and break off one end just prior to use. Since it is impractical to fire-polish the broken end, exercise CARE to prevent lacerations.
- (f) Continuously Adjustable Transformer, such as No. V-5 Variac, General Radio Co.
- (g) Variable Speed Stirrer, such as Eastern Industries, Model No. 3, New Haven, Conn.

DETERMINATION

1. Reduce a portion of the sample to a very fine powder if necessary, using a mortar and pestle. When specified, dry the sample at the time and temperature designated in the specification for the product. If the determination is to be made on the sample as is, transfer a portion to a clean piece of glassine paper for filling the melting point tube.
2. By tamping the open end in a small mound of sample, charge the tube with sufficient material to form a column in the bottom of the tube from 4.0 to 5.0 mm high when packed down as closely as possible by moderate reverse tapping on a solid surface. The filling of the tube can be aided by carefully scratching the side of the tube with a round file during the process. The remainder of the sample in the mortar or on the glassine paper should be discarded since it may contain slivers of glass.

DETERMINATION (cont'd)

3. Select the appropriate bath and make certain that it is filled to the proper height. The liquid should be filled so that it is within 1/2" of the top of the bath when at the melting point temperature. Adjust the thermometer to the proper immersion point, making certain that the thermometer is as close to the center of the bath as possible.
4. Start the stirrer and heat the bath rapidly to within approximately 15°C of the expected melting point. Regulate the rise with a Variac to 1°C per minute $\pm$ 10 seconds and maintain at this rate throughout the remainder of the test.
5. Insert the melting point tube exactly 10°C below the expected melting point so that the sample portion of the tube is adjacent to the bulb of the thermometer and clearly visible through the bath window.
6. Read the temperature of the melt at the following two points:
 - a. Meniscus Point - When the liquid phase of the sample forms a distinct meniscus at any level within the sample.
 - b. Clear Point (End of Melt) - When the last evidence of solid disappears and the sample is completely liquid.
7. Apply the thermometer calibration correction and report to 0.1°C the corrected temperatures as the melting point range.

DISCUSSION

The material must be in a fine powdered state to facilitate proper packing of the melting point tube. Granular or crystalline material will allow air voids, thus affecting the amount of material used. The volume of material is critical since amounts varying from the recommended 4.0-5.0 mm can affect the melting range by 0.2°C or more.

DISCUSSION (cont'd)

The position of the melting point tube and the thermometer are critical. They should be placed adjacent to each other and centered in the bath to minimize variations.

It is essential that the proper liquid level be maintained in order to facilitate efficient stirring of the bath. To counteract the greater convection currents in the highest temperature bath, this bath has been made 1" deeper than the other baths. The same liquid level should be maintained, namely 1/2" from the top.

Table 1 lists the significant differences in procedure between this method, the USP XVI, and the BP (1948) methods.

TABLE 1

COMPARISON OF THE MONSANTO MELTING POINT METHOD,
THE USP XVI, AND THE BP (1948) METHODS (1)

<u>Method</u>	<u>Visual Aid</u>	<u>Rate of Rise in Temp., °C per Minute</u>	<u>Alteration in Rate of Temp. Rise</u>	<u>Immersion Point in °C Below Ex- pected Melt</u>	<u>Preheating Time From Immersion to Expected Melt in Minutes</u>
Monsanto	Magnifier	1°	None	10	10
USP XVI	None	3°C	Change rate to 1°C per minute at 3° below the expected melt	30	12
BP (1948)	None	3°	None	10	3-1/3

DISCUSSION (cont'd)

Following are the five possible observable points during the transition from solid to liquid with notations of the points taken in each of the above three methods:

1. Sintering Point - The sample shows evidence of shrinking in volume and its appearance is altered.
2. Liquefaction Point - The sample shows pronounced shrinkage in volume, and the first distinct minute droplet of liquid forms at any level in the sample or on the wall of the capillary. This point closely approximates, or is identical with, the "Beginning of Melting" of the USP XVI method.
3. Flow Point - The sample forms droplets which readily begin to flow together and droplets adhering to the wall of the capillary begin to spread and flow. Shrinkage of the sample is virtually complete.
4. Meniscus Point - The liquid phase of the sample forms a distinct meniscus at any level within the sample. This point is identical with the "Meniscus Point" of this method and the BP (1948).
5. Clear Point (End of Melt) - The point where the last evidence of solid disappears and the sample is completely liquid. This point is identical with the "End of Melt" of this method and the USP XVI.

SCOPE

This method is intended for the determination of the melting point (range) of essentially pure organic compounds. It involves a capillary tube technique similar to a method proposed by Felker (1) for adoption as a standard procedure by the United States Pharmacopeia.

PRECISION AND ACCURACY

Results by this method, compared to those by the USP XVI and the BP (1948) methods will differ product for product since the transition from the solid to liquid phase varies as the complexity of the molecular structure varies. Products such as caffeine and acetophenetidin will respond as follows:

1. The meniscus point will be approximately 0.4°C higher than the USP XVI liquefaction point and 0.2°C lower than the BP (1948) meniscus point, owing primarily to the different rates of heating.
2. The end of melt will be essentially the same in this method and the USP XVI.

The precision of this method is dependent on the purity of the material being tested. Duplicate results on essentially pure material should agree within 0.3°C.

SAFETY

Mercury is a cumulative poison with a relatively high vapor pressure. In the event that a thermometer is broken, collect the mercury and place in a waste jar under water.

REFERENCES

Felker, W. S., Drug Standards. 20, 169 (1952).

/jw
7/22/71

RSV 0014987

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT 60
METHOD 602.245
May 31, 1971
Page 1 of 2

MATERIAL

Acrylamide

ANALYSIS REQUIRED

Solution appearance.

SCOPE

This method is intended to provide an indication of the completeness of solution of materials at specified concentrations in specific solvents. It involves comparison with an equivalent solution of an approved standard, or subjective description of the solution when no standard is available.

EQUIPMENT

- a) Clark Turbidimeter, Catalog No. T-11315, The Chemical Rubber Co., Cleveland, Ohio. (Lamp replaced with a 150 watt showcase lamp.)
- b) Appropriate comparison glassware such as test tubes, Nessler tubes, flasks, or bottles, matched as to volume and tint of glass; thoroughly cleaned, rinsed with distilled water, and dried.
- c) Solvent as designated in the specification for the product being tested. The solvent must be clear and colorless and should be filtered or distilled if necessary.

DETERMINATION

By Comparison with Standards

1. Read the specification for the product and, using a balance of 0.1 g. sensitivity, weigh the designated amount of sample and transfer to the appropriate glassware.

RSV 0014988

DETERMINATION (cont'd)

2. Add the specified volume of clear, colorless solvent and dissolve the sample by agitation at room temperature, or at the temperature designated in the specification.
3. Compare the color of the sample solution with the prescribed standard for color, referring to Method 38-D, E if an APHA color is designated, and noting any abnormality in hue.
4. Compare the turbidity of the sample solution under the turbidimeter in parallel with the prescribed standard, recording the appropriate degree of comparison as: better, slightly better, very slightly better, equivalent, very slightly more turbid, slightly more turbid, or more turbid, and note any difference in particle size, type or opalescence.

NOTE: When product minimum standards selected by production departments are designated, solutions of them are prepared in parallel with the sample, observing the same precautions.

By Description

1. Prepare the specified sample solution as directed in Method 112-A.
2. Describe the color of the solution.
3. Observe the solution under the turbidimeter and describe the degree of turbidity, opalescence, particles or sediment.

SAFETY

Take cognizance of any potential hazards in the solvents and samples, such as flammability, corrosiveness, lachrymatory or dermatological properties, or reactivity, and exercise adequate precautions.

/jw
7/22/71

RSV 0014989

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.246
May 31, 1971
Page 1 of 7

MATERIAL

Acrylamide

ANALYSIS REQUIRED

Moisture

SCOPE

This method is intended as a general directive for the determination of water, titratable with Karl Fischer reagent, in various organic and inorganic substances. It involves the solubilization or dispersion of a sample of the material to be analyzed in a suitable solvent and the titrimetric determination of the moisture present with Karl Fischer reagent. The end point may be detected colorimetrically in colorless or nearly colorless solutions, or electrometrically in colored solutions.

EQUIPMENT

- a) N/100 methanolic iodine - Dilute 10 ml. of N/10 iodine to 100 ml. with methanol.
- b) Standardized regular and low alcohol Karl Fischer reagents in automatic burettes protected with anhydrous silica gel traps. For details of preparation see Method 1 - Section VII.
- c) Auto-Aquatrator, Catalog No. S-29742, Precision Scientific Co., Chicago, Ill.
- d) Common laboratory apparatus and reagents.

APPLICATION

Solids are measured by weight. Liquids of known gravity may be measured from an appropriate volumetric pipette, multiplying volume by specific gravity to obtain weight.

RSV 0014990

<u>Expected Moisture</u>	<u>Sample Size</u>
0-0.5%	10 $\pm$ 0.1 g.
0.5-1.0%	5 $\pm$ 0.1 g.
1.0-5.0%	2 $\pm$ 0.01g.
5.0-12.0%	1 $\pm$ 0.01g.

DETERMINATIONA-Electrometric

1. Turn the line switch ON.
2. Set the selector to "DIRECT ADJUST".
3. Adjust the meter to full scale with the "DIRECT ADJUST" knob.
4. Set the selector to "DIRECT TITRATION".
5. Set the meter Low Set Point (left red needle) to 8 microamperes. (Where the moisture level is 0.05% or less, low set point must be set lower.)
6. Set the meter High Set Point (right red needle) to 15 microamperes.
7. Set the "End Point Hold" to 30 seconds.
8. Fill the K.F. burette.
9. To the titration vessel add a sufficient quantity of the solvent designated on the product specification sheet, to cover the electrode tip.
10. Place the titration vessel under the electrode, swing the support under the vessel and tighten the support until the vessel is snug against the cover. (The system must be air-tight.) Put the cover stopper in place.

DETERMINATION (cont'd)

11. Start the stirrer (stir slowly), apply pressure to the cover stopper and evacuate moist air from the system by pulling on the vacuum switch for 15 seconds.
12. Depress the "Start" button and allow the unit to titrate to the instrument end point.
13. While the solvent is being titrated, weigh-in a suitable container - an appropriate amount of the sample for moisture determination.
14. When the instrument indicates the end of the solvent titration, turn off the stirrer and remove the cover stopper.
15. Carefully add the sample to the titration vessel.
16. Replace the stopper, turn on the stirrer (stir slowly) and evacuate the system for 15 seconds.
17. Fill the K.F. burette, depress the start button, and while the unit is titrating to the instrument endpoint; reweigh the sample container to determine the weight of the sample used.
18. Read the burette and calculate.
$$\% \text{ Moisture} = \frac{\text{ml. titre} \times \text{K.F. factor} \times 100}{\text{sample weight}}$$

B-Visual

1. Measure about 25 ml. of the solvent designated in the specification for the product (or see Discussion) into a dry 125 ml. Erlenmeyer flask.
2. Titre with Karl Fischer reagent to an end point that matches the color of the N/100 methanolic iodine solution, sweeping moisture out of the flask by swirling the solution up the walls of the flask.

DETERMINATION (cont'd)

3. Zero the Karl Fischer burette.
4. Add the prescribed weight of sample to the blanked solvent, swirl to dissolve, and immediately titrate to the same endpoint.

NOTE: If the material to be analyzed is an insoluble porous solid, disperse finely divided sample in the solvent and allow to stand in a stoppered flask with occasional shaking for periods up to one hour at room temperature before titrating with Karl Fischer reagent.

5. Calculate:

$$\% \text{ Moisture} = \frac{\text{ml. titre} \times \text{K.F. factor} \times 100}{\text{sample weight}}$$

NOTE: In cases where pre-treatment of the sample is necessary to eliminate interaction of the material to be analyzed and Karl Fischer reagent, a separate solvent or reagent blank may be necessary. The calculation then becomes:

$$\% \text{ Moisture} = \frac{(\text{Sample ml.} - \text{Blank ml.}) \times \text{K.F. factor} \times 100}{\text{sample weight}}$$

PRECISION AND ACCURACY

- a) Sensitivity, precision, and accuracy depend on several factors, for example, concentration of the Karl Fischer reagent, titration technique, apparatus, quantity of water titrated, and nature of material being analyzed.
- b) The sensitivity is about 0.1 mg of water for visual titrations. Less than 0.02 mg can be measured by electrometric titration.

PRECISION AND ACCURACY (cont'd)

- c) The following is an example of the precision attained in an interlaboratory study on two samples of acetone containing 0.1 percent and 0.4 percent water and two samples of methyl ethyl ketone containing 0.05 percent and 0.17 percent water:



Repeatability - Two results (each the average of duplicate determinations) obtained by the same analyst should be considered suspect if they differ by more than 0.013 percent, absolute (95 percent confidence level). Duplicate determinations which agree within 0.008 percent are acceptable for averaging.

Reproducibility - Two results (each the average of duplicate determinations) obtained by analysts in different laboratories should be considered suspect if they differ by more than 0.028 percent, absolute (95 percent confidence level).



NOTE: The interlaboratory study was carried out by ASTM Committee D-1 on Paint, Varnish, Lacquer, and Related Product, Subcommittee V on Solvents, Plasticizers, and Chemical Intermediates. Seven laboratories participated with a single analyst performing duplicate determinations on each of two days, using two methods on the four samples described above. The Method of Test for Water in Lacquer Solvent and Diluents (Fischer Reagent Titration Method) (ASTM Designation: D1364)⁵ was the subject of the test program being compared with each laboratory's own version of a Karl Fischer method. As neither the means nor the variances of the two sets of data proved significantly different, all of the results were pooled to give estimates of the repeatability based on 55 degrees of freedom and reproducibility based on 47 degrees of freedom.

DISCUSSION

The stoichiometry of the reaction between water and Karl Fischer reagent is as follows:



and



Pyridine is used to combine with acidic products, which will cause a reversal in the following reaction:



It has the added advantage of combining with the sulfur dioxide, thereby reducing the vapor pressure of the latter.

Most substances are inert to Karl Fischer reagent and may be analyzed without special treatment. However, a number of compounds and classes of compounds react stoichiometrically with one or more of the components of Fischer reagent. Included among these are: - Ascorbic Acid, Hydrazine Salts, Substituted Hydrazine Salts, Mercaptans, Alkali Bicarbonates, Alkali Carbonates, Alkali Sulfites, Alkali Pyrosulfites, Boric Acid and Oxides, Carbonyl Compounds, Cupric Salts, Ferric Salts, Metal Hydroxides, Metal Oxides, Sodium Arsenate, Sodium Arsenite, Sodium Tetraborate, Sodium Thiosulfate, and Stannous Chloride.

In general these reactive compounds either may be rendered inert (the use of glacial acetic acid as a solvent eliminates amine and hydrazine interference) or may be estimated by an independent method and a stoichiometric correction applied to the water determination.

DISCUSSION

Although the cyanohydrin technique is the only general method for the determination of water in the presence of carbonyl compounds, low alcohol Karl Fischer reagent and/or special solvents may be used to inhibit interfering side reactions of many ketones and some aldehydes. By greatly reducing the methanol and increasing the pyridine in the regular Fischer solution, the tendency of carbonyls to form acetals and ketals by reacting with the methanol of the standard reagent is reduced. Appreciable quantities of lower alcohols in the sample will permit the interfering acetal or ketal reaction even with the use of modified reagent.

Details of other method variations are beyond the scope of this method and will be given in the methods of analysis for the specific compounds.

In general, methanol may be used as a solvent for the determination of water in inert materials. Amines stronger than benzyl amine ($K_B = 2.4 \times 10^{-5}$) and other basic compounds are preferably treated with or dissolved in an excess of glacial acetic acid or spent Fischer reagent, which in both cases combines with the amines. Some organic acids, formic, acetic, adipic, etc., tend to give slightly high results when methanol is used as a diluent, because of the ease with which they esterify. When equal parts of dry methanol and pyridine are used as a diluent, the inhibiting influence of the amine is sufficient to discourage any esterification of the acid at room temperature.

SAFETY

No unusual hazards are involved in this operation when cognizance of the inflammability and toxicity of the reagents and solvents is taken.

REFERENCE

1. "Aquametry", Mitchell and Smith, Interscience Publishers, N.Y., 1948.
2. ASTM E203-64, 31 pp. 567-578, 1965.

7/22/71

RSV 0014996

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.247
May 31, 1971
Page 1 of 2

MATERIAL

Acrylamide, Sodium lauryl sulfate

ANALYSIS REQUIRED

Sulfate (as SO_4^{2-})

SCOPE

This method is intended for the determination of soluble sulfates in concentrations around 0.1% in essentially water soluble materials. It involves gravimetric determination as barium sulfate.

EQUIPMENT

- a) Usual laboratory apparatus and reagents.
- b) 20 ml. fine porosity Sela filtering crucible, Catalog No. 8-227, Fisher Scientific Co.

DETERMINATION

1. Dissolve 10 grams of sample, weighed on a prescription balance, in 100ml. of water in a 600ml. beaker.
2. Neutralize the solution to methyl orange with 10% HCl and add 4 ml. in excess.

NOTE: If the solution is not completely clear at this point, filter by gravity through Whatman 41 paper into a 600 ml. beaker.

3. Dilute the solution to about 400 ml. and heat to boiling.
4. While boiling gently, add 25 ml. and heat to boiling.
5. Cover the beaker with a watch glass and digest the solution on a steam bath for one hour, or until the supernatant liquid is clear.

RSV 0014997

DETERMINATION (cont'd)

6. Filter through an ignited and tared fine porosity Sela crucible, transferring the precipitate quantitatively with a stream of water from a wash bottle, rinsing the beaker thoroughly, and washing down the inside walls of the crucible.
7. Ignite the crucible in a muffle furnace at 800-900°C for 20 minutes, cool in a desiccator for 30 minutes and accurately weigh.
8. Calculate the sulfate content as the specific sulfate designated in the specification for the product:

$$\% \text{ Specified Sulfate} = \frac{\text{net wt. of ppt.} \times \text{Factor}}{\text{sample weight}}$$

Table of Sulfate Factors

<u>Specified Sulfates</u>	<u>BaSO₄ Factor</u>
CaSO ₄	58.3
Fe ₂ (SO ₄) ₃	57.1
H ₂ SO ₄	42.0
K ₂ SO ₄	74.7
MgSO ₄	51.6
MnSO ₄	64.7
Na ₂ SO ₄	60.9

PRECISION AND ACCURACY

Duplicate determinations should agree to within 0.002%.

SAFETY

Exercise due care in the use of the high temperature muffle furnace. Handle crucible with long muffle tongs, preferably held with an asbestos glove.

/jw
7/22/71

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.248
May 31, 1971
Page 1 of 2

MATERIAL

Ammonia Persulfate

ANALYSIS REQUIRED

Active Oxygen

APPARATUS

- a) Analytical balance
- b) 250 ml Erlenmeyer flask
- c) Titrating burette
- d) Usual laboratory equipment

REAGENTS

- a) 1N H_2SO_4
- b) 0.5N ferrous ammonium sulfate solution
- c) 0.5N $KMnO_4$

PROCEDURE

The active oxygen content of the persulfates can be conveniently determined by the following method:

Weigh approximately 1 gram of the sample to the nearest milligram and dissolve it in 50 ml of approximately 1 normal H_2SO_4 in an Erlenmeyer flask. Add 40 ml of ferrous ammonium sulfate solution (0.5N). Swirl constantly while adding the iron solution. Let stand for one minute and titrate with 0.5N $KMnO_4$. Run a blank titration on 40 ml of ferrous ammonium sulfate solution, as used above, in 50 ml of the 1 normal H_2SO_4 .

RSV 0014999

CALCULATIONS

$$\text{Percent active oxygen} = \frac{(A-B)C \times 0.8}{D}$$

$$\text{Percent Ammonium Persulfate} = \frac{(A-B)C \times 11.4}{D}$$

A = ml KmnO_4 solution used for titrating the blank

B = ml KmnO_4 solution used for titrating sample

C = Normality of the KmnO_4 solution used

D = Weight of the sample in grams

/jw
7/22/71

RSV 0015000

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.249
May 31, 1971
Page 1 of 1

MATERIAL

Ferric Ammonium Sulfate
Sodium Formaldehyde Sulfoxylate
Tetrasodium ethylene diamine tetraacetate

ANALYSIS REQUIRED

Odor

SCOPE

This method is intended for the characterization of the odors of materials by olfactory perception.

DETERMINATION

1. Remove sufficient material to allow about two inches of head room in the sample bottle above the contents to permit adequate trapping of the odor evolved.
2. If the nature of the sample permits, agitate the contents of the bottle.
3. Let the bottle stand closed at least 30 minutes before performing the odor evaluation.
4. If unfamiliar with the odor of the product being tested, carefully open a previously approved sample and cautiously smell the air in the neck of the bottle.
5. Describe the sample odor with respect to the specification requirement.

/jw
7/22/71

RSV 0015001

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.250
May 31, 1971
Page 1 of 4

MATERIAL

Sodium Lauryl Sulfate
Tetrasodium ethylene diamine tetraacetate

ANALYSIS REQUIRED

pH

SCOPE

This method is intended to be a guide to proper techniques in determining the pH of aqueous solutions by pH meter with the glass electrode. The manufacturer's instruction manual should be consulted for proper operation of the pH meter used.

EQUIPMENT

1. pH meter, such as the Leeds and Northrup pH meter 7664 - Leeds and Northrup Co., 4901 Stention Avenue, Philadelphia 44, Pa., or the Beckman Model H2 pH meter, Scientific Instruments Division, Beckman Instruments, Inc., Fullerton, California.
2. Suitable calomel reference electrode.
3. Suitable glass electrode.
4. Brushless stirring motor such as Fultork Labmotor 14-502, Fisher Scientific Co.

STANDARDIZATION

1. Ascertain that the meter is connected to the proper line voltage, the electrodes are properly plugged into the proper electrode connections, and that the meter is set for reading pH.

RSV 0015002

STANDARDIZATION

2. With the meter in a "Zero" or "Neutral" position, remove the water bath in which the electrodes are immersed and check visually to see that they are clean and unbroken and that the salt bridge is filled. Oils, etc. that adhere to the electrodes may be removed with a suitable solvent and the electrodes rinsed with water.
3. Immerse the electrodes in a suitable standard buffer, start the stirrer and check the temperature of the buffer. Adjust the temperature compensator on the meter to this temperature.
4. Switch the meter from "Neutral" to "Read" position and adjust the standardization control so that the meter reads the exact pH of the buffer solution.
5. After two minutes, read the meter again. If it has drifted from the original setting, adjust the meter to the buffer pH value with fresh buffer solution and once more allow two minutes to check for drift. Continue until a stable reading is obtained.
6. Switch the meter from "Read" to "Neutral" position and turn off the stirrer. Remove the buffer solution from the electrodes, rinse them with water, and proceed to measurement of the sample solutions.

pH MEASUREMENT

1. After standardization of the meter, rinse the electrodes with a portion of the sample solution and immerse the electrodes in a fresh portion of the sample solution.
2. Turn on the stirrer and allow to stir several minutes; check the sample temperature and adjust the temperature compensator if necessary. Then switch the meter from "Neutral" to "Read" position.
3. Observe the pH reading. Wait two minutes and read again. If the second reading does not vary more than ± 0.05 pH units from the first reading, the pH value is acceptable.

pH MEASUREMENT (cont'd)

4. If the second reading varies more than ± 0.05 pH units from the first reading, shut off the stirrer, switch the meter to "Neutral" and repeat steps 2 and 3 with a fresh portion of the sample solution.
5. Continue as in steps 2, 3 and 4 until a reading is obtained that does not vary more than ± 0.05 pH units within two minutes.
6. Turn off stirrer, switch meter to "Neutral" position, remove sample solution from electrodes and clean them with rinses of water or a solvent that will effectively clean them.
7. Leave the electrodes immersed in distilled water.

DISCUSSION

The meter should be standardized with a buffer, the pH of which is within two pH units of the pH of the sample solution to be measured. Glass electrodes will occasionally lose response, so it is good practice to check the pH reading of a second buffer solution, bracketing the anticipated pH of the sample solution. If the buffers are known to be accurate, and the meter does not accurately relate them, then either the meter or the electrode system is not functioning properly, and an accurate pH value cannot be obtained.

It should be noted that pH meters are often used with non-aqueous systems that remove water from the bridge in the calomel reference cell making a high resistance junction with resultant change in reference potential and instability. The pH response of the glass electrode can also be impaired by dehydration by non-aqueous solutions. The electrodes may be restored by soaking in distilled water.

SAFETY

Line-operated pH meters should be grounded.

REFERENCES

1. Willard, Merritt, and Dean, "Instrumental Methods of Analysis", 3rd. ed., 1958, D. Van Nostrand Co., New York, N.Y.
2. Manov, G. G., and Acree, S. F., "The pH of Some Standard Buffer Solutions from 0 to 60°C, and the Calibration of Glass Electrode pH Meters", ASTM Bulletin No. 137, Dec. 1945.
3. Pharmacopeia of the United States XV, "Potentiometric Determination of pH", 1. 933.

/jw
7/21/71

RSV 0015005

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.251
July 13, 1971
Page 1 of 2

MATERIAL

Reactor Overhead Gases

ANALYSIS REQUIRED

Vinyl chloride

APPARATUS

Perkin-Elmer Vapor Fractometer, Model 154B or equivalent
Gas charging valve with 5 cc sample loop
1/4" O.D. stainless steel tubing

PROCEDURE

A. Column Preparation

1. Pack 16 ft of 1/4" stainless tubing with 33% dimethyl sulfolane on 45-60 mesh chromosorb.
2. Install column in the Fractometer. Adjust the temperature to 30°C and helium flow to 50 ml/min. Determine the optimum detector current and set at this value. Allow instrument to reach equilibrium at these conditions.

B. Analysis

1. Turn recorder on and set recorder range to 1. Recorder trace should show no drift for five minutes before Fractometer is ready for use.
2. Purge 5 cc sample loop with sample to be analyzed; shut off valve from sample container, and enter 5cc sample into Fractometer.
3. Record peaks for ethylene and vinyl chloride.

RSV 0015006

CALCULATIONS

1. Draw base lines for both peaks.
2. Use rule with 50 divisions per inch and measure peak height and width at one-half peak height for all peaks. Measurements should be made to the nearest one-tenth of a division.

3. Calculate areas for both peaks:

$$A = H \times w$$

4. Calculate mole percent VCM:

$$\text{mol \%} = \frac{\text{Area VCM}}{\text{Area VCM} + \text{Area Ethylene}}$$

/jw
7/22/71

RSV 0015007

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.252
July 13, 1971
Page 1 of 2

MATERIAL

Caustic

ANALYSIS REQUIRED

Assay

APPARATUS

- a) Hydrometers capable of reading gravity to 0.001 ranging from 1.010 to 1.525.
- b) Constant temperature both @ 20°C.

DETERMINATION

- a) Cool 250 ml of sample to 1-2°C below the specified temperature and fill a hydrometer jar to within 2-3 inches of the top.
- b) Place the specified hydrometer in the liquid and quickly adjust the sample temperature to the specified temperature.
- c) Allow the hydrometer to float to rest.
- d) Read the bottom of the meniscus and record.

CALCULATION

Convert specific gravity to percent NaOH by following graph.

/jw
7/22/71


RSV 0015008

METHOD 602.252

METHOD 602.252

7/13/71

Page 2 of 2



MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.253
July 13, 1971
Page 1 of 2

MATERIAL

Acrylamide - Soap Solution

ANALYSIS REQUIRED

% Acrylamide

APPARATUS

- a) Erlenmeyer flask, 250 ml heavy wall, fitted with a 24/40 Standard Taper joint (Ace Glass Co., No. 2671).
- b) Filling funnel, 50 ml cylindrical, open top with Standard Taper stopcock and 24/40 Standard Taper joint on bottom (Kontes Glass Co., No. K-63125). In the top of the funnel insert a rubber stopper fitted with a short glass tube for connection to vacuum line.

REAGENTS

- a) 0.1 N KBrO_3 - KBr solution (2.79 g potassium bromate and 10 g of potassium bromide per liter)
- b) 6 N H_2SO_4
- c) 20% KI solution
- d) Std. 0.1 N $\text{Na}_2\text{S}_2\text{O}_3$
- e) 1% Starch solution

PROCEDURE

Weigh accurately a sample containing 3-4 g Stream #4 of acrylamide into a 500 ml volumetric flask, dilute to volume with water and mix. Pipet a 25.00 ml aliquot of this solution and 10 ml of bromate-bromide reagent into the 250 ml iodine. Draw 5 ml of 6N sulfuric acid into the flask and seal with water. Allow the solution to stand in the dark for 20 min. and shake frequently. Add 15 ml of potassium iodide solution, mix and then break the vacuum. Titrate the solution immediately with standard thiosulfate using starch indicator. Run a blank determination.

RSV 0015010

METHOD 602.253

July 13, 1971

Page 2 of 2

CALCULATION

$$\% \text{ Acrylamide} = \frac{(\text{blank titer-sample titer}) (\text{N Na}_2\text{S}_2\text{O}_3) (3.554)}{\text{Weight of sample in aliquot}}$$

/jw
7/22/71

RSV 0015011

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.254
July 13, 1971
Page 1 of 2

MATERIAL

Acrylamide - Soap Solution

ANALYSIS REQUIRED

% Alkyl Sulfate

REAGENTS

p-toluidine hydrochloride solution 6.8%
Formula 30 Alcohol
Carbon tetrachloride (reagent grade)
Meta-cresol purple indicator (0.1% solution)
Phenol red indicator (0.1% solution)
0.1 N NaOH

APPARATUS

150 ml beaker
250 ml volumetric flask
25 ml pipette
500 ml separatory funnel
Analytical balance

SCOPE

This procedure employs the extraction of the organic sulfates by p-toluidine hydrochloride precipitation.

PROCEDURE

Using an analytical balance, transfer 25 g of Stream #4 sample to a 150 ml beaker. Dissolve in distilled water and transfer to a 250 ml volumetric flask (a small amount of F-30 Alcohol may be used to dissipate foam) and dilute to 250 ml mark. Mix well and withdraw a 25 ml aliquot portion (equivalent to 0.5 g of sample).

RSV 0015012

METHOD 602.254

July 13, 1971

Page 2 of 2

PROCEDURE (cont'd)

Transfer the aliquot portion to a 500 ml separatory funnel. Add 2 drops of phenol red indicator (red=alkaline, yellow=acid). If solution is alkaline, adjust to neutral with .1 N HCl. If solution is acid, adjust to neutral with .1 N NaOH. To neutral solution add 30 ml of Formula 30 Alcohol, 30 ml of 6.8% p-toluidine hydrochloride solution and 30 ml of carbon tetrachloride. Shake mixture for three minutes (occasionally releasing gas from separatory funnel). Allow to stand and separate. Cut lower (carbon tetrachloride) layer into 250 ml containing 50 ml of Formula 30 Alcohol and 3 cc H₂O adjusted to the alkaline (purple) end-point of meta cresol purple. Add 30 ml carbon tetrachloride to contents of separatory funnel and repeat extraction. Allow to stand and separate and again cut lower layer into same receiving flask.

Repeat extraction once more. These combined extractions in the alcohol medium contain all the organic sulfates and are titrated directly with standard 0.1 N NaOH to determine the alkyl sulfate content. The end point is taken as the first change in color from gray to purple that remains for 15 seconds.

CALCULATION

$$\frac{\text{ml NaOH} \times \text{N} \times 300 \times 100}{\text{Wt. of sample in aliquot} \times 1000} = \% \text{ Alkyl Sulfate as sodium lauryl sulfate}$$

/jw
7/22/71

RSV 0015013

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.255
July 13, 1971
Page 1 of 1

MATERIAL

Product Latex B, BH and E
Spent THF Solvent

ANALYSIS REQUIRED

Percent Total Solids

APPARATUS

Analytical balance capable of weighing to nearest 0.1 mg.
Aluminum weighing dish
Oven
Dessicator
200 mesh sample screen

PROCEDURE

Screen latex sample through 200 mesh screen. Tare duplicate aluminum weighing dishes to 0.1 mg. Transfer 3-4 g of sample to each dish and obtain gross weight to 0.1 mg. Place the dishes in an oven at 130°C for 3 hours. Cool dish in a dessicator and weigh to 0.1 mg. Report % solids. If two samples differ by more than 0.3%, repeat.

CALCULATIONS

$$\% \text{ Solids} = \frac{\text{net weight after drying}}{\text{net weight before drying}} \times 100$$

/jw
7/22/71

RSV 0015014

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.256
July 13, 1971
Page 1 of 1

MATERIAL

Product Latex B, BH and E

ANALYSIS REQUIRED

Latex Viscosity

APPARATUS

Brookfield viscometer, Model LVT with #2 spindle or Model
RVT with #2 spindle.
400 ml beaker
200 mesh sample screen (Stainless Steel)

PROCEDURE

Screen latex sample through 200 mesh screen. Sample temperature must be $25^{\circ} \pm 1^{\circ}\text{C}$. Pour 300 ml sample into a 400 ml beaker. Place a #2 spindle in the Brookfield viscometer model LVT and lower spindle into the sample to the level indicated by the indentation in the spindle shaft. Set the RPM selector at 6 and turn on the viscometer. After the spindle has completed 3 revolutions, depress the clutch and record the viscometer reading. While the clutch is still depressed turn the RPM selector to 60, turn on the viscometer, release the clutch and allow the spindle to complete 10 revolutions. Depress the clutch and record the viscometer reading.

CALCULATIONS

Convert the readings to centipoise.

Reading at 6 RPM x 50 = centipoise (LVT)

Reading at 60 RPM x 5 = centipoise (LVT)

Reading at 100 RPM (RVT)

/jw
7/22/71

RSV 0015015

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.257
July 13, 1971
Page 1 of 4

MATERIAL

Product Latex B, BH, and E

ANALYSIS REQUIRED

pH

APPARATUS

- a) pH Meter, such as the Leeds and Northrup pH meter 7664 - Leeds and Northrup Co., 4901 Stenton Avenue, Philadelphia 44, Pa., or the Beckman Model H2 pH meter, Scientific Instruments Division, Beckman Instruments, Inc., Fullerton, California.
- b) Suitable calomel reference electrode.
- c) Suitable glass electrode.
- d) Brushless stirring motor such as Fultork Labmotor 14-502, Fisher Scientific Co.

STANDARDIZATION

- a) Ascertain that the meter is connected to the proper line voltage, the electrodes are properly plugged into the proper electrode connections, and that the meter is set for reading pH.
- b) With the meter in a "Zero" or "Neutral" position, remove the water bath in which the electrodes are immersed and check visually to see that they are clean and unbroken and that the salt bridge is filled. Oils, etc. that adhere to the electrodes may be removed with a suitable solvent and the electrodes rinsed with water.
- c) Immerse the electrodes in a suitable standard buffer, start the stirrer and check the temperature of the buffer. Adjust the temperature compensator on the meter to this temperature.

RSV 0015016

STANDARDIZATION (cont'd)

- d) Switch the meter from "Neutral" to "Read" position and adjust the standardization control so that the meter reads the exact pH of the buffer solution.
- e) After two minutes, read the meter again. If it has drifted from the original setting, adjust the meter to the buffer pH value with fresh buffer solution and once more allow two minutes to check for drift. Continue until a stable reading is obtained.
- f) Switch the meter from "Read" to "Neutral" position and turn off the stirrer. Remove the buffer solution from the electrodes, rinse them with water, and proceed to measurement of the sample solutions.

pH MEASUREMENT

- a) Screen the latex sample through a 200 mesh stainless steel screen. After standardization of the meter, rinse the electrodes with a portion of the sample solution and immerse the electrodes in a fresh portion of the sample solution.
- b) Turn on the stirrer and allow to stir several minutes; check the sample temperature and adjust the temperature compensator if necessary. Then switch the meter from "Neutral" to "Read" position.
- c) Observe the pH reading. Wait two minutes and read again. If the second reading does not vary more than ± 0.05 pH units from the first reading, the pH value is acceptable.
- d) If the second reading varies more than ± 0.05 pH units from the first reading, shut off the stirrer, switch the meter to "Neutral" and repeat steps b and c with a fresh portion of the sample solution.
- e) Continue as in steps b, c and d until a reading is obtained that does not vary more than ± 0.05 pH units within two minutes.

pH MEASUREMENT (cont'd)

- f) Turn off stirrer, switch meter to "Neutral" position, remove sample solution from electrodes and clean them with rinses of water or a solvent that will effectively clean them.
- g) Leave the electrodes immersed in distilled water.

SCOPE

This method is intended to be a guide to proper techniques in determining the pH of aqueous solutions by pH meter with the glass electrode. The manufacturer's instruction manual should be consulted for proper operation of the pH meter used.

DISCUSSION

The meter should be standardized with a buffer, the pH of which is within two pH units of the pH of the sample solution to be measured. Glass electrodes will occasionally lose response, so it is good practice to check the pH reading of a second buffer solution, bracketing the anticipated pH of the sample solution. If the buffers are known to be accurate, and the meter does not accurately relate them, then either the meter or the electrode system is not functioning properly, and an accurate pH value cannot be obtained.

It should be noted that pH meters are often used with non-aqueous systems that remove water from the bridge in the calomel reference cell making a high resistance junction with resultant change in reference potential and instability. The pH response of the glass electrode can also be impaired by dehydration by non-aqueous solutions. The electrodes may be restored by soaking in distilled water.

SAFETY

Line-operated pH meters should be grounded.

REFERENCES

- a) Willard, Merritt, and Dean, "Instrumental Methods of Analysis", 3rd ed., 1958, D. Van Nostrand Co., New York, N.Y.
- b) Manov, G. G., and Acree, S. F., "The pH of Some Standard Buffer Solutions from 0 to 60°C, and the Calibration of Glass Electrode pH Meters", ASTM Bulletin No. 137, Dec. 1945.
- c) Pharmacopeia of the United States XV, "Potentiometric determination of pH", 1. 933.

/jw
7/22/71

RSV 0015019

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.258
July 13, 1971
Page 1 of 1

MATERIAL

Product Latex B, BH, and E

ANALYSIS REQUIRED

% Grit

APPARATUS

325 mesh stainless steel screen
Analytical balance capable of weighing to nearest 0.1 mg
Oven
Evaporating dish

PROCEDURE

Take 500 g of latex and filter through a 325 mesh screen. Wash the latex on the screen with a gentle stream of distilled or deionized water until no further latex can be washed through. Collect the retained grit by washing onto a weighed evaporating dish. Dry the grit at 120°C and determine the weight.

CALCULATIONS

$$\% \text{ grit} = \frac{\text{wt. grit}}{500 \text{ g}} \times 100$$

/jw
7/22/71

RSV 0015020

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.259
July 13, 1971
Page 1 of 5

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

Surface Tension

APPARATUS

Fisher tensiometer
Shallow sample container
200 mesh stainless steel screen

PROCEDURE

Screen sample through 200 mesh screen.
In making a determination of surface tension, careful preparation of the sample and the Surface Tensiometer must precede the actual manipulation of the instrument.

The sample should be placed in a glass beaker or cylindrical vessel with a diameter of at least 45 millimeters. For testing oil samples according to ASTM Method D-971, the glassware should be cleaned according to a definite procedure. Any residual oil from the previous sample is removed with petroleum naphtha or benzene followed by several washes with methyl ethyl ketone and water, then the glassware is immersed in a hot cleaning solution of chromic acid. The glassware also should be rinsed thoroughly with tap water, then with distilled water. It should be drained in an inverted position over a clean cloth, unless it is to be used immediately.

The platinum-iridium ring should be cleaned by rinsing it in petroleum naphtha or benzene, then by rinsing in methyl ethyl ketone. The ring should then be heated in the oxidizing portion of a gas flame.

RSV 0015021

PROCEDURE (cont'd)

The preceding instructions apply particularly to preparation for determinations in oil samples. A comparable degree of cleanliness should be maintained for other samples.

Since surface tension is dependent upon temperature, consideration must be given to this factor. For theoretical work, temperatures must be specified; however, 25°C is the temperature most commonly used. For control work, the operator should use the same temperature for each type of measurement.

Since the Surface Tensiometer may be employed as a manual or semi-automatic instrument, separate procedures are described below for both operating modes.

MANUAL DETERMINATIONSMeasuring Surface Tension

The cleaned platinum-iridium ring should first be attached to the hook at the end of the lever arm. The arrest mechanism should be holding the arm at this time.

The liquid to be measured is transferred to the clean glass vessel and placed on the sample table. The sample table is moved around until it is directly beneath the platinum-iridium ring. Raise the sample table until the ring is immersed in the test liquid. The ring should be in the liquid, beneath the surface so that the entire ring will be wetted. About 1/8 inch immersion is generally considered sufficient.

The torsion arm is now released and the instrument adjusted to a zero reading. Adjust the knob on the right side of the case until the index and its image are exactly in line with the reference mark on the mirror. Be careful to keep the ring in the liquid during this manipulation, raising or lowering the sample table if necessary by means of the knob adjustment underneath of the table. Now turn the knob beneath the main dial on the front of the case until the vernier reads zero on the outer scale of the dial.

MANUAL DETERMINATIONS (cont'd)

Lower the sample table until the ring is in the surface of the liquid, while at the same time adjusting the knob on the right side of the case to keep the index lined up with the reference mark on the mirror. The surface of the liquid will become distended but the index must be kept on the reference. Continue the two simultaneous adjustments until the distended film at the surface of the liquid breaks. The scale reading at the breaking point of the distended film is the apparent surface tension.

This reading is the one obtained after the sample has been adjusted for cross hairs mesh. Make only one repeat determination after this adjustment. Repeat determinations without cleaning the loop will give higher and higher readings due to the latex coating the wire.

SEMI-AUTOMATIC DETERMINATIONS

The semi-automatic operation of the FISHER Surface Tensiometer is basically similar to the manual mode of operation - the only essential difference being the use of small electric motor to twist the torsion wire and thus raise or lower the ring in the test solutions. No adjustments or complicated procedures are involved in changing from either mode of operation to the other, the use of a switch to activate the semi-automatic feature and a slightly different technique in regard to the level of the surface or interface are the sole differences. The motor is connected to the semi-automatic mechanism by means of a clutch which does not interfere with the manual mode.

The semi-automatic drive mechanism actuates the torsion arm and consequently the ring; the position of the sample table is not changed or controlled by the drive mechanism. As a consequence, the simultaneous adjustment of torsion arm and sample table is not feasible. Therefore, it is necessary to determine the level of the surface or interface at the instant the film is broken by means of a manual determination and to maintain this level of surface or interface at the instant the film is broken by means of a

SEMI-AUTOMATIC DETERMINATIONS (cont'd)

manual determination and to maintain this level of surface or interface for subsequent determinations of similar samples. Provided the subsequent samples are of the same character (approximately the same density and surface or interfacial tension) no appreciable error is introduced by presetting the level of the surface or interface. Since a major advantage of the Surface Tensiomat is its ability to make repetitive determinations on similar samples rapidly, this determination of a level is a simple and worthwhile operation.

A practical method for attaining a reproducible level for each sample is to use a sample container with a mark showing the level. More than one sample container can be so marked and used in succession. Also, it is usually practical to introduce reproducible volumes of sample into containers of reasonably uniform size. For example, Petri dishes with identical volumes of sample will produce levels which are for all practical purposes the same.

Once the level of the surface has been determined and the approximate reading for the type of sample ascertained, the semi-automatic mode of operation is very simple. The sample is placed on the sample table and the ring immersed in the proper liquid. The knob on the right side of the case is put in the UP position for surface tension measurements or interfacial tension measurements where the ring is pulled up through the interface, and in the DOWN position for interfacial tension determinations where the ring is forced down through the interface. When the surface or interface film ruptures, the level arm actuates a contact (No. 5 of Figure 1) which in turn actuates a relay and stops the drive motor. The value is then read from the proper scale on the dial.

A pilot light on the front of the case indicates when the Surface Tensiomat is running. When the film ruptures, the motor stops and the light goes out. The switch must be returned to the NEUTRAL position and the lever arm returned to the proper position prior to making the next determination.

METHOD 602.259
July 13, 1971
Page 5 of 5

SEMI-AUTOMATIC DETERMINATIONS (cont'd)

Experience has shown that readings obtained using semi-automatic operation will be slightly different from those obtained manually. However, with any given set of conditions, these variations can be determined by manual and semi-automatic experiments, and if necessary, the variations may then be taken into account by the user.

/jw
7/22/71

RSV 0015025

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.260
July 13, 1971
Page 1 of 5

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

Glass Transition Temperature (T_g)

APPARATUS AND MATERIALS

duPont 900 Differential Thermal Analyzer or its equivalent
Vacuum Pump
Liquid nitrogen and reservoir
Dry nitrogen gas
Bunsen burner

TEST SPECIMEN

The test specimens are prepared by air drying a film of latex on a glass plate.

PROCEDURE

A condensed operating procedure is described here. For detailed operating instructions, refer to duPont 900 DTA instruction manual.

a) Preparation for the Run

- 1) Cold junction thermocouples should be at room temperature.
- 2) Fill macro sample tube with the sample to depth of 5 mm. Brush the loaded sample tube through a small Bunsen flame two or three times, and pack the sample tightly by pushing the centering ceramic sleeve to the bottom of the tube. A thermocouple is then placed in the sample so that the junction touches the bottom of the tube. The ceramic sleeve should be positioned over the thermocouple to center it in the macro tube.

RSV 0015026

PROCEDURE (cont'd)

- 3) Fill control and reference tubes with quartz powder to depth of 5 mm. Insert thermocouples and centering sleeves in tubes.
 - 4) Wrap glass wool around thermocouple connections at the receptacle in cell assembly in order to shield from drafts during cooling.
 - 5) Place tubes in heating block, bend leads so junctions bottom in tubes, and put cell assembly in place.
 - 6) Place chart paper on recorder.
- b) Set the Recorder Controls
- 1) Switch POWER to RECORD (5 minutes warm-up).
 - 2) Set T ZERO SHIFT at +5.
 - 3) Set T SCALE at 10°C/in. (Scale on chart: -50°C to + 50°C)
 - 4) Turn Δ T ZERO SHIFT until pen is near top of chart.
 - 5) Set Δ T SCALE at 0.2°C/in.
 - 6) Set BASELINE SLOPE to 0.
- c) Set Sample Atmosphere Controls
- 1) Connect nitrogen gas supply to PURGE, liquid nitrogen to COOL and vacuum to VACUUM connections.
 - 2) Open VACUUM and COOL valves fully, this introduces liquid nitrogen into the cell.
 - 3) When pen reaches chart reading of -20°C, close COOL valve and adjust the VACUUM gauge to read approximately 15. Set GAS FLOW of purge nitrogen at 2 SCFH.

PROCEDURE (cont'd)

d) Set Temperature Program Controls

- 1) Push RESET.
- 2) Set TEMPERATURE RATE at 20°C/min.
- 3) Turn STARTING TEMPERATURE clockwise until heater voltage indicator reads minimum value. Then turn counterclockwise until voltage indicator just begins to move upscale.

e) Start the Run

- 1) Set PROGRAM MODE at HEAT.
- 2) Set POWER at RECORD.
- 3) The thermogram will not be recorded on the chart until the pen reaches on the right hand edge of the paper. The ΔT SCALE or ΔT ZERO SHIFT controls may be adjusted at any time during the run to keep strong endotherms or exotherms on the chart paper. Similarly, the BASELINE SLOPE control may be adjusted at any time during the run to make the baseline more nearly horizontal.


f) Conclude the Run

- 1) Turn PROGRAM MODE to HEATER OFF.
- 2) Turn POWER to STAND BY.
- 3) After each run the sample tube may be discarded. The sample thermocouple and the centering ceramic sleeve should be cleaned by burning over a small Bunsen flame to burn all organic compounds, they can be reused for twenty or more runs.

PROCEDURE (cont'd)

g) Definition of Transition Temperature

The literature of DTA contains many definitions for the location of significant transition temperatures. This procedure defines the transition temperature as a temperature of extrapolated onset as shown in the figure.



Unlike first-order thermodynamic transitions, the glass transitions of amorphous polymers may not exhibit a peak but only show a shift in baseline. Similarly, the T_g is defined as the temperature of the extrapolated onset in the thermogram.

h) Report

The report shall include the following:

- 1) Complete identification of the material tested.
- 2) Size of sample tube.
- 3) Reference material.
- 4) Heating, rate, $^{\circ}\text{C}/\text{min}$.
- 5) Nitrogen purge rate, SCFH.
- 6) Temperature scale and shift factor.
- 7) ΔT of sensitivity, $^{\circ}\text{C}/\text{in}$.

SCOPE

This method describes a test procedure for the determination of Tg of the polymer using a differential thermal analyzer.

SIGNIFICANCE

Differential thermal analysis (DTA) is a technique for studying the thermal behavior of materials as they undergo physical and chemical changes during heating or cooling. The glass transition temperature (Tg) of amorphous polymer is only one of many thermal fingerprints of polymer obtainable with DTA. Determination of Tg, which is one of the most important fundamental properties of polymers, by DTA represents one of the simplest and fastest methods.

/jw
7/23/71

RSV 0015030

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.261
July 13, 1971
Page 1 of 4

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

Intrinsic Viscosity

APPARATUS

1. Viscometers:
 - a. Cannon-Ubbelohde Dilution Type, Size 50, uncalibrated
 - b. Cannon-Fenske #50
2. Constant temperature water bath, $\pm 0.05^{\circ}\text{C}$
3. Volumetric pipettes: 10 ml, 20 ml, 100 ml
4. Aluminum weighing dishes
5. Forced air oven
6. Stopwatch
7. 1 Liter volumetric flask
8. Tetrahydrofuran (THF), Fisher T-397)

PROCEDURE

1. Measure the solids content of the latex sample and dilute it to 45% solids.
2. Weigh 2 grams of the 45% solids latex into a 200 ml jar.
3. Add 100 ml of pure THF (Tetrahydrofuran) to the jar, cap tightly, and shake. Heat the sample at $40-50^{\circ}\text{C}$ for 1 hour with intermittent shaking. The solution may still appear hazy.

RSV 0015031

PROCEDURE (cont'd)

4. Prepare the dilution solvent by adding 10.9 ml H₂O to a 1000 ml volumetric flask and making up to the mark with pure THF.
5. Obtain tare weights on two aluminum weighing dishes. Pipette 10 ml of the polymer solution into each dish and dry to constant weight at 40-50°C. This will take approximately one hour. Record the weight of residue.
6. Determine the concentration of the polymer solution.

Wt. residue x 10 equals the concentration expressed as grams per deciliter. The initial concentration should be just below 1 gm/deciliter.
7. With a small size coarse glass filter adapted to a 10 ml pipette, draw up 10 ml of the dilution solvent and deliver to a viscometer. The viscometer water bath should be maintained at 30°C.
8. Record the flux time of the dilution solvent. A viscometer should be chosen which gives a 100-150 sec. flux time for the dilution solvent. Make three determinations and record the average.
9. Replace the solvent in the viscometer with the polymer solution as in No. 7. Determine the flux time at this concentration.
- 10a. If a dilution type viscometer is used, simply add 10 ml of the dilution solvent to the polymer solution already in the viscometer and determine the flux time for this solution. Calculate the actual concentration of the solution.
- b. Make a second dilution by adding 20 ml of dilution solvent to the viscometer after Step 10a. Determine the flux time for this solution and calculate the exact concentration.
11. If a dilution type viscometer is not available, the two additional solutions can be made up separately by adding 10 ml and 30 ml of dilution solvent to each of two 10 ml samples of the original polymer solution.

PROCEDURE (cont'd)

12. For each of the three solution concentrations, calculate the reduced viscosity according to the following formula:


$$\text{Reduced viscosity} = \frac{\left(\frac{\text{Flux time of solution}}{\text{Flux time of solvent}} - 1 \right)}{\text{Actual concentration in gm/deciliter}}$$

13. Plot the reduced viscosity versus the actual concentration. The data should fall on a straight line which is extended to zero concentration.

The intrinsic viscosity is given by the value of the reduced viscosity at zero concentration and will generally be in the range of 0.6-0.8 deciliters/gm.
(See attached Figure.)

/jw
7/23/71

RSV 0015033



RSV 0015034

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.262
July 13, 1971
Page 1 of 5

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

Particle Size

APPARATUS

Amberlite XAD-2 resin
Analytical balance
Tumbler
Buchner funnel, monofilament nylon filter cloth
Aluminum weighing dish
Fisher Tensiomat

REAGENTS

Sodium lauryl sulfate (sapon WD)

PROCEDURE

- A. Wash Amberlite XAD-2 absorbent with distilled water. Remove excess water from washed Amberlite and store Amberlite in closed jar to keep moist.
- B. Weigh 150 g of screened (200 mesh) latex (40-50% T.S.) and 300 g wet Amberlite into 16 oz. wide mouth bottle.
- C. Cap bottle, shake to mix well, and keep mixed by slow tumbling motion for one hour.
- D. Filter out latex from mixture with Buchner funnel and monofilament nylon filter cloth under suction.
- E. Determine surface tension of latex. Repeat Steps B-D if surface tension is less than 70 dynes/cm.

RSV 0015035

PROCEDURE (cont'd)

- F. Determine total solids of the Amberlite treated latex stock:

Weigh out 1-2 g of stock in aluminum dish. Add enough water to thin out latex and barely cover dish bottom. Dry in oven at 120°C to constant weight (approx. 1 hour). Calculate % T.S. from dry weight.

- G. Prepare two latex samples in 4 oz. bottles with following compositions:

Sample 1: 60 g Amberlite treated latex stock
11 g 10% SLS (Sipon WD)
39 g Distilled water

110 g latex with 1% SLS

Sample 2: 60 g Amberlite treated latex stock
50 g Distilled water

110 g latex without SLS

Shake to mix well.

- H. Mix Samples 1 and 2 in 1 oz. bottles at various proportions as shown in Column I of Table I. Shake to mix well. The corresponding SLS concentration of each sample is shown in Column II of Table I.
- I. Determine surface tension of the 11 samples. A hypothetical set of surface tension values is given in Column III of Table I for illustrative purposes.
- J. Plot surface tension vs. SLS concentration of sample in semi-log graph paper as shown in Fig. 1. The resulting graph will typically show a break which separates the two rather distinct portions of the graph, i.e., the initial portion where the surface tension decreases rapidly with increasing SLS concentration, and the latter portion where surface tension practically remains at a constant low value.

PROCEDURE (cont'd)

- K. Compute the particle size using the following equation:

$$\text{Particle size } (\text{\AA}) = 0.24 \frac{\% \text{ T.S. of latex stock} \times (60/110)}{\text{Titration end point in } \frac{\text{g SLS}}{\text{g latex}}}$$

- L. When particle size is much larger than 1000 $\text{\AA}$ or much smaller than 500 $\text{\AA}$, the titration end point will move too close to the low or high end of the covered SLS concentration range, making a precise determination of the end point difficult.

When the end point is too low ($<0.002 \frac{\text{g SLS}}{\text{g latex}}$), repeat the determination using the same procedure except change the composition of Sample I of Step G as follows:

60.0 g Amberlite treated latex stock
5.5 g 10% SLS
<u>44.5 g Distilled water</u>
100.0 g latex with 0.5% SLS

and correct the computed particle size of Step K by a factor of 2.

When the end point is too high ($>0.008 \frac{\text{g SLS}}{\text{g latex}}$), repeat determination. Change composition of Sample I of Step G to:

60 g Amberlite treated latex stock
22 g 10% SLS
<u>28 g Distilled water</u>
110 g latex with 2% SLS

and then correct particle size computed in Step K by a factor of 0.5.

Table I


I	II	III
$\frac{\text{gr Sample I}}{\text{gr Sample II}}$	$\frac{\text{gr SLS}}{\text{gr latex}}$	Surface tension $\neq$ (dynes/cm)
0/20	0	71.5
2/18	0.001	64.8
4/14	0.002	56.7
6/14	0.003	49.6
8/12	0.004	43.5
10/10	0.005	38.2
12/8	0.006	34.5
14/6	0.007	33.3
16/4	0.008	33.1
18/2	0.009	33.1
20/0	0.010	33.0

$\neq$ Hypothetical set of values for illustrative purposes.

/jw
7/23/71

RSV 0015038

K-E SEMI-LOGARITHMIC 48 4873
2 CYCLES X 10 DIVISIONS
RUFFEL & SON CO.



RSV 0015039

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.263
July 17, 1971
Page 1 of 3

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

Instron tests for tensile stress, elongation, and 100% modulus.

APPARATUS

Instron Tester
Glass Plate and Razor

PROCEDURE

1. Draw down a wet film of screened (200 mesh) latex using an 8-9 mil bar.
2. Slowly air dry the film. If drying proceeds too fast, the film will crack and craze excessively. Drying can be slowed down by placing a flat box over the wet film immediately after it is drawn down. Drying will take 2-3 hours.
3. Remove the film from the glass plate and place it onto a Teflon covered plate.


If the film adheres too tightly to the glass plate, it may be necessary to slightly moisten it to remove it. Also, it is usually better to cut small sections of film with a razor and remove these from the glass plate rather than trying to remove the entire film intact.
4. Allow the film to air dry completely, then cure in a hot air circulating oven for 5 minutes at 150°C.
5. Accurately cut the cured film into pieces 1/2" x 2" for testing. Use a sharp razor for cutting the test pieces and avoid any edge defects. Also avoid using sections of the cured film with obvious defects, such as air bubbles or inclusions of foreign particles.

If a specimen cutting die is available as described in ASTM D882-GIT, it can be used.
6. Determine the thickness of each test strip, using any convenient calipering device which can measure and read out to ± 0.05 mil.


RSV 0015040

Procedure (cont'd)

7. Reinforce the ends of each test specimen, using Scotch tape. This allows more secure clamping of the specimen and helps to minimize slippage in the jaws. One-half inch of each end should be taped.



8. Place the specimen in the Instron with a jaw span set at 1" and strain the specimen at a rate of 1"/min. until it fails.
9. Discard samples which show slippage in the jaws or within the taped ends. Also discard the results if an obvious defect develops upon straining.
10. The output curve will look something like this:

CALCULATIONS

- a. Nominal tensile stress at break =
(psi)

$$\frac{\text{Load at break (A)}}{(\text{Sample width in inches}) \times (\text{Initial thickness in inches})}$$

- b. % Elongation at break =

$$\frac{\text{Extension in inches at break (B)}}{\text{Initial sample length (1")}} \times 100$$

Calculations (cont'd)

- c. Initial tensile modulus =
(psi) - from line C

$$\frac{\text{Nominal tensile stress at break (A)}}{\text{Extension E in inches/initial span (1")}}$$

- d. Secant modulus at 100% extension =
(psi) - from line D

$$\frac{\text{Load at elongation of 1" (F)}}{(\text{Sample width in inches}) \times (\text{Initial thickness in inches})}$$

:ff
7/19/71

RSV 0015042

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.264
July 17, 1971
Page 1 of 3

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

% Acrylamide

APPARATUS

Kjeldahl apparatus
2000 ml beaker
Laboratory agitator for 2000 ml beaker
Buchner funnel
Vacuum source
Vacuum desiccator

REAGENTS

Isopropanol

PROCEDURE

Twenty ml latex are slowly added to a stirred beaker containing 1200 ml isopropanol. The coagulated latex and isopropanol are then stirred for 5 minutes. The polymer is filtered out using two isopropanol rinses. The vacuum to the filter is released before the polymer packs into the filter paper. The polymer is shredded up and stirred in 1500 ml of H₂O for 2 hours. It is then filtered using five 300 ml portions of H₂O as rinses. The polymer is dried in a vacuum desiccator for 24 hours at room temperature or until completely dry.

Total nitrogen in the dry polymer is obtained by the standard Kjeldahl method. Total nitrogen is assumed to be present as acrylamide. Therefore,

$$\% \text{ acrylamide} = \frac{(71)}{(14)} (\%N)$$

The total acrylamide put into the latex based on solids content may be calculated.

$$\% \text{ Acryl.} = \frac{(\text{lbs. Termonomer Soln. Added})}{(\text{lbs. Latex Out}) (\% \text{ Solids})} (\% \text{ Acryl. in Soln}) (100) + 0.$$

RSV 0015043

NITROGEN DETERMINATION - KJELDAHLEQUIPMENT

1. Kjeldahl digestion rack - gas fired.
2. Kjeldahl digestion flask - 100 ml capacity.
3. Distillation unit - JFQ Plant standard apparatus drawing 2.
4. Kjeldahl digestion catalyst - thoroughly mix 150 g. of potassium sulfate anhydrous and 10 g. of yellow mercuric oxide powder.
5. Bromcresol green methyl red mixed indicator solution - mix 50 ml of bromcresol green solution (0.1% in 95% ethanol) and 10 ml of methyl red solution (0.1% in 95% ethanol).
6. Zinc metal, 20 mesh.
7. Phenolphthalein solution - 1.0% in Formula 30 alcohol.
8. Boric acid solution - 4% in distilled water.
9. Sodium thiosulfate - 44% aqueous solution.
10. Sodium hydroxide - 50% aqueous solution.
11. Volhard absorbers - JFQ Plant standard apparatus drawing 22.
12. Usual laboratory reagents and equipment.

DETERMINATION

1. Accurately weigh 0.20 to 0.25 g. of sample into a 100 ml Kjeldahl flask, add 15 g. of the digestion catalyst and 20 ml of concentrated sulfuric acid and swirl to mix. Carry through a blank digestion omitting the sample.
2. Add several glass beads or boiling chips and place the flask in a 110°C oven for over night charring.
3. Remove the sample from the oven and digest slowly, gradually increasing the heat until the mixture boils briskly. Continue to heat until the mixture becomes clear.
4. Continue the digestion for an additional hour.

Nitrogen Determination (cont'd)

5. Cool to room temperature and cautiously add approximately 75 ml of distilled water with mixing.
6. Transfer quantitatively with the aid of distilled water to a 500 ml round bottom flask containing 10 g. of zinc metal. Add several drops of phenolphthalein solution and attach the flask to the distillation assembly. Carry through the blank in a similar manner.
7. Add sufficient 4% boric acid solution to a Volhard absorber so that the bottom overflow exit is completely covered and add 4-5 drops of the mixed indicator.
8. Attach the absorber tightly to the condenser exit.
9. Add to the round bottom flask, via the distillation side arm, sufficient 50% sodium hydroxide to make the sample alkaline to phenolphthalein. Add 5 ml of 44% sodium thiosulfate in the same manner.
10. Make sure all connections are tight, apply heat to the flask and distill approximately 100 ml of distillate into the Volhard absorber.
11. At the end of the distillation, disconnect the Volhard absorber and titrate its contents with standardized N/10 HCl solution.
12. Calculate:

$$\frac{(\text{ml titer for blank} - \text{ml titer for sample}) \times 0.14\%}{\text{sample weight}} = \% \text{ Nitrogen}$$

:ff

7/19/71

RSV 0015045

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.266
July 17, 1971
Page 1 of 1

MATERIAL

Product Latex B and E

ANALYSIS REQUIRED

Film Clarity

APPARATUS

Glass plate
Drawdown bar

PROCEDURE

1. Draw down a latex film from screened latex (200 mesh) on a clean glass plate, using an 8 mil applicator bar. The exact thickness of the film is not critical.
2. Allow the wet film to slowly dry at room temperature. The film should be about 3 mils thick. If the film is too thick, it will show excessive crazing upon drying.
3. Observe the film and note the following points:
 - a. Continuity of the film formed.
 - b. Presence of grit.
 - c. Presence of incompatible deposits in the film.
 - d. Excessive roughness or graininess in the film.
 - e. Clarity and color of the film.

:ff
7/19/71

RSV 0015046

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.267
July 17, 1971
Page 1 of 6

MATERIAL

Product Latex BH

ANALYSIS REQUIRED

I.G.T. Test

APPARATUS

Brookfield LVT or RVT viscometer
Holder for wire wound coating rods
Analytical balance
Supply of paper and board substrates
Paper cutter
I.G.T. Printability Tester, consisting of:
 Printmaker equipped with pendulum and one inking disc
 Ink distributing apparatus with a set of two polymethane rollers.
I.G.T. Spring Drive Device to increase speed to the 0-650 ft/min range.
Four double width inking discs (20 mm), aluminum
I.G.T. ink pipette
Extra mouthpiece for I.G.T. ink pipette
Set of two glue-glycerine composition rollers
10 sets of 6-paper packings 25 mm wide for I.G.T. Tester
Mylar film - 1 mil thick (Herculene or Stabiline Drawing film from Keuffel and Esser).
Solvent for clean-up. 2-5 gal cans of Varn Wash V-120 from Vain Products Co., Flushing, N.Y.
Two, 1/2# tubes of #4, #5, #6 black tack graded inks from Interchemical Corp., I.P.I. Printing Ink Division, 4168 Meramec, St. Louis, Mo. 63105
Three, 50# bags of Lustra Coating Day from Freeport Kadin Co., 405 Lexington Ave., N.Y., N.Y 10017.
A binocular microscope with a 15-20X magnification.
One - Series 2000 standard Premier laboratory Dispersator, 1/2 HP single phase, universal type motor, powerstat, complete with stand and mounting bracket, B/M 2053 equipped with 1-1/2" Hi-Vis head and stub shaft assembly, with shaft 3/4" dia. x 11" long constructed of 316 stainless steel. Order from Premier Mill Corp., 224 Fifth Ave., New York, N.Y.
Pyrex radiant panel heater, cat. no. 604001, 250 volts, 3500 watts, Corning Glass Works, Corning, N.Y.

RSV 0015047

Apparatus (cont'd)

- One - ATC percentage input timer, type 304B-07-B-O-XX, 230 volts, 60 cycle, 20 amp., Automatic Timing & Controls, Inc., King of Prussia, Pa. 19406.
- 3 each of the following wire wound drawdown coating rods:
#8, #10, #12, R.D. Specialties Co., P.O. Box 6, Webster, N.Y. 14580.
- One - #162D case for surface mounting of type 304 ATCOTROL manual set percentage timer, general purpose sheet steel (Order with #15 above)

PROCEDURE FOR MEASURING PIGMENT BINDING
STRENGTH OF LATEX - IGT PICK STRENGTH

The procedure can be broken down into four parts: 1) preparation of the coating formulation, 2) coating and drying of the substrate, 3) printing of the samples on the tester, 4) determination of end-point.

1. Preparation of the Coating Formulation

The coating formulation is a simple mixture of latex plus coating clay diluted to a total solids content of 60%. The amount of latex used is such as to give a mixture of 18 parts dry latex solids per 100 parts dry clay solids.

A 70% solids clay slip is prepared as follows:

- a. Add 1.5 grams tetrasodium pyrophosphate (TSPP) to 215 grams of deionized or distilled water and stir until the TSPP is dissolved.
- b. Place the solution under the Premier Dispersator and slowly add in 500 grams of Lustra Coating Clay. It will be necessary to increase the speed of the agitator as the clay is being added.
- c. After all the clay has been added, disperse at the highest speed possible for 10 minutes.
- d. Determine the solids content of the clay slurry.

To make up the coating formulation, add 286 grams of the clay slurry (equal to 200 grams of oven dry clay) to a 250 ml beaker. Add latex to the beaker to give the equivalent of 36 grams of oven dry latex.

Procedure (cont'd)

The amount of latex to add can be found from the following table:

<u>% Solids of Latex</u>	<u>Grams of Latex to Add To Clay Slurry</u>
45	80.0
46	78.3
47	76.7
48	75.0
49 -	73.5
50	72.0

After the latex has been added, add water to make the solids content of the total coating formulation equal to 60%. The amount of water needed will depend on the solids content of the latex, assuming that the clay slurry has been made up properly at 70% solids. The following table can be used as a guide for the water addition:

<u>% Solids of Latex</u>	<u>Grams of Water to Add To Formulation</u>
45	27.0
46	28.7
47	30.3
48	32.0
49	33.5
50	35.0

The coating mixture should be stirred for 10 minutes to insure proper mixing of all components. Filter the coating through a conical 200 mesh stainless steel screen to remove any aggregates of dried material and to remove excess entrained air.

The viscosity of the coating formulation should be measured with a Brookfield viscometer, Model RVT, using a #2 spindle at 10 RPM, room temperature in the 250 ml beaker.

2. Coating and Drying of the Substrate

The substrate to be used is a bleached solid sulfate board (15 point Fold-Brite White Tag Bleached Board from Riegel Paper Corp., Riegelwood, N.C.). The coating is applied to the "wire" side of the sheet. The wire side of the sheet can be discerned by the faint screen-like markings on the

Procedure (cont'd)

sheet when viewed with glancing light. It is important that the proper side be coated since each side may give different results.

A dry coat weight of 10-12#/3300 ft.² is applied using the appropriate wire wound rod. This is equivalent to a dry coat weight of .00957-.01148 gm./in.².

The uncoated base sheets are supplied cut to 12" x 15" with the grain direction in the long dimension. The base sheet should be weighed and the base sheet weight calculated as follows:

$$\frac{\text{grams}}{12" \times 15"} = \text{base sheet weight in gm./in.}^2$$

Immediately after the coating is applied, the sheet is dried in the infrared oven for one minute at a temperature of 120°C. If the sheet is not properly dried, the results may be too low.

After drying, the sheet is allowed to come to equilibrium in a constant temperature and humidity room and trimmed to a size of 10" x 11". It is important to keep track of the grain direction. A good convention to follow is to keep the long dimension the grain direction.

Weigh the coated sheet and calculate the weight per unit area:

$$\frac{\text{weight of coated sheet in grams}}{10" \times 11"} = \text{Coated sheet wt. in gm./in.}^2$$

The coat weight is the difference between the coated sheet weight and base sheet weight.

3. Printing of the Samples

The samples to be tested should be allowed to equilibrate to constant temperature and humidity. Control of both room temperature and humidity is very important in obtaining reproducible results. T.A.P.P.I. standard conditions of 73°F. and 50% relative humidity are recommended if at all

Procedure (cont'd)

possible to achieve with tolerable temperature variations of $\pm 3.5^{\circ}\text{F}$ and relative humidity variations of $\pm 2\%$.

The samples should be cut as strips one inch wide by ten inches long; in this case with the grain direction being the short dimension. Avoid touching the coated surface or in any other way contaminating the coated surface.

The tester should be set up and inked with 2 ml of #5 ink according to the instructions supplied with the tester. The B scale spring drive is used to give a velocity range of from 0 to 630 fpm. A loading of 70 kg. should be used on the printing nip to insure good contact of inked roll and coated sample.

At least two strips from each sample should be printed. In addition, a control latex should be run with each set of samples to insure reproducibility of results from one day to the next.

After every four strips, the ink supply on the distribution system should be replenished with 0.16 ml of ink. The glue-glycerine roll should be used as the top roll of the ink distribution device. After two hours the ink distribution system should be cleaned up and a fresh supply of ink used.

4. Determination of the End-Point

The test is designed to measure the tack stresses on the surface of the sheet which cause rupture of the coating film. These stresses are a function of the ink viscosity and the speed of printing. In the I.G.T. test the ink viscosity is kept constant by the use of a special tack graded ink. The strip which is printed is accelerated through the printing nip so that the speed of the sample through the nip is a logarithmic function of distance along the nip.

The end point measured is that speed which just causes rupture of the coating film. In order to define this point the sample is removed from the tester and examined under a low power microscope. Starting with the lowest speed end (the clamped end) draw the sample under the microscope and look for areas

Procedure (cont'd)

where the coating has been torn away. These will show up as distinct white spots against the black ink printed surface. Do not confuse white spots caused by the ink not covering the surface with picked areas. A little practice will be needed to tell the difference between a non-inked area and a pick. As the strip is scanned, the severity of picking should increase. The end point is marked as the point at which picking appears to increase in severity. Occasionally, a few isolated will be observed on the early part of the sample. These should be ignored if the picking does not increase in severity with distance. There is a considerable judgment factor involved in consistently calling the end point so it is often advisable to have only one person perform the test.

The end point is marked on the strip and converted to a speed reading (f.p.m.) using the scale supplied with the instrument.

:ff
7/20/71

RSV 0015052

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.268
July 17, 1971
Page 1 of 2

MATERIAL

Aqueous Charge

ANALYSIS REQUIRED

SLS Concentration

APPARATUS

100 ml, 200 ml volumetric flasks
20 ml pipet
125 ml separatory funnel
50 ml pipet
Cotton
Beckman Spectrophotometer

REAGENTS

Methylene Blue Solution
Chloroform

PROCEDURE

Accurately weigh 1 g sample into a 200 ml volumetric flask.
Dissolve in distilled water and make up to volume.

Pipette 20 ml into a 125 ml separatory funnel.

Add 30 ml of distilled water and 25 ml of methylene blue
solution. Mix well by swirling.

Extract with 15 ml CHCl_3 , then 25 ml CHCl_3 collecting the
fractions in a 100 ml volumetric flask.

NOTE: Filter the fractions through a cotton plug that has
been previously wetted with CHCl_3 . Bring to volume with
 CHCl_3 .

Read at 652 m μ on the Beckman vs. CHCl_3 for blank.

Settings Required: Phototube in Load
Resistor #3
Slit width 0.04 to 0.1
Tungsten lamp

RSV 0015053

CALCULATIONS

$$\% \text{ SLS} = \frac{(A652)(0.38)}{\text{Sample Weight}}$$

SCOPE

This method is for small concentrations of S.L.S. (40.5%) in process streams such as the initial aqueous charge. It involves complexing the SLS with methylene blue; extraction with CHCl_3 , and subsequent U.V. reading at 652 m μ .

DISCUSSION

1. Samples should be read within 90 min.
2. Sample size may have to be varied with stream charge variations.
3. Methylene blue, N.F. - Matheson, Coleman and Bell.

Methylene Blue Solution - dissolve 0.1 g methylene blue, N.F. contained in a 100 ml volumetric flask in distilled water and dilute to 100 ml. Transfer 30 ml of this solution to a 1 liter volumetric flask, containing 500 ml of distilled water, 6.8 ml of concentrated sulfuric acid, and 50 g of sodium sulfate anhydrous, A.R. grade. Shake until solution is complete and dilute to volume with distilled water.

:ff
7/20/71

RSV 0015054

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.289-A (WC)
July 17, 1971
Page 1 of 4

MATERIAL

Vinyl Chloride

ANALYSIS REQUIRED

Determination of Iron in Vinyl Chloride

SAMPLE

100 ml glass pressure sample bomb

APPARATUS

Beckman Spectrophotometer
Dry Ice Chest
Two, 250 ml Erlenmeyer Flasks
100 ml Graduated Cylinder
10 ml Graduated Cylinder
100 ml Volumetric Flask

REAGENTS

- 10% Aqueous solution of Hydroxylamine hydrochloride
- 2 molar solution of C.P. sodium acetate (164 ± 0.01 g in 1 liter of distilled water)
- 0.1% solution of brom phenol blue (0.1 g dissolved in 100 ml of 50% ethyl alcohol)
- 0.12% aqueous solution of o-phenanthroline
- Alcohol water solution (25 parts 2B alcohol and 30 parts water)
- Concentrated C.P. HCl

PROCEDURE

Place 25 ml of the alcohol-water solution in a 250 ml erlenmeyer flask. From the sample bomb, previously cooled for two hours in the dry ice chest, measure out 100 ml of vinyl chloride into a chilled graduated cylinder. Add the vinyl chloride to the 250 ml erlenmeyer flask and allow to evaporate. Then swirl the flask to remove the last trace of vinyl chloride. A vented hood is needed for the above evaporation.

RSV 0015055

Procedure (cont'd)

Quantitatively transfer the residue to a 100 ml volumetric flask.

Acidify the solution in the volumetric flask by adding two drops of concentrated hydrochloric acid. Add two drops of brom phenol blue to the flask and then add the 2 molar sodium acetate solution dropwise until the solution turns a faint blue. Next, add four ml of 10% hydroxylamine hydrochloride solution, and then 10 ml of 0.12% o-phenanthroline. Dilute to the mark with distilled water; mix well, and allow the solution to stand for at least 20 minutes.

With the Beckman spectrophotometer measure the light transmission at 510 mu, with visible light.

CALCULATIONS

Refer to the attached graph and report the iron concentration to the nearest 0.05 ppm.

PREPARATION OF REAGENTS

10% aqueous hydroxylamine HCl - dissolve 100 g of the salt in 900 ml of distilled water. Be sure water is iron free. Store in P.E. bottle.

2 molar sodium acetate - dissolve 164 gms of C.P. salt in one liter of distilled water. Mix. Store in clean bottle.

.1% solution of brom phenol blue - dissolve .1 gm of the indicator in 50 ml of EtOH. Dilute to 100 ml with distilled water.

.12% of o-phenanthroline - dissolve .12 gms of o-phenanthroline in 100 ml distilled water.

Alcohol water solution - Mix 455 ml of 2B alcohol and distilled water up to one liter.

1000 ppm Fe Standard - dissolve .1000 gms of C.P. iron wire in 10 ml of 20% HCl (iron free). Cool and dilute to 100 ml with distilled water. Mix. Keep for stock reagent.

10 ppm Standard - dilute 1 ml of 1000 ppm standard to 100 ml with demineralized water. Mix. Keeps only one or two days.

STANDARDIZATION AND/OR CALIBRATION

To five 100 ml volume flasks add 0, 3, 5, 7, and 10 ml of the 10 ppm standard. These correspond to 0, 30, 50, 70 and 100 μ gms of Fe.

Follow the procedure outlined above from paragraph three on.

Determine absorbance in 20 mm cells against the zero standard above at 510 m μ .

Plot the absorbances so obtained against the corresponding μ gms of Fe. Draw in the best average line. Use this as a working curve or calculate a factor f from it such that f x absorbance = μ gms Fe.

SAFETY

Handle VCM with caution. Keep vapors away from sparks and flames. Guard against cold burns, handle with gloves. Use normal precautions with chemicals, glassware, and electrical apparatus.

ANALYTICAL TIME

1/2 hour

PRECISION OF RESULTS & SIGNIFICANT FIGURES TO REPORT

95% confidence limits are $\pm .1$ ppm at a level of 1 ppm or $\pm 10\%$ of the level found. Report results to nearest .05 ppm.

SCOPE

VCM sample is evaporated leaving any iron in a non-freezing mixture of water and alcohol. The iron is reduced to valence two and reacted with 1,10 phenanthroline to a red product analytical to iron at 510 m μ . These are known interferences.

REFERENCES

Plastics Division Control Lab. procedure 112-70-3.
Texas City C.L. Method 22-18-10, G.L. 4/22/54.

:ff

7/20/71

RSV 0015057

METHOD 602.289-A (WC)

July 17, 1971

Page 4 of 4

IRON IN VCM CALIBRATION CURVE

Wave Length: 510 mu
Cell Length: 20 mm
Standard: Blank (with reagents)

<u>ppm Fe</u>	<u>Optical Density</u>
0-0.1	0-0.05
0.2	0.05-0.09
0.3	0.090-0.112
0.4	0.112-0.150
0.5	0.150-0.19
0.6	0.19-0.23
0.7	0.23-0.27
0.8	0.27-0.31
0.9	0.31-0.34
1.0	0.34-0.40

RSV 0015058

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.331 (G.C.)
July 17, 1971
Page 1 of 3

MATERIAL

Product ethylene

ANALYSIS REQUIRED

Methanol

SAMPLE

200 ml stainless steel flow thru bomb

APPARATUS

Gas chromatograph equipped as follows:

1. Hydrogen flame ionization detector which should have a sensitivity of approximately 1 coulomb per mole for a C_4 hydrocarbon.
2. Amplifier and recorder with sensitivity sufficient to yield a full scale response to 1×10^{-11} amps of ion current and a noise level equivalent to less than 1×10^{-13} amps.
3. Gas sampling valve with a 500 sample loop.
4. Temperature controller for operation at 80°C.
5. Column consisting of 20 feet of 1/8 inch O.D. aluminum tubing packed with 30-60 mesh red chromosorb coated with 20% EAC-1-R- 296 polyester (available from Cambridge Industries Co., Inc., Cambridge, Mass.) Methylene chloride is the solvent used in coating the chromosorb with the polyester.

REAGENTS

1. Ethylene - Pure grade from Phillips Petroleum Co.
2. Methanol - Reagent grade
3. Hydrogen - Electrolytic grade
4. Helium - Commercial
5. Air - Free of water, organic compounds and particulate matter.

RSV 0015059

PROCEDURE

1. Install the column in the instrument and adjust oven temperature to 80°C and helium carrier flow to 45 cc per minute.
2. Adjust the hydrogen and air flow to the detector for maximum signal-to-noise ratio. Vary hydrogen flow first then air.
3. Allow instrument to "line out" at these conditions as indicated by a stable baseline on the recorder for five minutes or more using maximum sensitivity on the amplifier.
4. Using the gas sample valve, introduce 5cc of sample into the instrument.
5. Record the chromatogram. Methanol is eluted in approximately 7 minutes.

CALCULATIONS

1. Measure the area of the methanol peak.
2. Calculate percent methanol as follows:

$$\text{Peak area} \times \text{attenuation} \times \text{calibration factor} = \% \text{ MeOH}$$

CALIBRATION

1. Prepare synthetic mixtures of methanol in ethylene to cover concentration range of interest. The importance of using clean bombs for preparation of synthetic mixtures cannot be over emphasized.
2. Analyze the synthetic mixtures according to the above procedure.
3. Calculate the calibration factor.

SAFETY PRECAUTIONS

Normal safety precautions on handling hydrocarbon gases should be observed.

ANALYTICAL TIME

Average time required for this analysis is 0.25 hours.

PRECISION OF ANALYSIS

The precision of this method is ± 5 ppm in the range of 5 to 50 ppm methanol in ethylene.

SCOPE

This is a gas chromatographic procedure for the determination of methanol in the concentration range of 1-500 ppm in ethylene. The minimum detection limit is approximately 5 ppm of methanol. C_4 's and up to 1% C_5 's will not interfere with this analysis. Heating bomb above ambient temperature before injecting a sample into the chromatograph was found unnecessary.

REFERENCES

1. Research Method TC - 1180 (GC)
2. Research Notebook HC No. 1, pp. 24-33.

:ff
7/20/71

RSV 0015061

MONSANTO COMPANY
TEXAS CITY, TEXAS

METHOD 602.343 (W.C.)
August 29, 1961
Page 1 of 3

MATERIAL

Vinyl Chloride Monomer

ANALYSIS REQUIRED

Acidity as HCl

SAMPLE SIZE

Glass Bomb

APPARATUS

500 ml Iodine Flasks

5 ml Micro Burette .01 ml Graduations

Mettler Balance

REAGENTS

0.1% Thymol Blue in Methanol

C.P. Methanol

Ethylene Dichloride (Neutral)

Methanolic Caustic (.003 N)

Solvent Mixture Prepared from above.

PROCEDURE

1. Add 100 ml solvent mixture to each of two 500 ml iodine flasks. To the "blank" flask add 60 ml of neutral EDC.
2. Neutralize both flasks to an orange-yellow end point with .003 N MeOH-NaOH. Bring "sample" flask and "blank" flask to the same color. Check the color by looking down through the liquid from the top. Stopper flasks.
3. Place flasks in dry ice chest for 2-3 minutes. Remove and match colors again. Adjust, if necessary, with MeOH-NaOH to the same color again.

RSV 0015062

PROCEDURE (Continued)

4. Attach bomb to 1/8" SS coil in cooling bath (-35°C). Place sample flask under outlet of tubing. Open valve on bomb and drain VCM into sample flask. Tubing should be under the surface of the liquid. When sample stops bubbling, close valve and remove bomb. Purge a small amount of nitrogen through the SS tubing.
5. Remove flask, swirl to mix, and titrate to end point of the blank by matching the colors. Note volume of MeOH-NaOH used.

NOTE: Obtain sample weight by weighing bomb on Mettler balance before and after entering sample into titration flask.

CALCULATION

$$\text{ppm HCl} = \frac{V \text{ MeOH} - \text{NaOH} \times N \text{ MeOH-NaOH} \times 36,500}{\text{Sample Weight}}$$

PREPARATION OF REAGENTS

Titration solvent: 920 ml neutral EDC, 30 ml 0.1% thymol blue and 50 ml CP. Methanol per liter of solution.

Neutral EDC: Refer to Solutions Manual.

.003 MeOH-NaOH: Refer to Solutions Manual.

0.1% Thymol Blue in Methanol: Refer to Solutions Manual.

STANDARDIZATION AND CALIBRATION:

Refer to Solutions Manual for standardization of .03 MeOH-NaOH. This in turn is diluted 1:10 with methanol for .003 N MeOH-NaOH.

SAFETY PRECAUTIONS

When disposing of material in sample flask, be sure that all VCM has been evaporated from the solution before the EDC is put in solvent drain.

ANALYTICAL TIME

0.3 Hours.

RSV 0015063

METHOD 602.343 (W.C.)
August 29, 1961
Page 3 of 3

PRECISION AND SIGNIFICANT FIGURES

$\pm$ 0.1 ppm at 1 ppm level. Report to the nearest 0.1 ppm.

SCOPE

This is a non-aqueous acid-base titration capable of detecting about 10 micrograms of HCl.

REFERENCE

Research Notebook #2271, pp. 148746-53.

/jw
7/29/71

RSV 0015064

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.344 (WC)
May 14, 1964
Page 1 of 5

MATERIAL

Product vinyl chloride and 22E11

ANALYSIS REQUIRED

Determination of aldehydes

SAMPLE SIZE

Quart thermos bottle

APPARATUS

Beckman Model "B" Spectrophotometer

100 ml volumetric flask

Coca-cola bottle and cap

Bottle opener

Ice pick

250 ml Erlenmeyer flask

100 ml graduated cylinder

25 ml graduated cylinder

1 ml pipette

Triple beam balance

REAGENTS

Alpha - methyl indole solution

Ammonium hydroxide (NH_4OH), (SP conc.)

Acetaldehyde (CH_3CHO)

RSV 0015065

PROCEDURE

Select a clean, dry Coca-cola bottle free from large scratches. Pipette 1.0 ml of concentrated NH_4OH into the bottle and add 10 ml H_2O .

Pour 100 ml of vinyl chloride into a 100 ml graduate which has been cooled in dry ice for at least 15 minutes. Transfer to the Coca-cola bottle and cap the bottle immediately.

Shake the sealed bottle vigorously and place in the cabinet for one hour (CAUTION: Wear face shield while handling vinyl chloride). Puncture the cap with an ice pick and then remove cap. Allow the vinyl chloride to evaporate in the hood.

When all of the vinyl chloride has evaporated, wash the contents of the bottle into a 100 ml volumetric flask with several portions of distilled water. The volume of distilled water shall not exceed 50 ml. Add 20 ml of α -methylindole solution to the washings. Shake to mix thoroughly and allow the solution to stand for 10 minutes. Dilute to 100 ml with distilled water and again shake thoroughly.

Read the optical density of this solution on the Beckman "B" Spectrophotometer at 550 m μ using 20 mm cells and distilled water as the blank.

If the turbidity exceeds the range of the chart, repeat the analysis using a smaller volume of vinyl chloride.

CALCULATIONS

The results are read directly from the chart to the nearest part per million. Duplicate determinations should check within 2 ppm. See chart on last page of procedure.

PREPARATION OF REAGENT SOLUTIONS

Prepare fresh solutions daily by adding 60 ml concentrated CP HCl to 30 ml of distilled water in a 250 ml Erlenmeyer flask. Then dissolve in this 0.15 ± 0.05 grams of alpha methyl indole.

STANDARDIZATION AND/OR CALIBRATIONS

Chill a 5 ml Mohr pipette graduated in 0.1 ml increments to a temperature of 0°C by placing in an ice chest. Wrap pipette with a paper towel to avoid moisture contamination. After ten minutes, remove and pipette 5 ml of technical grade acetaldehyde. Introduce 1.25 ml of this chilled acetaldehyde into a 1000 ml chilled volumetric flask containing distilled water. Make up to volume with distilled water and shake thoroughly. The contents of this flask contains 1000 ppm acetaldehyde. Call this Solution A.

10 millileters of A diluted to 1 liter water = 10 ppm AcH

5 millileters of A diluted to 1 liter water = 5 ppm AcH

1 millileter of solution A diluted to 1 liter = 1 ppm AcH

Intermediate parts are obtained by diluting various parts of A to 1 liter. These aliquot portions are then run in a similar manner to that of the sample and the various scale readings recorded. Duplicate determinations should be run for a check.

These scale readings will relate to the concentration of the acetaldehyde present.

SAFETY PRECAUTIONS

When depressurizing the VCM in the Coca-cola bottle, be sure that the ice pick is removed slowly and the bottle is kept in a safety cage and wearer has on a face shield.

METHOD 602.344 (WC)

May 14, 1964

Page 4 of 5

ANALYTICAL TIME

0.2 hour

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES FOR REPORTING

Report to nearest 1-ppm; duplicate determinations should check each other $\pm$ 1 ppm at 1 ppm level.

SCOPE

Alpha methyl indole reacts with aldehydes to form a turbidity. A measure (spectrophotometer) of turbidity determines the amount of aldehydes present.

/jw
7/29/71

RSV 0015068

METHOD 602.344 (WC)
May 14, 1964
Page 5 of 5

ACETALDEHYDE IN SPECIFICATION VCM

20 mm Cells

Beckman "B" Spectrophotometer

<u>Optical Density</u>	<u>ppm AcH</u>
0 - 0.027	0
0.028 - 0.082	1
0.083 - 0.136	2
0.137 - 0.190	3
0.191 - 0.245	4
0.246 - 0.299	5
0.300 - 0.354	6
0.355 - 0.409	7
0.410 - 0.463	8
0.464 - 0.518	9
0.519 - 0.572	10

RSV 0015069

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.431
May 23, 1971
Page 1 of 1

MATERIAL

Acrylamide
Ferric Ammonium Sulfate
Sodium Formaldehyde Sulfoxylate
Tetrasodium Ethylene diamine Tetraacetate

ANALYSIS REQUIRED

Appearance of Solids

SAMPLE SIZE

8-ounce bottle

APPARATUS

None

REAGENTS

None

PROCEDURE

1. If the material consists of free flowing relatively stable non-hazardous particles, spread a small amount on a white card. Otherwise observe in the sample container.
2. Describe the color of the sample.
3. Describe the physical size, form, and uniformity of particles such as powder, crystals, flakes, briquettes, lumps, fused mass, etc., with respect to the specification requirement.
4. Note the presence of any discolored areas or particles and/or the presence of extraneous materials.

SCOPE

This method is intended to provide a general physical description of solid materials.

:ff
7/20/71

RSV 0015070

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.432
May 23, 1971
Page 1 of 1

MATERIAL

Sodium Lauryl Sulfate
Tetrahydrofuran

ANALYSIS REQUIRED

Appearance of Liquids

SAMPLE SIZE

8-ounce bottle

APPARATUS

Source of strong oblique illumination shaded from the observer's eye, such as Turbidimeter, Clark, Catalog No. T11315, The Chemical Rubber Co., Cleveland, Ohio.

REAGENTS

None

PROCEDURE

1. Rinse the outside of the sample container, a clear glass bottle, with methanol and water, wipe it with a towel and place it on a white card or opal glass background.
2. Describe the color of the sample.
3. Describe the clarity with respect to haze, sediment, globules or layers of immiscible liquids or other insolubles.
4. Swirl the sample, avoiding the introduction of air bubbles, and observe under the Clark Turbidimeter, describing any turbidity, opalescence, or sediment.
5. Describe the fluidity with respect to the specification requirement.

SCOPE

This method is intended to provide general physical description of liquid materials when viewed under a constant set of conditions, and to define any visibly apparent contamination.

:ff
7/20/71

RSV 0015071

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.433
May 23, 1971
Page 1 of 2

MATERIAL

Acrylamide

ANALYSIS REQUIRED

Assay

SAMPLE SIZE

8-ounce bottle

APPARATUS

1. Erlenmeyer flask, 250 ml heavy wall, fitted with a 24/40 Standard Taper joint (Ace Glass Co., No. 2671).
2. Filling funnel, 50 ml cylindrical, open top with Standard Taper stop-cock and 24/40 Standard Taper joint on bottom (Kontes Glass Co., No. K-63125). In the top of the funnel insert a rubber stopper fitted with a short glass tube for connection to vacuum line.

REAGENTS

1. 0.1 N KBrO_3 - KBr solution (2.79 g. potassium bromate and 10 g. of potassium bromide per liter).
2. 6 N H_2SO_4
3. 20% KI solution
4. Std. 0.1 N $\text{Na}_2\text{S}_2\text{O}_3$
5. 1% starch solution

PROCEDURE

Weigh accurately a sample containing 1.0 - 1.5 g. of acrylamide into a 500 ml volumetric flask, dilute to volume with water and mix. Pipet a 25.00 ml aliquot of this solution and 25.00 ml of bromate-bromide reagent into the 250 ml Erlenmeyer flask. Attach the filling funnel and rubber stopper, evacuate the flask, then close the stopcock and remove the stopper. Draw 5 ml of 6N sulfuric acid into the flask while preserving the vacuum. Allow the solution to stand in the dark for twenty minutes and shake frequently. Add 15 ml of potassium iodide solution, mix and then break the vacuum. Titrate the solution immediately with standard thiosulfate using starch indicator. Run a blank determination.

RSV 0015072

METHOD 602.433

May 23, 1971

Page 2 of 2

CALCULATION

$$\% \text{ Acrylamide} = \frac{(\text{blank titer} - \text{sample titer}) (N \text{ Na}_2\text{S}_2\text{O}_3) (3.554)}{\text{Weight of sample in aliquot}}$$

SCOPE

Acrylamide may be determined by bromination of the double bond. This method, which is used in the laboratories of the American Cyanamid Company, is applicable to solids and solutions.

:ff
7/20/71

RSV 0015073

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.434
May 23, 1971
Page 1 of 4

MATERIAL

Sodium Formaldehyde Sulfoxylate

ANALYSIS REQUIRED

Ash

SAMPLE SIZE

8-ounce bottle

APPARATUS

- a. Porcelain wide form crucible, such as Coors No. 2, Catalog No. 7-955, Fisher Scientific Company.
- b. Platinum dish, 25 ml. American Platinum Works
- c. Burner, such as Catalog No. 3-902, Fisher Scientific Company.
- d. Nichrome Triangles, such as Catalog No. 15-260, Fisher Scientific Company.
- e. Nichrome Triangles with Fused Silica Tubes, such as Catalog No. 15-280, Fisher Scientific Company.
- f. Desiccator, such as Catalog No. 8-631, Fisher Scientific Co.
- g. Desiccant, such as anhydrous 4-mesh calcium chloride.
- h. Furnace - muffle, such as Catalog No. 10-526, Fisher Scientific Co.

REAGENTS

Sulfuric Acid, T.S., Diluted - Cautiously add 57 ml of sulfuric acid (95-98%) to sufficient water to make 1000 ml.

PROCEDURE

NOTE: This procedure specifies a platinum crucible for general applicability. A porcelain crucible may be used if previous experience with the product being ashed indicates no attack on porcelain.

RSV 0015074

Procedure (cont'd)

1. Sand polish and mold the platinum crucible.
2. Heat the crucible on a silica sleeved Nichrome triangle for 5-10 minutes being careful not to allow the blue reducing flame of the burner to come in contact with the crucible.
3. Allow the crucible to pre-cool for 10 minutes on the silica covered triangle, place in a desiccator for 30 minutes if using platinum or 1 hour if using porcelain, and weigh accurately.
4. Add 2.00 ± 0.01 grams of sample and 10 ml of diluted sulfuric acid T.S. to the crucible.
5. Place the crucible on a porcelain top hot plate and heat gently until the sample is dry and thoroughly charred, then continue heating until the sample has been volatilized or nearly all the carbon has been oxidized.
6. Allow the crucible to cool and then carefully moisten the residue with 0.1 ml of H_2SO_4 .
7. Heat as in Step 5 until the remainder of the sample and any excess sulfuric acid have been volatilized.
8. Finish the ignition in a muffle adjusted to $800 \pm 25^\circ C$.
9. Pre-cool the crucible for 10 minutes on a silica covered triangle, cool in a desiccator for 30 minutes if using platinum or 1 hour if using porcelain, and weigh accurately.
10. Calculate:
$$\% \text{ Ash} = \frac{\text{net weight of residue} \times 100}{\text{sample weight}}$$

DISCUSSION

To obtain constant weight of the porcelain crucible within a reasonable time, it is imperative that it be pre-cooled on a porcelain or aluminum block for 10-15 minutes prior to placing in a desiccator and then further cooled in the desiccator for one hour. Elimination of the pre-cooling step necessitates a three hour or more cooling cycle in the desiccator. One hour cooling periods without the pre-cooling period could incur errors of 5-10 milligrams, depending on the size and number of crucibles and the size of the desiccator.

The initial combustion should be made with as low a temperature as possible. Only when the carbon has fully disappeared from the

Discussion (cont'd)

sides of the crucible should it be placed in a muffle furnace. A further consideration to be kept in mind is that if the flame is allowed to envelop the crucible fully, the atmosphere within the crucible is not air but largely a mixture of nitrogen and the products of combustion, of which water vapor is an important one. The chief effect of the water vapor is to raise the temperature necessary for complete dehydration. This is the reason long ignition periods are often required for the conversion of a metal salt to the oxide, whereas under proper conditions the conversion could be made in minutes.

Once the carbon is completely removed from the wall of the crucible, it is only a matter of expedience as to whether the muffle or the Meker burner is used to burn the carbonaceous matter from the residue (Ref. 1).

If the residue after ignition contains sodium or potassium carbonate, it may fuse and occlude carbon, making complete removal of the latter difficult. Increasing the sample weight only increases the tendency toward false higher answer.

Results will always be low in the presence of easily reducible oxides of volatile metals.

Due to the fluffiness of some ash residues, extreme care must always be exercised to prevent loss by air currents.

PRECISION

The following criteria should be used for judging the acceptability of results (95% confidence):

- a. Repeatability. Duplicate results by the same operator should be considered suspect if they differ by more than the following amount - 0.0 to 0.15 Ash, %; 0.003 Repeatability.
- b. Reproducibility. The results submitted by each of two laboratories should be considered suspect if the two results differ by more than the following amount - 0.0 to 0.15 Ash, %; 0.005 Reproducibility.

SAFETY

Handle hot crucibles with tongs.

Operate in a fume hood, and behind a safety shield, if there is any possibility of the material having explosive tendencies.

Safety (cont'd)

Make sure the ground glass flange of the desiccator is greased to prevent sticking.

REFERENCES

1. "Applied Inorganic Analysis", Hillebrand, Lundell, Bright, and Hoffman, John Wiley & Sons, London, 1953.
2. 1968 Book of ASTM Standards, part 17, page 195.
3. Food Chemicals Codex I.
4. USP XVIII.

SCOPE

This method is intended for the determination of Ash contents of less than 1.0% in pharmaceuticals, oils, etc. The scope is similar to that of the ASTM D-482-63 method for "Ash From Petroleum Products." The method does not necessarily give a true measure of the metallo-organic compounds or the inorganic compounds present because of possible losses by vaporization or reduction of the metal compounds present.

:ff
7/20/71

RSV 0015077

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.435
May 23, 1971
Page 1 of 1

MATERIAL

Ammonium Persulfate

ANALYSIS REQUIRED

Assay

SAMPLE SIZE

8-ounce bottle

APPARATUS

- a. Analytical balance
- b. 250 cc glass stoppered flask
- c. Titrating burets
- d. Usual laboratory equipment

REAGENTS

- a. Standard Acid Ferrous Sulfate solution (0.2 N)
- b. 0.1 N potassium permanganate

PROCEDURE

Add 0.5 gm of sample, accurately weighed, to 25.0 ml of Standard Acid Ferrous Sulfate solution (about 0.2 N) in a glass-stoppered flask. Stopper the flask, allow it to stand for 1 hour with frequent shaking and titrate the excess ferrous sulfate with 0.1 N potassium permanganate. Titrate 25.0 ml of the Standard Acid Ferrous Sulfate solution with 0.1 N potassium permanganate. From the volume of 0.1 N permanganate required, subtract the volume of 0.1 N permanganate required in the titration of the sample solution. Each ml of the difference, consumed by the sample, is equivalent to 0.01141 gm of $(\text{NH}_4)_2\text{S}_2\text{O}_8$.

:ff
7/20/71

RSV 0015078

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.436
May 23, 1971
Page 1 of 1

MATERIAL

Ferric Ammonium Sulfate

ANALYSIS REQUIRED

Assay

SAMPLE SIZE

8-ounce bottle

APPARATUS

- a. Analytical balance
- b. 250 cc glass stoppered flask
- c. Titrating burets
- d. Usual laboratory equipment

REAGENTS

- a. 0.1 N hydrochloric acid
- b. 0.1 N potassium iodide
- c. 0.1 N sodium thiosulfate
- d. Starch indicator

PROCEDURE

Dissolve about 1.5 gm of sample, accurately weighed, in 50 ml of water in a glass-stoppered flask. Add 3 ml of hydrochloric acid and 3 gm of potassium iodide and let stand for 30 minutes in the dark. Dilute with 100 ml of water and titrate the liberated iodine with 0.1 N sodium thiosulfate adding starch indicator solution toward the end of the titration. Correct for a blank. Each ml of 0.1 N sodium thiosulfate consumed is equivalent to 0.04822 gm of $\text{FeNH}_4(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$.

:ff
7/20/71

RSV 0015079

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.437
May 23, 1971
Page 1 of 1

MATERIAL

Sodium Formaldehyde Sulfoxylate

ANALYSIS REQUIRED

Assay

SAMPLE SIZE

8-ounce bottle

APPARATUS

- a. Analytical balance
- b. 1000 ml glass flask; 250 ml glass flask
- c. Titrating buret
- d. Usual laboratory equipment

REAGENTS

- a. 0.10 N iodine
- b. 10% formaldehyde solution
- c. Starch indicator

PROCEDURE

Weigh out 10 grams of SODIUM FORMALDEHYDE SULFOXYLATE, dissolve in 150 cc of warm water, add distilled water to make a total volume of 500 cc. Take a 4 cc aliquot in 100 cc distilled water and titrate with 10 N Iodine (adding first 5 cc of 10% Formaldehyde solution) and Starch Indicator.

CALCULATIONS

$$\frac{\text{cc .10N Iodine} \times .00385 \times 100}{\text{Wt. of sample is 0.08 grams}} = \% \text{ Sodium Formaldehyde Sulfoxylate}$$

:ff
7/20/71

RSV 0015080

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.438
May 23, 1971
Page 1 of 2

MATERIAL

Sodium Lauryl Sulfate

ANALYSIS REQUIRED

Assay

SAMPLE SIZE

8-ounce bottle

APPARATUS

150 ml beaker
250 ml volumetric flask
25 ml pipet
500 ml separatory funnel
Analytical balance

REAGENTS

p-toluidine hydrochloride solution, 6.8%
3-A Alcohol, Form #30
Carbon Tetrachloride (reagent grade)
Meta-cresol purple indicator (0.1% solution)
Phenol red indicator (0.1% solution)
0.1 N NaOH

PROCEDURE

Using an analytical balance, transfer 5.00 grams of SIPON WD to a 150 ml beaker. Dissolve in distilled water and transfer to a 250 ml volumetric flask (a small amount of 3-A Alcohol may be used to dissipate foam) and dilute to 250 ml mark. Mix well and withdraw a 25 ml aliquot portion (equivalent to 0.5 grams of sample).

Transfer the aliquot portion to a 500 ml separatory funnel. Add 2 drops of phenol red indicator (red-alkaline, yellow-acid). If solution is alkaline, adjust to neutral with .1N HCl. If solution is acid, adjust to neutral with .1N NaOH. To neutral solution add 30 ml of 3-A Alcohol, 30 ml of 6.8% p-toluidine hydrochloride solution and 30 ml of carbon tetrachloride. Shake mixture for

RSV 0015081

Procedure (cont'd)

three minutes (occasionally releasing gas from separatory funnel). Allow to stand and separate. Cut lower (carbon tetrachloride) layer into 250 ml wide-mouth Erlenmeyer flask, containing 50 ml of 3-A Alcohol and 3 cc H₂O adjusted to the alkaline (purple) end-point of meta cresol purple. Add 30 ml carbon tetrachloride to contents of separatory funnel and repeat extraction. Allow to stand and separate and again cut lower layer into same receiving flask.

Repeat extraction once more. These combined extractions in the alcohol medium contain all the organic sulfates and are titrated directly with standard 0.1N NaOH to determine the alkyl sulfate content. The end-point is taken as the first change in color from gray to purple that remains for 15 seconds.

CALCULATION
$$\frac{\text{ml NaOH} \times N \times 300 \times 100}{\text{Wt. of sample in the aliquot} \times 1000} = \% \text{ Alkyl Sulfate as sodium lauryl sulfate}$$
SCOPE

This procedure employs the extraction of the organic sulfates by p-toluidine hydrochloride precipitation.

:ff
7/20/71

RSV 0015082

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.439
May 23, 1971
Page 1 of 2

MATERIAL

Tetrasodium Ethylenediamine Tetraacetate

ANALYSIS REQUIRED

Chelation value

SAMPLE SIZE

8 ounce bottle

APPARATUS

- a. Analytical balance
- b. 250 ml Erlenmeyer flask
- c. Titrating burette
- e. Usual laboratory equipment

REAGENTS

- a) Ammonium oxalate, saturated solution. Dissolve approximately 60 grams of ACS reagent grade (calcium-free) ammonium oxalate $(\text{NH}_4)_2\text{C}_2\text{O}_4 \cdot \text{H}_2\text{O}$ in a liter of water and allow to cool to room temperature.
- b) Calcium chloride, standard solution. Dissolve 0.5 molar of primary standard calcium carbonate (Mallincrodt No. 4071) in 300 ml of distilled water by slowly adding 86 ml of concentrated hydrochloric acid to completely dissolve the calcium carbonate. Heat to boiling to remove carbon dioxide and completely dissolve the sample. Cool to room temperature and dilute this solution to exactly one liter. One ml of this solution contains the equivalent of 50 mg of calcium carbonate.
- c) Ammonium hydroxide, 1:1 solution. Add one volume of concentrated ammonium hydroxide to one volume of distilled water.

RSV 0015083

PROCEDURE

Weigh, to the nearest milligram, a 10 gram sample of chelating agent into a clean 250 ml Erlenmeyer flask, and add 85 ml of distilled water and 5 ml of saturated ammonium oxalate solution. Titrate the sample to the first faint permanent turbidity with standard calcium chloride. The pH of the solution should be checked by using pH paper after the endpoint has been reached; if the pH is below 11, a dilute solution of ammonium hydroxide should be added to raise the pH above 11. After the pH has been adjusted, the titration should be completed.

CALCULATION

$$\frac{\text{ml CaCl}_2 \text{ solution} \times 50}{\text{gms of sample}} = \text{mg CaCO}_3 \text{ per gm of chelating agent}$$

SCOPE

The total chelation value is determined by titrating with a solution of calcium salt in the presence of oxalate ion as an indicator at a pH of 11 or slightly higher. The pH is important because the chelating ability decreases when the pH of the solution drops below 11 and consequently low results are obtained.

/jw
7/23/71

RSV 0015084

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.440
May 23, 1971
Page 1 of 2

MATERIAL

Ferric Ammonium Sulfate, Sodium Lauryl Sulfate

ANALYSIS REQUIRED

Chlorides

SAMPLE SIZE

8 ounce bottle

APPARATUS

- a) Usual laboratory apparatus
- b) Potentiometer sensitive to potential changes of 5 mv. or less, coupled to the system:

Ag::Solution::2M MgSO_4 : Hg_2SO_4 :Hg

REAGENTS

- a) 10% H_2SO_4
- b) Methyl orange indicator.
- c) 0.01 N AgNO_3

PROCEDURE

1. Dissolve 5 grams of sample, weighed on a prescription balance, in 40 ml of water in a 150 ml beaker.
2. Slowly neutralize to methyl orange with 10% H_2SO_4 and add 0.5 ml in excess.
3. Cool the solution to 0-5°C in an ice bath and stir mechanically under the silver electrode system.

NOTE: If a precipitate forms, as with some organic salts, add just sufficient acetone to redissolve.

RSV 0015085

PROCEDURE (cont'd)

4. Titrate fairly rapidly with N/100 AgNO₃ until the potential begins to decrease markedly with small increments of titrant.
5. Continue the titration in 0.1 increments, recording the potential after each, then underscore the total volume of titrant at the point corresponding to the largest potential change per increment of titrant.
6. Calculate the chloride content as the specific chloride designated in the specification for the product:

$$\% \text{ Specified chloride} = \frac{\text{ml titre} \times \text{Factor}}{\text{sample weight}}$$

Table of Chloride Factors

<u>Specified Chloride</u>	<u>Factor for N/100 AgNO₃</u>
Cl ⁻	0.0355
CaCl ₂	0.0555
HCl	0.0365
MgCl ₂	0.0476
MnCl ₂	0.0629
NaCl	0.0585
NH ₄ Cl	0.0535

PRECISION AND ACCURACY

Duplicate determinations should agree to within 0.02%.

SAFETY

No unusual hazards are involved in this operation.

SCOPE

This method is intended for the determination of 0.01 to 0.5% ionizable chloride in water soluble materials. It involves a potentiometric titration with silver nitrate solution.

7/23/71

RSV 0015086

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.441
May 23, 1971
Page 1 of 2

MATERIAL

Sodium Lauryl Sulfate

ANALYSIS REQUIRED

Cloud point Determination

SAMPLE SIZE

8 ounce bottle

APPARATUS

- a) ASTM cloud and pour point kerosene cooling baths adjusted to 30°F, 0°F and -30°F with dry ice.
- b) Test jars, gaskets and ASTM Cloud and Pour Point thermometers.
- c) Usual laboratory equipment.

REAGENTS

None

PROCEDURE

1. Bring the fluid to be tested to a temperature at least 25°F above the specification cloud point in the bath adjusted to the next lower temperature.
2. If moisture is present, remove by filtration thru dry filter paper until the fluid is perfectly clear. The temperature at which the filtration is made must be at least 25°F above the cloud point.
3. Fill a test jar to a point slightly above the mark with a portion of the sample.

RSV 0015087

PROCEDURE (cont'd)

4. Close the jar tightly with a stopper carrying a cloud and pour point test thermometer in a vertical position in the center of the jar, with the bulb touching the bottom of the jar.
5. Adjust a gasket around the jar one inch from the bottom and place in the jacket of the bath adjusted to the first temperature below the temperature of the sample.
6. As the temperature falls at each reading that is a multiple of 2, remove the test jar from the jacket, quickly without disturbing the sample, observe against a good light for the appearance of a cloudy haze, and replace the jar in the bath within 3 seconds.
7. If the sample does not develop a cloud when it has cooled to within 20°F of the bath temperature, place the jar in the next lower bath and continue observing at 2°F intervals.
8. Record as the cloud point the temperature at which the first distinct cloudiness or haze appears in the bottom of the test jar.

SAFETY

Kerosene is a flammable liquid. Avoid open flames or excessive heat. Take cognizance of any dermatological, lachrymatory, or corrosive properties of samples and avoid skin contact and inhalation of vapors.

SCOPE

This method is intended for the determination of the cloud point of fluids which are transparent in layers 1-1/2" in thickness. It has the same general scope as ASTM D97-57 as applied to the determination of cloud points in petroleum products.

/jw
7/23/71

RSV 0015088

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.442
May 24, 1971
Page 1 of 5

MATERIAL

Tetrahydrofuran
Tetrasodium ethylenediamine tetraacetate

ANALYSIS REQUIRED

Color of Liquids, Melts, and Solutions (APHA)

APPARATUS AND REAGENTS

- a) Matched short form Nessler tubes of 50 ml capacity and tall form Nessler tubes of 100 ml capacity.
- b) 0.05 M Periodic Acid - Dissolve 11.4 g. of $\text{HIO}_4 \cdot 2\text{H}_2\text{O}$ in water and dilute to 1000 ml.
- c) Color Standards - Prepare APHA 500 standard by dissolving 1.0 g of cobaltous chloride ($\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$) and 1.245 g of potassium chloroplatinate (K_2PtCl_6) in water, adding 100 ml of conc. HCl, and diluting to 1000 ml with water.

NOTE: The absorbance of the 500 platinum-cobalt stock solution must fall within the limits given below. The measurements shall be made in a cell having a 10 mm light path using a Beckman Model B Spectrophotometer at a sensitivity setting of one, or in another spectrophotometer having equivalent resolution and sensitivity. Reagent water in a matched cell shall be used as the reference solution.

<u>Wavelength, mμ</u>	<u>Absorbance</u>
430	0.110 to 0.120
455	0.130 to 0.145
480	0.105 to 0.120
510	0.055 to 0.065

RSV 0015089

APPARATUS AND REAGENTS (cont'd)

From a burette add the amount of APHA 500 standard indicated in the following table to the respective tube, add 0.5 ml of 0.05 M periodic acid and dilute to the mark with double distilled water. If the marks on a series of tubes vary in depth, a final adjustment of all standards to the depth of the shortest should be made. All tubes should be kept covered when not in use and wrapped in black paper so that no light enters from the side.

<u>APHA No.</u>	<u>Short (50 ml.) Tubes</u> <u>APHA 500</u>	<u>Tall (100 ml.) Tubes</u> <u>APHA 500</u>
5	--	1.0 ml.
6	--	1.2 ml.
7	--	1.4 ml.
8	--	1.6 ml.
9	00	1.8 ml.
10	1.0 ml.	2.0 ml.
11	--	2.2 ml.
12	--	2.4 ml.
13	--	2.6 ml.
14	--	2.8 ml.
15	1.5 ml.	3.0 ml.
16	--	3.2 ml.
18	--	3.6 ml.
20	2.0 ml.	4.0 ml.
25	2.5 ml.	5.0 ml.
30	3.0 ml.	6.0 ml.
35	3.5 ml.	7.0 ml.
40	4.0 ml.	8.0 ml.
45	4.5 ml.	9.0 ml.
50	5.0 ml.	10.0 ml.
60	6.0 ml.	12.0 ml.
70	7.0 ml.	14.0 ml.
80	8.0 ml.	16.0 ml.
90	9.0 ml.	18.0 ml.
100	10.0 ml.	20.0 ml.
125	12.5 ml.	
150	15.0 ml.	
175	17.5 ml.	
200	20.0 ml.	
250	25.0 ml.	
300	30.0 ml.	
350	35.0 ml.	
400	40.0 ml.	
450	45.0 ml.	

PROCEDURE

1. Fill a (black paper-wrapped) Nessler tube, using a 50 ml low form tube for off-white fluids or a 100 ml tall form tube for water-white fluids, to the depth of the standards with a portion of the sample, melt, or specified solution.

NOTE: If the sample has any visible turbidity, filter it. (Note this on the analysis ticket).

2. Compare with the APHA standards viewing vertically against a white glass plate.
3. Record the number of the APHA standard which most nearly matches the sample, making a note of any abnormality in tint or clarity of the sample. In the event the color lies midway between two standards, report the darker of the two.

SAFETY

Take cognizance of any potential hazards in the material being tested, such as the temperature melts, corrosiveness, lachrymatory or dermatological properties, and exercise appropriate precautions.

PRECISION

Results should not differ from the mean by more than the following amounts:

Platinum-Cobalt Color	Repeatability	Reproducibility
	One Operator & Apparatus	Different Operators & Apparatus
Under No. 50	2.5	5
No's. 50 to 100	5	10
Over No. 100	10	15

DISCUSSION

The function of the periodic acid is to stabilize the colors of the standards. It will not, however, protect against contamination. Hence, care must be exercised against splashing samples into the standards. In any event, the standards should be renewed monthly, unless obvious discrepancies dictate more frequent renewal.

The chart on page four shows some comparisons of color intensities as designated in various systems. In the comparison of color intensities of the ASTM D-1500 and the NPA systems, the asterisk (*) indicates the colors are not matched but they are close approximations. See also, references 3, 4.

SCOPE

This method is intended for the estimation of the color of liquids in terms of platinum-cobalt (Hazen) standards established by the American Public Health Association (ref. 1). It involves visual comparisons of the sample with prepared standards in matched Nessler tubes.

REFERENCES

1. "Standards Methods for the Examination of Water, Sewage, and Industrial Wastes", American Public Health Association, Inc., New York, 1955.
2. 1968 ASTM Standards, part 20, pages 555 and 817.
3. Federal Standard Stock Catalog, Section IV, Part 5, TT-P-141b, Method 458.1.
4. Quality Assurance Register Directory/Handbook and Specifications Catalog, Hitchcock Publishing Company, Hitchcock Building, Wheaton, Illinois 60187.

/jw
7/26/71

RSV 0015092

COLOR OF LIQUIDS, MELTS, AND SOLUTIONS (APHA)

COLOR INTENSITY COMPARISONS

Unit No.	APHA	Barrett	ASTM D-1500	NPA	Sapthali	Gardner 1933	Gardner/Neale 1931	Hellige S-P	Print	DuPont	Fac	Purkin	K ₂ Cr ₂ O ₇ H ₂ SO ₄	Lut/Lead Y/R	Photometric Index (A. O. C. S. I) 440/550 nm	Reine Standard	Unit No.
1	5	1	0.5	0.5	0.5	1	1	1	1	1	1	1	1	1	1	1	1
2	10	2	1.0	1.0	1.0	2	2	2	2	2	2	2	2	2	2	2	2
3	15	3	1.5	1.5	1.5	3	3	3	3	3	3	3	3	3	3	3	3
4	20	4	2.0	2.0	2.0	4	4	4	4	4	4	4	4	4	4	4	4
5	25	5	2.5	2.5	2.5	5	5	5	5	5	5	5	5	5	5	5	5
6	30	6	3.0	3.0	3.0	6	6	6	6	6	6	6	6	6	6	6	6
7	35	7	3.5	3.5	3.5	7	7	7	7	7	7	7	7	7	7	7	7
8	40	8	4.0	4.0	4.0	8	8	8	8	8	8	8	8	8	8	8	8
9	45	9	4.5	4.5	4.5	9	9	9	9	9	9	9	9	9	9	9	9
10	50	10	5.0	5.0	5.0	10	10	10	10	10	10	10	10	10	10	10	10
11	55	11	5.5	5.5	5.5	11	11	11	11	11	11	11	11	11	11	11	11
12	60	12	6.0	6.0	6.0	12	12	12	12	12	12	12	12	12	12	12	12
13	65	13	6.5	6.5	6.5	13	13	13	13	13	13	13	13	13	13	13	13
14	70	14	7.0	7.0	7.0	14	14	14	14	14	14	14	14	14	14	14	14
15	75	15	7.5	7.5	7.5	15	15	15	15	15	15	15	15	15	15	15	15
16	80	16	8.0	8.0	8.0	16	16	16	16	16	16	16	16	16	16	16	16
17	85	17	8.5	8.5	8.5	17	17	17	17	17	17	17	17	17	17	17	17
18	90	18	9.0	9.0	9.0	18	18	18	18	18	18	18	18	18	18	18	18
19	95	19	9.5	9.5	9.5	19	19	19	19	19	19	19	19	19	19	19	19
20	100	20	10.0	10.0	10.0	20	20	20	20	20	20	20	20	20	20	20	20
21	105	21	10.5	10.5	10.5	21	21	21	21	21	21	21	21	21	21	21	21
22	110	22	11.0	11.0	11.0	22	22	22	22	22	22	22	22	22	22	22	22
23	115	23	11.5	11.5	11.5	23	23	23	23	23	23	23	23	23	23	23	23
24	120	24	12.0	12.0	12.0	24	24	24	24	24	24	24	24	24	24	24	24
25	125	25	12.5	12.5	12.5	25	25	25	25	25	25	25	25	25	25	25	25
26	130	26	13.0	13.0	13.0	26	26	26	26	26	26	26	26	26	26	26	26
27	135	27	13.5	13.5	13.5	27	27	27	27	27	27	27	27	27	27	27	27
28	140	28	14.0	14.0	14.0	28	28	28	28	28	28	28	28	28	28	28	28
29	145	29	14.5	14.5	14.5	29	29	29	29	29	29	29	29	29	29	29	29
30	150	30	15.0	15.0	15.0	30	30	30	30	30	30	30	30	30	30	30	30
31	155	31	15.5	15.5	15.5	31	31	31	31	31	31	31	31	31	31	31	31
32	160	32	16.0	16.0	16.0	32	32	32	32	32	32	32	32	32	32	32	32
33	165	33	16.5	16.5	16.5	33	33	33	33	33	33	33	33	33	33	33	33
34	170	34	17.0	17.0	17.0	34	34	34	34	34	34	34	34	34	34	34	34
35	175	35	17.5	17.5	17.5	35	35	35	35	35	35	35	35	35	35	35	35
36	180	36	18.0	18.0	18.0	36	36	36	36	36	36	36	36	36	36	36	36
37	185	37	18.5	18.5	18.5	37	37	37	37	37	37	37	37	37	37	37	37
38	190	38	19.0	19.0	19.0	38	38	38	38	38	38	38	38	38	38	38	38
39	195	39	19.5	19.5	19.5	39	39	39	39	39	39	39	39	39	39	39	39
40	200	40	20.0	20.0	20.0	40	40	40	40	40	40	40	40	40	40	40	40
41	205	41	20.5	20.5	20.5	41	41	41	41	41	41	41	41	41	41	41	41
42	210	42	21.0	21.0	21.0	42	42	42	42	42	42	42	42	42	42	42	42
43	215	43	21.5	21.5	21.5	43	43	43	43	43	43	43	43	43	43	43	43
44	220	44	22.0	22.0	22.0	44	44	44	44	44	44	44	44	44	44	44	44
45	225	45	22.5	22.5	22.5	45	45	45	45	45	45	45	45	45	45	45	45
46	230	46	23.0	23.0	23.0	46	46	46	46	46	46	46	46	46	46	46	46
47	235	47	23.5	23.5	23.5	47	47	47	47	47	47	47	47	47	47	47	47
48	240	48	24.0	24.0	24.0	48	48	48	48	48	48	48	48	48	48	48	48
49	245	49	24.5	24.5	24.5	49	49	49	49	49	49	49	49	49	49	49	49
50	250	50	25.0	25.0	25.0	50	50	50	50	50	50	50	50	50	50	50	50
51	255	51	25.5	25.5	25.5	51	51	51	51	51	51	51	51	51	51	51	51
52	260	52	26.0	26.0	26.0	52	52	52	52	52	52	52	52	52	52	52	52
53	265	53	26.5	26.5	26.5	53	53	53	53	53	53	53	53	53	53	53	53
54	270	54	27.0	27.0	27.0	54	54	54	54	54	54	54	54	54	54	54	54
55	275	55	27.5	27.5	27.5	55	55	55	55	55	55	55	55	55	55	55	55
56	280	56	28.0	28.0	28.0	56	56	56	56	56	56	56	56	56	56	56	56
57	285	57	28.5	28.5	28.5	57	57	57	57	57	57	57	57	57	57	57	57
58	290	58	29.0	29.0	29.0	58	58	58	58	58	58	58	58	58	58	58	58
59	295	59	29.5	29.5	29.5	59	59	59	59	59	59	59	59	59	59	59	59
60	300	60	30.0	30.0	30.0	60	60	60	60	60	60	60	60	60	60	60	60
61	305	61	30.5	30.5	30.5	61	61	61	61	61	61	61	61	61	61	61	61
62	310	62	31.0	31.0	31.0	62	62	62	62	62	62	62	62	62	62	62	62
63	315	63	31.5	31.5	31.5	63	63	63	63	63	63	63	63	63	63	63	63
64	320	64	32.0	32.0	32.0	64	64	64	64	64	64	64	64	64	64	64	64
65	325	65	32.5	32.5	32.5	65	65	65	65	65	65	65	65	65	65	65	65
66	330	66	33.0	33.0	33.0	66	66	66	66	66	66	66	66	66	66	66	66
67	335	67	33.5	33.5	33.5	67	67	67	67	67	67	67	67	67	67	67	67
68	340	68	34.0	34.0	34.0	68	68	68	68	68	68	68	68	68	68	68	68
69	345	69	34.5	34.5	34.5	69	69	69	69	69	69	69	69	69	69	69	69
70	350	70	35.0	35.0	35.0	70	70	70	70	70	70	70	70	70	70	70	70
71	355	71	35.5	35.5	35.5	71	71	71	71	71	71	71	71	71	71	71	71
72	360	72	36.0	36.0	36.0	72	72	72	72	72	72	72	72	72	72	72	72
73	365	73	36.5	36.5	36.5	73	73	73	73	73	73	73	73	73	73	73	73
74	370	74	37.0	37.0	37.0	74	74	74	74	74	74	74	74	74	74	74	74
75	375	75	37.5	37.5	37.5	75	75	75	75	75	75	75	75	75	75	75	75
76	380	76	38.0	38.0	38.0	76	76	76	76	76	76	76	76	76	76	76	76
77	385	77	38.5	38.5	38.5	77	77	77	77	77	77	77	77	77	77	77	77
78	390	78	39.0	39.0	39.0	78	78	78	78	78	78	78	78	78	78	78	78
79	395	79	39.5	39.5	39.5	79	79	79	79	79	79	79	79	79	79	79	79
80	400	80	40.0	40.0	40.0	80	80	80	80	80	80	80	80	80	80	80	80
81	405	81	40.5	40.5	40.5	81	81	81	81	81	81	81	81	81	81	81	81
82	410	82	41.0	41.0	41.0	82	82	82	82	82	82	82	82	82	82	82	82
83	415	83	41.5	41.5	41.5	83	83	83	83	83	83	83	83	83	83	83	83
84	420	84	42.0	42.0	42.0	84	84	84	84	84	84	84	84	84	84	84	84
85	425	85	42.5	42.5	42.5	85	85	85	85	85	85	85	85	85	85	85	85
86	430	86	43.0	43.0	43.0	86	86	86	86	86	86	86	86	86	86	86	86
87	435	87	43.5	43.5	43.5	87	87	87	87	87	87	87	87	87	87	87	87
88	440	88	44.0	44.0	44.0	88	88	88	88	88	88	88	88	88	88	88	88
89	445	89	44.5	44.5	44.5												

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.443
May 31, 1971
Page 1 of 2

MATERIAL

Tetrahydrofuran

ANALYSIS REQUIRED

Distilling range

APPARATUS

- a. 100 ml Engler flask
- b. + 50°C to +70°C, grade AA, 76 mm immersion thermometer graduated in 0.1°C
- c. Asbestos ring with 1 1/4 inch diameter
- d. 100 ml graduated cylinder
- e. 9/16" O.D. No. 20GA seamless brass tubing 2 feed long for condenser
- f. Bunsen burner, ring stand and steam bath
- g. Usually laboratory equipment
- h. Apparatus setup pictured in Figure 1




Figure 1
BOILING RANGE APPARATUS

RSV 0015094

PROCEDURE AND CALCULATION

Prior to conducting the boiling range test, use Method 602.446 to check the THF sample for THF-peroxide. If more than trace amounts of peroxide are present, add a few caustic pellets to the distillation flask to be used in the determination and proceed with the following test.*

Measure a 100 ml sample of THF at approximately 20°C, into a standard 100 ml Engler flask. Insert a +50°C to +70°C, grade AA, 76 mm immersion thermometer graduated in 0.1°C. (Wm. Hiergesell & Sons Catalog No. 6536) into the neck of the Engler flask so that the top of the mercury in the expansion bulb is level with the inside of the bottom of the flask vapor outlet tube. Support the flask by an asbestos ring or a "Transite"² plate having a 1 1/4 inch diameter hole and assemble the distillation apparatus as shown in Figure 1. Place a 100 ml graduated cylinder at the outlet of the condenser tube with the tube extending at least one inch into, but not below the 100 ml mark of, the cylinder. Cover the top of the graduated cylinder closely with blotting paper cut to fit the condenser tube tightly.

Heat the contents of the flask over an enclosed steam bath, taking care that the water in the steam bath is actively boiling before the flask is placed in position.

Place a "Transite" board guard on all sides around the legs of the steam bath except that facing away from the condenser assembly.

Record as the initial boiling point the temperature at which the first drop of condensate falls from the end of the condenser tube. Adjust the receiving graduate so that the end of the condenser tube touches the side of the cylinder and continue distilling. Maintain the rate of distillation between 3 ml and 5 ml per minute and record the temperature at which 95 ml has been collected, halt the distillation.

Record the barometric pressure and room temperature at the time of distillation and correct the recorded boiling range temperatures to 760 mm of mercury according to the following formula:

temperature correction (°C) = $(760 - \text{barometric reading reduced to } 0^\circ\text{C}) \times 0.041$

* Based on ASTM D 1078-58 test method.

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.444
May 31, 1971
Page 1 of 1

MATERIAL

Tetrahydrofuran

ANALYSIS REQUIRED

Specific Gravity

APPARATUS

Westphal balance

PROCEDURE

Determine the specific gravity at 20°C by means of a Westphal balance which has been standardized with distilled water at 4°C.

/jw
7/26/71

RSV 0015096

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.445
May 31, 1971
Page 1 of 1

MATERIAL

Tetrahydrofuran

ANALYSIS REQUIRED

Moisture

APPARATUS

Usual laboratory equipment including 100 ml titrating buret.

REAGENTS

Anhydrous methanol
Karl Fischer reagent

PROCEDURE AND CALCULATIONS

Titrate 50 ml of anhydrous methanol to a permanent deep-red color with Karl Fischer reagent. Add a 10 ml sample of tetrahydrofuran and titrate immediately to the same endpoint with Karl Fischer reagent, using an "Aquameter" or the visual endpoint. Calculate percent moisture by use of the following equation.

$$\% \text{ moisture} = \frac{(\text{ml reagent}) (\text{water equivalent}) \times 100}{(\text{specific gravity}) \times 10}$$

/jw
7/26/71

RSV 0015097

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.446
May 31, 1971
Page 1 of 2

MATERIAL

Tetrahydrofuran

ANALYSIS REQUIRED

Peroxide

APPARATUS

- a. 2 x 250 ml Erlenmeyer flasks and stoppers
- b. 25 ml graduated cylinders
- c. Titrating buret and equipment

REAGENTS

- a. 1:4 concentrated sulfuric acid to water
- b. 10% potassium iodide solution
- c. 0.02 N sodium thiosulfate

PROCEDURE AND CALCULATION

Unstabilized tetrahydrofuran upon exposure to air may slowly form THF hydroperoxide, also designated as THF-peroxide. A simple, accurate analysis for THF-peroxide may be made as follows:

1. Place 100 ml of distilled water in each of two Erlenmeyer flasks.
2. To each flask add 25 ml of a 1:4 mixture of concentrated sulfuric acid and water, and then add 25 ml of a 10% potassium iodide solution.
3. Add a 25 ml sample of tetrahydrofuran to one flask, stopper both flasks, shake, and place in the dark for 15 - 20 minutes.
4. Titrate the contents of each flask with 0.02 N sodium thiosulfate until the solutions become colorless.

RSV 0015098

METHOD 602.446
May 31, 1971
Page 2 of 2

PROCEDURE AND CALCULATION

5. Calculate the THF-peroxide content of the THF sample by
by use of the following formula:

$$\begin{aligned} \% \text{ THF peroxide} &= (\text{ml thio in sample} - \text{ml thio in blank}) \times \\ &\quad N \text{ of thio} \times 0.23 \end{aligned}$$

/jw
7/26/71

RSV 0015099

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.447
May 31, 1971
Page 1 of 3

MATERIAL

Tetrahydrofuran

ANALYSIS REQUIRED

Carbonyl and hydroxyl compounds

APPARATUS

- a. Vapor fractometer equipped with one millivolt recorder
- b. 6 ft column filled with 15% carbowax 1500 on 80-100 mesh chromosorb W
- c. 60 ml serum bottle
- d. 0.05 ml Hamilton syringe fitted with Chaney adapter
- e. Usual laboratory equipment

REAGENTS

- a. Helium
- b. Acetone
- c. Isopropyl alcohol
- d. n-propyl alcohol
- e. n-butyl alcohol

PROCEDURE AND CALCULATIONS

The amounts of acetone, isopropyl alcohol, n-propyl alcohol, and n-butyl alcohol in tetrahydrofuran may be determined by gas chromatography, using a 6-ft column filled with 15% "Carbowax" 1500⁴ on 80-100 mesh "Chromosorb" W.² Apparatus suitable for this determination includes a vapor "Fractometer" (equipped with a 1-millivolt recorder)³ or an equivalent gas chromatograph. The column may be purchased as a complete unit (Perkin-Elmer type "K" column or equivalent) or prepared in accordance with instructions outlined below.

To prepare the column, dissolve 3.5 g of "Carbowax" 1500 in 75 ml of methylene chloride and mix with 20 g of 80-100 mesh "Chromosorb" W to make a slurry. Evaporate the methylene chloride on a steam plate with occasional stirring to obtain the dried adsorbent. Fill a 6 ft (1/4" o.d.)

RSV 0015100

PROCEDURE AND CALCULATIONS (cont'd)

thin-walled tube of stainless steel with the dried adsorbent, tapping the tube continually while filling to insure uniform packing. Plug the ends of the tube with glass wool and insert into the "Fractometer", bending the tube as required to fit.

Prepare a standard sample by adding 50 ml of THF known to contain only minute quantities of other components to a tared serum bottle (60 ml) fitted with a rubber stopper. Into the THF standard, measure 0.01 ml acetone, 0.01 ml isopropyl alcohol, 0.01 ml n-propyl alcohol, and 0.02 ml n-butyl alcohol. Calculate the percentage by weight of each constituent added to the THF standard.

Before operating the "Fractometer", adjust the detector to 8.0 volts and sensitivity 1, and allow the "Fractometer" to attain a temperature of 73°C. Apply a helium gas pressure of 15 psi at a flowmeter reading of 5.5 (70 ml/min of helium), and inject 0.01 ml of the THF standard into the "Fractometer" by means of a Hamilton syringe (0.05 ml) fitted with a Chaney adapter. As the carrier gas passes through the "Fractometer", the components of the sample appear at the detector after various lengths of time according to their partition coefficients. The times of elution of the components in the sample are plotted automatically against recorder response on the recorder chart (Figure 2). Deflections in recorder response are indicated on the chart as peaks characteristic of each component.




Figure 2
TYPICAL CHROMATOGRAPH FOR THF

RSV 0015101

PROCEDURE AND CALCULATIONS (cont'd)

Repeat the chromatographic procedure, injecting a 0.01 ml sample of the original THF into the "Fractometer". Measure the heights of the acetone, isopropyl alcohol, n-propyl alcohol, and n-butyl alcohol peaks and compare with the THF standard. Calculate the percentage by weight of each component in the standard in accordance with the following equation:

$$\% \text{ in std} = \frac{\% \text{ added} \times \text{peak height (std)}}{\text{peak height (std)} - \text{peak height (sample)}}$$

Chromatograph the samples of THF to be analyzed under the conditions specified above, and measure the peak heights corresponding to each component. Calculate the percentage by weight of component in the unknown by means of the following equation:

$$\% \text{ in sample} = \frac{\% \text{ (std)} \times \text{peak height (sample)}}{\text{peak height (std)}}$$

/jw
7/26/71

RSV 0015102

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.448
May 31, 1971
Page 1 of 6

MATERIAL

Tetrasodium ethylene diamine tetraacetate

ANALYSIS REQUIRED

Heavy metals (Pb)

APPARATUS

- A. Matched 100 ml low form Nessler tubes.
- B. pHydron test papers, 3.0 to 5.5 (Micro Essential Laboratory) or glass electrode pH indicator.

REAGENTS

- A. 6% acetic acid
- B. Saturated hydrogen sulfide water
- C. Concentrated ammonium hydroxide
- D. 10% hydrochloric acid
- E. Redistilled acetone
- F. Lead stock solution (100 mg/l) - Dissolve exactly 0.1598 gram of lead nitrate in 100 ml of water containing 1 ml of conc. nitric acid and dilute to 1 liter in a volumetric flask.
- G. Standard lead solution (10 mg/l) - Accurately pipette 10.00 ml of the 100 mg/l stock solution into a 100 ml volumetric flask and dilute to volume with water. Prepare fresh daily.
- J. 0.1% p-Nitrophenol indicator solution

DETERMINATION

- 1. Ascertain into which the following classifications the product being tested falls and treat as indicated:
 - a. Soluble Products: Weigh 2.0 grams (see note on page 2) into a 100 ml Nessler tube and treat as directed in the table.

RSV 0015103

DETERMINATION (cont'd)

If the product under test is not listed but is soluble in water or acetone, dissolve in 25 ml of acetone or water; if acid or basic adjust to pH 3.5 with NH_4OH or 10% HCl , or if neutral add 4 ml of 6% HOAc and dilute to 50 ml with the alternate solvent to obtain an approximate 50% aqueous-acetone solution. (Some products may tend to separate on dilution. This difficulty can sometimes be overcome by a judicious selection of acetone-water ratio, taking care to duplicate this ratio in preparation of the standard (Step 4). Alternately, the product may be treated as described below under b, c, d or e.)

NOTE: With strong buffers such as the glycerophosphates, start with 1.0 gram of sample and 50 ml of water, omit acetone from both sample and standard (Step 4) and adjust to pH 3.5 by glass electrode.

<u>Product</u>	<u>Dissolve In</u>	<u>Add</u>	<u>Dilute to 50 ml with</u>
Aspirin	25 ml acetone	water	water
Aminoacetic Acid (Glycine)	25 ml water	4.0 ml 6% HOAc	acetone
Amobarbital	25 ml acetone	water	water
Benzoic Acid	25 ml acetone	water	water
Caffeine	20 ml hot H_2O	4.0 ml 6% HOAc	acetone
Caffeine Citrate	25 ml water	0.5 ml NH_4OH	acetone
Calcium Carbonate	20 ml 10% HCl	NH_4OH to pH 3.5	acetone
Calcium Chloride	20 ml water	4.0 ml 6% HOAc	acetone
Calcium Cyclamate	25 ml water	10% HCl to pH 3.5	acetone
Calcium Saccharin	25 ml water	10% HCl to pH 3.5	acetone
Chloramphenicol	25 ml acetone	4.0 ml 6% HOAc	water
Chlorothenyl- pyramine Citrate (Chlo- rothen Citrate)	25 ml water	1.0 ml 10% HCl	acetone

<u>Product</u>	<u>Dissolve In</u>	<u>Add</u>	<u>Dilute to</u> <u>50 ml with</u>
Maleic			
Anhydride	25 ml water	NH ₄ OH to pH 3.5	acetone
Phenacetin	25 ml acetone	4.0 ml 6% HOAc	water
Saccharin			
Insoluble	25 ml acetone	0.8 ml NH ₄ OH	water
Salicylamide	25 ml acetone	4.0 ml 6% HOAc	water
Salicylic Acid	25 ml acetone	1.0 ml NH ₄ OH	water
Sodium Acetate	20 ml water	4.25 ml 10% HCl	acetone
Sodium Benzoate	25 ml water	3.5 ml 10% HCl	acetone
Sodium Carbonate	15 ml water	10.0 ml 10% HCl	acetone
Sodium Cyclamate	25 ml water	10% HCl to pH 3.5	acetone
Sodium Saccharin	25 ml water	0.25 ml 10% HCl	acetone
Sodium			
Salicylate	25 ml water	1.0 ml 10% HCl	acetone
Sodium Seco-			
barbital	25 ml water	2.5 10% HCl	acetone
Sulfanilamide	25 ml acetone	4.0 ml 6% HOAc	water
Sulfathiazole	25 ml acetone	4.0 ml 6% HOAc	water
Sulfuric Acid	25 ml water	NH ₄ OH to pH 3.5	acetone
Thenylpyramine			
Fumarate (methyl-			
pyriline Fu-			
marate)	25 ml water	4.0 ml 6% HOAc	acetone
Thenylpyramine			
Hydrochloride			
(Methapyriline			
HCl)	25 ml water	4.0 ml 6% HOAc	acetone
Urea	25 ml water	4.0 ml 6% HOAc	acetone

- b. Insoluble Powders - (such as Acramer, Phthalyl Sulfathiazole, Succinyl Sulfathiazole, Theobromine):
Moisten 2.0 grams with 10 ml of 10% HCl in a 100 ml beaker, hold on the steam bath for 5 minutes, cool to about 20°C, let stand until the supernatant liquid is clear, filter by vacuum and wash the residue with a small amount of water.

To the filtrate add 1 drop of 0.1% p-nitrophenol, titrate to a faint yellow endpoint with NH₄OH, add 4.0 ml of 6% HOAc and 25 ml of acetone, and dilute to 50 ml with water.

- c. Oils (such as Methyl Salicylate): Extract 2.0 grams with 10 ml of 10% HCl in a 60 ml separatory funnel and filter the aqueous layer through a small circle of wet filter paper.

To the extract add 1 drop of 0.1% p-nitrophenol, titrate to a faint yellow endpoint with NH_4OH , add 4.0 ml of 6% HOAc and 25 ml of acetone, and dilute to 50 ml with water.

- d. Colored Products (such as MHA or Kalcitor Acid): Moisten 2.0 grams with 4 ml sulfuric acid in a 250 ml beaker. Hold the beaker with tongs and swirl until any initial reaction subsides. With constant swirling, heat to boiling over a moderate flame to carbonize. Remove the flame and add 30% hydrogen peroxide ((Superoxol) DROP-WISE from a 15 ml pipette, allowing the reaction to subside somewhat between drops. When there is noticeable diminution in the vigor of the reaction upon addition of hydrogen peroxide, heat cautiously to fumes, remove the heat and resume addition of hydrogen peroxide, cool the beaker in ice, add 10 ml of water and 1 drop 0.1% p-nitrophenol. Neutralize to a faint yellow endpoint with ammonium hydroxide, add 4 ml of 6% acetic acid and dilute to 50 ml with water.
- e. Volatile Products (such as Acetic Anhydride): Weigh 2.0 grams into a clean evaporating dish and take to dryness on the steam bath. Dissolve the residue in 4 ml of 6% acetic acid, add 25 ml acetone and dilute to 50 ml with water.

2. Add (an additional) 25 ml of acetone to the prepared solution.
3. Verify the pH by a spot test on pHyrion paper, adjusting to 3.5 with NH_4OH or 10% HCl if necessary.

DETERMINATION (cont'd)

4. Prepare a standard by pipetting exactly 1.00 ml of 10 mg/l standard lead solution for each 5 ppm of the heavy metals specification, into a second Nessler tube containing 20 ml of water, 4 ml of 6% acetic acid, and the same quantity of acetone as is in the sample solution. (NOTE: if a 1.00 gram sample was taken, use 0.5 ml of standard lead solution for each 5 ppm of the heavy metals specification).
5. Dilute the standard and sample solutions to 75-80 ml with water.
6. Add 20 ml of saturated hydrogen sulfide water to the standard and to the sample solution, dilute both solutions to 100 ml with water, mix well and note the time.
7. Exactly 10 minutes after mixing the solutions compare longitudinally over a white, well lighted background.
8. Record the sample as containing less, equal to, or more than the ppm of heavy metals represented by the comparison standard, depending on whether the sample solution is lighter, equivalent to, or darker than the standard solution.

SAFETY

Prepare hydrogen sulfide water in an efficient fume hood and avoid inhalation during use. Conduct wet digestions (Determination, step 1, part d) in an efficient fume hood and behind a safety shield.

DISCUSSION

This method is essentially equivalent to the USP limit test, with no provision for compensating for inherent sample color. It is not directly applicable to ferric and manganese glycerophosphates owing to complex precipitation and oxidation of the sulfides.

DISCUSSION (cont'd)

Although some indication is obtained for most metals of the copper and arsenic groups, the method is most sensitive to lead. Thus, for example, the absorptivity of an equivalent amount of copper sulfide is about half that for lead. Since the test is intended primarily for the detection of toxic metals, the inclusion of iron is avoided by conducting the test in acetate solution below pH 5.

Glassware used in this test should be rinsed with concentrated hydrochloric acid and distilled water before use. Even reagent grade acetone should be redistilled. Oxidizing agents, such as traces of chlorine in the diluent water used, will cause the precipitation of elemental sulfur and invalidate the results.

REFERENCES

1. Stewart, F. N. and Strode, C. W., J. Am. Pharm. Assoc., Sci. Ed., 41, 242 (1952)
2. U.S.P. SVII.

SCOPE

This method is intended for the detection of apparent heavy metal contamination beyond a designated amount in fine chemicals. It involves formation of the sulfides at pH 3.5 and visual comparison with a standard.

/jw
7/26/71

RSV 0015108

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.449
May 31, 1971
Page 1 of 4

MATERIAL

Acrylamide
Sodium Lauryl Sulfate

ANALYSIS REQUIRED

Iron

APPARATUS

- a. Usual laboratory apparatus and reagents including matched 50 ml short form Nessler tubes.
- b. Colorimeter such as the Bausch and Lomb Spectronic 20, with matched 1" test tube cells.

REAGENTS

- a. 15% Ammonium Thiocyanate - Dissolve 75 g of ammonium thiocyanate A.R. in water and dilute to 500 ml.
- b. 100 mg/l Iron - Clean primary standard iron wire with emery paper and cleansing tissue, holding with tweezers. Accurately weigh 0.1000 g into a 1000 ml volumetric flask, add 10 ml HCl, heat on a steam bath until dissolved, add 1.0 g. ammonium persulfate, and dilute to volume with distilled water at room temperature.
- c. 5 mg/l Iron - Pipette 25 ml of 100 mg/l iron into a 500 ml volumetric flask and dilute to volume with water.
- d. Acetone, A.R. Grade

CALIBRATION

1. Pipette 0, 5, 10, 20, 25, and 30 ml aliquots of 5 mg/l iron solution into separate 50 ml glass stoppered cylinders, add 2 ml of conc. HCl and dilute each solution to about 35 ml with distilled water. (After treatment and dilution to volume, these solutions will represent concentrations of 0, 0.5, 1.0, 1.5, 2.0, 2.5, and 3.0 mg/l.)

RSV 0015109

CALIBRATION (cont'd)

2. To each solution add 0.05 g of ammonium persulfate, 5.0 ml of 15% ammonium thiocyanate, and 5 ml of acetone; then dilute to volume with water and mix well.
3. Using 1" test tube cells in the Spectronic 20 colorimeter and a wavelength setting of 475, immediately measure the absorbance of the 1.5 mg/l solution against the 0, 0.5, and 1.0 mg/l solutions; measure the absorbances of the 2.0, 2.5, and 3.0 mg/l solutions against the 1.5 mg/l solution. -
4. Plot the absorbances against concentrations on coordinate paper to obtain a "V" shaped curve with the apex at 1.5 mg/l. Divide (1.5 - respective solution concentration) by the absorbances of the first three solutions, divide (respective solution concentration - 1.5) by the absorbances of the last three solutions, and average the results as calibration constant K.

NOTE: For JFQ Lab Spectronic 20 Colorimeter, Serial No. RD 8233, calibration of 12-2-55, $K = 4.7$.

DETERMINATION

A-1 - Spectrophotometrically

1. Prepare 35 ml of clear, colorless, oxidized extract or solution of the sample or ash containing the equivalent of 2 ml of free conc. HCl in a 50 ml glass stoppered cylinder, as directed for the product being tested.

NOTE: The optimum sample size in grams is between 25 and 125 divided by the expected iron content in ppm.
2. Add 0.5 g of ammonium persulfate, 5.0 ml of 15% ammonium thiocyanate, and 5 ml acetone; then dilute to 50 ml with water and mix well.
3. Concurrently prepare a separate mixture of 15.0 ml of 5 mg/l iron solution, 2 ml of conc. HCl, 0.05g of ammonium persulfate, 5.0 ml of 15% ammonium thiocyanate, and 5 ml acetone; then dilute to 50 ml with water and mix well.

RSV 0015110

DETERMINATION

4. Immediately measure the absorbance of the darker of the two solutions against the lighter solution in the Spectronic 20 colorimeter, using 1" test tube cells and a wavelength setting of 475.

5. Calculate:

$$\text{ppm Fe} = \frac{50 (1.5 + KA)}{\text{sample wt. in solution}}$$

$$(\% - \text{ppm}/10,000)$$

where:

K = the calibration constant, which has a positive value when the sample solution is darkest and a negative value when the 1.5 mg standard solution is darkest.

and

A = the absorbance found

A-2 - Visually

1. Prepare 35 ml of clear, colorless oxidized extract or solution of the sample or ash containing the equivalent of 2 ml of free conc. HCl in a 50 ml Nessler tube, as directed for the product being tested.
2. Add 0.05 g of ammonium persulfate and 5.0 ml of 15% ammonium thiocyanate, and 5 ml of acetone.
3. To a blank of 35 ml of water and 2 ml of conc. HCl in a second Nessler tube add 0.05 g of ammonium persulfate, 5.0 ml of 15% ammonium thiocyanate and 5 ml of acetone.
4. Using a Mohr pipette, titrate the blank with 5 mg/l iron solution, observing the tubes transversely over a white background, until the color of the blank matches that of the sample.

DETERMINATION (cont'd)

5. Calculate:

$$\text{ppm Fe} = \frac{\text{ml iron solution} \times 5}{\text{sample wt. in solution}}$$

(% = ppm/10,000)

DISCUSSION

In concentrations between approximately 0.1 and 100 ppm of iron in various materials the precision and accuracy of this method may vary from 5 to 10% of the amount present in the absence of interferences.

Fluoride, pyrophosphate and oxalate interfere, as may colored ions, metals such as mercury or zinc which react with thiocyanate, or substances present in sufficient amount to exhibit color before the addition of thiocyanate. High salt concentration may inhibit color formation, and strong light may cause fading. The color is not stable for more than 30 minutes.

SCOPE

This method is intended for the determination of small amounts of iron in aqueous solutions of soluble colorless materials, ashes or extracts. It involves formation of the ferric thiocyanate complex and colorimetric estimation of the concentration.

SAFETY

With proper regard to observance of handling instructions on bottles of ammonium persulfate, no unusual hazards are involved in this operation.

REFERENCE

Snell and Snell, "Colorimetric Methods of Analysis", Vol. II, D. Van Nostrand Company, Inc., New York, 1954.

/jw
7/26/71

RSV 0015112

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.626 (M.C.)
April 26, 1962
Page 1 of 7

MATERIAL

Vinyl Chloride Monomer (VCM)

ANALYSIS

Heavy-end Impurities:

MVA

Ethyl Chloride

2-CPY

Vinylidene Chloride (1, 1-DCY)

t-1,2-DCY

CB

SAMPLE

Liquid Sample in VCM Bomb

APPARATUS

P. E. Model 154-C Vapor Fractometer equipped with gas inlet valve, or equivalent, valved for multi-column operation (See attached drawing.)

Hydrogen Flame Ionization Detector - capable of detecting 2×10^{-10} moles of chlorinated compounds.

Meeco Vaporizing Valve, Type 1S, or equivalent, for flashing liquid VCM into gas sample loop of chromatograph.

Flow-meter for determining rate of purge on sample loop.

Recorder - L. & N. Speedomax Type G, or equivalent, with 1 second pen and 5 mv full-scale deflection. Recorder should be modified to switch to either the thermistor detector or the flame detector.

RSV 0015113

REAGENTS AND SUPPLIES

Aluminum Tubing - 3/16 inch O.D.
Helium - 2 regulated supplies.
Chromosorb - red, 30/60 mesh.
Narcoil 40 - dinonyl phthalate from National Research Corp.
Hydrogen - electrolytic grade.
Air - well-filtered, clean supply.
Ethyl chloride - commercial grade.
2-Chloropropene - Columbia Organic Chemical Co., purified
by gas chromatography.
2-Chloro-1,3-butadiene - Monomer-Polymer, Inc., purified
by gas chromatography.
1,1-Dichloroethylene - Monomer-Polymer, Inc., purified
by gas chromatography.
Vinylacetylene - E. I. DuPont Co., purified by gas
chromatography.
Vinyl Chloride - Free of, or very low concentration of,
impurities.
Acetone - for column packing preparation.

PROCEDURE

A. Column Preparation

1. Mix 15 grams of Narcoil 40 with approximately 100 ml of acetone.
2. Pour mixture over 85 grams of 30/60 mesh chromosorb and add more acetone if necessary to make a slurry.
3. Mix continuously while evaporating the acetone. If available, use rotary-film evaporator with baffled flask for mixing. Apply heat from hot air gun and vacuum to flask to aid in evaporation.

RSV 0015114

PROCEDURE (cont'd)

4. When the acetone is evaporated, run screen to remove fines and pack one 20-foot length and one 35-foot length of 3/16-inch tubing with the 30/60 mesh fraction.
5. Install the columns in the instrument as shown on the attached drawing.
6. Adjust the two helium pressures and restrictor valve to give the following flows:
 - a. Column I: 75 ml/minute.
 - b. Column II (with restrictor in flow path): 85 ml/min.
 - c. Columns I & II (in backflush position): 85 ml/min.
7. Adjust temperature to 50° C. and allow instrument and columns to reach equilibrium.
8. Light the hydrogen burner and adjust the hydrogen and air flows so as to yield maximum signal-to-noise ratio when MVA passes through it.
9. Turn on amplifier and recorder and set for maximum sensitivity. The instrument should show no drift for at least five minutes before it is ready for use.

B. Analysis

1. Purge the 5 cc sample loop for at least one minute with sample flashed through the Meco Vaporizing Valve, then inject the sample into the flowing carrier stream with the columns in Position A and the recorder turned to the thermal conductivity detector. The attenuator switch of the detector should be set on an attenuation of 2. The valve between the thermal detector and the flame detector should be in the vent position to prevent VCM from going to the flame detector.
2. Allow the sample to flow in Column I as shown in Position A for approximately three minutes (valve should be switched when recorder pen reaches 50% of full scale coming down on VCM peak). This allows the components lighter than VCM and most of the VCM to be vented from Column I.

PROCEDURE (cont'd)

3. Just before MVA is eluted from Column I, switch the valve to Position B. This allows everything behind the VCM peak to be backflushed onto Column II where they are separated and then flow to the detector. Sometime after backflushing and prior to the backflushed VCM being eluted from Column II, switch the recorder to the flame detector and the valve between the thermal detector and the flame detector to the position such that the flow is through the flame detector.
4. Record all peaks eluted from Column II. A typical chromatogram is attached.

CALCULATIONS

The concentration of each component is calculated as follows:

1. Measure the peak height and width at one-half the peak height using a scale of 50 divisions to 1 inch.
2. Calculate the area of each peak according to the following equation:

area =

peak height x peak width at one-half peak height x attenuation

3. Calculate the concentration of each component as follows:


Concentration (ppm) = peak area x factor

PREPARATION OF REAGENT SOLUTIONS

None necessary.

CALIBRATION



1. Prepare synthetics of MVA, ethyl chloride, 2-CPY, 1,1-DCY, trans-1,2-DCY, and CB in VCM to cover the concentration range of interest. At least two synthetics should be prepared in each concentration range. Synthetics should be prepared in clean, stainless steel bombs of any convenient size. The pressure should be sufficient to permit extended use of the mixtures. The blending system used must be completely free of leaks and all kinds of foreign materials.




- V₁ - WALL 8-PORT VALVE
- V₂ - CIRCLE SEAL DOUBLE 4-PORT VALVE
- V₃ - 6-WAY GAS CHANGING VALVE
- V₄ - NEEDLE VALVE
- C₁ - COLUMN I
- C₂ - COLUMN II
- R_{1,2} - HELIUM REGULATORS
- G_{1,2} - PRESSURE GAUGES, 0-600 IF
- V₅ - CIRCLE SEAL 4-PORT VALVE

METHOD 602.626 (G.O.)
April 26, 1962
Page 6 of 7

Figure II



METRO 602.626 (A.C.)
 April 26, 1962
 Page 7 of 7



RSV 0015119

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.628 (G.C.)
June 18, 1962
Page 1 of 7

MATERIAL

Vinyl Chloride Monomer (VCM).

ANALYSIS

Light-end Impurities:

C₁ - C₄ Hydrocarbons
Methyl Chloride

SAMPLE

Liquid sample in VCM bomb.

APPARATUS

Gas chromatograph equipped with valve(s) for multi-column operation and two regulated helium supplies. (See Note)

Hydrogen flame ionization detector, Barber-Colman Model A-4504 or equivalent: sensitivity should be approximately one coulomb per mole.

Amplifier and Recorder: Over-all sensitivity should be sufficient to yield full-scale deflection for 1×10^{-11} amp with a noise level of less than 1×10^{-13} amp. Barber-Colman Electrometer for above flame detector is recommended.

Integrator: Perkin-Elmer Printing Integrator, or equivalent.

Tubing: 3/16" O.D. Aluminum.

Pressure Regulators: Two for helium, one for hydrogen, and one for air.

Meeco Vaporizing Valve, Type 1S, or equivalent, for flashing liquid VCM into gas sample loop of chromatograph.

REAGENTS

Helium - Two regulated supplies.

Chromosorb - red, 30/60 mesh.

Trioctyl Phosphate: Union Carbide Chemical Company.

RSV 0015120

REAGENTS (cont'd)

Tributyl Phosphate: Fisher Scientific.

C₁ - C₄ Hydrocarbons: Phillips Research Grade. •

Methyl Chloride: The Matheson Company.

Vinyl Chloride Monomer: Free of, or very low concentration
of, impurities.

Hydrogen: Electrolytic Grade.

Air: Purified source, such as cylinder breathing air.

Acetone: For column packing preparation.

PROCEDURE

A. Column and Instrument Preparation

1. Prepare two columns as follows:
 - a. Column I: 20 feet of 15% w/w Trioctyl Phosphate on 30/60 mesh red chromosorb packed in 3/16" O.D. tubing.
 - b. Column II: 45 feet of 15% w/w Tributyl Phosphate on 30/60 mesh red chromosorb packed in 3/16" O.D. tubing.
2. Install the columns in the gas chromatograph as shown in Position A and adjust the helium pressure on Column I to give a flow of 100 cc/minute. When the valve is switched to Position B, the flow through Columns I and II should be 100 cc/minute. The columns are operated at room temperature.
3. Light the hydrogen burner and adjust the hydrogen and air flows so as to yield maximum signal-to-noise ratio when a C₃ or C₄ hydrocarbon passes through it.
4. Turn on amplifier and recorder and set at maximum sensitivity. The instrument should show no drift for at least five minutes before it is ready for use.
5. Prepare the integrator for use.

PROCEDURE (cont'd)

b. Analysis

The analysis of light-ends may be made using the following procedure:

1. Line-out instrument in Position B, then switch to Position A for entering sample.
2. Purge the 5 cc sample loop for at least one minute with sample flashed through the Meeco Vaporizing Valve, then inject the sample into the flowing carrier stream with the columns in Position A.
3. Fifteen (15) seconds after entering the sample, switch the valve to diagram Position B. This allows the entire sample to be backflushed onto Column II where the components are separated and then flow to the detector. Using the integrator, record areas for all peaks eluted ahead of VCM.
4. When VCM begins to be eluted, switch the valve to diagram Position A. This allows the VCM to be vented instead of passing through the detector.

CALCULATIONS

The concentration of each light-end component is calculated as follows:

$$\text{ppm} = \text{attenuation} \times \text{peak area} \times \text{factor}$$

PREPARATION OF REAGENT SOLUTIONS

None necessary.

CALIBRATION

1. Prepare synthetics of the light-end impurities, hydrocarbons and methyl chloride, in VCM to cover the concentration range of interest. At least two synthetics should be prepared. Synthetics should be prepared in clean stainless steel bombs of any convenient size. The pressure should be sufficient to permit extended use of the mixtures. The blending system must be clean and completely free of leaks.

CALIBRATION (cont'd)

2. Analyse synthetic mixtures as described under the Procedure for Analysis.
3. Calculate factors for each of the light-end impurities according to the following equation:

$$\text{factor} = \frac{\text{concentration of component in synthetic}}{\text{integrated area} \times \text{attenuation}}$$

SAFETY PRECUATIONS

VCM is a toxic, volatile, and highly flammable material and should be treated as such. Sample used to purge sample loop should be vented outside of building.

ANALYTICAL TIME

The average analytical time for this analysis should be 0.5 hour.

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES FOR REPORTING RESULTS

95% confidence limits using this method is $\pm 10\%$ of the amount present.

Report results to nearest 0.1 ppm.

SCOPE

This method is designed for light-end impurities in VCM; however, it may be used for the analysis of light-end impurities in volatile pure grade hydrocarbons as long as the impurities are hydrocarbons in the $C_1 - C_4$ range.

REFERENCE

Monsanto Research Notebook - No. 2330, pp. 154660-154665.

NOTE:

Valving recommended for this procedure is as follows:

- a. Backflushing and column switching may be done with a Wall 8-Port Valve purchased from Tycam Engineering and Manufacturing Co., Inc., Houston, Texas. A drawing (may be a rough sketch) should be furnished to indicate number and location of channels in teflon slide.


REVISION 602.628 (G.C.)
June 15, 1962
Page 5 of 7

NOTE (cont'd)

- b. A double 4-Port Circle Seal Valve may be used for the control of air to the air actuator of the Well Valve as well as switching the flow of the second helium supply when backflushing. (See attached drawing.)



/jw
7/29/71

RSV 0015124



METHOD 602.628 (G.C.)
June 18, 1952
Page 6 of 7

FIGURE I



Sample: Vinyl Chloride

Column 1: 20 feet of Trioctyl Phosphate
on 30/60 mesh red Chromasorb
packed in 1/4" O.D. tubing.

II: 45 feet of Tributyl Phosphate
on 30/50 mesh red Chromasorb
packed in 3/16" O.D. tubing.

Temperature: Ambient.
Carrier Flow: 100 cc/min.
Detector: Hydrogen Flame.
Sample Size: 5 cc Gas.

MONSANTO COMPANY
TEXAS CITY, TEXAS

DEPT. 60
METHOD 602.639 (MC)
May 28, 1964
Page 1 of 3

MATERIAL

Product vinyl chloride monomer

ANALYSIS REQUIRED

Phenol

SAMPLE SIZE

Quart thermos jug

APPARATUS

250 ml Erlenmeyer flask

50 ml beaker

2 ml pipette

100 ml graduate

Beckman Model "B" Spectrophotometer

50 mm absorption cells (matched pair)

REAGENTS

Millons reagent -

Prepare in well ventilated hood.

PROCEDURE

Evaporate 100 mls of sample in a 250 ml Erlenmeyer flask. When sample is completely evaporated, add 25 ml distilled H₂O. Pipette 2 ml of Millons Reagent into flask. Boil for 1 minute on hot plate. Prepare a blank in the same manner, omit the sample. When solution has cooled, filter and read sample against blank on Beckman Model "B". Use 59 mm cells, sensitivity of 3, and wavelength of 440 mμ.

RSV 0015127

CALCULATIONS

Multiply O.D. by 10 to obtain ppm phenol. Report to nearest 0.1 ppm.

Note:

Use of sensitivity of 3 is not normal for our procedures. Use setting of 3 to adjust instrument on blank. This will allow the use of much smaller slit and will give some improvement in reproducibility of results.

PREPARATION OF REAGENT SOLUTIONS

To 100 grams of fuming nitric acid, add 100 grams of mercury. Allow sufficient time for mercury to react and go completely in solution. Add 20 mls of distilled H_2O while stirring mixture. Place in a glass stoppered container.

STANDARDIZATION AND/OR CALIBRATIONS

See Solutions Manual.

SAFETY PRECAUTIONS

Handle mercury with care, avoid breathing vapors and dust spilled mercury with sulfur. Mercury reacts vigorously with nitric acid liberating large quantities of toxic NO_2 . Perform the preparation of this reagent in the hood with good ventilation.

ANALYTICAL TIME

.2 hour

METHOD 602.639 (WC)

May 28, 1964

Page 3 of 3

PRECISION AND SIGNIFICANT FIGURES

95% confidence limit ± 1 ppm at 5 ppm level.
Report to nearest 1 ppm level.

SCOPE

Mercuric nitrate reacts with the phenol in the sample to produce a color complex.

/jw
7/29/71

RSV 0015129

MONSANTO CHEMICAL COMPANY
ORGANIC DIVISION
TEXAS CITY, TEXAS

DEPT. 33, 60
METHOD 602.657 (W.C.)
February 4, 1964
Page 1 of 4

MATERIAL

Caustic receipts

ANALYSIS REQUIRED

Iron

SAMPLE

1 pint

APPARATUS

Beckman "B" spectrophotometer

50 mm matched cells

5 ml pipettes

100 ml volumetric flasks

Analytical balances

Phydrion paper 4.0-6.0 range

Small weighing bottle

REAGENTS AND SOLUTIONS

1. 0.15% 1, 10 - phenanthroline (aqueous)
2. 10% hydroxylamine hydrochloride (aqueous)
3. Iron standard 100 μ gm/ml (aqueous)

RSV 0015130

4. A. R. grade ammonium hydroxide
5. .5N NaOH
6. .5N HCl
7. 1-1 HCl solution (iron free)

PROCEDURE

1. Fill a small weighing bottle with sample and weigh.
2. Place 2 gm of sample (40-50 drops) in a 100 ml beaker and dilute with 10-15 ml of water. Reweigh weighing bottle to obtain sample weight.
3. Add 10 ml of 1-1 HCl to the beaker and boil down until approximately 5 ml remains.
4. Pipette 5 ml of 10% NH_4OH , HCl and 5 ml of .15% O-phenanthroline to the beaker. Then pipette 4 ml of concentrated NH_4OH into the mixture.
5. Adjust the pH to 6-7 using phydron paper with .5N HCl or .5N NaOH.
6. Transfer contents of the beaker to a 100 ml volumetric, rinsing the beaker several times to insure removal of all the sample. Adjust the volume with distilled water and mix.
7. Prepare a reagent blank.
8. Wait 20 minutes and then determine the absorbance of the sample at 510 wavelength in a 50 mm cells against a reagent blank. Obtain net absorbance.

CALCULATIONS

$$\text{PPM Fe} = \frac{\text{abs.} \times 526}{\text{Weight of sample} \times 5}$$

The 5 in the calculation denotes the cell length in CM.

PREPARATION OF REAGENTS

1. 0.15% 1, 10 phenanthroline - dissolve .75 μ gm of reagent grade material in 5 ml of glacial acetic acid. Dilute to 500 ml with distilled water.
2. 10% hydroxylamine HCl - dissolve 50 μ gms of AR grade material in 500 ml of distilled water.
3. Iron standard 100 μ gm/ml - weigh .1000 gms of standard iron wire. Dissolve in 10 ml HCl, dilute to 1 liter.
4. Iron standard 20 μ gm/ml - dilute 20 ml of the standard 100 μ gm/ml to 100 ml iron free distilled water.
5. Iron standard 5 μ gm/ml - dilute 5 ml of the 100 μ gm/ml standard to 100 ml with iron free distilled water.

STANDARDIZATION AND/OR CALIBRATIONS

Pipette 5, 10, 15, and 20 ml of the 5 μ gm/ml standard to each of 4-100 ml volumetric flasks. These will contain 25, 50, 75 and 100 μ gm/ml of iron. Similarly pipette 10, 15, 20 and 25 ml of the 20 μ gm/ml standard to each of 4-100 volumetric flasks. These will contain 200, 300, 400 and 500 μ gm of Fe. Prepare a reagent blank. Analyze following the method in the procedure. Plot the absorbances found against the corresponding μ gm Fe. Draw in the best straight line through the points obtained. Re-check any points that are over 10% from the line value.

From this line, find the μ gm of Fe corresponding to unit absorbance.

The factor 526 μ gm Fe/unit absorbance was determined by reading the absorbance at 500 μ gm from the curve and dividing 500 by the unit absorbance at this point (.95 abs.)

PRECISION AND SIGNIFICANT FIGURES

95% confidence limits are $\pm$ 0.1 PPM at a level of 1.05 PPM. Report results to two significant figures.

ANALYTICAL TIME

0.5 hours

SAFETY

Use normal precautions in handling laboratory chemicals, glassware and electrical equipment.

Strong caustic can cause severe burns if spilled on the skin. If such a spill does occur, wash with large amounts of water and report to first aid.

SCOPE

Iron is reduced to the divalent state with hydroxylamine at a pH 4 to 7. Fe(II) forms a red complex with 1.10 orthophenthroline. This red color is analytical W.R.T. iron.

REFERENCES

J.F.Q. General Method 70-F

/jw
7/26/71

RSV 0015133

MONSANTO CHEMICAL COMPANY
ORGANIC DIVISION
TEXAS CITY, TEXAS

DEPT. 33.60
METHOD 602.658
February 4, 1964
Page 1 of 3

MATERIAL

Caustic Receipts

ANALYSIS REQUIRED

Sodium hydroxide and sodium carbonate

SAMPLE

1 pint bottle

APPARATUS

50 ml acid burette

Small weighing bottle with dropper attached

250 ml Erlenmeyer flask

Analytical balance

REAGENTS

0.5 N standard HCl

Distilled Water

Phenolphthalein Indicator

Methyl Orange Indicator

PROCEDURE

1. Fill a small weighing bottle with the sample and weigh. This sample bottle can be used for all subsequent analyses of the caustic.
2. Place approximately 50 ml of distilled water in a 250 ml Erlenmeyer flask and add 2 gms of sample from the weighing bottle (approx. 40-50 drops). Reweigh to obtain the weight of sample.

RSV 0015134

3. Add 2-3 drops of phenolphthalein indicator and titrate with 0.5 N HCl until the solution is just colorless. Record this volume of HCl as V_1 .
4. Refill HCl burette to zero. Add 2-3 drops of methyl orange indicator to the titration flask and titrate with 0.5 HCl until the solution is red orange. Record this volume of HCl as V_2 .

CALCULATIONS

$$\frac{(V_1 - V_2) (N_{\text{HCl}}) 4.0}{\text{Weight of Sample}} = \% \text{ NaOH}$$

$$\frac{(2V_2) (N_{\text{HCl}}) 5.3}{\text{Weight of Sample}} = \% \text{ Na}_2\text{CO}_3$$

PREPARATION OF REAGENTS

Prepare indicator solutions and normal acid solution outlined in Reagent Manual.

STANDARDIZATION AND/OR CALIBRATIONS

Standardize HCl solution according to Reagent Manual.

SAFETY

Use care in handling strong caustic solutions, as caustic can cause severe burns if spilled on the skin. If a spill should occur, flush the area with large amounts of water and report to first aid. Use normal precautions with chemicals, glassware, and electrical equipment.

ANALYTICAL TIME



0.5 hours

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES


No precision has been determined, but duplicate analysis should agree within 0.5%.

SCOPE

In titrating to the phenolphthalein end point, the following two reactions occur:



The subsequent titration to the methyl orange end point completes the titration of the NaHCO_3 by the following reaction.



REFERENCES

Method 22-10-1, Revised 3/17/53

GK/mlw

RSV 0015136

MONSANTO CHEMICAL COMPANY
ORGANIC DIVISION
TEXAS CITY, TEXAS

DEPT. 33,60
METHOD 602.659
February 4, 1964
Page 1 of 2

MATERIAL

Caustic Receipts

ANALYSIS

Chlorides

SAMPLE

1 pint

APPARATUS

Small weighing bottles

250 ml beaker

50 ml burette

pH meter

Magnetic stirrer with stirrer bar

Silver electrode

Hg-HgSO₄ - K₂SO₄ electrode

REAGENTS

Concentrated nitric acid

0.1 N silver nitrate solution

PROCEDURE

1. Weigh 2 gms of sample (40-50 drops) and place in a 250 ml beaker. Add 50 ml of distilled water and 5 ml of concentrated nitric acid.
2. Place on magnetic stirrer and immerse the electrodes in the solution.
3. While solution is being agitated, titrate with 0.1 AgNO₃ to the first break, which will be about 100 MV.

RSV 0015137

CALCULATIONS

$$\frac{(\text{Volume of AgNO}_3) (\text{Normality}) (3.55)}{\text{Weight of Sample}} = \% \text{ Cl}^-$$

Weight of Sample

PREPARATION OF REAGENTS

Prepare 0.1 AgNO₃ according to Reagents Manual

STANDARD AND/OR CALIBRATIONS

Standard AgNO₃ by procedure in Reagents Manual

SAFETY

Handle all chemicals, glassware, and electrical equipment with normal precautions.

Caustic can cause severe burns to the skin. Handle with caution and if a spill should occur, wash the area with large amounts of water and report to first aid.

ANALYTICAL TIME

0.5 hours

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES

No precision has been determined, but should be accurate to $\pm .01\%$. Report 2 significant figures.

SCOPE

Chlorides in an acid solution precipitate as AgCl. An excess of chlorides in AgNO₃ causes a sharp change in potential, which is taken as the end point of the titration.

REFERENCES

Method 20-29-50, Rev. 7/27/58 (Modified)

GK/mlw

RSV 0015138

MONSANTO CHEMICAL COMPANY
ORGANIC DIVISION
TEXAS CITY, TEXAS

DEPT. 33.60
METHOD 602.660 (W.C.)
February 4, 1964
Page 1 of 4

MATERIAL

Caustic Receipts

ANALYSIS REQUIRED

Sulfates

SAMPLE

1 pint

APPARATUS

100 ml beaker

Small weighing bottle

50 ml stoppered graduated cylinders

Pyrex hydrosol filter holders, Cat. #XX1002500

Millipore Filter Corp., Bedford, Massachusetts

Millipore filter discs., 25 mm dia. type AA 0.8

Scoop, glass, Cat. #1H-12. Laboratory Equipment Company

St. Joseph, Michigan

Beckman Model B

50 mm cells (matched)

REAGENTS

1. Barium chloride 20-30 mesh parr turbidimetric grade, J. T. Baker Chemical Company.
2. Dispersing reagent
3. Concentrated HCl
4. Phenolphthalein indicator

RSV 0015139

PROCEDURE

1. Weigh 2 gms of sample (40-50 drops from a weighing bottle) into a 100 ml beaker. Add 10 ml of distilled water and 1-2 drops of phenolphthalein indicator to the beaker. Carry a reagent blank along with the sample.
2. Drop-wise, add concentrated HCl until the solution is just colorless, then filter the sample through a millipore filter type AA. Wash the beaker several times with 5 ml portions of distilled water, passing the wash water through the filter after each washing. Transfer the contents of the filter flask to a 50 ml stoppered graduate cylinder, again washing the filter flask several times with 5 ml portions of distilled water.
3. Add 5 ml of dispersing reagent to the graduate and dilute to 50 ml with distilled water. Mix well. Then add 1 scoop (0.3 gm) of BaCl_2 and shake for 1 minute.
4. Allow the graduate to stand for five minutes, and then transfer to a 50 mm cell and read against the blank at 420 m μ on the Beckman "B."

CALCULATIONS

From the attached calibration chart, determine the amount of SO_4 present.

$$\frac{\mu\text{gm of SO}_4}{\text{Sample weight}} = \text{PPM SO}_4$$

PREPARATION OF REAGENTS

1. Dispersing reagent - Mix 25 ml concentrated HCl, 75 ml of deeminac water, 300 ml Cp glycerine, and 600 ml of 95% ethanol.

STANDARDIZATION AND/OR CALIBRATION

1. Measure .148 gm of AR Na_2SO_4 (Sodium sulfate) into a liter volumetric flask and dilute to a liter. This gives a 100 $\mu\text{gm/ml}$ of SO_4 solution.

2. Pipette 10, 30, 50, 70 ml of the 100 μ gm/ml standard solution into a 100 ml volumetric flask and dilute to volume. This gives a 10, 30, 50, 70 μ gm/ml solution.
3. Use 1 ml of each solution and carry through the standard procedure for SO_4 determination. Plot the absorbance for each sample against the μ gm SO_4 present. Draw a straight line connecting these values for an absorbance versus μ gm chart.

SAFETY

Use normal laboratory safety precautions in handling chemicals, glassware and electrical equipment.

Caustic can cause a severe burn if spilled on the skin. If a spill should occur, wash with large amounts of water and report to first aid.

ANALYTICAL TIME

0.5 hours

PRECISION OF ANALYSIS AND SIGNIFICANT FIGURES

Report to the nearest PPM SO_4 . 95% confidence limits are ± 5 PPM at a 50 PPM level.

SCOPE

The SO_4 present reacts with Ba to produce a turbidity, which is read as absorbance on the Beckman "B."

REFERENCES

Method 602.503 (W.C.) (Modified)

Micrograms of SO_4 , versus absorbance on
Beckman Model "B" spectrophotometer at
420 millimicrons with a 50 mm cell.

METHOD 602.660 (W.C.)
February 4, 1964
Page 4 of 4

Absorbance	$\mu\text{gm SO}_4$	Absorbance	$\mu\text{gm SO}_4$
.002	3	.049	102
.003	6	.054	117
.004	9	.059	126
.005	12	.062	135
.007	15	.068	144
.009	18	.075	153
.010	21	.080	162
.012	24	.085	171
.013	27	.090	180
.014	30	.095	189
.015	33	.101	198
.017	36	.106	207
.019	39	.110	216
.020	42	.115	225
.021	45	.120	234
.022	48	.125	243
.023	51	.130	252
.024	54	.135	261
.026	57	.140	270
.027	60	.146	279
.028	63	.155	288
.030	66	.160	297
.031	69	.166	306
.032	72	.172	315
.034	75	.180	324
.035	78	.186	333
.036	81	.194	342
.038	84	.199	351
.040	87	.209	366
.042	90	.219	381
.043	93	.230	396
.044	96	.240	411
.045	99	.251	426

GK/mlw

RSV 0015142

008

RSV 0015143

TANK CALIBRATION TABLE

TANK NUMBER: 60N2 Page 1 of 2
 TANK NAME: Caustic-Finished Latex Mix Tank MMS/ff
 TANK CONTENTS: Caustic & Hydrolyzed Latex Rev. 2/25/72
 METHOD OF GAUGING: LI 214
 TANK DIMENSIONS: 4' OD x 6'6"
 TYPE: Vertical
 HEADS: Bottom Flat; Top Dished
 EQUIPPED WITH: An agitator
 TRANSMITTER RANGE: 82.5 inches H₂O
 SP. GR.: 1.145 GALLONS/INCH: 7.83

Instrument Reading(%)	Gallons		Instrument Reading(%)	Gallons	
0	47		37	256	
1	53		38	261	
2	58		39	267	
3	64		40	273	
4	70		41	278	
5	75		42	284	
6	81		43	290	
7	86		44	295	
8	92		45	301	
9	98		46	307	
10	103		47	312	
11	109		48	318	
12	115		49	323	
13	120		50	329	
14	126		51	335	
15	132		52	340	
16	137		53	346	
17	143		54	352	
18	149		55	357	
19	154		56	363	
20	160		57	369	
21	165		58	374	
22	171		59	380	
23	177		60	386	
24	182		61	391	
25	188		62	397	
26	194		63	402	
27	199		64	408	
28	205		65	414	
29	211		66	419	
30	216		67	425	
31	222		68	431	
32	228		69	436	
33	233		70	442	
34	239		71	448	
35	244		72	453	
36	250		73	459	

TANK CALIBRATION TABLE

TANK NUMBER: 60N2			Page 2 of 2		
TANK NAME: Caustic-Finished Latex Mix Tank			MMS/ff		
TANK CONTENTS: Caustic & Hydrolyzed Latex			Rev. 2/25/72		
METHOD OF GAUGING: LI 214					
TANK DIMENSIONS: 4' OD x 6'6"					
TYPE: Vertical					
HEADS: Bottom Flat; Top Dished					
EQUIPPED WITH: An agitator					
TRANSMITTER RANGE: 82.5 inches H ₂ O					
SP.GR.: 1.145			GALLONS/INCH: 7.83		

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
74	464	5.6			
75	470				
76	476				
77	481				
78	487				
79	493				
80	498				
81	504				
82	510				
83	515				
84	521				
85	527				
86	532				
87	538				
88	543				
89	549				
90	555				
91	560				
92	566				
93	572				
94	577				
95	583				
96	589				
97	594				
98	600				
99	605				
100	611				

RSV 0015145

TANK CALIBRATION TABLE

TANK NUMBER: 60S1-1 & 2
 TANK NAME: Blowdown Separator #1 & #2
 TANK CONTENTS: Latex
 METHOD OF GAUGING: LI 201-1&2
 TANK DIMENSIONS: 11' dia. x 15' SS
 TYPE: Vertical
 HEADS: Flanged & Dished
 EQPT. WITH: Impingement Plate
 TRANSMITTER RANGE: 16.5 to 80.5 inches H₂O
 SP. GR.: 1.10 GALLONS/INCH: 59

Page 1 of 1
 MMS/ff
 Rev. 2/25/72

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	0		69	1890	34.4
37.5-Top of Dish	807	807	70	1924	
38	824	17	71	1959	
39	858	34.4	72	1993	
40	892		73	2028	
41	927		74	2062	
42	961		75	2096	
43	996		76	2131	
44	1030		77	2165	
45	1064		78	2200	
46	1099		79	2234	
47	1133		80	2268	
48	1168		81	2303	
49	1202		82	2337	
50	1236		83	2372	
51	1271		84	2406	
52	1305		85	2440	
53	1340		86	2475	
54	1374		87	2509	
55	1408		88	2544	
56	1443		89	2578	
57	1477		90	2612	
58	1512		91	2647	
59	1546		92	2681	
60	1580		93	2716	
61	1615		94	2750	
62	1650		95	2784	
63	1684		96	2819	
64	1718		97	2853	
65	1752		98	2888	
66	1787		99	2922	
67	1821		100	2956	
68	1856				

RSV 0015146

TANK CALIBRATION TABLE

TANK NUMBER: 60S2			Page 1 of 1		
TANK NAME: THF Recovery Column Still Pot			MMS/ff		
TANK CONTENTS: THF/Water/Latex Solids			Rev. 2/16/72		
METHOD OF GAUGING: LIC 305					
TANK DIMENSIONS: 6'4" I.D. x 7'8" SS					
TYPE: Vertical					
HEADS: Bottom double w/one inverted; Top F&D					
EQPT. WITH: Steam Jacket and Agitator					
TRANSMITTER TYPE: Displacer					
SP. GR.: 1.0 GALLONS/INCH: 17.6					
Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	88	—	52	748	25.4
2	113	25.4	54	773	
4	138		56	799	
6	164		58	824	
8	189		60	850	
10	215		62	875	
12	240		64	900	
14	265		66	926	
16	291		68	951	
18	316		70	977	
20	342		72	1002	
22	367		74	1027	
24	392		76	1053	
26	418		78	1078	
28	443		80	1104	
30	469		82	1129	
32	494		84	1154	
34	519		86	1180	
36	545		88	1205	
38	570		90	1231	
40	596		92	1256	
42	621		94	1281	
44	646		96	1307	
46	672		98	1332	
48	697		100	1358	
50	723				

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60S3
 TANK NAME: THF Separator
 TANK CONTENTS: THF
 METHOD OF GAUGING: LIC 304
 TANK DIMENSIONS: 30" OD x 39"
 TYPE: Vertical
 HEADS: Flanged and Dished
 EQPT. WITH:

Page 1 of 2
 JSM/ff
 7/4/71

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	21.7	.43	35	36.8	.43
1	22.2		36	37.2	
2	22.6		37	37.6	
3	23.0		38	38.1	
4	23.4		39	38.5	
5	23.9		40	38.9	
6	24.3		41	49.4	
7	24.7		42	39.8	
8	25.2		43	40.2	
9	25.6		44	40.6	
10	26.0		45	41.1	
11	26.4		46	41.5	
12	26.9		47	41.9	
13	27.3		48	42.4	
14	27.7		49	42.8	
15	28.1		50	43.2	
16	28.6		51	43.6	
17	29.0		52	44.1	
18	29.5		53	44.5	
19	29.9		54	44.9	
20	30.3		55	45.4	
21	30.8		56	45.8	
22	31.2		57	46.2	
23	31.6		58	46.7	
24	32.0		59	47.1	
25	32.5		60	47.5	
26	32.9		61	47.9	
27	33.3		62	48.4	
28	33.8		63	48.8	
29	34.2		64	49.2	
30	34.6		65	49.7	
31	35.0		66	50.1	
32	35.5		67	50.5	
33	35.9		68	51.0	
34	36.3		69	51.4	

RSV 0015148

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60S3 (continued)Page 2 of 2
JSM/ff
7/4/71

<u>Instrument Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	51.8	.43	86	58.7	.43
71	52.3	↓	87	59.1	↓
72	52.7		88	59.6	
73	53.1		89	60.0	
74	53.5		90	60.4	
75	54.0		91	60.9	
76	54.4		92	61.3	
77	54.8		93	61.7	
78	55.3		94	62.1	
79	55.7		95	62.6	
80	56.1		96	63.0	
81	56.6		97	63.4	
82	57.0		98	63.9	
83	57.4		99	64.3	
84	57.8		100	64.7	
85	58.3				

RSV 0015149

TANK CALIBRATION TABLE

TANK NUMBER: 60T1-1 & 2
 TANK NAME: Vinyl Chloride Storage
 TANK CONTENTS: Vinyl Chloride
 METHOD OF GAUGING: LI 111-1&2
 TANK DIMENSIONS: 10' dia x 40'
 TYPE: Horizontal
 HEADS: Elliptical
 TRANSMITTER RANGE: 98 Inches H₂O
 SP. GR.: .919

Page 1 of 2
 MMS/ff
 Rev. 2/25/72

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	453		35	8218	280
1	572	124	36	8499	281
2	712	135	37	8782	283
3	857	145	38	9066	284
4	1011	154	39	9353	287
5	1171	160	40	9640	287
6	1336	165	41	9928	288
7	1506	170	42	10216	288
8	1686	180	43	10504	288
9	1876	190	44	10793	289
10	2071	195	45	11083	290
11	2271	200	46	11373	290
12	2476	205	47	11663	290
13	2686	210	48	11954	291
14	2901	215	49	12245	291
15	3122	221	50	12540	295
16	3345	223	51	12833	291
17	3572	227	52	13124	291
18	3801	229	53	13414	290
19	4034	233	54	13704	290
20	4274	240	55	13994	291
21	4517	243	56	14285	291
22	4763	246	57	14576	291
23	5012	149	58	14867	291
24	5263	254	59	15157	290
25	5520	257	60	15447	290
26	5780	260	61	15735	288
27	6042	262	62	16021	286
28	6306	264	63	16305	284
29	6573	267	64	16587	282
30	6840	267	65	16867	280
31	7113	273	66	17147	280
32	7387	274	67	17426	279
33	7662	275	68	17704	278
34	7938	276	69	17979	275

TANK CALIBRATION TABLE

TANK NUMBER: <u>60T1-1 & 2</u> (continued)			Page <u>2</u> of <u>2</u> MMS/ff Rev. 2/25/72		
<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	18254	275	86	22270	224
71	18527	273	87	22490	220
72	18798	271	88	22706	216
73	19067	269	89	22918	212
74	19334	267	90	23126	208
75	19598	264	91	23326	200
76	19858	260	92	23520	194
77	20114	256	93	23708	188
78	20368	254	94	23890	182
79	20618	250	95	24065	175
80	20866	248	96	24235	170
81	21110	244	97	24397	162
82	21350	240	98	24551	154
83	21586	236	99	24698	147
84	21818	232	100	24837	139
85	22046	228			

TANK CALIBRATION TABLE

TANK NO.: 60T2-1&2
 TANK NAME: THF Storage Tank 1 & 2
 TANK CONTENTS: Tetrahydrofuran
 METHOD OF GAUGING: LI 301-1&2
 TANK DIMENSIONS: 94.8" O.D. x 30' Str. Side
 TYPE: Horizontal Cylindrical
 HEADS: ASME Flanged & Dished
 TRANSMITTER RANGE: 75 inches H₂O
 SPECIFIC GRAVITY: 0.89

Page 1 of 2
 MMS/ff
 6/21/72

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	495	0	37	4533	129
1	570	75	38	4662	129
2	647	77	39	4792	130
3	720	80	40	4922	130
4	802	82	41	5053	131
5	898	86	42	5184	131
6	987	89	43	5315	131
7	1078	91	44	5447	132
8	1161	93	45	5579	132
9	1266	95	46	5711	132
10	1363	97	47	5843	132
11	1462	99	48	5975	132
12	1564	102	49	6107	132
13	1668	104	50	6239	132
14	1773	105	51	6370	131
15	1880	107	52	6500	130
16	1989	109	53	6629	129
17	2099	110	54	6757	128
18	2210	111	55	6885	128
19	2322	112	56	7013	128
20	2436	114	57	7141	128
21	2552	116	58	7269	128
22	2669	117	59	7396	127
23	2787	118	60	7523	127
24	2906	119	61	7649	126
25	3026	120	62	7774	125
26	3147	121	63	7898	124
27	3269	122	64	7921	123
28	3393	124	65	8193	122
29	3518	125	66	8264	121
30	3643	125	67	8384	120
31	3769	126	68	8503	119
32	3895	126	69	8621	118
33	4022	127	70	8737	116
34	4149	127	71	8852	115
35	4276	127	72	8966	114
36	4404	128	73	9079	113

6/26/72


RSV 0015152

TANK CALIBRATION TABLE

TANK NO.: 60T2-1&2Page 2 of 2
MMS/ff
6/21/72

<u>Instrument Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
74	9191	112
75	9300	109
76	9408	108
77	9515	107
78	9620	105
79	9723	103
80	9824	101
81	9924	100
82	10022	98
83	10119	97
84	10215	96
85	10300	95
86	10390	90
87	10475	85
88	10557	82
89	10636	79
90	10711	75
91	10783	72
92	10853	70
93	10921	68
94	10987	66
95	11050	63
96	11107	57
97	11158	51
98	11200	42
99	11233	33
100	11258	25

<u>Instrument Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
-------------------------------	----------------	----------



6/26/72

RSV 0015153

TANK CALIBRATION TABLE

TANK NUMBER:	60T3		Page 1 of 1
TANK NAME:	SLS Storage Tank		MMS /ff
TANK CONTENTS:	Sodium Lauryl Sulfate		Rev. 2/25/72
METHOD OF GAUGING:	LI 102		
TANK DIMENSIONS:	12' OD x 15' SS		
TYPE:	Vertical		
HEADS:	Top - Conical; Bottom - Flat		
EQPT. WITH:	Side-mounted agitator and heating coil		
TRANSMITTER RANGE:	175 inches H ₂ O		
SP. GR.:	1.031	GALLONS/INCH:	70.46

Instrument	Reading (%)	Gallons	Δ	Instrument	Reading (%)	Gallons	Δ
0		846	119	35		5031	119
1		965		36		5151	119
2		1085		38		5390	239
3		1205		40		5630	
4		1325		42		5869	
5		1444		44		6109	
6		1564		46		6348	
7		1684		48		6588	
8		1804		50		6828	
9		1924		52		7067	
10		2043		54		7307	
11		2163		56		7546	
12		2283		58		7779	
13		2403		60		8018	
14		2523		62		8258	
15		2642		64		8498	
16		2762		66		8737	
17		2882		68		8977	
18		3002		70		9216	
19		3121		72		9456	
20		3234		74		9695	
21		3354		76		9935	
22		3474		78		10174	
23		3594		80		10414	
24		3713		82		10654	
25		3833		84		10893	
26		3953		86		11133	
27		4073		88		11372	
28		4192		90		11541	
29		4312		92		11851	
30		4432		94		12091	
31		4552		96		12324	
32		4672		98		12563	
33		4791		100		12803	
34		4911					

TANK CALIBRATION TABLE

TANK NUMBER: 60T5
 TANK NAME: Ammonia Make-up Tank
 TANK CONTENTS: NH_4OH (22.2%)
 METHOD OF GAUGING: LG 121
 TANK DIMENSIONS: 4' dia x 6' SS
 TYPE: Horizontal
 HEADS: Flanged & Dished
 EQPT. WITH: Cooling Coil

Page 1 of 1
 JSM/ff
 7/4/71

Height		Gallons	Δ	Height		Gallons	Δ
0'	0"	45	12	2'	0"	425	16
	1	57	12		1	441	16
	2	69	12		2	457	15
	3	82	13		3	472	15
	4	95	13		4	487	15
	5	110	15		5	502	15
	6	125	15		6	517	15
	7	140	15		7	532	14
	8	155	15		8	546	13
	9	170	15		9	559	13
	10	185	15		10	572	12
	11	200	15		11	584	12
1'	0"	216	16	3'	0"	596	11
	1	232	16		1	607	9
	2	249	17		2	616	8
	3	267	18		3	624	7
	4	285	18		4	631	6
	5	303	19		5	637	6
	6	322	19		6	641	4
	7	341	19				
	8	359	18				
	9	377	16				
	10	393	16				
	11	409	16				

RSV 0015155

TANK CALIBRATION TABLE

TANK NUMBER: 60T6
 TANK NAME: 20% Caustic Storage Tank
 TANK CONTENTS: 20% Sodium Hydroxide
 METHOD OF GAUGING: LG 210
 TANK DIMENSIONS: 6'6" OD x 15'
 TYPE: Horizontal
 HEADS: Flanged & Dished
 EQPT. WITH:

Page 1 of 1
 JSM/ff
 7/4/71

Height	Gallons	△	Height	Gallons	△
0' 0"	135	35	3' 0"	2243	66
1	170	37	1	2309	66
2	207	39	2	2375	65
3	246	40	3	2440	65
4	286	44	4	2505	65
5	330	48	5	2570	65
6	378	48	6	2635	65
7	426	50	7	2700	65
8	476	50	8	2765	64
9	526	51	9	2829	64
10	577	52	10	2893	64
11	629	54	11	2957	64
1' 0"	683	56	4' 0"	3021	63
1	739	57	1	3084	60
2	796	59	2	3144	59
3	855	59	3	3203	59
4	914	60	4	3262	57
5	974	63	5	3319	56
6	1037	64	6	3375	54
7	1101	64	7	3429	52
8	1165	64	8	3481	51
9	1220	64	9	3532	50
10	1293	65	10	3582	50
11	1358	65	11	3632	48
2' 0"	1423	65	5' 0"	3680	48
1	1488	65	1	3728	44
2	1553	65	2	3772	40
3	1618	65	3	3812	39
4	1683	66	4	3851	37
5	1749	66	5	3888	35
6	1815	71	6	3923	
7	1886	71			
8	1957	72			
9	2029	72			
10	2101	71			
11	2172	71			

RSV 0015156

TANK CALIBRATION TABLE

TANK NUMBER: 60T7
 TANK NAME: Condensate Storage Tank
 TANK CONTENTS: Condensate
 METHOD OF GAUGING: LIC 101
 TANK DIMENSIONS: 12' ID x 19'11" SS
 TYPE: Vertical
 HEADS: Top - Conical; Bottom - Flat
 TRANSMITTER RANGE: 154 inches H₂O
 SP. GR.: 1.00 GALLONS/INCHES: 70.5

Page 1 of 1
 MMS /ff
 Rev. 2/25/72

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
0	846	543
5	1388	
10	1956	
15	2473	
20	3016	
25	3558	
30	4101	
35	4643	
40	5186	
45	5728	
50	6271	

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
55	6814	543
60	7356	
65	7899	
70	8441	
75	8984	
80	9526	
85	10069	
90	10611	
95	11154	
100	11696	

TANK CALIBRATION TABLE

TANK NUMBER: 60TB
 TANK NAME: Aqueous Make-Tank
 TANK CONTENTS: Aqueous Solution
 METHOD OF GAUGING: None
 TANK DIMENSIONS: 90" OD x 8'6" SS
 TYPE: Vertical
 HEADS: Flanged & Dished
 EQPT. WITH: An agitator

Page 1 of 1
 JSM/ff
 7/4/71

<u>Liquid Level</u> <u>% or Inches</u>	<u>Gallons</u>	<u>Δ</u>
0	256	
5	396	140
10	536	
15	676	
20	816	
25	956	
30	1096	
35	1236	
40	1376	
45	1516	
50	1656	

<u>Liquid Level</u> <u>% or Inches</u>	<u>Gallons</u>	<u>Δ</u>
55	1796	140
60	1936	
65	2076	
70	2216	
75	2356	
80	2496	
85	2636	
90	2776	
95	2916	
100	3056	

TANK CALIBRATION TABLE

TANK NUMBER:	60T9	Page <u>1</u> of <u>1</u>
TANK NAME:	Aqueous Storage Tank	MMS/ff
TANK CONTENTS:	Aqueous Mix	Rev. 2/25/72
METHOD OF GAUGING:	LI 110	
TANK DIMENSIONS:	11' dia. x 15' SS	
TYPE:	Vertical	
HEADS:	Top - Conical; Bottom - Flat	
EQPT. WITH:	Internal Coils	
TRANSMITTER RANGE:	157 inches H ₂ O	
SP. GR.: 1.00	GALLONS/INCH: 59.2	

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	710		52	5541	184
2	894	184	54	5731	
4	1083		56	5914	
6	1267		58	6104	
8	1456		60	6287	
10	1640		62	6471	
12	1823		64	6661	
14	2013		66	6844	
16	2196		68	7033	
18	2386		70	7217	
20	2569		72	7400	
22	2753		74	7589	
24	2942		76	7773	
26	3126		78	7962	
28	3309		80	8146	
30	3499		82	8329	
32	3682		84	8519	
34	3872		86	8702	
36	4055		88	8892	
38	4245		90	9075	
40	4428		92	9259	
42	4612		94	9448	
44	4801		96	9632	
46	4985		98	9821	
48	5174		100	10005	
50	5358				

TANK CALIBRATION TABLE

TANK NUMBER: 60T10
 TANK NAME: APS Make-Up Tank
 TANK CONTENTS: APS Solution
 METHOD OF GAUGING: None
 TANK DIMENSIONS: 4'6" OD x 5'3"
 TYPE: Vertical
 HEADS: Flanged and Dished
 EQPT. WITH: An agitator and Heating Coils

Page 1 of 1
 JSM/ff
 7/4/71

<u>Liquid Level</u> <u>%</u>	<u>Gallons</u>	<u>Δ</u>	<u>Liquid Level</u> <u>%</u>	<u>Gallons</u>	<u>Δ</u>
0	55	31	55	396	31
5	86	31	60	427	31
10	117	31	65	458	31
15	148	31	70	489	31
20	179	31	75	520	31
25	210	31	80	551	31
30	241	31	85	582	31
35	272	31	90	613	31
40	303	31	95	644	31
45	334	31	100	675	31
50	365	31			

TANK CALIBRATION TABLE

TANK NUMBER: 60T11				Page 1 of 1			
TANK NAME: APS Feed Tank				JSM/ff			
TANK CONTENTS: APS Solution				7/4/71			
METHOD OF GAUGING: LG 118							
TANK DIMENSIONS: 36" OD x 3'0"							
TYPE: Vertical							
HEADS: Flanged & Dished							

LG Height		Gallons	Δ	LG Height		Gallons	Δ
Ft.	In.			Ft.	In.		
0	0	34	4.4	2	0	140	4.4
	1	39	↓		1	144	↓
	2	43			2	149	
	3	47			3	153	
	4	52			4	157	
	5	56			5	162	
	6	61			6	166	
	7	65			7	171	
	8	69			8	175	
	9	74			9	179	
	10	78			10	184	
	11	83			11	188	
1	0	87	↓	3	0	193	
	1	91					
	2	96					
	3	100					
	4	105					
	5	109					
	6	113					
	7	118					
	8	122					
	9	127					
	10	131					
	11	135					

TANK CALIBRATION TABLE

TANK NUMBER: 60T12
 TANK NAME: SFS Make-Up Tank
 TANK CONTENTS: SFS Solution
 METHOD OF GALGING: None
 TANK DIMENSIONS: 42" OD x 4' 0" SS
 TYPE: Vertical
 HEADS: Flanged & Dished
 EQPT. WITH: Agitator and Heating Coils

Page 1 of 1
 JSM/ff
 7/4/71

Liquid Level %	Gallons	Δ	Liquid Level %	Gallons	Δ
0	26	14.4	55	184	14.4
5	40	↓	60	199	↓
10	55		65	213	
15	69		70	228	
20	84		75	242	
25	98		80	256	
30	112		85	271	
35	127		90	285	
40	141		95	300	
45	156		100	314	
50	170				

TANK CALIBRATION TABLE

TANK NUMBER: 60T13
 TANK NAME: SFS Feed Tank
 TANK CONTENTS: SFS Solution
 METHOD OF GAUGING: LG 114
 TANK DIMENSIONS: 36" OD x 3'0"
 TYPE: Vertical
 HEADS: Flanged & Dished

Page 1 of 1
 JSM/ff
 7/4/71

LG Height		Gallons	Δ
Ft.	In.		
0	0	34	4.4
	1	39	↓
	2	43	
	3	47	
	4	52	
	5	56	
	6	61	
	7	65	
	8	69	
	9	74	
	10	78	
	11	83	
1	0	87	↓
	1	91	
	2	96	
	3	100	
	4	105	
	5	109	
	6	113	
	7	118	
	8	122	
	9	127	
	10	131	
	11	135	

LG Height		Gallons	Δ
Ft.	In.		
2	0	140	4.4
	1	144	↓
	2	149	
	3	153	
	4	157	
	5	162	
	6	166	
	7	171	
	8	175	
	9	179	
	10	184	
	11	188	
3	0	193	

TANK CALIBRATION TABLE

TANK NUMBER:	<u>60T14</u>	Page <u>1</u> of <u>1</u>
TANK NAME:	Acrylamide Make-Up Tank	JSM/ff
TANK CONTENTS:	Acrylamide/SLS Solution	7/4/71
METHOD OF GAUGING:	None	
TANK DIMENSIONS:	60" x 63"	
TYPE:	Vertical	
HEADS:	Flanged & Dished	
EQPT. WITH:	Agitator and Heating Coil	

Liquid Level %	Gallons	Δ	Liquid Level %	Gallons	Δ
0	131	82	55	1033	82
5	213	↓	60	1115	↓
10	295		65	1197	
15	377		70	1279	
20	459		75	1361	
25	541		80	1443	
30	623		85	1525	
35	705		90	1607	
40	787		95	1689	
45	869		100	1771	
50	951				

TANK CALIBRATION TABLE

TANK NUMBER: 60T15 Page 1 of 1
 TANK NAME: Acrylamide Feed Tank MMS /ff
 TANK CONTENTS: Acrylamide/SLS Solution Rev. 2/25/72
 METHOD OF GAUGING: LI 122
 TANK DIMENSIONS: 5' ID x 6'2"
 TYPE: Vertical
 HEADS: 5' radius
 EQPT. WITH: External Jacket (not used).
 TRANSMITTER TYPE: Displacer
 SP. GR.: 1.04 GALLONS/INCH: 12.24

<u>Instrument</u> <u>Reading(%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading(%)</u>	<u>Gallons</u>	<u>Δ</u>
0	7		50	428	18
2	15	8	52	446	
4	28	13	54	464	
6	42	14	56	481	
8	59	17	58	499	
10	76	18	60	577	
12	93		62	534	
14	111		64	552	
16	129		66	569	
18	147		68	588	
20	164		70	605	
22	181		72	622	
24	200		74	640	
26	217		76	657	
28	235		78	676	
30	252		80	693	
32	269		82	710	
34	288		84	728	
36	305		86	746	
38	323		88	764	
40	340		90	781	
42	357		92	798	
44	376		94	816	
46	392		96	834	
48	411		98	852	
50	428		100	869	

TANK CALIBRATION TABLE

TANK NUMBER: 60T16-1 & 2, 60T17			Page <u>1</u> of <u>2</u>		
TANK NAME: Hydrolyzer #1 & 2, Filter FD. TK			AKS/ff		
TANK CONTENTS: Product Latex			Rev. 2/11/72		
METHOD OF GAUGING: LI 203-1&2 (0-110" WC)					
TANK DIMENSIONS: 7.5' dia. x 8'9"					
TYPE: Vertical					
HEADS: Flanged & Dished					
EQPT. WITH: Agitator and Heating Coils (#1 only)					

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	411	27.3	35	1367	27.3
1	438		36	1394	
2	466		37	1421	
3	493		38	1448	
4	520		39	1476	
5	548		40	1503	
6	575		41	1530	
7	602		42	1558	
8	629		43	1585	
9	657		44	1612	
10	684		45	1640	
11	711		46	1667	
12	739		47	1694	
13	766		48	1721	
14	793		49	1749	
15	821		50	1776	
16	848		51	1803	
17	875		52	1831	
18	902		53	1858	
19	930		54	1885	
20	957		55	1913	
21	984		56	1940	
22	1012		57	1967	
23	1040		58	1994	
24	1066		59	2022	
25	1094		60	2049	
26	1121		61	2076	
27	1148		62	2104	
28	1175		63	2131	
29	1203		64	2158	
30	1230		65	2186	
31	1257		66	2213	
32	1285		67	2240	
33	1312		68	2267	
34	1339		69	2295	

TANK CALIBRATION TABLE

TANK NUMBER: 60T16-1 & 2, 60T17 (continued)Page 2 of 2
AKS/ff
Rev. 2/11/72

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	2322	27.3	86	2759	27.3
71	2349	↓	87	2786	↓
72	2377		88	2813	
73	2404		89	2841	
74	2431		90	2868	
75	2459		91	2895	
76	2486		92	2923	
77	2513		93	2950	
78	2540		94	2977	
79	2568		95	3005	
80	2595		96	3032	
81	2622		97	3059	
82	2650		98	3086	
83	2677		99	3114	
84	2704		100	3141	
85	2732				

TANK CALIBRATION TABLE

TANK NUMBER: 60T18			Page 1 of 2		
TANK NAME: Latex Receiver			AKS/ff		
TANK CONTENTS: Filtered Latex			Rev. 2/11/72		
METHOD OF GAUGING: LIC 206 (0-50" WC)					
TANK DIMENSIONS: 3'6" OD x 4'3"					
TYPE: Vertical					
HEADS: Flanged and Dished					
EQPT. WITH:					
Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	62.0	2.7	35	156.5	2.7
1	64.7		36	159.2	2.7
2	67.4		37	161.9	
3	70.1		38	164.6	
4	72.8		39	167.3	
5	75.5		40	170.0	
6	78.2		41	172.7	
7	80.9		42	175.4	
8	83.6		43	178.1	
9	86.3		44	180.8	
10	89.0		45	183.5	
11	91.7		46	186.2	
12	94.4		47	188.9	
13	97.1		48	191.6	
14	99.8		49	194.3	
15	102.5		50	197.0	
16	105.2		51	199.7	
17	107.9		52	202.4	
18	110.6		53	205.1	
19	113.3		54	207.8	
20	116.0		55	210.5	
21	118.7		56	213.2	
22	121.4		57	215.9	
23	124.1		58	218.6	
24	126.8		59	221.3	
25	129.5		60	224.0	
26	132.2		61	226.7	
27	134.9		62	229.4	
28	137.6		63	232.1	
29	140.3		64	234.8	
30	143.0		65	237.5	
31	145.7		66	240.2	
32	148.4		67	242.9	
33	151.1		68	245.6	
34	153.8		69	248.3	

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: <u>60T18</u> (continued)			Page <u>2</u> of <u>2</u> AKS /ff Rev. 2/11/72		
<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	251.0	2.7	86	294.2	2.7
71	253.7	↓	87	296.9	↓
72	256.4		88	299.6	
73	259.1		89	302.3	
74	261.8		90	305.0	
75	264.5		91	307.7	
76	267.2		92	310.4	
77	269.9		93	313.1	
78	272.6		94	315.8	
79	275.3		95	318.5	
80	278.0		96	321.2	
81	280.7		97	323.9	
82	283.4		98	326.6	
83	286.1		99	329.3	
84	288.8		100	332.0	
85	291.5				

RSV 0015169

TANK CALIBRATION TABLE

TANK NUMBER: 60T19
 TANK NAME: Latex B Storage Tank
 TANK CONTENTS: Latex B
 METHOD OF GAUGING: LI 208 (0-366.5"WC)
 TANK DIMENSIONS: 13' dia. x 28' SS
 TYPE: Vertical
 HEADS: Top - Conical; Bottom - Flat
 EQPT. WITH: An eductor

Page 1 of 2
 AKS/ff
 Rev. 2/14/72

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
0	910	273	35	10465	273
1	1183		36	10738	
2	1456		37	11011	
3	1729		38	11284	
4	2002		39	11557	
5	2275		40	11830	
6	2548		41	12103	
7	2821		42	12376	
8	3094		43	12649	
9	3367		44	12922	
10	3640		45	13195	
11	3913		46	13468	
12	4186		47	13741	
13	4459		48	14014	
14	4732		49	14287	
15	5005		50	14560	
16	5278		51	14833	
17	5551		52	15106	
18	5824		53	15379	
19	6097		54	15652	
20	6370		55	15925	
21	6643		56	16198	
22	6916		57	16471	
23	7189		58	16744	
24	7462		59	17017	
25	7735		60	17290	
26	8008		61	17563	
27	8281		62	17836	
28	8554		63	18109	
29	8827		64	18382	
30	9100		65	18655	
31	9373		66	18928	
32	9646		67	19201	
33	9919		68	19474	
34	10192		69	19747	

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: <u>60T19</u> (continued)			Page <u>2</u> of <u>2</u> AKS /ff Rev. 2/14/72		
<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	20020	273	86	24388	273
71	20293	↓	87	24661	↓
72	20566		88	24934	
73	20839		89	25207	
74	21112		90	25480	
75	21385		91	25753	
76	21658		92	26026	
77	21931		93	26299	
78	22204		94	26572	
79	22477		95	26845	
80	22750		96	27118	
81	23023		97	27391	
82	23296		98	27664	
83	23569		99	27937	
84	23842		100	28210	
85	24115				

RSV 0015171

TANK CALIBRATION TABLE

TANK NUMBER: 60T20			Page 1 of 2		
TANK NAME: Latex BH Storage Tank			RTH /ff		
TANK CONTENTS: Latex BH Product			Rev. 1/25/72		
METHOD OF GAUGING: LI 207 (0-218.5" WC)					
TANK DIMENSIONS: 17' ID x 18'					
TYPE: Vertical					
HEADS: Bottom Flat; Top Conical					
EQPT. WITH: An Agitator					

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	707		35	10473	279
1	987	279	36	10752	
2	1266		37	11031	
3	1545		38	11310	
4	1824		39	11589	
5	2103		40	11868	
6	2382		41	12147	
7	2661		42	12426	
8	2940		43	12705	
9	3219		44	12984	
10	3498		45	13263	
11	3777		46	13542	
12	4056		47	13821	
13	4335		48	14100	
14	4614		49	14379	
15	4893		50	14658	
16	5172		51	14937	
17	5451		52	15216	
18	5730		53	15495	
19	6009		54	15774	
20	6288		55	16053	
21	6567		56	16332	
22	6846		57	16611	
23	7125		58	16890	
24	7404		59	17169	
25	7683		60	17448	
26	7962		61	17727	
27	8341		62	18006	
28	8520		63	18285	
29	8799		64	18564	
30	9078		65	18843	
31	9357		66	19122	
32	9636		67	19401	
33	9915		68	19680	
34	10194		69	19959	

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: <u>60T20</u> (continued)			Page <u>2</u> of <u>2</u> RTH/ff Rev. 1/25/72		
<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	20238	279	86	24702	279
71	20517	↓	87	24981	↓
72	20796		88	25260	
73	21075		89	25539	
74	21354		90	25818	
75	21633		91	26097	
76	21912		92	26376	
77	22191		93	26655	
78	22470		94	26934	
79	22749		95	27213	
80	23028		96	27492	
81	23307		97	27771	
82	23586		98	28050	
83	23865		99	28329	
84	24144		100	28608	
85	24423				

RSV 0015173

TANK CALIBRATION TABLE

TANK NUMBER:	60T21	Page <u>1</u> of <u>2</u>
TANK NAME:	Misc. Product Storage	AKS/ff
TANK CONTENTS:	Latex Product	Rev. 2/14/72
METHOD OF GAUGING:	LI 209 (0-323.9" WC)	
TANK DIMENSIONS:	13' ID x 25'	
TYPE:	Vertical	
HEADS:	Top - Conical; Bottom - Flat	
EQPT. WITH:	An agitator	

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
0	703	241.3 ↓	35	9149	241.3 ↓
1	944		36	9390	
2	1186		37	9631	
3	1427		38	9872	
4	1668		39	10114	
5	1910		40	10355	
6	2151		41	10596	
7	2392		42	10838	
8	2633		43	11079	
9	2875		44	11320	
10	3116		45	11562	
11	3357		46	11803	
12	3599		47	12044	
13	3840		48	12285	
14	4081		49	12527	
15	4323		50	12768	
16	4564		51	13009	
17	4805		52	13251	
18	5046		53	13492	
19	5288		54	13733	
20	5529		55	13975	
21	5770		56	14216	
22	6012		57	14457	
23	6253		58	14698	
24	6494		59	14940	
25	6736		60	15181	
26	6977		61	15422	
27	7218		62	15664	
28	7459		63	15905	
29	7701		64	16146	
30	7942		65	16388	
31	8183		66	16629	
32	8425		67	16870	
33	8666		68	17111	
34	8907		69	17353	

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60T21 (continued)

Page 2 of 2
AKS /ff
Rev. 2/14/72

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
70	17594	241.3	86	21455	241.3
71	17835	↓	87	21692	↓
72	18077		88	21937	
73	18318		89	22179	
74	18559		90	22420	
75	18801		91	22661	
76	19042		92	22903	
77	19283		93	23144	
78	19524		94	23385	
79	19766		95	23627	
80	20007		96	23868	
81	20248		97	24109	
82	20490		98	24350	
83	20731		99	24592	
84	20972		100	24833	
85	21214				

RSV 0015175

TANK CALIBRATION TABLE

TANK NUMBER:	60T22	Page <u>1</u> of <u>2</u>
TANK NAME:	Waste Latex & Vent Tank	AKS/ff
TANK CONTENTS:	Off-Spec. Latex	Rev. 2/11/72
METHOD OF GAUGING:	LT 123 (0-304" WC)	
TANK DIMENSIONS:	9'0" ID x 27' SS	
TYPE:	Vertical	
HEADS:	8'6" Radius	
EQUIPPED WITH:	Impingement Plate	

Instrument Reading(%)	Gallons H ₂ O -	Δ gal/%	Instrument Reading(%)	Gallons H ₂ O	Δ gal/%
0	795	120.5	35	5013	120.5
1	916		36	5133	
2	1036		37	5254	
3	1157		38	5374	
4	1277		39	5495	
5	1398		40	5615	
6	1518		41	5736	
7	1639		42	5856	
8	1759		43	5977	
9	1880		44	6097	
10	2000		45	6218	
11	2121		46	6338	
12	2241		47	6459	
13	2362		48	6579	
14	2482		49	6700	
15	2603		50	6820	
16	2723		51	6941	
17	2844		52	7061	
18	2964		53	7182	
19	3085		54	7302	
20	3205		55	7423	
21	3326		56	7543	
22	3446		57	7664	
23	3567		58	7784	
24	3687		59	7905	
25	3808		60	8025	
26	3928		61	8146	
27	4049		62	8266	
28	4169		63	8387	
29	4290		64	8507	
30	4410		65	8628	
31	4531		66	8748	
32	4651		67	8869	
33	4772		68	8989	
34	4892		69	9110	

TANK CALIBRATION TABLE

TANK NUMBER: 60T22 (continued)Page 2 of 2

AKS/ff

Rev. 2/11/72

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons H₂O</u>	<u>Δgal/%</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons H₂O</u>	<u>Δgal/%</u>
70	9230	120.5	86	11158	120.5
71	9351		87	11279	
72	9471		88	11399	
73	9592		89	11520	
74	9712		90	11640	
75	9833		91	11761	
76	9953		92	11881	
77	10074		93	12002	
78	10194		94	12122	
79	10315		95	12243	
80	10435		96	12363	
81	10556		97	12484	
82	10676		98	12604	
83	10797		99	12725	
84	10917		100	12845	
85	11038				

RSV 0015177

TANK CALIBRATION TABLE

TANK NUMBER: 60T24
 TANK NAME: Waste Tank
 TANK CONTENTS: Water/Latex Solids
 METHOD OF GAUGING: LI 303 (0-100" WC)
 TANK DIMENSIONS: 8'6" dia. x 8'8" SS
 TYPE: Vertical
 HEADS: 8'6" Radius
 EQPT. WITH:

Page 1 of 2
 AKS/ff
 Rev. 2/14/72

Instrument Reading (%)	Gallons H ₂ O	Δ	Instrument Reading (%)	Gallons H ₂ O	Δ
0	460	35.4	35	1699	35.4
1	495		36	1734	
2	531		37	1770	
3	566		38	1805	
4	602		39	1841	
5	637		40	1876	
6	672		41	1911	
7	708		42	1947	
8	743		43	1982	
9	779		44	2018	
10	814		45	2053	
11	849		46	2088	
12	885		47	2124	
13	920		48	2159	
14	956		49	2195	
15	991		50	2230	
16	1026		51	2265	
17	1062		52	2301	
18	1097		53	2336	
19	1133		54	2372	
20	1168		55	2407	
21	1203		56	2442	
22	1239		57	2478	
23	1274		58	2513	
24	1310		59	2549	
25	1345		60	2584	
26	1380		61	2619	
27	1416		62	2655	
28	1451		63	2690	
29	1487		64	2726	
30	1522		65	2761	
31	1557		66	2796	
32	1593		67	2832	
33	1628		68	2867	
34	1664		69	2903	

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60T24 (continued)

Page 2 of 2
AKS/ff
Rev. 2/14/72

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons H₂O</u>	<u>Δ</u>
70	2938	35.4	86	3504	35.4
71	2973	↓	87	3540	↓
72	3009		88	3575	
73	3044		89	3611	
74	3080		90	3646	
75	3115		91	3681	
76	3150		92	3717	
77	3186		93	3752	
78	3221		94	3788	
79	3257		95	3823	
80	3292		96	3858	
81	3327		97	3894	
82	3363		98	3929	
83	3398		99	3965	
84	3434		100	4000	
85	3469				

RSV 0015179

TANK CALIBRATION TABLE

TANK NUMBER: 60T25
 TANK NAME: Seal Pot
 TANK CONTENTS: Filtered Water
 METHOD OF GAUGING: LG 403
 TANK DIMENSIONS: 4' dia. x 10' SS
 TYPE: Horizontal
 HEADS: 4' Radius
 EQPT. WITH:

Page 1 of 1
 JSM/ff
 7/4/71

Height		Gallons	Δ
Ft.	In.		
0	0	37	32
	2	69	37
	4	106	41
	6	147	44
	8	191	45
	10	236	48
1	0	284	52
	2	336	52
	4	388	54
	6	442	54
	8	496	54
	10	550	54

Height		Gallons	Δ
Ft.	In.		
2	0	604	52
	2	656	52
	4	708	48
	6	756	45
	8	804	44
	10	845	41
3	0	886	

TANK CALIBRATION TABLE

TANK NUMBER: 60T26			Page 1 of 2		
TANK NAME: Chilled Water Tank			JSM/ff		
TANK CONTENTS: Chilled Condensate			7/4/71		
METHOD OF GAUGING: LI 501 (0-114" WC)					
TANK DIMENSIONS: 10' dia. x 40'					
TYPE: Horizontal					
HEADS: Elliptical					
EQPT. WITH:					

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	448	133	35	8871	305
1	581	144	36	9176	306
2	725	152	37	9482	308
3	877	160	38	9790	308
4	1037	169	39	10098	309
5	1206	179	40	10407	311
6	1385	189	41	10718	313
7	1574	196	42	11031	316
8	1770	201	43	11347	316
9	1971	206	44	11663	316
10	2177	218	45	11979	316
11	2395	222	46	12295	316
12	2617	227	47	12611	317
13	2844	231	48	12928	317
14	3075	235	49	13245	316
15	3310	242	50	13561	316
16	3552	244	51	13877	316
17	3796	249	52	14193	314
18	4045	256	53	14507	312
19	4301	264	54	14819	312
20	4565	270	55	15131	311
21	4835	273	56	15442	309
22	5108	274	57	15751	309
23	5382	276	58	16060	308
24	5658	277	59	16368	306
25	5935	280	60	16674	302
26	6215	281	61	16976	300
27	6496	285	62	17276	299
28	6781	286	63	17575	298
29	7067	296	64	17873	296
30	7363	297	65	18169	294
31	7660	300	66	18463	288
32	7960	302	67	18751	285
33	8262	304	68	19036	282
34	8566	305	69	19318	280

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60T26 (continued)Page 2 of 2
JSM/ff
7/4/71

<u>Instrument Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	19598	277	86	23542	198
71	19875	274	87	23740	193
72	20149	272	88	23933	185
73	20421	269	89	24118	178
74	20690	266	90	24296	167
75	20956	263	91	24463	156
76	21219	257	92	24619	146
77	21476	252	93	24765	137
78	21728	247	94	24902	130
79	21975	241	95	25032	127
80	22216	236	96	25159	100
81	22452	230	97	25259	89
82	22682	224	98	25348	73
83	22906	219	99	25421	39
84	23125	212	100	25460	
85	23337	205			

RSV 0015182

TANK CALIBRATION TABLE

TANK NUMBER: 60T36
 TANK NAME: Catch Tank
 TANK CONTENTS: Filtered Water
 METHOD OF GAUGING: None
 TANK DIMENSIONS: 8' I.D. x 14'0"
 TYPE: Horizontal
 HEADS: Dished Heads
 EQPT. WITH:

Page 1 of 2
 JSM/ff
 7/4/71

Height				Height			
Ft.	In.	Gallons	Δ	Ft.	In.	Gallons	Δ
0	0			3	0	1925	74
	1	16	16		1	1999	74
	2	32	17		2	2073	75
	3	49	29		3	2148	75
	4	78	30		4	2223	75
	5	108	31		5	2298	76
	6	139	38		6	2374	76
	7	177	39		7	2450	77
	8	216	39		8	2527	77
	9	255	44		9	2604	78
	10	299	44		10	2682	78
	11	343	48		11	2760	78
1	0	391	50	4	0	2838	78
	1	441	50		1	2916	78
	2	491	53		2	2994	78
	3	544	53		3	3072	77
	4	597	54		4	3149	77
	5	651	55		5	3226	76
	6	706	61		6	3302	76
	7	767	61		7	3378	75
	8	828	63		8	3453	75
	9	891	63		9	3528	75
	10	954	63		10	3603	74
	11	1017	65		11	3677	74
2	0	1082	65	5	0	3751	74
	1	1147	67		1	3825	73
	2	1214	68		2	3898	73
	3	1282	69		3	3971	72
	4	1351	69		4	4043	72
	5	1420	70		5	4115	71
	6	1490	71		6	4186	70
	7	1561	72		7	4256	69
	8	1633	72		8	4325	69
	9	1705	73		9	4394	68
	10	1778	73		10	4462	67
	11	1851	74		11	4529	65

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: <u>60T36</u> (continued)				Page <u>2</u> of <u>2</u> JSM/ff 7/4/71			
Height				Height			
<u>Ft.</u>	<u>In.</u>	<u>Gallons</u>	<u>Δ</u>	<u>Ft.</u>	<u>In.</u>	<u>Gallons</u>	<u>Δ</u>
6	0	4594	65	7	0	5285	48
	1	4659	63		1	5333	44
	2	4722	63		2	5377	44
	3	4785	63		3	5421	39
	4	4848	61		4	5460	39
	5	4909	61		5	5499	38
	6	4970	55		6	5537	31
	7	5025	54		7	5568	30
	8	5079	53		8	5598	29
	9	5132	53		9	5627	17
	10	5185	50		10	5644	16
	11	5235	50		11	5660	16
				8	0	5676	

RSV 0015184

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60T38
 TANK NAME: Chemlime Storage Tank
 TANK CONTENTS: Chemlime
 METHOD OF GAUGING: LI 402
 TANK DIMENSIONS: 12' OD x 14' 0"
 TYPE: Vertical
 HEADS: Top - Conical; Bottom - Flat
 EQPT. WITH: Agitator and Internal Coils

Page 1 of 2
 JSM/ff
 7/4/71

<u>Instrument</u> <u>Reading(%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading(%)</u>	<u>Gallons</u>	<u>Δ</u>
0	846	110	35	4696	110
1	956		36	4806	
2	1066		37	4916	
3	1176		38	5026	
4	1286		39	5136	
5	1396		40	5246	
6	1506		41	5356	
7	1616		42	5466	
8	1726		43	5576	
9	1836		44	5686	
10	1946		45	5796	
11	2056		46	5906	
12	2166		47	6016	
13	2276		48	6126	
14	2386		49	6236	
15	2496		50	6346	
16	2606		51	6456	
17	2716		52	6566	
18	2826		53	6676	
19	2936		54	6786	
20	3046		55	6896	
21	3156		56	7006	
22	3266		57	7116	
23	3376		58	7225	
24	3486		59	7335	
25	3596		60	7445	
26	3706		61	7555	
27	3816		62	7665	
28	3926		63	7775	
29	4036		64	7885	
30	4146		65	7995	
31	4256		66	8105	
32	4366		67	8215	
33	4476		68	8325	
34	4586		69	8435	

RSV 0015185

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60T38 (continued)Page 2 of 2
JSM/ff
7/4/71

<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u> <u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
70	8545	110	86	10305	110
71	8655	↓	87	10415	↓
72	8765		88	10525	
73	8875		89	10635	
74	8985		90	10745	
75	9095		91	10855	
76	9205		92	10965	
77	9315		93	11075	
78	9425		94	11185	
79	9535		95	11295	
80	9645		96	11405	
81	9755		97	11515	
82	9865		98	11625	
83	9975		99	11735	
84	10085		100	11845	
85	10195				

RSV 0015186

60-801.00

TANK CALIBRATION TABLE

TANK NUMBER: 60T40			Page <u>1</u> of <u>2</u>		
TANK NAME: Latex E Storage Tank			AKS/ff		
TANK CONTENTS: Latex E Product			Rev. 1/25/72		
METHOD OF GAUGING: LI 216 (0-230.2"WC)					
TANK DIMENSIONS: 20'0" ID x 18'					
TYPE: Vertical					
HEADS: Bottom - Flat; Top - Conical					
EQPT. WITH: An Eductor					

<u>Instrument</u>	<u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>	<u>Instrument</u>	<u>Reading (%)</u>	<u>Gallons</u>	<u>Δ</u>
0		2153	405.9	35		16360	405.9
1		2559		36		16765	
2		2965		37		17171	
3		3371		38		17577	
4		3777		39		17983	
5		4183		40		18389	
6		4588		41		18795	
7		4994		42		19201	
8		5400		43		19607	
9		5806		44		20013	
10		6212		45		20419	
11		6618		46		20824	
12		7024		47		21230	
13		7430		48		21636	
14		7836		49		22042	
15		8242		50		22448	
16		8647		51		22854	
17		9053		52		23260	
18		9459		53		23666	
19		9865		54		24072	
20		10271		55		24478	
21		10677		56		24883	
22		11083		57		25289	
23		11489		58		25695	
24		11895		59		26101	
25		12301		60		26507	
26		12706		61		26913	
27		13112		62		27319	
28		13518		63		27725	
29		13924		64		28131	
30		14330		65		28537	
31		14736		66		28942	
32		15142		67		29348	
33		15548		68		29754	
34		15548		69		29754	

RSV 0015187

Monsanto

FROM (NAME & LOCATION)

Texas City, Texas

DATE

CC

SUBJECT

Additional Material for Dept.(s) 60
Operating Manual/Specification Manual

REFERENCE

TO

M E Gibbs - SL

Please insert the attached material into its correct
location in your Department(s) _____
Operating Manual/Specification Manual, Copy No. _____.

Freida Fredericks
Standards Section

The attached material is
COMPANY-CONFIDENTIAL

ANY MATERIAL THAT IS BEING SUPERSEDED BY THE ATTACHED
SHOULD BE DESTROYED OR RETURNED TO THE STANDARDS SECTION

RSV 0015188

TANK CALIBRATION TABLE

TANK NUMBER: 60T41
 TANK NAME: LATEX STORAGE TANK
 TANK CONTENTS: LATEX PRODUCT (Sp.Gr. 1.11)
 METHOD OF GAUGING: LI (0-360" W.C.)
 TANK DIMENSIONS: 16' 0" OD x 30'
 TYPE: VERTICAL
 HEADS: BOTTOM-FLAT; TOP-CONICAL
 EQPT. WITH: BAFFLES FOR FUTURE AGITATOR

Page 1 of 2
 JSM/ff
 9/7/73

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
0	1,467	404.9	35	15,639	404.9
1	1,872		36	16,044	
2	2,277		37	25,357	
3	2,682		38	16,854	
4	3,087		39	17,259	
5	3,492		40	17,664	
6	3,897		41	18,069	
7	4,302		42	18,474	
8	4,706		43	18,879	
9	5,111		44	19,284	
10	5,516		45	19,689	
11	5,921		46	20,094	
12	6,326		47	20,498	
13	6,731		48	20,903	
14	7,136		49	21,308	
15	7,540		50	21,714	
16	7,946		51	22,118	
17	8,351		52	22,523	
18	8,756		53	22,928	
19	9,161		54	23,333	
20	9,566		55	23,738	
21	9,970		56	24,143	
22	10,375		57	24,548	
23	10,780		58	24,953	
24	11,185		59	25,357	
25	11,590		60	25,762	
26	11,995		61	26,167	
27	12,400		62	26,572	
28	12,805		63	26,977	
29	13,210		64	27,382	
30	13,615		65	27,787	
31	14,020		66	28,192	
32	14,425		67	22,597	
33	14,830		68	29,002	
34	15,234		69	29,407	

TANK CALIBRATION TABLE

TANK NUMBER: 60T41 (continued)Page 2 of 2
JSM/ff
9/7/73

Instrument Reading (%)	Gallons	Δ	Instrument Reading (%)	Gallons	Δ
70	29,812	404.9	86	36,290	404.9
71	30,217	↓	87	36,695	↓
72	30,621		88	37,100	
73	31,026		89	37,505	
74	31,431		90	37,910	
75	31,836		91	38,315	
76	32,241		92	38,920	
77	32,646		93	39,125	
78	33,051		94	39,530	
79	33,459		95	39,935	
80	33,861		96	40,340	
81	34,266		97	40,745	
82	34,671		98	41,149	
83	35,676		99	41,554	
84	35,481		100	41,959	
85	35,885				

MONSANTO COMPANY
MONSANTO INDUSTRIAL CHEMICALS DIVISION
TEXAS CITY PLANT

STANDARD MANUFACTURING PROCESS
FOR
ETHYLENE VINYL CHLORIDE
(Dept. 60)

Prepared by: M. E. Gibbs
R. T. Hammann
G. A. Hart
M. M. Siegel
A. Shadensack
J. Morgan

Approved by: J. L. Rasmussen
J. L. Rasmussen
Manufacturing Superintendent

R. W. Flint
R. W. Flint
Guest Superintendent

M. E. Gibbs
M. E. Gibbs
Research Group Leader

J. R. Savage
J. R. Savage
Manufacturing Manager

Copy 1 of 20

Issued:

OCT 1 1972

CONFIDENTIAL INFORMATION

This process is the property of Monsanto Company and the recipient is responsible for its safekeeping. It contains confidential information of Monsanto Company which must not be reproduced, revealed to unauthorized persons or sent outside the Company without proper authorization. Either retain in a secure file or return to Standards Department.

RSV 0015191

DISTRIBUTION

<u>COPY NO.</u>	<u>DISTRIBUTION TO:</u>
1.	Technical Information Center St. Louis Research
2.	Technical Information Center St. Louis Research
3.	Vital Records Center General Offices
4.	Business Group Research Manager Specialty Products Group Administration
5.	Business Group Research Group Leader Specialty Products Group Administration
6.	Manufacturing Manager Specialty Products Group Administration
7.	General Superintendent MICC Manufacturing, Texas City
8.	Manufacturing Superintendent Texas City
9.	Manufacturing Supervisor Texas City
10.	Manufacturing Foreman Texas City
11.	Process Technology Superintendent MICC Manufacturing, Texas City
12.	Central Files, Texas City
13.	Library, Texas City
14.	Safety Review Committee
15.	Safety Review Committee
16.	Safety Review Committee
17-20	Extra


RSV 0015192

DISTRIBUTION

Date: _____

COPY NO.

DISTRIBUTION TO:

1.  Technical Information Center - St. Louis Research
2. Technical Information Center - St. Louis Research
3. Vital Record Center - General Offices - St. Louis
4. W. R. Richard - St. Louis
5. M. E. Gibbs - St. Louis
6. J. R. Savage - St. Louis
7. R. W. Flint - Texas City
8. J. L. Rasmussen - Texas City
9. R. T. Hammann - Texas City
10. G. A. Hart - Texas City
11. W. L. Smull - Texas City
12. Central Files - Texas City
13. Library - Texas City
14. J. F. Quinn - St. Louis
15. M. A. Terpstra - St. Louis
16. R. W. Radue - St. Louis
- 17-20 Extra Copies

CONFIDENTIAL

RSV 0015193

TABLE OF CONTENTS

TITLE PAGE

DISTRIBUTION

TABLE OF CONTENTS

PRE-SECTION* - AMENDMENTS
SECTION I - SYNOPSIS OF PROCESS

SECTION II - FLOW SHEETS

- A. Makeup System, Raw Material Storage
- B. Autoclave Feed System
- C. Autoclave System
- D. Blowdown, Hydrolysis and Filter Feed
- E. Filtration and Product Storage
- F. THF Storage and Recovery
- G. Vent System and Incineration
- H. Chilled Water and Hot Water Systems

SECTION III - PROCESS IN DETAIL

- A. Raw Materials
- B. Solution Makeup
- C. Autoclave
- D. Hydrolysis
- E. Filtration
- F. Storage and Packout

Table of Contents (cont'd)

SECTION IV - COMMENTS ON PROCESS

A. Raw Materials and Feed Makeup

1. Acrylamide Solution
2. Ammonium Persulfate Solution
3. Aqueous Solution
4. Sodium Formaldehyde Sulfoxylate Sol'n.
5. Sodium Lauryl Sulfate

B. Polymerization

1. Catalyst
2. Polymerization Reaction
 - a. Pressure
 - b. Temperature
 - c. Heat Generation and Removal
 - d. Initial Batch Charging
 - e. Reaction Start
 - f. Rate and Timing of Acrylamide Addition
 - g. Catalyst Addition Rates
 - h. Total Amount of VCM Addition
3. Blowdown
4. Agitation

C. Hydrolysis

1. Temperature
2. Caustic Concentration
3. Amount of Caustic
4. Percent Total Solids

D. Screening and Latex Transfer

1. Temperature
2. Pressure
3. Grit Content
4. Latex Transfer

Table of Contents (cont'd)

SECTION IV - COMMENTS ON PROCESS (cont'd)

E. Auxiliary Systems

1. Incinerator
2. Chilled Water
3. Reactor Sequencer
4. Ethylene
5. Solvent System

SECTION V - PROCESS CONTROL PROCEDURES

- A. Introduction
- B. Safety and Housekeeping
- C. Quality
- D. Yields and Rates
- E. Performance - Key Variables

APPENDIX A - MATERIAL BALANCES

1. Material Balance for Latex B Production
2. Material Balance for Latex BH Production
3. Material Balance for Latex E Production

APPENDIX B - SPECIFICATIONS, FINISHED GOODS AND RAW MATERIALS

1. Raw Materials
 - a. Acrylamide
 - b. Ammonia
 - c. Ammonium Persulfate
 - d. Butylated Hydroxy Toluene
 - e. Ethylene
 - f. Ferric Ammonium Sulfate
 - g. Methanol
 - h. Mono Methyl Ether Hydroquinone
 - i. Sodium Formaldehyde Sulfoxylate
 - j. Caustic
 - k. Sodium Lauryl Sulfate
 - l. Tetrasodium Ethylenediamine Tetraacetate
 - m. Tetrahydrofuran
 - n. Vinyl Chloride
2. Finished Product
 - a. Latex B
 - b. Latex BH
 - c. Latex E

APPENDIX C - PHYSICAL CONSTANTS

APPENDIX D - PACKAGING AND SHIPPING

APPENDIX E - SAFETY DATA

1. Rules and Precautions
 - a. Personnel
 - b. Equipment
 - c. Operating Practices

II. Toxicity

APPENDIX F	-	<u>ENVIRONMENTAL CONTROL AND WASTE DISCHARGE</u>
	I.	Incinerator Operation
	II.	Comments on Incineration
	III.	Water Emissions
	IV.	Emission Summary for Latexes B & BH
	V.	Emission Summary for Latex E
APPENDIX G	-	<u>EQUIPMENT LIST</u>
APPENDIX H	-	<u>UTILITY SUMMARY</u>
APPENDIX I	-	<u>PERSONNEL</u>
APPENDIX J	-	<u>TIME CYCLES</u>
APPENDIX K	-	<u>QUALITY CONSIDERATIONS</u>
APPENDIX L	-	<u>FORM FOR TENTATIVE AMENDMENT TO THE STANDARD MANUFACTURING PROCESS</u>

Monsanto


NAME & LOCATION G. A. Danner - Texas City

DATE March 12, 1974

SUBJECT EVC Tentative Amendment B

REFERENCE

TO : Ethylene Vinyl Chloride SMP Holders #1



Attached is a copy of Amendment B to the Ethylene Vinyl Chloride SMP. Please add this to your copy of the SMP.

G. A. Danner jp
G. A. Danner

/jp
Attach.

Date: January 24, 1974

Tentative Amendment No. B

STANDARD MANUFACTURING PROCESS FOR
ETHYLENE VINYL CHLORIDE

Product Code No.'s
840.61, 840.62, 840.63, 840.64

ISSUANCE OF CURRENT PROCESS DOCUMENT 1972

TITLE OF AMENDMENT: Autoclave Heatup for Hydrolyzing EVC-BH

PRESENT PRACTICE:

At the end of the reaction period when producing EVC-BH, the batch is heated from 55 to 60°C by the use of steam on the autoclave coils so that the product in the hydrolyzer will be greater than 45°C.

SUGGESTED CHANGE:

When the VCM addition is 90% complete, the steam valve on the autoclave will be manually blocked in and the temperature controller will be increased from 55 to 60°C letting the heat of reaction increase the temp. Due to a period of low VCM addition because of pressure control on the autoclave, the reaction cycle will be increased from 4 1/3 to 4 2/3 hours per batch.

JUSTIFICATION FOR CHANGE:

There is some indication that heating the latex using steam on the coils may have contributed to fouling of the autoclave coils. This suggested approach will allow the department to make hydrolyzed EVC-BH without using steam on the autoclave coils.

SAFETY:

No new property or personnel hazards will be introduced by the proposed change. All temperature and pressure shutdown devices will continue to be in service.

QUALITY:

No quality problems are expected from this change.

EXTENT OF DEMONSTRATION:

Permission is requested to run for three months. If successful this amendment will become part of the Standard Manufacturing Process.

Originated by:

R. T. Hammann
R. T. Hammann

Endorsed by:

J. L. Rasmussen
J. L. Rasmussen

APPROVED BY:

R. W. Flint
R. W. Flint

Supt., Mfg.

M. E. Gibbs
M. E. Gibbs

Research Group Leader

W. R. Richard
W. R. Richard

Research Manager

J. R. Savage
J. R. Savage

Mfg. Manager

mm


Monsanto

(NAME & LOCATION)

J. S. Morgan - Texas City

DATE February 23, 1973 CC
SUBJECT EVC Tentative Amendment A
REFERENCE
TO : Ethylene Vinyl Chloride SMP Holders #2

Amendment A to the Ethylene Vinyl Chloride SMP was successfully demonstrated in October and November, 1972. Attached is a copy of the amendment to be included in your SMP.



J. S. Morgan

/jp
Attach.

RSV 0015202

Tentative Amendment A
to the
Ethylene Vinyl Chloride
Tentative Process

DATE: March 20, 1972

TENTATIVE AMENDMENT: A

TENTATIVE PROCESS REPORT FOR PRODUCTION OF POLVIN PIGMENT BINDER,
MONSAT, AND MONFLEX

January 7, 1970. Report No. P-1556.

TITLE: Eliminate Second Water Wash

PRESENT PRACTICE:

The EVC reactors are washed with 55°C steam condensate twice after batch blowdown to remove any residual latex and grit particles that will foul the cooling coils. Every fifteenth reaction batch also includes a solvent (tetrahydrofuran - THF) wash between the first and second condensate washes to dissolve any coagulated latex that has built up on the cooling coils. The wash cycle consists of the following steps:

- | | |
|--|------------------|
| 1. Fill reactor (2,150 gallons) with condensate (or THF) | - 20 min. |
| 2. Agitate reactor contents | - 15 min. |
| 3. Pump out reactor contents to sewer (or THF storage) | - <u>10 min.</u> |
| TOTAL | - 45 min. |

PROPOSED CHANGE:

This amendment proposes the elimination of the second condensate wash - retaining the first condensate wash on each batch and the THF wash every 15 batches. Each THF wash will be followed by a second condensate wash to remove residual THF from the reactor.

JUSTIFICATION

The primary justification for eliminating the second water wash is the reduction of 2,150 gallons H₂O/batch (15.9 M gallons/day at standard production rate) to the sewer. This represents a

JUSTIFICATION (cont'd.)

raw material savings of \$5 M/year and \$16 M in proposed secondary waste treating capital. These figures are for standard product rate - 23.7 M lbs./yr. (January, 1972 Monthly Report - M. E. Gibbs). The wash elimination will result in an individual batch cycle reduction of 45 minutes (8 1/2%). No savings due to increased unit capacity can be claimed at this time. →

The first water wash contains approximately 20 lbs. (1,100 ppm) of dissolved polymer solids and approximately 1/2 lb. of grit. The color of this stream is chalky white. The second water wash contains approximately 2 lbs. (110 ppm) polymer solids and has been observed to be clear with some grit (settleable) present. Each wash is pumped to the sewer which flows through the Primary Waste Treatment Area (startup 6/1/72) and the proposed secondary treatment facility and eventually to Galveston Bay.

Any latex left in the reactor prior to batch start will result in the production of some larger than normal particle sizes. The 2 lbs. of latex from the second wash will be left in the reactor prior to the next batch start. The amount of latex remaining, however, is slight and will probably not noticeably affect latex particle size.

BACKGROUND

The low solids content of the second water wash has been observed on every batch since the unit startup in November, 1971. The Tentative Process gives approximately 100 ppm as the solids concentration of the second wash. It has been observed the amount of solids in the washes can be reduced by optimizing the 2 low pressure (100 psig) reactor blowdowns.

SAFETY

Safety will not be affected.

QUALITY

Product quality will be monitored through regular analysis throughout the demonstration. The solids content of the first wash will also be monitored to determine if the solids concentration increases. Reactor heat transfer data will be taken to assure that the THF wash frequency is not increased to greater than 1/15 ba due to coil fouling. Particle size will be obtained every fifth batch during the demonstration.

EXTENT OF DEMONSTRATION

Permission is requested to run for one month - in order to get at least 15 batches on one reactor - with the second water wash eliminated. If in this period no "off-spec" product is produced and the THF wash frequency is not increased to greater than 1/15 batches (U overall should remain > 100 BTU/hr. ft.²°F at max. heat load), the Tentative Amendment will be considered permanent.

Originated by: Alan K. Shadensack
A. K. Shadensack

Endorsed by: R. T. Hammann
R. T. Hammann

Approved by:

W. L. Smull /W. L. Smull M. E. Gibbs /M. E. Gibbs
J. L. Rasmussen /J. L. Rasmussen W. R. Richard /W. R. Richard
R. W. Flint /R. W. Flint

Safety Review Committee

R. W. Radue ^{5/1/72} /~~R. W. Radue~~
F. Quinn ^{5/1/75} /F. Quinn
M. A. Terpstra ^{5/1/71} /M. A. Terpstra
J. R. Savage /J. R. Savage

/jp
4/7/72

Monsanto

NAME & LOCATION: J. S. Morgan - Texas City

DATE : February 20, 1973
SUBJECT : EVC Tentative Amendment A
REFERENCE : Plant Demonstration
TO : Ethylene Vinyl Chloride SMP Holders

Tentative Amendment A was successfully demonstrated October 28 - November 8, 1972. Objective was to eliminate the second autoclave water wash from the normal batch cycle. The results are shown below.

Objectives and Results
For Plant Demonstration

<u>Performance Criteria</u>	<u>Objective</u>	<u>Results</u>
THF Wash Frequency	15	15
Heat Transfer Coefficient (BTU/hr. F ² °F)	>100	>178
Product Particle Size (A°)	650-900	685-827
Off-Spec Product Produced (lbs.)	0	0
Condensate Savings @ 23.7 M lbs./yr. (\$/yr.)	5 M	5 M
Secondary Treatment Capital Savings (\$)	16 M	16 M
Batch Cycle Reduction (%)	8.5	8.5

A tabulation of the complete test results is attached.

Conclusions

The second water wash can be eliminated from the batch cycle without effecting product quality or the frequency of the washes. The \$5 M/yr. condensate savings (at 23.7 M lbs./yr.) and the 8 1/2% reduction in batch cycle time will be achieved.

This amendment does not include:

1. Elimination of the two water washes during the THF wash sequence.
2. New products with significantly different heat transfer requirements and fouling characteristics as compared to Latex B or E.

Disposition of Tentative Amendment D

The amendment is accepted. No modifications will be required to permanently implement the new procedure. The sequencer currently

RSV 0015206

halts after the first water wash and must be manually advanced to the subsequent step. Under Amendment A the sequencer is advanced to Step 25 and held until the next batch is started.

A handwritten signature in black ink, appearing to read "J. S. Morgan", with a stylized, cursive script.

J. S. Morgan

/jp
Attach.

<u>Date</u>	<u>Batch Number</u>	<u>Heat Trans. Coefficient*</u>	<u>Particle Size</u>	<u>Since THF Wash</u>
10/28/72	223	238		1
10/29/72	225	236	795	2
10/29/72	227	211		3
10/30/72	229	-205		4
10/30/72	230	213		5
10/31/72	233	206		6
10/31/72	235	185	827	7
11/1/72	237	188		8
11/2/72	239	182		9
11/3/72	241	196		10
11/4/72	243	192	822	11
11/5/72	246	181		12
11/6/72	248	169		13
11/7/72	250	182	685	14
11/8/72	252	178		15

* Assuming Coil Area = 468 ft.²
H₂O Flow = 340 GPM

JSM/jp
2/20/73

RSV 0015208

RSV 0015209

I. SYNOPSIS OF PROCESS

I. SYNOPSIS OF PROCESS


The EVC department located at the Texas City plant was designed to produce 19 M pounds of EVC latex (weight on basis of dry polymer) per year via a Monsanto developed batch process. Three different latices are made, denoted as Latex B, Latex BH and Latex E. Capacity for Latex B and Latex E has been demonstrated at 22-23 M pounds per year. Capacity of BH is as yet not defined but is somewhat lower at present.

The basic raw materials include the three monomers: vinyl chloride, ethylene, acrylamide, and a surfacant, sodium lauryl sulfate. Ethylene is obtained via pipeline from the Polymers and Petrochemical Division facilities at both the Texas City and Chocolate Bayou plants. Vinyl chloride is shipped by rail car from Dow. Sodium lauryl sulfate is received in bulk and acrylamide in 50 pound bags.


The manufacture of EVC latex is centered in the polymerization step and accomplished in 2150 gallon, 2500 psi, stainless steel reactors. The reactors are equipped with high intensity agitation and internal coils which provide the ability to remove high heat loads. The reaction is a semi-batch type whereby initial charges of aqueous solution, vinyl chloride and ethylene are made batchwise and continuous charges of vinyl chloride and other aqueous streams are made on a programmed basis during ensuing reaction. Two basic reactions occur simultaneously: one, the redox catalyst reaction including reduction of ferric to ferrous ion by sulfoxylate and the coupling oxidation reaction regenerating ferric ion and decomposing persulfate ion to a free radical; and two, the subsequent free radical emulsion polymerization forming molecular weight chains of about 35,000 from the vinyl chloride, ethylene and acrylamide monomers. Molecular weight data and simplified equations follow:

Persulfate Redox Catalyst System

Reduction




Oxidation




Polymerization Reaction

Initiation




Chain Growth



The reaction is carried out at 55°C and 1750 - 1850 psi. Vinyl chloride is added to maintain constant pressure, as a pressure fall is associated with reaction of less dense monomers to more dense polymer. Surfactant addition maintains appropriate charge distribution on polymer particle surface to prevent coagulation. Acrylamide is added directly in proportion to degree of reaction. Two catalyst streams generate necessary free radicals. When reaction mass reaches a solids content of about 50%, reaction is terminated and blowdown to atmospheric pressure during which a gas/liquid phase separation occurs. The unreacted monomers are vented to an incineration system and the residual liquid latex is sent forward for further processing.

Primary difference in Latex B and Latex E is polymer composition which is determined during polymerization by feed rates and compositions. B type latices contain about 76% VCM, 21% ethylene and 3% acrylamide. E type latices contain more ethylene (29%) and less vinyl chloride (68%). No by-product, as such is made but some solid coagulum must be filtered out of the liquid latex before storage as final product. This is accomplished by a special type of vibrating filter screen.

Latex BH is produced by a further chemical processing step. B latex is hydrolyzed with caustic already diluted by previously hydrolyzed latex such that the caustic concentration is not too severe for virgin latex. High caustic concentrations can result in undesirable amounts of coagulation. The amide group contributed to the polymer of acrylamide is partially converted to carboxyl and imide as follows:



Hydrolysis is a relatively slow reaction requiring temperature, time and agitation to reach completion in a satisfactory batch cycle.

Following filtration, the latex is stored in bulk storage tanks. Loadout of the product is made from the storage facilities into tankcar, tank trucks, or into 55 gallon fibre drums.

A tetrahydrofuran solvent system is provided to solvent wash reactors at a given frequency (currently after every 15 batches). Between solvent washes heat transfer coefficients decrease due to a solid film deposition on the cooling coils. When the coefficients reach a minimum value, solvent washing is instituted. A solvent recovery system is also provided to remove water and polymer solids from spent solvent. The recovery system consists primarily of a 26 foot packed, distillation tower.

EVC yields of major components are as follows:

<u>Materials</u>	<u>Theory</u> <u>Usage</u>		<u>Design</u> <u>Usage</u>	
	<u>Lbs/Cwt B</u>	<u>Lbs/Cwt E</u>	<u>Lbs/Cwt B</u>	<u>Lbs/Cwt E</u>
VCM	72.4	65.2	82.7	79.5
Ethylene	21.0	28.0	36.1	42.8
Acrylamide	3.15	2.62	3.58	3.07
Sodium Lauryl Sulfate	7.29	7.11	8.35	8.06

TABLE I

ABBREVIATIONS AND MOLECULAR WEIGHTS


Compound	Formula	Abbreviation	Molecular Weight
Acrylamide	$\text{CH}_2 = \text{CHCONH}_2$	A	71.08
Ammonium hydroxide	NH_4OH	NH_4OH	35.05
Ammonium persulfate	$(\text{NH}_4)_2\text{S}_2\text{O}_8$	APS	228.21
Ethylene	$\text{H}_2\text{C}=\text{CH}_2$	E	28.05
Ferric Ammonium Sulfate	$\text{Fe}(\text{NH}_4)(\text{SO}_4)_2 \cdot 12 \text{H}_2\text{O}$	FAS	482.21
Sodium formaldehyde sulfoxylate	$\text{NaHSO}_2 \cdot \text{CH}_2\text{O} \cdot 2 \text{H}_2\text{O}$	SFS	154.12
Sodium hydroxide	NaOH	NaOH	40.01
Tetrasodium Ethylene-diamine tetraacetate	$\begin{array}{c} \text{NaOOCCH}_2 \diagdown \quad \text{CH}_2\text{COONa} \\ \quad \quad \quad \text{N}-\text{CH}_2-\text{CH}_2-\text{N} \\ \text{NaOOCCH}_2 \diagup \quad \quad \quad \diagdown \text{CH}_2\text{COONa} \end{array}$	Na_4EDTA	380.18
Sodium lauryl sulfate	$\text{C}_{12}\text{H}_{25}\text{OSO}_3\text{Na}$	SLS	288.49
Vinyl Chloride	$\text{H}_2\text{C}=\text{CHCL}$	VCl	62.50
Water	HOH	H_2O	18.02
Tetrahydrofuran	$\begin{array}{c} \text{H}_2\text{C} - \text{CH}_2 \\ \quad \quad \\ \text{H}_2\text{C} \quad \quad \text{CH}_2 \\ \quad \quad \quad \diagup \quad \diagdown \\ \quad \quad \quad \text{O} \end{array}$	FHF	72.10
Monomethyl ether hydroquinone	$\text{CH}_3\text{OC}_6\text{H}_4\text{OH}$	MEHQ	124.14


I-4


II. FLOW SHEETS


RSV 0015216

II-1










FLOW DIAGRAM # 8
HOT WATER AND
CHILLED WATER SYSTEMS
3-10-72 J.M.V.



RSV 0015224

III. PROCESS IN DETAIL

RSV 0015225

III. PROCESS IN DETAIL

RSV 0015226

III. PROCESS IN DETAILA. RAW MATERIALS SECTION

The following raw materials are used in the production of EVC Latexes:

1. Vinyl Chloride Monomer (VCM)

- a. Product: B, BH, E
- b. Received: Received from Dow in 20,000 gallon tank cars
- c. Storage: - Stored in two (2) 24,000 gallon, horizontal, carbon steel storage tanks under 90 psig pressure

2. Ethylene (C₂H₄)

- a. Product: B, BH, E
- b. Received: Received from a plant header at 350 psig
- c. Storage: Stored in two (2) - 372 cubic foot bottles under 3,000 psig of pressure

3. Acrylamide

- a. Product: B, BH, E
- b. Received: Received from American Cyanamid in 50 lb bags
- c. Storage: The bags of acrylamide are stored in a solids storage shed until ready for use.

4. Sodium Lauryl Sulfate (SLS)

- a. Product: B, BH, E
- b. Received: Received from Proctor and Gamble in 4,500 gal. trucks
- c. Storage: Stored in a 12,000 gallon, 304 SS storage tank.

5. Ammonium Persulfate (APS)

- a. Product: B, BH, E
- b. Received: Received from FMC Corp. in 250 lb fibre drums.
- c. Storage: The drums of APS are stored in a solids storage shed until ready for use.

A. RAW MATERIALS SECTION (cont'd)6. Ferric Ammonium Sulfate (FAS)

- a. Product: B, BH, E
- b. Received: Received from Mallinckrodt in 25 lb fibre drums
- c. Storage: The drums of FAS are stored in a solids storage shed until ready for use.

7. Sodium formaldehyde Sulfoxalate (SFS)

- a. Product: B, BH, E
- b. Received: Received from Rohm and Haas in 325 lb. fibre drums.
- c. Storage: The drums of SFS are stored in a solids storage shed until ready for use.

8. Tetrasodium Ethylenediamine Tetraacetate (EDTA)

- a. Product: B, BH, E
- b. Received: Received from Dow in 100 lb fibre drums
- c. Storage: The drums of EDTA are stored in a solids storage shed until ready for use.

9. Ammonia

- a. Product: B, BH, E
- b. Received: Received from a plant header at 50 psig
- c. Storage: Stored as a solution with water at 22.2%.

10. Caustic

- a. Product: BH
- b. Received: Received from a $\approx$ 25% plant header.
- c. Storage: Stored at 20% concn. in a 4,000 gallon, horizontal, C.S. storage tank.

11. Tetrahydrofuran (THF)

- a. Product: B, BH, E
- b. Received: Received from DuPont in 500 lb drums and 6,000 gallon tank trucks
- c. Storage: Stored in two (2) 11,200 gallon, horizontal storage tanks.