

Microbics Microtox Toxicity Test. Sample: FC-94-X (Li
salt of PFOS)

MICROBICS MICROTOX® TOXICITY TEST

TEST SUBSTANCE

Identity: Perfluorooctanesulfonate, Lithium salt; may also be referred to as PFOS Li salt, FC-94, or FC-94-X. (1-Octanesulfonic acid, 1,1,2,2,3,3,4,4,5,5,6,6,7,7,8,8,8-heptadecafluoro-, lithium salt, CAS # 29457-72-5)

Remarks: The test sample is a mixture of the test substance in water (approximately 24.5% test substance and 75.5% water). No calculations were made to adjust for the actual concentration of the test substance in the test sample.

METHOD

Method: Microbics Microtox® "BASIC" Procedure

Test type: Static acute

GLP: No

Year Completed: 1994

Species: *Photobacterium phosphoreum*

Analytical monitoring: pH and osmotic condition

Statistical methods: EC10 and EC50 values calculated by statistical linear regression program provided for Microtox®

Test organism source: Microbics, Reagent Lot # AMO25.

Test conditions

Dilution water: Millipore Milli-Q™ water

Stock and test solutions preparation: A primary stock solution was prepared in dilution water at 2000 mg/L. The pH of the stock solution was then adjusted from 5.9 to 7.0 using 0.02N NaOH and an osmotic adjustment was made using 200 mg NaCl dissolved in 10 mL of stock solution. Appearance of test solutions was noted as "clear and colorless". All test solutions were made by proportional dilution.

Exposure vessels: 3 mL glass cuvettes.

Number of replicates: 2

Number of concentrations: Four concentrations and controls

Element basis: Percent light loss

RESULTS

Nominal concentrations: 125, 250, 500, 1000 mg/L

Element value: 15-minute EC₁₀ = 564 (367-866) mg/L

15-minute EC₅₀ = >1000 mg/L

30-minute EC₁₀ = 487 (320-738) mg/L

30-minute EC₅₀ = >1000 mg/L

Remarks: Values reported are for the test sample. No calculations were made to adjust for the concentration of the test substance in the test sample.

CONCLUSIONS

The test sample containing 24.5% Perfluorooctanesulfonate, Lithium salt exhibited a 30-minute EC₅₀ for *Photobacterium phosphoreum* of >1000 mg/L.

Submitter: 3M Company, Environmental Laboratory, P.O. Box 33331, St. Paul, Minnesota, 55133

DATA QUALITY

Reliability: Klimisch ranking 2. This study, while well conducted, lacks analytical data for: determination of the test substance concentration in the test solutions, and determination of the sample purity.

REFERENCES

This study was conducted by the 3M Company, Environmental Laboratory, Lab Request number M1018, 2/07/94.

OTHER

Last changed: 5/2/00

ENVIRONMENTAL LABORATORY

LR M1018
JBD 2/7/94

MICROBICS MICROTOX® TOXICITY TEST

SAMPLE: FC-94-X

PERCENT LIGHT LOSS OF *P. phosphoreum* AFTER:

<u>Test Conc.</u> <u>(mg/L)</u>	<u>5-Min. Exposure</u>	<u>15-Min. Exposure</u>	<u>30-Min. Exposure</u>
125	0	7	8
250	0	5	7
500	0	9	10
1000	2	15	18
EC50, mg/L (95% C.I.)	> 1000 (-)	> 1000 (-)	> 1000 (-)

Sample Description

Preparation of Stock Solution
For Microtox TestingDate: 2/7/94

FC-94-X

Analyst: JBD

1) Preparation of stock solution:

- * Diluted 500 mg (0.5g) of test substance in a 250 mL volumetric flask with Milli-Q™ water.

$$\text{Conc.} = 500 \text{ mg} / 0.25 \text{ L} = 2000 \text{ mg/L}$$

2) pH adjustment of stock solution:

- * Initial pH = 5.9

- * pH adjusted to 7.0 with 0.02N NaOH

3) Osmotic adjustment of stock solution:

- * 200 mg NaCl dissolved in 10 mL of stock solution (= 2% NaCl).

4) Appearance of test solutions: clear & colorless

30 min. EC_{50} : > 1000 mg/L

ויוזר שט

TEST DATE: 2-7-94

Investigator: JBD

Approved by: _____

Procedure: BASIC

Osmotic Adjustment: NaCl

Dilution Factor : 2

Test Time: 30 minutes

Concentration Units: mg/L

NUMBER	IO/IT	CONC.	GAMMA
1	97.94/ 77.53	125.0000	0.07276
2	97.79/ 73.86	125.0000	0.09498
3	89.54/ 73.12	250.0000	0.04017#
4	76.83/ 74.15	250.0000	0.10924#
5	96.54/ 73.54	500.0000	0.11509#
6	90.53/ 70.23	500.0000	0.09495#
7	93.75/ 64.43	1000.0000	0.23577#
8	91.26/ 64.43	1000.0000	0.20221#

micron
Reagent
Lot # AM035

CONTROL IT/IO = 81.28/ 96.69 83.78/ 97.62**

RATIOS = 0.5406 0.8532**

CORRECTION FACTOR = 0.6474

[illegible]
$$.35 \text{SLOPE} = 0.8607$$

G
A
M
M
A

[illegible]

EC303048.3767 (95% CONFIDENCE RANGE: 600.9886 TO 715462.1895)

```
# Used for calculations
```

30 min.

STATISTICAL REPORT FOR DATA FILE M1018.K30

ESTIMATING EQUATION: $\text{LOG } C = 0.8351 \times \text{LOG } b + 3.4841$

NUMBER OF DATA PAIRS USED: 6

RESIDUAL VARIANCE: 0.02548

95% CONFIDENCE FACTOR: 5.07227 FOR EC50

CORRELATION COEFFICIENT: $R = 0.84781$

COEFFICIENT OF DETERMINATION: $R^2 = 0.71879$

EFFECTIVE CONCENTRATIONS (95% CONFIDENCE RANGE):

EC50 = 3048.3767 mg/L (600.9386 TO 15462.1895)***

EC20 = 957.7798 mg/L (474.9458 TO 1931.4672)

EC10 = 486.5695 mg/L (320.5834 TO 738.4969)

EC01 = 65.6808 mg/L (10.7405 TO 401.6549)***

*** EXTRAPOLATED RESULTS

.....

INITIAL DATA

M101G.K30

NUMBER	CONCENTRATION	GAMMA	% EFFECT
1	125.0000	0.073	6.78
2	125.0000	0.095	8.67
3	250.0000	0.040	3.86
4	250.0000	0.109	9.85
5	500.0000	0.115	10.32
6	500.0000	0.095	8.67
7	1000.0000	0.236	19.09
8	1000.0000	0.202	16.62

15 min.

m1018

STATISTICAL REPORT FOR DATA FILE M1018.K15

ESTIMATING EQUATION: $\text{LOG } C = 0.7754 \times \text{LOG } b + 3.4909$

NUMBER OF DATA PAIRS USED: 6

RESIDUAL VARIANCE: 0.02548

95% CONFIDENCE FACTOR: 5.14018 FOR EC50

CORRELATION COEFFICIENT: $R = 0.84780$

COEFFICIENT OF DETERMINATION: $R^2 = 0.71876$

EFFECTIVE CONCENTRATIONS (95% CONFIDENCE RANGE):

EC50 = 3096.7422 mg/L (602.4584 TO 15917.8008)***

EC20 = 1056.9373 mg/L (488.4282 TO 2287.1658)***

EC10 = 563.5806 mg/L (366.8522 TO 865.8067)

EC01 = 67.7866 mg/L (18.3214 TO 420.6264)***

*** EXTRAPOLATED RESULTS

rr11018

M101E.K15

% EFFECT

1	125.0000	0.065	3.11	
2	125.0000	0.082	8.06	7.08
3	250.0000	0.028	2.74	
4	250.0000	0.087	8.02	5.38
5	500.0000	0.099	9.02	
6	500.0000	0.095	8.70	8.66
7	1000.0000	0.189	15.88	
8	1000.0000	0.170	14.53	15.21

MICROTOX DATA REPORT

FILE NAME: M1018.K5

TEST DATE: 2-7-94

Investigator: JRD

TEST TIME: _____

Approved by: _____

Sample Description: Microtox Testing of FC-94-X

Procedure: BASIC

Osmotic Adjustment: NaCl

Initial Concentration : 1000 mg/L

Dilution Factor : 2

Test Time: 5 minutes

Concentration Units: mg/L

NUMBER	IO/IT	CONC.	GAMMA
1	97.94/ 85.74	125.0000	-0.01908*
2	97.79/ 83.68	125.0000	0.00352*
3	89.54/ 80.52	250.0000	-0.04508*
4	96.63/ 82.60	250.0000	0.00667*
5	96.54/ 84.41	500.0000	-0.01767*
6	90.53/ 80.48	500.0000	-0.03404*
7	93.75/ 77.91	1000.0000	0.03332
8	91.26/ 78.10	1000.0000	0.00342*

Microtox
Reagent:
Lot # Amc25

CONTROL IT/IO = 82.95/ 96.69 83.91/ 97.62**

RATIOS = 0.8579 0.8596**

CORRECTION FACTOR = 0.8597

.....NO TOXICITY AT CONCENTRATIONS TESTED

* Invalid data or controls

100

M101B.K5

NUMBER

CONCENTRATION

GAMMA

% EFFECT

1	125.0000	-0.019 *	-1.95	
2	125.0000	0.004 *	0.35	-0.80
3	250.0000	-0.045 *	-4.72	
4	250.0000	0.007 *	0.66	-2.03
5	500.0000	-0.018 *	-1.82	
6	500.0000	-0.034 *	-3.52	-2.67
7	1000.0000	0.033	3.22	
8	1000.0000	0.003 *	0.34	1.78
* 7 INVALID DATA REMOVED FOR CALCULATION				