

INTER-OFFICE MEMO **TENNECO CHEMICALS, INC.**

TO A. Gentilini - AT Burlington Lab. DATE August 28, 1969
FROM H. B. Carr AT Burlington COPY TO J. C. Fisher
SUBJECT E. L. Kanyok
R. J. Likar
E. R. Parker
Personnel File

Your suggestion, Al, on cleaning the Brabender pot is excellent and is very much appreciated. Earl will institute this suggestion as an operating procedure and I'm sure that it will make a necessary but bothersome job more efficient and easy for our Quality Control Laboratory.

H. B. Carr
H. B. Carr

HBC/rjl

Attachment

8-27-69

To: H.B. Carr O.C. Manager
Via: E.L. Kanyok O.C. Supervisor
From: A. Gentilini

Subject: Cleaning of Babender Pot.

- 1- After initially discovering that the use of Thermolite #31 was a good agent for periodic cleaning of the Babender Pot, it wasn't quite the absolute solution.
- 2- I ventured further, and discovered that by adding approximately 1 gram of wood chips to the regular cleaning compound it did a fairly thorough job of cleaning and returned the reading to zero.
3. The wood is very soft and leaves no scars and scratches and leaves the metal shiny and clean.

4- We could buy a pound or two of soft cedar and have the pot operator use in their routine cleaning of the pot.

I hope this is helpful
Respect

A. Gentilini