

NIOSH
MANUAL OF ANALYTICAL METHODS

Physical and Chemical Analysis Branch
Division of Laboratories and Criteria Development
1014 Broadway, Cincinnati, Ohio 45202

U.S. DEPARTMENT OF HEALTH, EDUCATION, AND WELFARE
Public Health Service
Center for Disease Control
National Institute for Occupational Safety and Health
1974

BFG25863

23522001

Mr. John V. Crable and Dr. David G. Taylor, Chief and Assistant Chief respectively, of the Physical and Chemical Analysis Branch were responsible for the development and preparation of the manual.

BFG25864

23522002

TABLE OF CONTENTS

<u>P & CAM No.</u>	<u>Pages</u>	<u>Title</u>
101	1-6	Lead in Blood and Urine
102	1-9	Lead in Blood and Urine
103	1-5	Mercury in Urine
105	1-5	Mercury in Blood
106	1-7	Colorimetric Method for Free Silica
107	1-6	Antimony in Urine
108	1-12	Nitrogen Dioxide in Air
109	1-7	Free Silica (Quartz, Cristobalite, Tridymite) in Atmospheric Dust
110	1-6	Quartz in Coal Dust by Infrared Spectroscopy
112	1-5	Carbon Monoxide in Air
114	1-5	Fluoride in Urine
116	1-6	Cyanide in Air
117	1-6	Fluoride and Hydrogen Fluoride in Air
118	1-6	Acrolein in Air
121	1-6	Beryllium in Air
123	1-8	Beryllium in Air
125	1-9	Formaldehyde in Air
126	1-13	Hydrogen Sulfide in Air
127	1-11	Organic Solvents in Air
139	1-8	Arsenic in Urine and Air
140	1-7	Arsenic in Urine and Air
141	1-8	2,4-Toluenediisocyanate (TDI) in Air
142	1-7	p,p-Diphenylmethanediisocyanate (MDI) in Air
146	1-7	Sulfur Dioxide in Air
151	1-5	Cadmium in Air
152	1-5	Total Particulate Chromium in Air
153	1-6	Ozone in Air
155	1-5	Lead in Air
158	1-9	Parathion in Air
159	1-5	Oil Mist in Air
160	1-10	Sulfur Dioxide in Air
163	1-7	Sulfur Dioxide in Air
165	1-7	Mercury in Urine

TABLE OF CONTENTS (Continued)

<u>P & CAM No.</u>	<u>Pages</u>	<u>Title</u>
167	1-8	Mercury in Blood
168	1-11	Aromatic Amines in Air
169	1-6	Chromic Acid Mist in Air
173	1-8	General Procedure for Metals
175	1-15	Mercury in Air
Appendix A		Description, Installation, and Operation of Two Stage Thermal Desorption Unit for the Mercury Determination

INDEX OF METHODS ARRANGED ACCORDING TO ANALYTE

<u>Analyte</u>	<u>Matrix</u>	<u>Procedure</u>		<u>P & CAM No.</u>
		<u>Sampling</u>	<u>Analysis</u>	
Acrolein	Air	Impinger	Colorimetric	118
Antimony	Urine		Colorimetric	107
Aromatic Amines	Air	Silica gel Adsorption	GC	168
Arsenic	Air, Urine	Membrane Filter	AA	139
Arsenic	Air, Urine	Membrane Filter	Colorimetric	140
Beryllium	Air, Settled Dust	Membrane Filter	AA	121
	Ore, Swipe Samples			
Beryllium	Air, Settled Dust	Membrane Filter	Emission Spectroscopy	123
	Ore, Swipe Samples			
Cadmium	Air	Membrane Filter	AA	151
Carbon Monoxide	Air	Grab-bag	IR	112
Chromic Acid Mist	Air	Membrane Filter	Colorimetric	169
Chromium (Particulate)	Air	Membrane Filter	AA	152
Cyanide	Air	Impinger	Ion-Specific Electrode	116
Diphenylmethane- diisocyanate	Air	Impinger	Colorimetric	142
Fluoride	Urine		Ion-Specific Electrode	114
Fluorides and Hydrogen Fluoride	Air	Impinger	Ion-Specific Electrode	117
Formaldehyde	Air	Impinger	Colorimetric	125

<u>Analyte</u>	<u>Matrix</u>	<u>Procedure</u>		<u>P & CAM No.</u>
		<u>Sampling</u>	<u>Analysis</u>	
Hydrogen Sulfide	Air	Impinger	Colorimetric	126
Lead	Air	Membrane Filter	AA	155
Lead	Blood, Urine	Vacutainers	AA	101
Lead	Blood, Urine	Vacutainers	Colorimetric	102
Mercury (Particulate, Metallic, Organic Vapor)	Air	Three Section Solid Phase Sampler	Flameless AA	175
Mercury	Blood	Vacutainers	Flameless AA	105
Mercury	Blood	Vacutainers	Flameless AA	167
Mercury	Urine		Flameless AA	103
Mercury	Urine		Flameless AA	165
Metals (General Procedure)	Industrial Hygiene Samples	Membrane Filter	AA	173
Nitrogen Dioxide	Air	Impinger	Colorimetric	108
Oil Mist	Air	Membrane Filter	Fluorescence	159
Organic Solvents	Air	Charcoal Adsorption	GC	127
Ozone	Air	Impinger	Colorimetric	153
Parathion	Air	Impinger	GC	158
Quartz	Coal Dust	Membrane Filter	Infrared	110
Silica (Quartz)	Air, Dust	Membrane Filter	Colorimetric	106
	Biological Tissue			
Silica (Quartz, Cristobalite, Tridymite)	Air	Membrane Filter	X-Ray Diffraction	109
Sulfur Dioxide	Air	Impinger	Titration	146
Sulfur Dioxide	Air	Impinger	Colorimetric	160

<u>Analyte</u>	<u>Matrix</u>	<u>Procedure</u>		<u>P & CAM No.</u>
		<u>Sampling</u>	<u>Analysis</u>	
Sulfur Dioxide	Air	Impinger	Titration	163
2,4-Toluene-diisocyanate	Air	Impinger	Colorimetric	141