



"Roger J. Zipfel" <ZIPFEL@WWPS-A1.EMAIL.DUPONT.COM> on 09/11/99
02:56:00 AM

To: Roger J Zipfel/AE/DuPont
cc:
Subject: [C-8 ANALYTICAL/18-Feb-1994]First draft C-8 sampling + analys

Date: Fri, 18 Feb 1994 13:46:00 EDT
From: Rik Brandenburg <BRANDERA@EUOIS2-A1.EMAIL.DUPONT.COM>
Subject: First draft C-8 sampling + analysis procedure
To: Gerard Rietbroek <RIETBRG@EUOIS2-A1.EMAIL.DUPONT.COM>, Willem Huisman
<HUISMAW@EUOIS2-A1.EMAIL.DUPONT.COM>, "Richard W. Quarles,, Jr."
<QUARLES@wwps-a1.email.dupont.com>, YASUSHI NISHIKAWA
<NISHIKAY@wwps-a1.email.dupont.com>, AKITO ABE
<ABEA@YOVX00-A1.EMAIL.DUPONT.COM>, IKUKO SAWATAKI
<SAWATAI@YOVX00-A1.EMAIL.DUPONT.COM>, Chris J Primus
<PRIMUSCJ@wwps-a1.email.dupont.com>, ZIPFEL@wwps-a1.email.dupont.com
Cc: YOKOYAK@YOVX00-A1.EMAIL.DUPONT.COM
MIME-version: 1.0
Content-type: MULTIPART/MIXED; BOUNDARY="Boundary_(ID_5IRvr06z4ItTd4eMFHp7gQ)"
Delivery-date: Fri, 18 Feb 1994 08:33:00 EDT
Posting-date: Fri, 18 Feb 1994 14:27:00 EDT
Importance: normal
Al-type: MAIL

To: see distribution list
From: R.A. Brandenburg

As promised I would come back in February with a proposal for uniform sampling and analysis procedures for C-8. The reason for this is that so far, most C-8 results have been generated in a non-routine way. The choice of sample-taking equipment, sample preparation, and actual analytical technique was often based on individual preference. Although the analytical techniques can easily be verified by using standard solutions, this is not the case with the sample collection, concentration and preparation. Therefore special emphasis should be given to the way of sampling and preparation.

We want to be able to compare results from different sites, and -most importantly- to compare last years' results with next years' result. This is only possible if we use the same methods all the time.

Attached is a matrix with 4 columns. The first column lists the various samples we may have to deal with. Column "C" indicates the proposed sampling Collection /Concentration method. Column "P" is for the sample Preparation, and column "A" is for the Analysis.

RJZ013779

EID081896

This is just a first shot so if you have other ideas or additions, please don't hesitate to let me know. Also I should like to know if Parkersburg and MDF have value for this global approach, and are willing to adjust their procedures where necessary to achieve some degree of global uniformity. Please answer before April 1.

best regards,

Rik Brandenburg

Date: Tue, 25 Jan 1994 15:08:00 EDT
 MIME-version: 1.0
 Content-type: text/plain
 Posting-date: Fri, 18 Feb 1994 00:00:00 EDT
 Importance: normal
 Al-type: DOCUMENT

DRAFT PROCEDURE FOR C-8 SAMPLING AND ANALYSIS

I. POSSIBLE SAMPLES

	COLLECTION/ CONCENTR.	PREPA- RATION	ANALYSIS METHOD
explanation of numbers, see page 2, code letter :	"C"	"P"	"A"
C-8 POWDER IMPURITIES	C3.1	P5	A3
C-8 POWDER ASSAY	C3.1	P9	A6
C-8 SOLUTIONS FROM VENDOR	C2.1	P0	A6,7
C-8 SOLUTIONS FROM RECOVERY SYSTEMS	C2.1	P0	A6,7
C8 DUST IN AIR (PERSONNEL EXPOSURE SAMPLES)	C1.3	P0	A5
C8 VAPOR IN AIR (PERSONNEL EXPOSURE SAMPLES)	C1.1	P4	A3
C8 VAPOR + DUST IN STACKS, EXHAUSTS	C1.2	P5	A3
C8 IN RIVER WATER (PPB)	C2.1	P0	A4
C8 IN GROUND WATER (PPB)	C2.1	P0	A4
C8 IN SCRUBBER EFFLUENTS	C2.1	P5	A3
C8 IN SUPERNATE	C2.1	P0	A1
C8 IN FEP WASTE WATER	C2.1	P2	A1
C8 IN PTFE WASTE WATER	C2.1	P2	A1
C8 IN VITON WASTE WATER	C2.1	P2	A1
C8 IN PFA WASTE WATER	C2.1	P2	A1
C8 IN WASTE WATER FROM FLUOROPOLYMERS	C2.1	P2	A1
C8 IN WATER FROM BIOTREATMENT	C2.1	P2	A1
C8 IN TEFLON DISPERSION	C2.2	P3	A1
C8 IN WET TEFLON POWDER	C3.1	P6	A1
C8 IN TEFLON POWDER FINISHED PRODUCT	C3.1	P3	A1
C8 IN WET VITON BEFORE BUSS KNEADER	C3.2	P ?	A ?
C8 IN VITON SHEET	C3.3	P ?	A ?
C8 IN FEP BEFORE TD DRYER	C3.1	P3	A1

RJZ013780

C8 IN FEP FINISHED PRODUCT	C3.1	P 2	A 2
C8 IN WAX	C3.4	P10	A1
C8 IN LANDFILL PTFE	C3.1	P3	A1
C8 IN DRYER PAPER	C3.5	P3	A1
C8 IN BIOTREATMENT SLUDGE	C3.6	P5+11	A3

II. EXPLANATION OF NUMBERS IN COLUMNS "C" ,"P" AND "A"

COLUMN "C": POSSIBLE SAMPLING COLLECTION AND CONCENTRATION METHODS

C1.: GAS:

- C1.1. On tenax
- C1.2. In a water scrubber or impinger
- C1.3. on a dust filter specially for air samples
- C1.4. in a plastic bag

C2. LIQUID

- C2.1. Just a water sample
- C2.2. dispersion or emulsion
- C2.3. absorption on active carbon
- C2.4. absorption on an ion-exchanger

C3. SOLID

- C3.1. wet or dry powders, pellets or cubes
- C3.2. wet Viton before buss kneader
- C3.3. Viton sheet
- C3.4. wax
- C3.5. paper
- C3.6. biotreatment sludge

COLUMN "P": POSSIBLE SAMPLE PREPARATION METHODS

- P0. none
- P1. dissolve in water
- P2. shake with water + CHCl₃
- P3. reflux with water + CHCl₃
- P4. desorption from tenax
- P5. esterification and extraction
- P6. desorption from active carbon
- P7. desorption from ion exchanger
- P8. vacuum-drying of wet solids (polymer, paper, sludge, wax)
- P9. dissolve in anhydrous acetic acid
- P10. dissolve in CHCl₃
- P11. filter solids after extraction

COLUMN "A": POSSIBLE METHODS OF ANALYSIS

- A1. methylene blue and spectrometry (10 ppm - 1%)
- A2. CTMAB titration (very inaccurate)
- A3. GC (ppm level).

RIZ013781

- A4. CH2M HILL method (ppb level)
- A5. weigh dust particles on filter

In addition: For pure C-8 or commercial solutions:

- A6. perchloric acid titration
- A7. density tests
- A8. refractive index
- A9. acid/base titration

RJZ013782

EID081899