

APR 11 1977
EST. WGT. 2402 D 1061 LBS.

D 1061 Jan

() BT-417-D (L) AJC 2-15-77
() BT-434-D (C) DCG 3/11/77
() () MTW 3/15/77
() ()
TYPED BY
NE 3/19/77

STD. HRS. TIME M

FIN. AREA: BT-434-D ONLY
SPLIT OFF ONE-DRUM ($\pm 150\#$)
To BT-417-D, Lot

Do NOT USE N-8710A Lot #1 or #2

NON-BLEEDING

PMENT	NAME:	USE LBS.	N CODE	ACCT. LBS.
SV 1	MONASTRAL® BLUE GF			
	Run 2-3" (± 200 -G) water into V-202			
	Use Station # 3, preset water meter per chart for water addition to full totes, then using mode 4, pump in 5 full totes ($\pm 2600\# 100\%$)		8710-A	
	GALS. OF WATER			
	SEQUENCER STEPS			
	Lot 6 11 16 #100%			
	298B 100 175 20 497 26.8 1855		8710-A	5435
	291A 100 175 20 509 28.9 1750			
	301B 100 175 20 478 27.0 1762			
	298C 100 175 20 446 26.2 1701			
	303A 100 175 20 504 28.9 1750			
	Total 24.32 # 100%			
	13 Weigh & add Assume 50%	14	7-0	7
	14 Stir 5 min			
	15 Record: pH 7.5 (7.0-8.0)			
	16 If necessary add to pH 7.0-8.0 Assume 50%		7-0	
	17 Record: pH 7.5 M		30-0	
	18 Heat to 120-125°F in 30-35 min			
	Record: Temp 122°F			
	19 Stop agitator. Record: Vol 56			
SC 2	20 Calc Est Wgt:			
	24.32 #100% x 0.9876 = 24.02 # Est. Wgt.			
	21 Record: Est. Wgt. on Make & Process cards			
	START 11:00 AM			
	FIN 6:15 PM			
	OPTR 114			