

MAY -0-81- 40001
EST. WGT. 2289 LBS.

STD. HRS. TIME M

FIN. AREA: BT-484-D ONLY
SPLIT OFF ONE-DRUM (+150#)
To BT-417-D, Lot

PROCESSES	APPROVED BY
() BT-417-D	(L) AJC 2-15-77
() BT-434-D	(S) DCG 3/11/77
()	() MTW 3/15/77
()	()
TYPED BY	
WE 3/19/77	

CLEAN UP

AS YOU GO ALONG

NON-BLEEDING

PMENT	NAME:	USE LBS.	N CODE	ACCT. LBS.
SV 1	MONASTRAL® BLUE GF			
	1 Run 2-3" (± 200-G) water into V-204			
	2 Use Station # 2, preset water meter per chart for water addition to full totes, then using mode 4, pump in 5 full totes (± 2600# 100%)		8710-A	
	GALS. OF WATER			
	SEQUENCER STEPS			
	Lot 6 11 16 #100%			
	29713 150 10 80 499 28.5 1730		8710-A	2289
	297 150 10 90 478 27.3 1732			
	29513 150 10 90 389 27.5 1416			
	295A 160 10 90 454 25.2 1820			
	298A 150 20 100 469 26.8 1750			
	Total 2289 # 100%			
	3 Weigh & add Assume 50% 14 7-B 7			
	4 Stir 5 min			
	5 Record: pH 9.0 (7.0-8.0)			
	6 If necessary adj to pH 7.0-8.0 Assume 50% 7-B 1			
	6 Record: pH 7.8 M-3 6389 2 30-II			
	7 Heat to 120-125°F in 30-35 min			
	Record: Temp 125 of			
	8 Stop agitator. Record: Vol 68			
ISC 2	9 Calc Est Wgt: 2289 #100% x 0.9876 = 2261 # Est. Wgt.			
	10 Record: Est. Wgt. on Make & Process cards			
	START 3PM	FIN	OPTR 787	

917