

PROCESSES APPROVED BY
() XBT-579-D ()

SEP 24 1982

K21105

[illegible]

EST. WT. ~~100~~ * LBS.

TYPED BY: JES

STANDARD HRS.

TIME 2-4M

PREPARATION NOTES: USE N-8554 LOT B-004 876#s, ONE TOTE

- CLEAN AS BT-284-D

EQUIP MENT	NAME: MONASTRAL BLUE EBF	1 PAGE	USE LBS.	N-CODE	ACCT LBS.
SV 1	RUN 2-3* (+/-200G) WATER INTO V- <u>202</u> USE STATION # <u>2</u> , PRESET WATER METER PER CHART FOR WATER ADDITION TO FULL TOTES, THEN USING MODE 3, PUMP IN 2 FULL TOTES (+/-1000 ± 100%) GALS OF WATER SEQUENCER STEPS				
	LOT	6	11	16	±100%
	<u>8-004</u>	<u>180</u>	<u>50</u>	<u>100</u>	<u>876</u>

MISC. 1-8 CALC EST WGT:
2-1 $\div 100\% \times 0.980 = 858 \div$ EST. WT.
1-9 RECORD EST. WT. ON MAKE & PROCESS CARDS

START	FIN	OPTR
4:00 AM	6:45 AM	101