

E1589

EST. WGT. 1094 LBS.

STD. HRS. _____ TIME _____ M

FIN. AREA: BT-434-D ONLY
SPLIT OFF ONE-DRUM ($\pm 150\#$)
To BT-417-D, Lot _____

() BT-417-D	(1) AJC 2-15-77
() BT-434-D	(3) DCG 3/11/77
()	() MTW 3/15/77
()	()
TYPED BY	
NE 3/19/77	

NON-BLEEDING

PMENT	NAME: MONASTRAL® BLUE GF.	USE LBS.	N CODE	ACCT. LBS.
SV	1 Run 2-3" (± 200 -G) water into V-			
1	2 Use Station # _____, preset water meter per chart for water addition to full totes, then using mode 4, pump in 5 full totes ($\pm 2600\#$ 100%)		8710-A	
	GALS. OF WATER			
	SEQUENCER STEPS			
	Lot 6 11 16 #100% %		8710-A	
	Total # 100%			
	3 Weigh & add Assume 50%	14	7-0	7
	4 Stir 5 min			
	5 Record: pH (7.0-8.0)			
	6 If necessary adj to pH 7.0-8.0 Assume 50%		7-0	
	& Record: pH M %		30-0	
	7 Heat to 120-125°F in 30-35 min			
	Record: Temp of			
	8 Stop agitator. Record: Vol			
sc	9 Calc Est Wgt:			
2	#100% x 0.9876 = # Est. Wgt.			
	10 Record: Est. Wgt. on Make & Process cards			
	START FIN OPTR			

