

D-19661 .5

APPENDI O. 3
FPG - WINCHESTER

ANALYSES OF AIR SAMPLES
FOR Methyl Ethyl Ketone

METHOD OF COLLECTION: ODZ-PS-Perbonal
sampling pump with activated charcoal absorption tubes
METHOD OF ANALYSIS: Gas Chromatography

DATE ANALYZED: 5/6/76
ANALYST: J. Desouky
TLV: 590 mg/M³

Sample No.		Date	Time	Total (M ³) Air Vol.	Location	Conc. mg/M ³	OPERATOR
Field	Lab.						
1	630-G	4/23/76	9:49am 10:20am	.024	Blok finish printer No. 2	3.13	Richard Renner
2	631-G	"	9:43am 10:23am	.03	Strip finish hand printer No. 2	2.17	Calvin Fox
3	632-G	"	10:01am 10:39am	.03	Strip finish printer machine No. 10	59.23	Stephen Michael

(*) Concentration in excess of TLV