

THOMPSON VITA-METER CORPORATION
6402 Cedar Avenue
Cleveland 3, Ohio

MIXING PROCEDURE FOR

"VITOL"

1. Mix $1/8$ part by weight of emulsifying base (Penola 2210 or equivalent) with $3/8$ parts by weight of kerosene. Stir well. Penola specific gravity is 1.07 and kerosene is about 0.780.
2. Mix $\frac{1}{2}$ part by volume of inhibitor from #1 above into 2 parts by volume of water. Stir or agitate until uniform white emulsion is obtained.
3. Separately mix 85 parts by volume of methanol with enough tetraethyl lead fluid to give 3.53 cc of tetraethyl lead per gallon of pure methanol.
4. To the 85 parts of leaded methanol add 13 parts of water slowly with agitation or stirring.
5. Slowly add inhibitor from #2 with agitation.

NOTE: Important - Total hardness of water must be less than 100 parts per million.

S/ C. H. Van Hartesveldt

CHVH:sg

RE 0020730

N22004