

AR 226 - 3432

36

000368

TABLE IIGC PRECISION

<u>SAMPLE</u>	<u>GC COUNTS/THOUSANDS</u>
5 ppm FC-143	232.5, 241.0, 234.1
1.0 ppm FC-143	63.2, 62.3, 65.2, 61.6, 59.0
0.5 ppm	33.4, 34.3, 38.8, 39.5
0.25 ppm	17.5, 15.5, 19.1, 18.4

BLOOD SAMPLESILC with PPD/S. Stafford

<u>SAMPLE</u>	<u>SPIKED LEVEL (ppm F)</u>	<u>MEASURED LEVEL (ppm F)</u>
1	.50	.56
2	1.5	1.9
3	Blank	.01
4	.027	0.026
5	0.17	0.145
6	1.0	1.2
7	Blank	.009
8	.054	.047
9	.34	.314
10	.054	.059