

AR226-1597

ANALYTICAL RESULTS: FC-143* IN BLOOD
POLYMER PRODUCTS - CHAMBERS WORKS

July-August 1981					
<u>Production</u>	<u>Building</u>	<u>First</u>	<u>Technicians</u>		<u>Method</u>
<u>Building</u>	<u>Mechanic</u>	<u>Line</u>	<u>Solution</u>	<u>Emulsion</u>	<u>Development</u>
<u>Operators</u>		<u>Supervision</u>	<u>Prep**</u>	<u>Prep</u>	<u>Technician</u>
FC-143 as ppm Fluorine in Blood					
0.014	0.059	0.013	(0.013)	0.020	N.D.
0.062		0.016	(0.015)	N.D.	
0.077		N.D.		<0.007	
0.140		0.032		<0.007	
0.150		0.015	(0.007)	0.010	
0.025		0.033	(0.007)		
0.084			(0.011)		
0.110			(0.016)		
0.039					
0.027					
0.250			(0.096)		
0.150			(0.120)		
0.014	-	0	0.007	0	- Low
0.250	0.059	0.033	0.120	0.020	0 High
0.094	0.059	0.018	0.036	0.009	0 Mean

Overall Range: 0-0.250, Overall Mean: 0.049

*Ammonium Perfluorooctanoate

**Duplicate numbers reported for four individuals.

First number is pre-campaign, second number post-campaign.

JCL/dmb
8/27/81

EID110665

BLD000186

000128