Ethanolamines/Glycol Ethers
Tier 3 Procedure

Section Number: 614
Issue Number: 1.3
Original Date: 12-17-91 **Author:** Paula Bourgeois
Revision Date: 12-29-95
Title: How To Use A Volumetric Flask
References : Glossary Of Equipment Identification (Text # 114)

==========

614. How To Use A Volumetric Flask

1.0 Objective:

1.1 To ensure correct use and accurate measurement.

2.0 Responsibility:

2.1 The operations supervisor is responsible for instructing or designating an instructor for training on all lab equipment.

2.2 All persons whether permanent or rotating are responsible for correctly using all lab equipment.

3.0 Procedure:

3.1 To fill a volumetric flask to the correct amount, pour liquid to be measured into the flask until it reaches the neck. Then with small amounts of the liquid, fill the neck until the bottom of the meniscus is on the indicator mark (usually a ring around the neck).

3.1.1 Meniscus - the bottom of the bubble at the surface of the liquid.

4.0 Revisions:

Revision No. 1	07-20-93	Yearly Review
Revision No. 2	12-21-94	Yearly Review
Revision No. 3	12-29-95	Yearly Review

Paper Copies Of This Document Are Uncontrolled & Not Subject To Update.
The Controlled Document Resides In File: Volumet.doc

DO A 059869
CONFIDENTIAL